

Total Control for Full Contour Confidence

CEREC®



MILLING

Basic Starter Kit A-D (Type B) – Ivoclar Vivadent

The IPS Empress CAD for CEREC Basic Starter Kit A-D comprises a selection of block shapes in the most popular A-D shades, as well as the necessary working accessories especially for chairside for intraoral application, the new finishers and polishers of the OptraFine Assortment, and a selection of OptraGate and OptraDam.

In order to enable shade determination with light tooth shades, the new Bleach module with the 4 Bleach shade BL1 to BL4 will also be integrated. The Basic Kit is supplied in the new materials cabinet.

Shades	Size	Qty
A1	I12	5 pcs.
A2	I12	5 pcs.
A3	I12	5 pcs.
B1	I12	5 pcs.
A1	C14	5 pcs.
A2	C14	5 pcs.
A3	C14	5 pcs.
B1	C14	5 pcs.
BL2	C14	5 pcs.
A1	C14 L	5 pcs.
A2	C14 L	5 pcs.
A3	C14 L	5 pcs.

Contains 1 - 50ml IPS Contrast Spray chairside, 1-18 pcs. OptraFine Assortment, 4 Small OptraGate, 4 Regular OptraGate, 4 Small OptraDam, 4 Regular OptraDam, 1 Bleach Shade Guide
 Each.....715-2903 1@ \$1,592.30



inLab Basic Starter Kit (Type B) – Ivoclar Vivadent

The IPS Empress CAD inLab Basic Starter Kit A-D comprises a selection of block shapes in the most popular A-D shades, as well as the necessary working accessories especially for labside applications: The Contrast Spray labside in the large delivery form, especially for extraoral application, as well as the Object Fix firing paste (for the cut-back technique). In order to enable shade determination with light tooth shades, the new Bleach shade guide with the 4 Bleach shade BL1 to BL4 will also be integrated.

The Basic Kit is supplied in the new materials cabinet and can be combined as desired with other Ivoclar Vivadent all-ceramic assortments, such as IPS e.max. The IPS Empress CAD for InLab Basic Kit contains no layering materials. Therefore, it is supplemented in the laboratory with the IPS Empress Veneer Kit (layering materials) for the cut-back technique, which is separately available.

Shades	Size	Qty
A1	I12	5 pcs.
A2	I12	5 pcs.
A3	I12	5 pcs.
B1	I12	5 pcs.
A1	C14	5 pcs.
A2	C14	5 pcs.
A3	C14	5 pcs.
B1	C14	5 pcs.
BL2	C14	5 pcs.
A1	C14 L	5 pcs.
A2	C14 L	5 pcs.
A3	C14 L	5 pcs.

Contains 1 - 200ml IPS Contrast Spray labside, 1 - 12gm IPS Object Fix, 1 Bleach Shade Guide
 Each715-2929 1@ \$1,592.30



IPS Empress CAD (Type B) – Ivoclar Vivadent

Leucite-Reinforced Glass Ceramic Block for CAD/CAM Dentistry

IPS Empress® CAD is the aesthetic leucite-reinforced glass ceramic for the efficient fabrication of single-unit restorations through CAD/CAM technology. By optimizing the manufacturing processes, a new material was created with more homogeneous leucite crystal distribution providing a glass ceramic with an increased flexural strength of 160 MPa.



High Translucency

Features a high degree of translucency, making the HT blocks ideally suited for the fabrication of smaller restorations (e.g., inlays and onlays). Restorations fabricated with the HT blocks demonstrate true-to-nature chameleon effect and outstanding adaptation to the residual tooth structure. The blocks simulate the life-like fluorescence in order to permit optimum integration into the natural residual tooth structure.

Chromascop

Shades	Size	Qty	Item #	
100	V12	5	715-2960	1@ \$108.95
100	I8	5	715-2978	3@ \$104.25
200	I8	5	715-2986	
300	I8	5	715-2994	
100	I10	5	715-3000	
200	I10	5	715-3018	
300	I10	5	715-3026	
100	I12	5	715-3034	
200	I12	5	715-3042	
300	I12	5	715-3059	

A-D

A1	V12	5	715-3067	1@ \$108.95
A2	V12	5	715-3075	3@ \$104.25
A3	V12	5	715-3083	
A1	I8	5	715-3091	
A2	I8	5	715-3109	
A3	I8	5	715-3117	
A3.5	I8	5	715-3125	
B1	I8	5	715-3133	
B2	I8	5	715-3141	
B3	I8	5	715-3158	
C2	I8	5	715-3166	
D3	I8	5	715-3174	
A1	I10	5	715-3182	
A2	I10	5	715-3190	
A3	I10	5	715-3208	
A3.5	I10	5	715-3216	
B1	I10	5	715-3224	
B2	I10	5	715-3232	
B3	I10	5	715-3240	
C2	I10	5	715-3257	
D3	I10	5	715-3265	
A1	I12	5	715-3273	
A2	I12	5	715-3281	
A3	I12	5	715-3299	
A3.5	I12	5	715-3307	
B1	I12	5	715-3315	
B2	I12	5	715-3323	
B3	I12	5	715-3331	
C2	I12	5	715-3349	
D3	I12	5	715-3356	



Low Translucency

Features a lower degree of translucency, making them ideal for most full coverage restorations (e.g., full crowns, large partial coverage restorations and veneers). The LT blocks are less translucent, which makes them appear higher in value, therefore resulting in a brighter restoration.

Chromascop

Shades	Size	Qty	Item #	
100	V12	5	715-3612	1@ \$108.95
100	I10	5	715-3620	3@ \$104.25
200	I10	5	715-3638	
300	I10	5	715-3646	
100	I12	5	715-3653	
200	I12	5	715-3661	
300	I12	5	715-3679	
100	C14	5	715-3687	
200	C14	5	715-3695	
300	C14	5	715-3703	

A-D

A1	V12	5	715-3711	1@ \$108.95
A2	V12	5	715-3729	3@ \$104.25
A3	V12	5	715-3737	
A1	I10	5	715-3745	
A2	I10	5	715-3752	
A3	I10	5	715-3760	
A3.5	I10	5	715-3778	

A-D (continued)

Shades	Size	Qty	Item #	
B1	I10	5	715-3786	1@ \$108.95
B2	I10	5	715-3794	3@ \$104.25
B3	I10	5	715-3802	
C2	I10	5	715-3810	
D3	I10	5	715-3828	
BL1	I10	5	715-4016	
BL2	I10	5	715-4024	
BL3	I10	5	715-4032	
BL4	I10	5	715-4040	
A1	I12	5	715-3836	
A2	I12	5	715-3844	
A3	I12	5	715-3851	
A3.5	I12	5	715-3869	
B1	I12	5	715-3877	
B2	I12	5	715-3885	
B3	I12	5	715-3893	
C2	I12	5	715-3901	
D3	I12	5	715-3919	
BL1	I12	5	715-4057	
BL2	I12	5	715-4065	
BL3	I12	5	715-4073	
BL4	I12	5	715-4081	
A1	C14	5	715-3927	
A2	C14	5	715-3935	
A3	C14	5	715-3943	
A3.5	C14	5	715-3950	
B1	C14	5	715-3968	
B2	C14	5	715-3976	
B3	C14	5	715-3984	
C2	C14	5	715-3992	
D3	C14	5	715-4008	
BL1	C14	5	715-4099	
BL2	C14	5	715-4107	
BL3	C14	5	715-4115	
BL4	C14	5	715-4123	



Multi

With multiple translucencies and a balanced blend of chroma from gingival to incisal, IPS Empress CAD Multi offers dental professionals a fusion of colour that provides the natural gradients to replicate dentin and enamel in a single block.

Shades	Size	Qty	Item #	
Multi A1	I12	5	715-4511	1@ \$139.95
Multi A2	I12	5	715-4529	3@ \$133.65
Multi A3	I12	5	715-4537	
Multi A3.5	I12	5	715-4545	
Multi B1	I12	5	715-4552	
Multi BL1	I12	5	714-3902	
Multi BL3	I12	5	714-3910	
Multi A1	C14	5	715-4560	
Multi A2	C14	5	715-4578	
Multi A3	C14	5	715-4586	
Multi A3.5	C14	5	715-4594	
Multi B1	C14	5	715-4602	
Multi BL1	C14	5	714-3928	
Multi BL3	C14	5	714-3936	
Multi A1	C14L	5	715-4818	1@ \$144.95
Multi A2	C14L	5	715-4826	3@ \$139.00
Multi A3	C14L	5	715-4834	
Multi A3.5	C14L	5	715-4842	
Multi B1	C14L	5	715-4859	
Multi BL1	C14L	5	714-3944	1@ \$144.95
Multi BL3	C14L	5	714-3951	1@ \$139.00

MILLING (CONT.)



IPS e.max CAD for CEREC Basic Kit LT – Ivoclar Vivadent

The IPS e.max CAD for CEREC Basic Kit LT comprises blocks as well as the necessary working accessories for chairside applications using the CEREC System (Sirona). The Basic Kit is supplied in a materials cabinet and can be expanded as desired with other IPS e.max Kits.

Contents

- 6x 5 pcs. IPS e.max CAD for CEREC and inLab Blocks C14, Shades: LT BL2, LT A1, LT A2, LT A3, LT A3.5, LT B1
- 7x 3g IPS e.max CAD Crystall./Shade, Shades: 0, 1, 2, 3, 4, I1, I2
- 7x 1g IPS e.max CAD Crystall./Stains, Shades: white, cream, sunset, copper, olive, khaki, mahogany
- 1x 3g IPS e.max CAD Crystall./Glaze Paste
- 1x 270ml IPS e.max CAD Crystall./Glaze Spray
- 1x 15ml IPS e.max CAD Crystall./Glaze Liquid
- 1x 5g IPS e.max CAD Crystall./Add-On
- 1x 5g IPS e.max CAD Crystall./Add-On Liquid
- 1x 50ml IPS Contrast Spray Chairside
- 1x IPS e.max CAD Crystallization tray
- 1x 10ml IPS object fix flow
- 1x 12ml IPS object fix putty
- 1x IPS Ceramic Ätzgel kit
- 1x IPS e.max Press/CAD LT Shade Guide
- Div. Zubehör

IPS e.max CAD for CEREC Basic Kit LT A-D

Each 717-4899 1@ \$1,399.00
6@ \$1,328.75



IPS e.max CAD – Ivoclar Vivadent

As soon as you see it you will know it is unique. Blue in colour, IPS e.max CAD is delivered in a precrystallized form that is ideal for milling. After milling, IPS e.max CAD is crystallized in a standard porcelain oven where it undergoes a dramatic change and the final physical properties are established. IPS e.max CAD is highly aesthetic and offers 360 MPa of strength.

IPS e.max CAD LT A-D I12/5

Shade	Size	Qty	Item #	
A1	I12	5	717-4519	1@ \$159.95
A2	I12	5	717-4527	3@ \$154.95
A3	I12	5	717-4535	
A3.5	I12	5	717-4543	
B1	I12	5	717-4550	
B2	I12	5	717-4568	
B3	I12	5	717-4576	
C2	I12	5	717-4584	
D3	I12	5	717-4592	

IPS e.max CAD LT A-D C14/5

Shade	Size	Qty	Item #	
A1	C14	5	717-4600	1@ \$159.95
A2	C14	5	717-4618	3@ \$154.95
A3	C14	5	717-4626	
A3.5	C14	5	717-4634	
B1	C14	5	717-4642	
B2	C14	5	717-4659	
B3	C14	5	717-4667	
C2	C14	5	717-4675	
D3	C14	5	717-4683	

IPS e.max CAD LT I12/5

Shade	Size	Qty	Item #	
BL1	I12	5	717-4691	1@ \$159.95
BL2	I12	5	717-4709	3@ \$154.95
BL3	I12	5	717-4717	
BL4	I12	5	717-4725	

IPS e.max CAD LT C14/5

Shade	Size	Qty	Item #	
BL1	C14	5	717-4733	1@ \$159.95
BL2	C14	5	717-4741	3@ \$154.95
BL3	C14	5	717-4758	
BL4	C14	5	717-4766	

IPS e.max CAD Accessories – Ivoclar Vivadent

IPS e.max CAD Crystall. Shade Refill – 3g

Shade - 0	717-5011	1@ \$50.00
Shade - 1	717-5029	
Shade - 2	717-5037	
Shade - 3	717-5045	
Shade - 4	717-5052	
Shade - 11	717-5060	
Shade - 12	717-5078	

IPS e.max CAD Crystall./Stains Refill – 1g

Stains - white	717-5169	1@ \$25.00
Stains - cream	717-5177	
Stains - sunset	717-5185	
Stains - copper	717-5193	
Stains - olive	717-5201	
Stains - khaki	717-5219	
Stains - mahogany	717-5227	

IPS e.max CAD Crystall./Glaze Paste – 3g

Each	717-5300	1@ \$50.00
------	----------	------------

IPS e.max CAD Crystall./Glaze Spray – 270ml

Each	717-5318	1@ \$129.95
------	----------	-------------

IPS e.max CAD Crystallization Tray

Each	717-5326	1@ \$120.00
------	----------	-------------

IPS e.max CAD Crystallization Pins

Pkg	717-5334	1@ \$76.00
-----	----------	------------

**inEos
Dental Digitizer**

The fastest scanner on the market

inEos makes even more of every productive minute with your inLab System. It's the best reason yet to integrate CAD/CAM technology into your dental laboratory.

inEos for inLab — The latest productivity enhancer from Sirona



The Dental Compan





Mark II VITABLOCS, 3D Master Shades - Vident

Mark II/3D Refills:

Fine particle nethaline feldspar ceramic. Indicated for anterior crowns, veneers and inlays and onlays.


Shade	Size	Qty	Type A	Type B	
1M1C	I8	10		707-3737	1@\$230.00
1M2C	I8	10		707-3745	3@\$223.00
2M1C	I8	10		707-3752	
2M2C	I8	10		707-3760	
2M3C	I8	10		707-3778	
3M1C	I8	10		707-3786	
3M2C	I8	10		707-3794	
3M3C	I8	10		707-3802	
4M2C	I8	10		707-3810	
1M1C	I10	10		707-3828	
1M2C	I10	10		707-3836	
2M1C	I10	10		707-3844	
2M2C	I10	10		707-3851	
2M3C	I10	10		707-3869	
3M1C	I10	10		707-3877	
3M2C	I10	10		707-3885	
3M3C	I10	10		707-3893	
4M2C	I10	10		707-3901	
0M1C	I12	10		707-3729	
1M1C	I12	10		707-3919	
1M2C	I12	10		707-3927	
2M1C	I12	10		707-3935	
2M2C	I12	10		707-3943	
2M3C	I12	10		707-3950	
3M1C	I12	10		707-3968	
3M2C	I12	10		707-3976	
3M3C	I12	10		707-3984	
4M2C	I12	10		707-3992	
0M1C	I14	10		707-4214	1@\$136.00
1M1C	I14	05		707-4222	3@\$131.00
1M2C	I14	05		707-4230	
2M1C	I14	05		707-4248	
2M2C	I14	05		707-4255	
2M3C	I14	05		707-4263	
3M1C	I14	05		707-4271	
3M2C	I14	05		707-4289	
3M3C	I14	05		707-4297	
4M2C	I14	05		707-4305	




Vita Classical VITABLOCS, 3D Master Shades - Vident

Fine particle nethaline feldspar ceramic. Indicated for anterior crowns, veneers and inlays and onlays.

Shade	Size	Qty	Type A	Type B	
A1	I8	10		707-3141	1@\$230.00
A2	I8	10		707-3158	3@\$223.00
A3	I8	10		707-3166	
A1	I10	10		707-3174	
A2	I10	10		707-3182	
A3	I10	10		707-3190	
A1	I12	10		707-3208	
A2	I12	10		707-3216	
A3	I12	10		707-3224	
A2	V5-12	10	707-9882		
A3	V5-12	10	707-9890		
A3.5	V5-12	10	707-9908		
B3	V5-12	10	707-9916		
A1	I14	05		707-3315	1@\$136.00
A2	I14	05		707-3323	1@\$132.00
A3	I14	05		707-3331	
A3.5	I14	05			



Type A
Mandrel
For use in CEREC 3/inLab only



Type B
Mandrel
For use in CEREC 3/
inLab/CEREC® MC XC/inLab MC XL



VITABLOCS Aesthetic Line - Vident

Translucent version of Mark II ceramic. Indicated for anterior crowns and veneers.

Shade	Size	Qty	Type A	
1M1C	V7	10	707-0485	1@ \$236.00

“I have been working with the CEREC® 3D for almost 2 years. I purchased it 3 months after opening a new office. All of my single restorations and veneers are done with CEREC. I’m happy, the patients are happy and the CEREC system continues to produce high-quality restorations that save me time and money.”

DR. OREST FRANGOPOL, *Rancho Santa Margarita, Calif.*

MILLING (CONT.)

Accu-Spray – Vident

Reflective spray scan powder specifically designed for lab use with a non-clinical propellant. Accu-Spray enhances the reflective properties of a prepared dental model fabricated with dental stone or scan stone. 200ml can.

Each707-1913 1@ \$250.00

Accu-Wax – Vident

A liquid version of scan powder that dries to a reflective powder prior to scanning. 10ml bottle specifically for lab use. Accu-Liquid is perfect for scanning implant abutments.

Each707-1939 1@ \$60.50

Accu-Liquid – Vident

This highly reflective wax provides increased laser-scanning accuracy. Perfect for blocking out undercuts, voids, defects and recontouring die stone prior to scanning.

Each707-1921 1@ \$141.00



Triluxe – Vident

Three Levels of Translucency and Saturation in One Block for CEREC 3

TriLuxe milling blocks offer three different colour saturations and three different translucencies integrated into one block, which makes it possible to copy the optical characteristics of a natural tooth. TriLuxe is separated into three softly graduated levels: the middle layer (body) represents a regular intensity, the top layer (enamel) offers a lower, less intense chroma with increased translucency and the lower layer (neck) offers high-colour saturation and a low amount of translucency. TriLuxe blocks are indicated for crowns in both the anterior and posterior regions:

Tri 12 (10mm x 12mm x 15mm) and Tri 14 (12mm x 14mm x 18mm).

Shade	Size	Qty	Type A	Type B	
1M2C	Tri 14/14	05	706-9651	707-4461	1@ \$195.00
2M2C	Tri 14/14	05	706-9669	707-4479	3@ \$189.00
3M2C	Tri 14/14	05	706-9677	707-4487	
1M2C	Tri 14	05	706-9362	707-4545	1@ \$168.00
2M2C	Tri 14	05		707-4552	3@ \$163.00
3M2C	Tri 14	05	706-9388	707-4560	
1M2C	Tri 12	10		707-4628	1@ \$317.00
2M2C	Tri 12	10		707-4636	3@ \$307.00
3M2C	Tri 12	10		707-4644	



Triluxe Forte Blocks (type B Only) – Vident

VITABLOCS TriLuxe forte are the next step in aesthetic CAD/CAM milling and are based on the outstanding aesthetics achieved with VITABLOCS TriLuxe for CEREC, which have proven their clinical efficacy in countless successful restorations since 2003. Improved are the finer nuances of the shade transition from the enamel to the neck layer. There is also more intense colour saturation and increased fluorescence in the cervical area. These new features produce a convincing and natural shade effect, and an improved chameleon blend, to the existing dentition. The result is an excellent aesthetic result without the need to individualize and characterize the shade. All-new TriLuxe forte blocks, available in sizes:

TF-12 (10mm x 12mm x 15mm),
TF-14 (12mm x 14mm x 18mm) and
TF-14/14 (14mm x 14mm x 18mm).

Shade	Size	Qty	Type B	
1M2C	TF 14/14	05	707-5633	1@ \$196.00
2M2C	TF 14/14	05	707-5641	3@ \$186.70
3M2C	TF 14/14	05	707-5658	
1M2C	TF 14	05	707-5716	1@ \$184.00
2M2C	TF 14	05	707-5724	3@ \$178.00
3M2C	TF 14	05	707-5732	
1M2C	TF 12	10	707-5781	1@ \$340.00
2M2C	TF 12	10	707-5799	3@ \$323.80
3M2C	TF 12	10	707-5807	

Vita VM9 Aesthetic Kit for VITABLOCS – Vident

Each707-1970 1@ \$1,299.00



CAD-TEMP Blocks – Vident

The VITA CAD-Temp composite blocks consist of a fiber-free, homogeneous, high-molecular and cross-linked acrylate-polymer with microfiller. CAD-Temp is used for the fabrication of multi-unit, fully or partially anatomical long-term temporary bridge restorations with up to 2 pontics, using the inLab or CEREC systems of Sirona. The bridge block can also be used for temporary crowns. The light-curing veneering composite VITA VM LC is used for building up or individualizing the shade of partially anatomical restorations. VITA CAD-Temp acrylate polymer blocks are available in the 4 shades 0M1T, 1M2T, 2M2T, 3M2T and in the size CT-40 (15.5mm x 19mm x 39mm).

(for MC XL)

Shade	Item#	
0M1T - 2 pack707-4677	1@ \$117.00
1M2T - 2 pack707-4685	6@ \$104.00
2M2T - 2 pack707-4693	
3M2T - 2 pack707-4701	
0M1T - 10 pack707-4347	1@ \$440.00
1M2T - 10 pack707-4354	6@ \$430.00
2M2T - 10 pack707-4362	
3M2T - 10 pack707-4370	

Vita Extreme Aesthetic Characterization Kit – Vident

Each707-1988 1@ \$965.00

“I have had the CEREC for a little over three years now. Dentistry became very enjoyable again after acquiring the CEREC. First, my patients appreciate the fact that I am concerned about their time. By using CEREC, I can minimize the number of appointments necessary. Secondly my practice revenue has increased by \$100,000 per year each year I have had CEREC. My dentistry has never been better! I know my clinical skills have improved. The restorations are fantastic.”

DR. BRIAN C. MCDOWELL, *Fitchburg, Mass.*




CEREC Full Contour Block Refills (Type B) – Sirona

Fine-structure feldspar ceramic blocks that are easy to polish and are similar to natural tooth enamel in terms of abrasion properties. Available in Translucent (T), Medium (M), Opaque (O) and Polychromatic (PC). Can be used in combination with the classic VITA A-D shade system as well as the VITA 3D Master shade system.

Size	Shade	Qty	Item #		Size	Shade	Qty	Item #	
8	S1-M	8	662-3011	1@ \$201.55	14	S0-M	8	662-3714	1@ \$241.50
8	S2-T	8	662-3026		14	S1-M	8	662-3722	
8	S2-M	8	662-3037		14	S2-T	8	662-3730	
8	S2-O	8	662-3045		14	S2-M	8	662-3748	
8	S3-T	8	662-3052		14	S2-O	8	662-3755	
8	S3-M	8	662-3060		14	S3-T	8	662-3763	
8	S3-O	8	662-3078		14	S3-M	8	662-3771	
8	S4-T	8	662-3086		14	S3-O	8	662-3789	
8	S4-M	8	662-3094		14	S4-T	8	662-3797	
8	S4-O	8	662-3102		14	S4-M	8	662-3805	
8	S5-M	8	662-3110		14	S4-O	8	662-3813	
10	S1-M	8	662-3128	1@ \$201.55	14	S5-M	8	662-3821	
10	S2-T	8	662-3136		12	S2-PC	8	662-4019	1@ \$299.95
10	S2-M	8	662-3144		12	S3-PC	8	662-4027	
10	S2-O	8	662-3151		12	S4-PC	8	662-4035	
10	S3-T	8	662-3169		14	S2-PC	8	662-4118	1@ \$318.50
10	S3-M	8	662-3177		14	S3-PC	8	662-4126	
10	S3-O	8	662-3185		14	S4-PC	8	662-4134	
10	S4-T	8	662-3193		14/14	S2-PC	8	662-4217	1@ \$339.50
10	S4-M	8	662-3201		14/14	S3-PC	8	662-4225	
10	S4-O	8	662-3219		14/14	S4-PC	8	662-4233	
10	S5-M	8	662-3227						
12	S0-M	8	662-3235	1@ \$201.55					
12	S1-M	8	662-3243						
12	S2-T	8	662-3250						
12	S2-M	8	662-3268						
12	S2-O	8	662-3276						
12	S3-T	8	662-3284						
12	S3-M	8	662-3292						
12	S3-O	8	662-3300						
12	S4-T	8	662-3318						
12	S4-M	8	662-3326						
12	S4-O	8	662-3334						
12	S5-M	8	662-3342						



Type A
Mandrel
For use in
CEREC 3/inLab only



Type B
Mandrel
For use in CEREC 3/
inLab/CEREC®
MC XC/inLab MC XL



Paradigm MZ100 – 3M ESPE

MZ100 Composite Blocks for CEREC

Made from 3M Z100 restorative, the Paradigm MZ100 Block is a strong, wear-resistant, highly aesthetic composite milling block for use in fabricating anterior and posterior CEREC restorations. Features:

- Easy finishing and polishing
- Extended cutting instrument life
- No glazing or firing required
- Repairable in the mouth
- No hydrofluoric acid etch

Paradigm MZ100 – Refills TYPE B

Size 10, 5 blocks/pkg		
Shade A1	502-6513	1@ \$151.45
Shade A2	502-6521	
Shade A3	502-6539	
Shade A3.5	502-6547	
Shade B3	502-6554	
Shade Enamel	502-6562	
Size 14, 5 blocks/pkg		
Shade A1	502-6760	1@ \$174.10
Shade A2	502-6778	
Shade A3	502-6786	
Shade A3.5	502-6794	
Shade B3	502-6802	
Shade Enamel	502-6810	

Paradigm C – 3M ESPE

Lucite-reinforced glass ceramic block for use with the CEREC systems. Brilliant fluorescence and translucency create a remarkable chameleon effect. Designed for inlays, onlays, crowns and veneers. 5 blocks/pkg

Shade	Size	Type B		Shade	Size	Type B	
A1	I8	502-5986	1@ \$125.85	A1	I12	502-6109	
A2	I8	502-5994		A2	I12	502-6117	
A3	I8	502-6000		A3	I12	502-6125	
A3.5	I8	502-6018		A3.5	I12	502-6133	
B3	I8	502-6026		B3	I12	502-6141	
Bleach	I8	502-6034		Bleach	I12	502-6158	
A1	I10	502-6042		A1	I14	502-6265	1@ \$125.85
A2	I10	502-6059		A2	I14	502-6273	
A3	I10	502-6067		A3	I14	502-6281	
A3.5	I10	502-6075		A3.5	I14	502-6299	
B3	I10	502-6083		B3	I14	502-6307	
Bleach	I10	502-6091		Bleach	I14	502-6315	



MILLING (CONT.)

Type A



Mandrel
For use in
CEREC 3/inLab only

Type B











Mandrel
For use in CEREC 3/
inLab/CEREC® MC XC/inLab
MC XL





Beauty that radiates from the core



inLab® Blocs – Vita

inLab® Blocs – Vita

 <p>CS-11 10/Pkg Spinell Type A 707-9494 1@ \$415.00 Type B 707-2978 3@ \$385.00</p>	 <p>BZ-33 2/Pkg Zirconia Type B 707-2911 1@ \$246.00 6@ \$240.00</p>
 <p>CA-12 10/Pkg Alumina Type B 707-2960 1@ \$288.00 6@ \$280.00</p>	 <p>BZ-40 2/Pkg Zirconia Type A 706-9891 1@ \$313.00 Type B 707-2929 1@ \$313.00 6@ \$305.00</p>
 <p>CZ-12 10/Pkg Zirconia Type A 707-2457 1@ \$349.00 Type B 707-2945 1@ \$343.00 6@ \$336.00</p>	 <p>AL-20 In-Ceram 2000 4/Pkg Type A 707-1012 1@ \$147.00 6@ \$143.00 24/Pkg Type A 707-1020 1@ \$815.00 6@ \$795.00 Type B 707-2879 1@ \$815.00 6@ \$795.00</p>
 <p>CZ-18 5/Pkg Zirconia Type A 706-9916 1@ \$230.00 Type B 707-2937 1@ \$381.00 6@ \$374.00</p>	 <p>AL-40 In-Ceram 2000 2/Pkg Type A 707-1038 1@ \$256.00 3@ \$250.00 10/Pkg Type A 707-1046 1@ \$1,058.00 3@ \$1,035.00 Type B 707-2887 1@ \$256.00 3@ \$250.00 Type B 707-2895 1@ \$1,058.00 3@ \$1,035.00</p>

 <p>YZ-20/15 Yttria Stabilized Zirconia 5/Pkg Type A 707-2119 1@ \$169.00 Type B 707-3000 6@ \$164.00 25/Pkg Type A 707-2127 1@ \$965.00 Type B 707-3018 6@ \$760.00</p>
 <p>YZ-20/19 Yttria Stabilized Zirconia 4/Pkg Type A 707-2010 1@ \$170.00 Type B 707-3026 24/Pkg Type A 707-2028 1@ \$930.00 Type B 707-3034</p>
 <p>YZ-40/15 Yttria Stabilized Zirconia 2/Pkg Type A 707-2135 1@ \$220.00 Type B 707-3042 10/Pkg Type A 707-2143 1@ \$1,030.00 Type B 707-3059</p>
 <p>YZ-40/19 Yttria Stabilized Zirconia 2/Pkg Type B 707-3067 1@ \$258.00 6@ \$253.00 10/Pkg Type A 707-2044 1@ \$1,295.00 Type B 707-3075 6@ \$1,185.00</p>



CAD-Waxx – Vita

High molecular translucent acrylpolymer for the construction of PFM metal frameworks for crowns and bridges. Burns away cleanly using the lost-wax technique.

CW-40/15		10/Pkg	
2/Pkg		Type A	
Type A	707-2051	1@ \$22.00	Type A
Type B	707-2994		Type B
			707-2069
			1@ \$91.00
			707-2986



55mm Flip Block – Vita

YZ-55, 1/Pkg, Yttria Stabilized Zirconia	
Type A	706-9636
Type B	707-3083
	1@ \$197.00

IPS. e.max ZirCAD Basic Kit (Type B) – Ivoclar Vivadent



IPS e.max ZirCAD blocks are made of yttrium oxide-stabilized zirconium oxide. In the partially sintered, chalk-like state, IPS e.max ZirCAD can be easily milled in the inLab (Sirona) unit. The milling process is conducted by enlarging the framework by 20%.

During the subsequent controlled sintering process in the Sintramat high-temperature furnace, the frameworks shrink to their final size and are imparted with the high fracture strength and high fracture toughness. With its outstanding strength, IPS e.max ZirCAD thus also fulfills clinical requirements related to masticatory forces, particularly in the posterior region.

A specially designed fluorescent zircon liner provides an optimum bond and is used to apply the basic shade on the frameworks.

The highlights:

- Extended range of indications thanks to the outstanding strength of 900 MPa and high-fracture toughness
- Optimum bond and application of the basic shade with the specially developed IPS e.max Ceram ZirLiner
- Highly aesthetic fluorapatite glass-ceramic ingot IPS e.max ZirPress or efficient use in the press-on technique
- Conventional cementation

C15	C15/L	B40	B40L
5/Pkg	5/Pkg	3/Pkg	3/Pkg
IPS e.max ZirCAD	IPS e.max ZirCAD	IPS e.max ZirCAD	IPS e.max ZirCAD
714-3704	714-3779	714-3712	714-3803
1@ \$126.50	1@ \$143.75	1@ \$210.45	1@ \$231.15
714-3696	714-3787	714-3795	714-3811
25/Pkg	25/Pkg	9/Pkg	9/Pkg
1@ \$603.75	1@ \$690.00	1@ \$610.65	1@ \$672.75

Kit
Each714-3761 1@ \$1,149.95

E.MAX CAD Basic Kit (Type B) – Ivoclar Vivadent



Lithium disilicate glass-ceramic block for the CAD/CAM technique. IPS e.max CAD blocks feature a characteristic and remarkable shading ranging from white and blue to blue-grey, which is a result of the composition of the glass-ceramic.

Thanks to an innovative manufacturing process, IPS e.max CAD features an impressive homogeneity and can be easily processed using the inLab or Everest system. The strength of the material in the pre-crystal blue phase can be compared with that of other glass ceramics currently available on the market.

After the blocks have been milled, the material undergoes a crystallization process in the Programat P300, P500 or EP 600 Combi furnace. The crystallization process results in a change of structure, which imparts the glass-ceramic with a strength of 360 MPa and the actual optical properties, such as shade, translucency and brightness.

The tooth-shaded IPS e.max CAD frameworks form the ideal basis for the fabrication of aesthetic crowns in the anterior and premolar region, which are completed with the optimally coordinated veneering ceramic IPS e.max Ceram.

The highlights:

- Efficient, economical processing due to easy grinding properties
- High final strength of 360 MPa
- Easy and fast crystallization (850°C, 30 minutes)
- Crystallization without shrinkage and complex infiltration processes
- Lifelike and translucent aesthetics

C14
5/Pkg
IPS e.max CAD
Shade 0 714-3613
Shade 1 714-3621
Shade 2 714-3639
Shade 3 714-3647
Shade 4 714-3654
1@ \$144.95

Kit
Each714-3753 1@ \$1,149.95

MILLING (CONT.)



InCoris ZL (Type B) – Sirona

Coping/Framework Blocks

Size	Shade	Qty	Item #		Size	Shade	Qty	Item #	
20/19	F0.5	6	662-2013	1@ \$319.90	40/19	F0.5	3	662-2450	1@ \$469.70
20/19	F1	6	662-2021	6@ \$303.10	40/19	F1	3	662-2468	6@ \$444.50
20/19	F2	6	662-2039		40/19	F2	3	662-2476	
20/19	F3	6	662-2047		40/19	F3	3	662-2484	
20/19	F4.5	6	662-2054		40/19	F4.5	3	662-2492	
20/19	F0.5	24	662-2120	1@ \$1,180.20	40/19	F0.5	12	662-2567	1@ \$1,779.40
20/19	F1	24	662-2138	6@ \$1,115.80	40/19	F1	12	662-2575	6@ \$1,684.20
20/19	F2	24	662-2146		40/19	F2	12	662-2583	
20/19	F3	24	662-2153		40/19	F3	12	662-2591	
20/19	F4.5	24	662-2161		40/19	F4.5	12	662-2609	
40/15	F0.5	3	662-2237	1@ \$426.30	55/19	F0.5	2	662-2674	1@ \$490.00
40/15	F1	3	662-2245	6@ \$403.20	55/19	F1	2	662-2682	6@ \$463.40
40/15	F2	3	662-2252		55/19	F2	2	662-2690	
40/15	F3	3	662-2260		55/19	F3	2	662-2708	
40/15	F4.5	3	662-2278		55/19	F4.5	2	662-2716	
40/15	F0.5	12	662-2344	1@ \$1,600.20					
40/15	F1	12	662-2351	6@ \$1,513.40					
40/15	F2	12	662-2369						
40/15	F3	12	662-2377						
40/15	F4.5	12	662-2385						



InCoris AL (Type B) – Sirona

Coping/Framework Blocks

Size	Shade	Qty	Item #	
20/19	F0.7	6	662-2773	1@ \$284.20 6@ \$264.60
20/19	F0.7	24	662-2781	1@ \$950.60 6@ \$886.20
40/15	F0.7	3	662-2799	1@ \$276.15 6@ \$261.10
40/15	F0.7	12	662-2807	1@ \$989.80 6@ \$936.60
40/19	F0.7	3	662-2815	1@ \$330.05 6@ \$308.00
40/19	F0.7	12	662-2823	1@ \$1,253.00 6@ \$1,216.60

Thank you
> for making us #1



PATTERSON
DENTAL/DENTAIRE



CAESY – Winner



Sirona CEREC® 3D – Winner



EagleSoft – Winner



Schick® CDR – Winner

**Your single-source
technology solution!**

Call 800.873.7683
to learn more about our
award-winning technology!

POWDER



Powder Meister – Dr David Lawler

New patented nozzle design with on/off control allows directional control of powder flow.

- Clog-free operation
- New – built-in air pressure regulator allows for easy powder control
- Air line filter removes moisture from the air line
- New – smaller bottle easily fits all hand sizes
- Used worldwide by more CEREC® users than any other powdering device

Each.....040-7007 1@ \$872.40

Powder Meister Holder (New)

Stainless steel container keeps the All-New Powder Meister upright and easy to reach.

Each.....461-5092 1@ \$118.45

Powder Meister Autoclavable Replacement Nozzle

Each.....461-5209 1@ \$104.90



PowderPerfekt™ – Our Glass LTD

- Ergonomically correct to reduce risk of repetitive stress injury and enhance productivity
- Newly designed applicator handle allows precise metering of the flow
- The air source comes from a filtered dry air syringe with single button control
- Thumb-activated button on the syringe easily controls the flow of powder
- Syringe design allows powder to be placed in the most difficult access cases
- With the provided handy Access® tip installed on the air syringe, it becomes a multi-purpose instrument useful for blowing clean dry air
- Tips and nozzles are autoclavable

Each040-2990 1@ \$859.85

Parts Kit

Each547-2709 1@ \$120.35

Nozzle, 3/Pkg

Each547-2717 1@ \$109.15

Regulator With Mini-Gauge, Each

Each547-2725 1@ \$97.95



CEREC® Spray Head – Vident

Used in conjunction with CEREC® propellant and powder. Propels CEREC® powder through powder application bottle onto tooth surface.

Each707-9940 1@ \$65.25

CEREC® Propellant – Vident

Special propellant for application of CEREC® powder. Produces a uniform film of powder on the preparation.

250ml Can707-9932 1@ \$104.90



CEREC® Powder – Vident

A white contrast powder used to aid in optical impression – surface of preparation becomes opaque non-reflective. Easily removable with water. Applied with CEREC® propellant.

12gm Bottle707-9924 1@ \$34.30



CEREC® Liquid – Vident

Adhesive for CEREC® powder.

6ml Bottle707-9957 1@ \$48.60



IPS Contrast Spray (Chairside and Labside) – Ivoclar Vivadent

A spray is applied before the scanning procedure of the preparation or the model in order to adjust the optical properties of dentin and enamel and/or the plaster model and thus enable an accurate scan.

With the atomizing nozzle, an optimum covering layer with a detailed representation of the edges is achieved easily and efficiently with only a short spray discharge.

50ml.....715-5005 1@ \$216.75

200ml.....715-5013 1@ \$273.95



Powder Pro for CEREC – ADI

The Powder Pro is a true breakthrough in precision imaging with the CEREC® 3D Optical Camera. Its unique design allows for the finest of CAD/CAM restoration.

- Unique ergonomical “Power Pen” offers unparalleled precision powder application
- Stationary Powder Chamber - ensures a uniform powder layer creating the finest of restorations
- On-board air valve allows for perfect titration of powder flow
- Rheostat control gives you the ability for spot application
- 30-second setup right next to your high-speed allows for instant accessibility
- Handpiece and tips are autoclavable

Each035-2047 1@ \$1,220.00

Powder Pro Tips, 10/Pkg

Each Package100-8804 1@ \$54.55



Scan Film – Vident

An alternative to powder and propellants, Scan Film is applied as a liquid and allowed to dry into a homogeneous powder surface. 10 micron film thickness.

5ml Bottle708-0070 1@ \$253.00

Scan Film Single Unit Dose

0.1ml/60 Pack708-0088 1@ \$277.00

POWDER (CONT.)



Scan Spray – Vident

Approximately 150 applications per container.
 Each708-0096 1@ \$233.00



Aesthetic - Base Gold – Dentona

Die Stone Suitable for CAD/CAM

- High-contrast colour that creates a harmonic overall picture
- Minimum expansion values - models have dimensional stability
- A pleasant vanilla scent

5kg
 Each279-3495 1@ \$38.95
 25kg
 Each279-3529 1@ \$99.95



Dentatec – Sirona

Cleaner and lubricant. 1,000ml. of water additive to enhance milling operations.
 Each636-6223 1@ \$102.55

ACCESSORIES



R-SI-Line Metal-Bite – R-Dental

Antagonist Registration for CEREC

For the powderless optical 3D data registration of occlusal zones of antagonists (CAD/CAM/CIM - Impression Technology). The material is dark grey, extremely hard and fast (hardening time in the mouth about 40 seconds) bite registration material based on vinyl polysiloxane (A-silicone). Additionally it has high dimension stability, high standing, is inflexible and very suitable for cutting and contouring – not crumbly. Safety cartridge system. Contains 2 - 50ml cartridges and 12 mixing cannulas

Box.....592-7108 1@ \$67.05
 Mixing Cannulas (Green)
 Box of 50.....592-7181 1@ \$48.30

Diamond Die Scan Stone – Hi Tech

Sirona-approved grey model stone. This gypsum product is approved for scanning with the inLab/inEos, internally or externally without the need for powder application.

25lb
 Each362-9003 1@ \$90.95



Replacement Filter – Sirona

Replacement filter for CEREC® MC XL and inLab MC XL milling units
 Each84050248 1@ \$38.50



Dentatec – Sirona

Cleaner and lubricant, 250ml water additive to enhance milling operations.
 Each636-6025 1@ \$64.05

Camera Support – Sirona

Protects lens from scratches. Camera Support is the most effective way to protect the camera – rather than supporting the camera lens by resting it on the distal tooth, which can scratch the lens.



6/Pkg.....636-6371 1@ \$15.05

Diamond Wheel Disc – Sirona

Package of 3 discs.
 Pkg.....707-9973 1@ \$234.00



CEREC® Blue Model – Sirona

Each636-6280









Calibration Test Stubs – Vident

Each636-6033 1@ \$5.80

MILLING BURS





CEREC MC XL AND INLAB MC XL MILLING INSTRUMENTS – SIRONA



	Qty.	Patterson #	Sirona #	Retail Price	
Cylinder Pointed Bur - 12mm	6/pkg	636-6967	61 43 866	\$249.80	
Cylinder Pointed Bur - 12mm "S"	6/pkg	636-6991	62 08 685	\$219.50	
Cylinder Pointed Bur - 20mm	6/pkg	636-6983	61 43 874	\$283.80	
Step Bur - 12mm "S"	6/pkg	636-7015	62 08 677	\$219.50	
Step Bur - 12mm	6/pkg	636-6959	61 43 841	\$249.80	
Step Bur - 20mm	6/pkg	636-6975	61 43 858	\$283.80	

CEREC 3 COMPACT MILLING INSTRUMENTS – SIRONA



	Qty.	Patterson #	Sirona #	Retail Price	
Cylinder Bur 1.2mm diameter	10/pkg	636-5746	54 66 151	\$219.45	
Cylinder Bur 1.6mm diameter	6/pkg	636-5803	54 66 193	\$230.30	
Cylinder Pointed Bur	6/pkg	636-6363	59 45 535	\$231.70	
Step Bur - 10mm	6/pkg	636-6926	60 89 010	\$199.55	

INLAB COMPACT MILLING INSTRUMENTS – SIRONA



	Qty.	Patterson #	Sirona #	Retail Price	
Cylinder Bur 1.2mm diameter	10/pkg	636-5746	54 66 151	\$219.45	
Cylinder Bur 1.6mm diameter	6/pkg	636-5803	54 66 193	\$230.30	
Cylinder Pointed Bur	6/pkg	636-6363	59 45 535	\$231.70	
Cone Bur - 10mm	6/pkg	636-6298	58 85 103	\$236.60	
Cone Bur - 12mm	6/pkg	636-6355	59 35 668	\$280.00	
Cone Bur - 14mm (Flip Block)	6/pkg	636-6397	59 99 771	\$182.70	
Step Bur - 10mm	6/pkg	636-6926	60 89 010	\$199.55	
Step Bur - 12mm	6/pkg	636-6918	60 52 265	\$253.40	
Step Bur - 14mm	6/pkg	636-6934	60 89 028	\$245.00	

ABRASIVES

LS-7564 Cerec® Finishing Set - Axis

Developed by Dr. James Klim

Each 115-1943 1@ \$122.55



FIG. NO. Size 1/10mm Length (mm)	P3035-3	P30035-3	P0800	H7901-MW	H7404-MW	F368 023 4.5	SF368 023 4.0
GRIT	CeraGlaze® Polishers		PDO™	MW Series		Fine	Superfine
SHANK	RA (21)			FG (31)			
REORDER #	121-7868	121-7801	121-8528	115-0911	115-0937	545-3295	545-6686

LS-910 CeraGlaze® Porcelain Adjusting and Polishing Lab Set - Axis

Each 115-4921 1@ \$250.20




FIG. NO. Size 1/10mm Length (mm)	P301	P3001	P30001
TYPE	Pre-polishing	Refining	Reglaze/High Shine
SHANK	HP (11)		
REORDER #	121-6142	121-6126	121-6118

LS-911 CeraGlaze® Extraoral Porcelain Adjusting and Polishing Set - Axis

Each 115-4913 1@ \$292.70




FIG. NO. Size 1/10mm Length (mm)	P344	P343	P3044	P304	P30004	P30043
TYPE	Pre-polishing		Refining		Reglaze/High Shine	
SHANK	HP (11)					
REORDER #	121-8015	121-8056	116-0118	116-0134	116-0142	116-0167

LS-506 CeraGlaze® Ultimate Porcelain Polishing Set (Intraoral) – Axis

Each.....121-8577 1@ \$181.25



FIG. NO.	P332	P335	P336	P3032	P3035	P3036	P30032	P30035	P30036
Size 1/10mm Length (mm)									
GRIT	Green Pre-polishing			Blue Refining			Yellow Reglazing		
SHANK	RA (21)								
REORDER #	121-7512	121-7520	121-7538	121-7546	121-7553	121-7561	116-0068	116-0084	121-7579

LS-904 CeraGlaze® Porcelain Polishing Set (Extraoral) – Axis

Each.....116-0209 1@ \$261.25



FIG. NO.	P0301	P3044	P3042	P3043	P30044	P30042	P30043
Size 1/10mm Length (mm)	220						
GRIT	White Pre-polish	CeraGlaze® Polishers					
SHANK	HP (11)						
REORDER #	115-8054	116-0118	116-0126	116-0134	116-0142	116-0159	116-0167

LS-7527 CEREC Prep N' Glaze Logic Set – Axis

Developed by Dr. Mark R. Morin and Dr. Douglas M. Voiers

Each.....122-6125 1@ \$187.05



FIG. NO.	807	889	862	379	837	80	P3032	P30032	P30033	P30035
Size 1/10mm Length (mm)	018 5.0	027 7.0	012 8.0	023 4.5	014M 8.0	014M 8.0				
GRIT	Coarse	Coarse	Fine	Fine	Coarse	Coarse	CeraGlaze® Polishers			
SHANK	FG (31)			FG SS (33)			RA (21)			
REORDER #	121-0194	121-2497	545-5605	545-3386	122-3106	120-9063	121-7546	116-0068	116-0076	116-0084

LS-7542 CEREC Art and Passion Intraoral Set – Axis

By Dr. Rich Masek and Dr. Steve Tsotsos

Each.....122-6083 1@ \$169.65



FIG. NO.	837KRT	830L	368	833	836KR	845KR	845KR	P332	P3032	P30032
Size 1/10mm Length (mm)	014 8.0	016 5.0	018 4.5	018 5.9	012 6.0	018 4.0	018 4.0			
GRIT	Turbo	Coarse	Fine	Fine	Fine	Fine	Fine	CeraGlaze® Polishers		
SHANK	FG (31)							RA (21)		
REORDER #	120-5699	121-0251	545-3261	121-0905	545-4046	121-0988	115-8245	121-7512	121-7546	116-0068

ABRASIVES (CONT.)

**LS-7543 CEREC Art and Passion
Extraoral Set – Axis**

By Dr. Rich Masek and Dr. Steve Tsotsos
Each122-6091 1@ \$319.95



FIG. NO.	9805 G8001 P344 P3044 P30044 P3043 P30043
Size 1/10mm	037 050
Length (mm)	2.5 13.0
GRIT	Medium All Ceramic CeraGlaze® Polishes
SHANK	HP (11)
REORDER #	120-7315 115-5035 121-8015 116-0118 116-0142 116-0134 116-0167

The CeraMaster Finishing and Polishing Kit (CA) - SHOFU

Designed to finish and polish porcelain restorations without the need for reglazing. The kit contains the exceptional performing Dura-Green Stones for fast contouring and finishing and DuraWhite Stones for fine finishing. The kit also contains CeraMaster Coarse and Reality 5 Star Award winning CeraMaster diamond impregnated polishing points designed to finish, polish and superpolish porcelain to the point of a high glaze.

Each634-4097 1@ \$162.40



Dura-Green Stones Dura-White CeraMaster Coarse CeraMaster

SHAPE	CN1	FL2	RE1	CN1	FL2	RE1	Knife	Cup	Mini-Pt.	Knife	Cup	Mini-Pt.
ISO No.	024	025	062	024	025	062	090	065	030	065	030	125
PN	0052	0060	0086	0221	0223	0229	122C	123C	124C	122	123	124

All-Ceramic Finishing and Polishing Kit (HP) - SHOFU

Designed to finish and polish porcelain restorations prior to insertion. The kit contains Dura-Green Stones, Dura-Green Fine Stones (Shofu's NEW stones are excellent for adding surface characteristics and adjusting details), CeraMaster Coarse and CeraMaster.

Each634-4105 1@ \$157.50



Dura-Green Stones Dura-Green Fine CeraMaster Coarse CeraMaster

SHAPE	WH6	IC9	PC1	TC1	FL3	PC1	Bullet	KN7	WH6	Bullet	KN7	WH6
ISO No.	130	135	030	020	030	030	040	125	120	090	120	220
PN	0048	0022	0030	039F	011F	030F	126C	127C	128C	0126	0127	0128

2 Striper® CEREC® 3D Kit – Premier Dental

Prep and Polish Kit

Each585-8584 1@ \$203.70




FIG. NO.	317	515	2003	260	2014	F1S	285 MF2
Abrasive Length (mm)	8.0	5.0	4.5	5.9	6.0	4.0	4.0
GRIT	Coarse	Coarse	Coarse	Fine	Extra Fine	Very Fine	20µ Polish
SHANK	FG						
REORDER #	572-4794	571-8416	572-7243	571-8473	572-7615	571-4951	572-4323



Diamond Twist – Premier Dental

Designed specifically for ceramic crowns, but also works on all porcelains and composites.

Diamond Twist SCL Kit

Unique Two-Step CEREC® Polishing System used primarily immediately after the crown is milled. A simple Two-Step procedure that uses the unique Fibra Points for smoothing the restoration and a special formula Diamond Twist polishing paste, used with 100% soft wool points to achieve an exceptional glaze appearance.

Kit509-9510 1@ \$177.75

Diamond Twist SCL Kit Refills

Diamond Twist Polishing Paste

3gm Syringe509-9536 1@ \$60.15

Felt Points

W-803 Wheel, 3/4" x 1/4" Soft

5/Pack509-9577 1@ \$7.00

W-802 Wheel, 3/4" x 1/8" Soft

12/Pack509-9585 1@ \$12.95

P-800 Point, Hard

Each509-9593 1@ \$11.35

B-801 Barrel, Soft

Each509-9601 1@ \$11.35

Fibra Points

FB3M Barrel

Each509-9619 1@ \$15.40

FP2M Taper

Each509-9627 1@ \$13.80

FW40M Wheel, Thin

6/Pack509-9635 1@ \$41.65

FW50M Wheel, Thick

4/Pack509-9643 1@ \$41.30

Mandrels

M804 Stainless Steel

Each509-9668 1@ \$4.70

M600 SurGrip Stainless Steel

Each509-9676 1@ \$14.35



Diamond Twist

Used chairside for Int contact adjustments polishing brushes for Micro-Cloth Disks for flat areas.

Kit509-9502 1@ \$96.55

Diamond Twist SCO Refills

Diamond Twist Polishing Paste

6gm Twist Tube509-9544 1@ \$82.25

Polishing Brush

BR70, RA 25/Pack.....509-9908 1@ \$39.10

Micro-Cloth Disks and Mandrel

MC 500, 1/2" disk

25/Pack509-9700 1@ \$27.30

MC 312, 5/16" disk

25/Pack509-9718 1@ \$29.35

230 RA Mandrels

4/Pack509-9650 1@ \$15.85



Occlusal Reduction Dr. Kometas

Benefits of Ideal Occlusion

- Minimize chairside adjustment
 - Ensure proper occlusal reduction
- Easy to use – first start the depth grooves in the central fossa and extend to all fissures. Allow the self-limiting feature to touch the occlusal surface and then push the bar laterally. Then use your preferred axial reduction bur to remove any remaining occlusal tooth.



Occlusal Reduction Bur Kit

Kit includes a bur block and 6 burs (1.2mm, 1.5mm, 1.8mm, 2.0mm, 2.2mm and 2.4mm).

Kit509-9726 1@ \$111.65

Reduction Bur Refills

5/Pack509-9742 2.0mm509-9775

1.5mm509-9759 2.2mm509-9783

1.8mm509-9767 2.4mm509-9791

Profin® – Dentatus

CEREC Specific Kits Contain the Recommended by CEREC® Instrument

Contains: one Profin instrument Lamineer tips, 100 polishing tip tip stand and a 12-minute videotape on CD-ROM.

Doriot

Fits Star and Midwest straight nose cones.

Kit765-6606 1@ \$1,026.20

E Type

Fits all other low-speed systems including

KaVo non-optics.

Kit765-6614 1@ \$1,120.00

CEREC Profin Lamineer Tips

Six Profin Lamineer tips recommended by CEREC instructors for use with CEREC® all-ceramic restorations.

Kit765-6598 1@ \$340.20



Occlubrush Polishers – Kerr

The Occlubrush is a fine-surface polishing brush system. Pre-polish and finish restorations to a high gloss with this instrument. Special impregnated fibers in each handle a polishing instrument, eliminating the need for paste.



Occlubrush Assortment Introductory Kit

Contains 5 each (Regular Cup, Small Cup, Point)

Kit430-6304 1@ \$109.65

Occlubrush Refills, 10/Pkg

Regular Cup430-6320 1@ \$75.50

Small Cup430-6338 3@ \$73.55

Point430-6346

STAINING & GLAZING



IPS Empress Universal Shade/Stains Kit – Ivoclar Vivadent

The IPS Empress Universal Shade/Stains Kit contain a selection of the most popular Dentin and intensive Stains. These stains can be used to apply surface staining and characterization. The Kit also contains a smooth and fine-grained glazing paste for final glazing of the IPS Empress restorations.

IPS Empress® Universal Shade/Stains Kit A-D
 Contains 1 each IPS Empress Universal Shade 3g A1, A2/A3/A3.5, B1, B2/B3/B4, C1/C2, D2/D3, 1 – IPS Empress Universal Shade Incisal 3g I1 and I2, 1 each IPS Empress Universal Stains 1g white, orange, mahogany, khaki and maroon, 1 – IPS Empress Universal Glaze 3g, 1 – IPS Empress Universal Glaze and Stain Liquid 15ml, 1 – Tweezer, 1 – Brush 00.
 Each 715-5039 **1@ \$440.00**

IPS Empress® Universal Shade/Stains Kit Chromascop
 Contains 1 each IPS Empress Universal Shade 3g 110/120, 130, 140/210, 220/230, 310, 410/420, 1 each IPS Empress Universal Shade Incisal 3g I1 and I2, 1 each IPS Empress Universal Stains 1g white, orange, mahogany, khaki and maroon, 1 – IPS Empress Universal Glaze Paste 3g, 1 – IPS Empress Universal Glaze and Stain Liquid 15ml, 1 – Tweezer, 1– Brush 00.
 Each 715-5021 **1@ \$440.00**

IPS Empress® Universal Shades
 The IPS Empress Universal Shades are available in 15 Chromascop and 9 A-D shades. They are used for the exact reproduction of the Chromascop and A-D shades.

Shade	Size	Item	
110/120	3g	715-5062	1@ \$47.95
140/210	3g	715-5070	
220/230	3g	715-5088	
410/420	3g	715-5096	
430/440	3g	715-5104	
130	3g	715-5112	
240	3g	715-5120	
310	3g	715-5138	
320	3g	715-5146	
330	3g	715-5153	
340	3g	715-5161	
510	3g	715-5179	
520	3g	715-5187	
530	3g	715-5195	
540	3g	715-5203	
A1	3g	715-5211	
A2/A3/A3.5	3g	715-5229	
A4	3g	715-5237	
B1	3g	715-5245	
B2/B3/B4	3g	715-5252	
C1/C2	3g	715-5260	
C3/C4	3g	715-5278	
D2/D3	3g	715-5286	
D4	3g	715-5294	

IPS Empress® Universal Stains
 These intensive stains are available in 14 different shades. They are used to imitate natural characteristics. With the three complementary shades basic blue, basic red and basic yellow, the shades may be individually adjusted.

Shade	Size	Item	
White	1g	715-5419	1@ \$25.00
Black	1g	715-5427	
Mahogany	1g	715-5437	
Khaki	1g	715-5443	
Orange	1g	715-5450	
Grey	1g	715-5468	
Vanilla	1g	715-5476	
Crackliner	1g	715-5484	
Olive	1g	715-5492	
Yellow	1g	715-5500	
Basic Blue	1g	715-5518	
Basic Yellow	1g	715-5526	
Basic Red	1g	715-5534	
Maroon	1g	715-5542	

IPS Empress® Paste
 Glaze 3g 715-5708 **1@ \$47.95**

IPS Empress® Liquid Stain
 Glaze 15ml 715-5716 **1@ \$21.25**



Sinfony Magic Set – 3M ESPE
For Characterization of Indirect Composites

Sinfony Magic Refills **1@ \$95.80**

I1 Ocean Blue (49410)	503-8054
I2 Atlantis (49420)	503-8062
I3 Chestnut (49430)	503-8070
I4 Havana (49440)	503-8088
I5 Orange (49450)	503-8096
I6 Khaki (49460)	503-8104
I7 Vanilla (49470)	503-8112
I8 Honey Yellow (49480)	503-8120
I9 Gingiva (49490)	503-8138
I10 Red (49491)	503-8146
I11 Terra Cotta (49492)	503-8153

Sinfony Dispenser Handpiece
 Each 503-8211 **1@ \$59.00**

Sinfony Transparent Opal Material
 T1 Neutral 503-8203 **1@ \$76.25**



Tetric Color – Ivoclar Vivadent


Light-cured restorative tinting and opaquing system derived from the Tetric family of products. Contains 7 - 1gm syringes (shades white, light yellow, bluish grey, ochre, medium brown, dark brown and black), 1 Vivapad, 40 (0.4mm) tips

Each 708-9998 **1@ \$165.95**

Tetric Color Refills
 Contains 1 - 1gm syringe and tips

White.....	708-7497	1@ \$28.25
Light Yellow	708-7505	
Bluish Grey	708-7513	
Ochre.....	708-7521	
Medium Brown.....	708-7539	
Dark Brown.....	708-7547	
Black	708-7554	

SHADE GUIDES



Classic Vita Lumin® Shade Guide – Vident

Simple, accurate selection. Corresponds to all Vita porcelain systems.
Each707-9601 1@ \$131.00



Vita Shade Guide for CEREC® Blocks – Vident

This shade guide contains 13 removable tabs indicating all Vita shades for CEREC® Blocks, including Vita Classic, 3D Master, Esthetic and bleach shade. (G55-900)
Each707-0089 1@ \$87.75



VITA System 3D-Master® Tooth Guide Kit – Vident

Dental professionals around the world have discovered precise and repeatable shade selection and improved communication, using the VITA 3D-Master® Tooth Guide. Based on true colour science, shades can now be selected with incredible accuracy using a fast, efficient, three-step process applying Value (lightness), Chroma (saturation) and Hue (Yellow-Red shift). VITA offers 29 shade tabs, including 3 bleached tab examples, in one autoclavable guide. System 3D-Master is supported by VITA tooth-coloured restorative and prosthetic materials used for the veneering of metal ceramics, densely sintered zirconia and alumina, as well as indirect laboratory composites, denture teeth and direct micro-hybrid composites.
Each708-0815 1@ \$114.00

VITA Akzent Stain Kit – Vident

Developed by Claude Sieber, these stains enable the ceramist or dentist to create more life-like three-dimensional restorations. The use of the VITA Akzent Finishing Agent eliminates the shadow caused by translucent glazes, resulting in a semi-translucent glaze finish for excellent light transmission.

Includes 20 - 3gm Akzent stains, plus 1 each of Akzent finishing agent (5gm), Akzent glaze (5gm), Akzent stain liquid (20ml), Akzent mixing palette, glass spatula, sable brush, Akzent shade sample disc, instruction and conversion chart.

Kit706-9008 1@ \$1,250.00



Firing Stands For Vident

Contains: 3 firing stands
Kit706-9909 1@ \$43.00



Fibrous Firing Pads – Vident

Set of 2707-0147 1@ \$14.60

VITA Akzent Stain Kit Refills – Vident

Akzent Finishing Agent, 5gm
Each706-9453 1@ \$77.00

Akzent Glaze, 5gm
Each706-9602 1@ \$60.50

Akzent Stain Liquid, 20ml
Each706-9610 1@ \$35.25

VITA Akzent Stain Kit Refills – Vident

Stain Refills, 1 Each, 3gm 1@ \$53.50

- Birch.....706-9115 Tabac.....706-9214
- Mellow Redwood.....706-9222
- Yellow.....706-9123 Shak.....706-9230
- Sun Kiss.....706-9131 Jet Black.....706-9248
- Inca Gold.....706-9149 Brown Stone..706-9255
- Gold Earth.....706-9156 Seal Brown....706-9263
- Curry.....706-9164 Niagara.....706-9271
- Adobe.....706-9172 Fumo 01.....706-9289
- Almond.....706-9180 Fumo 02.....706-9297
- Burnt Clay....706-9198 Fumo 03.....706-9305
- Muted Clay ..706-9206

CEREC® Firing Tray W – Vident

Each708-0658 1@ \$47.40



CEREC® Brush IC-4 – Vident

Each708-0633 1@ \$89.95

Tweezers For Crowns – Ivoclar Vivadent

Each708-1730 1@ \$89.95

VITA Shading Pastes – Vident

Fluorescent Shading

These fine grained staining pastes are designed for use on porcelain and designed for use on porcelain. Ready to use, they do not require mixing a powder and a liquid.

Kit contains: 8 - 4gm shading pastes (shades are listed at right under refills), 1 - 7gm Vita glaze paste, 1 - 15ml Vita stain liquid, 1 shade sample disc, 1 mixing palette, 1 agate mixing spatula, stain brush with tips and 1 flat squirrel hair brush.

Kit708-0674 1@ \$960.00



VITA Stain Kit Refills – Vident

Stain Refills 4gm Each 1@ \$76.00

- White.....708-0773
- Orange.....708-0781
- Ochre.....708-0799
- Khaki.....708-0831
- Olive.....708-0849
- Medium Brown.....708-0856
- Reddish Brown.....708-0864
- Mauve.....708-0872


VITA Stain Kit Refills – Vident

Misc. Refills

VITA Glaze Paste, 7gm
Each708-0948 1@ \$127.00

VITA Stain Liquid, 15ml
Each708-0955 1@ \$52.25

SHADE GUIDES (CONT.)



VITABlocs® Guide 3D-Master® – Vident

Improves the shade selection by matching Vitabloc® shades faster and now more precise than ever. The new Vitabloc shade tabs offer a realistic shade and translucency that's relevant to the final restoration appearance. All the shade tabs are fabricated from the same ceramic materials as used in the original Vitablocs. The new VITABLOCKS Guide is based on the original VITA System 3D-Master System, reproducing the three natural tooth shade dimensions of Lightness, Chroma and Hue perfectly. VITA future-oriented for more than 80 years.

Each707-5906 1@ \$114.00

HYDROFLUORIC ACID

CEREC® Ceramics Etch – Vident

A 5 percent hydrofluoric acid gel for ceramic restorations. Cannot be used on metal. Do not use on any surface requiring tooth structure.

6ml Bottle707-9965



Etch Ease – Densol

The Etch Ease device helps to prepare a porcelain restoration for bonding or cementation. Etch Ease connects to the suction line and, using suction, allows the porcelain restoration to be held in place during the etching process. Rinsing is a breeze and the product provides a safe way to handle porcelain restorations prior to bonding. Using the Etch Ease, the technician is protected and the restoration is kept free of contaminants.

3/Pkg
6ml Bottle258-0801 1@ \$215.40

**Porcelain Etch Gel
9.6 Percent Hydrofluoric Acid**

For etching porcelain veneers and textured porcelain surfaces. Does not stain composites or thixotropic, washes off easily. Etches porcelain in only 30 seconds.

Kit590-1905 1@ \$30.10
Porcelain Etch Gel Refill
3ml Syringe590-1913 1@ \$30.10



GLAZING FURNACES



Centurion Q50 – Dentsply/Neytech

- 50-program memory
- Quartz sleeve muffle
- Continuous step technology
- Enhanced temperature stability
- Closed muffle drying
- Programmable muffle movement
- Fast Cooling
- Bright, backlit LCD readout
- Icon graphic controls
- Auto calibration, auto purge, auto night mode
- Program override
- Large 10 cm (4") work platform
- Internal surge/spike protection
- Program advance key
- °C to °F conversion
- Full two-year warranty

Each042-4010 1@ \$4,871.95



Programat CS Furnace – Ivoclar Vivadent

Ivoclar Vivadent has developed a ceramic and crystallization furnace for the dentist. The Programat CS is an easy-to-operate, small furnace for the dental office. It is particularly suitable for glaze and crystallization firing of IPS e.max CAD LT blocks. The furnace features 20 programs and is very easy to use. Because of its small dimensions, the furnace occupies very little space in the dental office. Thanks to the vacuum built-in function, more complex firing cycles with high-quality results can be conducted.


- Easy operation
- Automatic temperature calibration
- State-of-the-art design
- 3 pre-programmed standard programs
- 17 individually created programs
- QTK (Quartz Tube Kanthal) heating technology
- Fast, turbo opening system (5 seconds)
- Compact ergonomic design
- 2-year warranty

Each048-8858 1@ \$5,495.00

We offer a wide variety of CEREC® marketing tools. Please contact your Patterson representative for more info.

MARKETING TOOLS


PATTERSON OFFICE SUPPLIES



CEREC® Before-and-After Business Card

- Use professional business cards to remind patients of your CEREC® services
- Personalization: up to five lines
- Printed in black ink, 3.5" W x 2" H


500/Pkg		
Package	052-5097	1@ \$79.50
		2@ \$49.50
		4@ \$40.25



CEREC® Recall Cards

- Perfect for reminding patients of an upcoming appointment, that it is time to make an appointment or to reactivate patients
- Use to inform patient base of your new CEREC®
- Personalization: up to five lines, printed in black ink
- Your choice of free messages suggested
- Size A: 4.25" x 6"
- Size D: 4.25" x 5.5" detached (laser)
- 100/Pkg


Standard - Size A.....	052-5220	1@ \$55.95
		2@ \$25.95
		4@ \$20.95
Custom - Size A	052-5246	1@ \$61.95
		2@ \$31.95
		4@ \$24.95
Laser - Size D.....	052-5873	1@ \$20.95
		2@ \$19.95
		4@ \$18.95



CEREC® Personalized Magnet

- Show patients what CEREC® can do for them with these before-and-after CEREC® magnets
- Personalize with up to five lines, typestyle as shown in black ink, 3.5" x 2"

500/Pkg		
Package	080-9293	1@ \$319.55
		2@ \$279.65
		5@ \$247.65



CEREC® Appointment Cards

- Build your professional image and use for appointment reminders
- Available on several paper stocks, printed in many ink colours, 3.5" x 2"

500/Pkg		
Package	069-1584	1@ \$79.50
		2@ \$49.50
		4@ \$40.25

CEREC® Poster

- Promote your CEREC® capabilities by playing this poster in your practice
- Matching frames available
- Printed on durable paper for long life
- Size: 18" x 24"


Each	060-4371	1@ \$17.95
		2@ \$17.45
		3@ \$16.95



CEREC® Patient Education

"The Choice is Simple"


100/Pkg.....	636-6405	1@ \$48.95
--------------	----------	------------



CEREC® Notepads

- Great for jotting down notes or a quick professional impression
- Imprinted with up to four lines of personalization in black ink
- Size: 5.5" x 8.5"
- White, 20lb bond paper

100 Sheets/Pad;		
5 Pads/Pkg		
Package	052-5097	1@ \$79.50



CD/DVD

This priceless marketing tool will help your marketing efforts with customizable ad templates, PR tips with press releases, radio script and a Flash movie for your Web site.

Each	011-6236	1@ \$24.95
------------	----------	------------