

# Evacuation Products

## SALIVA EJECTORS & ACCESSORIES

### PATTERSON

#### Saliva Ejectors

Disposable, available in two colours, each 6" long with standard diameter. Non-removable tip means it won't come free during procedure. Slotted design protects patient, won't pull or suck tissue. Non-crimping construction holds its shape when bent - no rebound.



<b>100/Pkg</b>		
Clear w/ White Tip ....089-4089	1@ \$4.10	
White w/ White Tip...089-4097	10@ \$3.80	
Clear w/ Blue Tip .....089-4105		

#### Formaject Saliva Ejector - Otis

Disposable  
Formable white plastic tubes with coloured cushion-tips.



<b>100/Package</b>	1@ \$27.90
Pink Regular .....	3@ \$27.25
Blue Hi-Flow.....	

<b>Formaject Tips Only</b>	1@ \$17.70
<b>100/Package</b>	3@ \$16.80

Pink Regular .....	547-7989
Blue Hi-Flow.....	547-7997

<b>500/Package</b>	1@ \$74.55
Pink Regular .....	547-8474
Blue Hi-Flow.....	547-8482

#### Hyoformic Saliva Ejectors - Orsing

Disposable  
Contains:  
100 ejectors  
and 5 adaptors.



<b>Pkg</b> .....	589-7749	1@ \$21.00
		20@ \$19.90

<b>Hyoformic Adapter</b>	
5/Pkg.....	589-8002 1@ \$3.65



#### Denta Pops - Primotec

Saliva Ejectors, Disposable

Mouth prop/tongue holder/suction device in one disposable unit. Patient holds bite block comfortably between teeth and assists in dental aspiration and tongue retraction. Universal left or right.

#### Intro Pack

Contains: 20 disposable dry field aspirators, 42" high vac adapter and instructional CD.  
Each .....569-3007 1@ \$30.10

#### 50/Pack

Contains: 50 disposable dry field aspirators and 42" high vac adapter.

Multicolour .....569-3023 1@ \$63.90

#### 100/Pack

Contains: 100 disposable dry field aspirators and 42" high vac adapter.

Multicolour .....569-3056 1@ \$112.45  
White.....569-3064

#### Flexo® Saliva Ejectors - Bosworth



Personal, hygienic saliva ejectors that will adapt to any mouth and maintain their shape. The specially designed soft, smooth tips provide patient comfort and will not pull tissue. Will not collapse when bent, maintaining a non-restricted flow of suction.

#### 100/Pkg

Clear Removable Tip .224-8375 1@ \$9.45  
White Glued Tip .....224-8409

#### Colour Tips

Blue.....224-8425 1@ \$9.45  
Green.....224-8441  
Red.....224-8458  
Yellow.....224-8433

#### Lingua-Fix - Zirc

Disposable Saliva Ejector

- Simple operation - easy to place and remove
- Chin-lock fixes Lingua-Fix in place
- Reduces high-risk waste
- Comfortable - no sharp edges or corners
- Large suction surface avoids blockage and is more comfortable



10/Pkg.....771-2615 1@ \$13.20

50/Pkg.....771-2607 1@ \$50.60

#### Saliva Ejectors - BDM

Disposable

100/pkg 1@ \$7.95  
Clear .....152-0246 3@ \$7.25

#### Monoject Saliva Ejectors - Covidien

Disposable

Blue tip, 5 1/2" length. 1@ \$35.60  
250/Pkg.....632-8306 12@ \$33.45

#### Saliva Ejectors - Crosstex

Disposable



Optimum suction without aspirating tissue. Soft, smooth non removable tip to assure maximum patient safety. Available in three colour combinations.

#### 100/Pkg

Clear w/ Blue Tip .....235-4371 1@ \$6.15  
Clear w/ White Tip ....235-4389 10@ \$5.80  
White w/ White Tip...235-4629 1@ \$6.70  
Clear w/ Clear Tip .....235-4637 10@ \$5.55



#### Bullfrog Mini Vac Saliva Ejector Assembly - Hager Worldwide

Solid Stainless Steel

- 360° swivel quick-disconnect
- Fits all sizes of saliva tubing
- Easily disassembled for service and cleaning
- Convenient lever for easy handling
- Twice the volume flow of competitive units
- Autoclavable

With Swivel.....68070611 1@ \$88.90

#### Bullfrog Saliva Ejector Tips - Hager Worldwide

Cold-Sterilizable

Black .....159-2674 1@ \$7.60



## HVE TIPS & ACCESSORIES

### PATTERSON

#### HVE Aspirator Tips

Virgin-grade plastic with no impurities – these High Volume Evacuation aspirator tips pass all FDA product requirements and testing. Made in the U.S.A. for high-quality control and consistent fit.

<b>Disposable, 50/Pkg</b>	
Non-Slotted (#D1) .....089-4287	1@ \$3.65 6@ \$3.30
Slotted (#D2).....089-4246	1@ \$5.50 6@ \$5.10

<b>Autoclavable, 1/Pkg</b>	
HVE Autoclavable .....089-4204	1@ \$4.00 6@ \$3.40

#### Combination HVE Aspirator Tips

Multi-use combination tips have both vented and non-vented ends on the same evacuator. "S" tip design with no sharp edges means greater patient comfort. Made in the U.S.A., this virgin-grade plastic is FDA approved with no impurities.

<b>100/Pkg</b>	
Package .....089-4220	1@ \$6.75 3@ \$5.80

#### Aspir-A-Tips

<b>Disposable</b>	
For use with Aspir-A-Tips and Tip-A-Dilly Holders	
200/Box.....083-0802	1@ \$6.70 6@ \$5.70

<b>Holder only</b>	
Each .....083-0810	1@ 41.25

#### Bullfrog Saliva Ejector Adaptor – Hager Worldwide



To HVE Valve  
 • Aluminum (Autoclavable)  
 • Fits all saliva ejectors  
 Each .....68070649 1@ \$16.50

#### Bullfrog Adapter – Hager Worldwide



HVE to Saliva Tubing  
 • Aluminum (Autoclavable)  
 • Fits all sizes of saliva tubing  
 Each .....68070718 1@ \$19.80

#### Bullfrog Evacuation Tip Holder – Hager Worldwide

Autoclavable, Grey  
 1/Pkg.....68070691 1@ \$10.35

#### E-VAC-U Saliva Valve – Bullfrog

Autoclavable  
 Stainless steel, no screen required to do internal designs. 360° swivel quick-disconnect. Fits all sizes of saliva tubing. Convenient lever for easy handling.

Saliva Valve Assembly with Swivel  
 Each .....68070570 1@ \$78.05

Combination Kit  
 Contains: 3-lever saliva valves, 1½" tubing adaptor swivel.  
 Each .....68070631 1@ \$172.80



#### Plastic Saliva Ejector Screens – Zirc

Convenient, inexpensive, disposable. Assists in reducing infectious cross-contamination. Easily inserted and removed.  
 25/Pkg.....771-2482 1@ \$10.00  
 250/Pkg White .....771-9016 1@ \$61.00

#### Rubber Saliva Ejector – Coltène/Whaledent

Tip Only  
 Black .....381-2484 1@ \$45.95

#### Saliva Ejector Assembly Without Tubing – A-dec



Autoclavable  
 Each .....68070566 1@ \$83.40

#### Saliva Ejector Basket Screen – Coltène/Whaledent

Each .....381-2591 1@ \$22.60

#### Saliva Ejector Tips – Hager Worldwide



Mirror Suction  
 Mouth Mirror with Mouth Suction

- Autoclavable
- Acts as tongue guard
- Acts as mouth retractor
- Provides excellent aspiration
- Is a front surface mirror

#### Mirror Suction, "Mini" Kit

Includes: 2 Handles w/Mirrors, 1 Hose, 1 - 11mm Adapter and 1 Cleaning Brush.  
 Each .....777-7402 1@ \$158.90



#### Mirror Suction Full Kit

Includes: 5 Handles w/Mirrors, 1 Hose, 1 - 11mm Adapter and 1 Cleaning Brush.  
 Each .....777-0183 1@ \$291.90



#### Mirror Suction Parts

Mirror Suction Mirror and Handle.  
 Red #254059 .....777-0225 1@ \$66.40  
 Green #254060 .....777-0233  
 Yellow #254061 .....777-0241  
 Blue #254062 .....777-0258  
 Grey #254063 .....777-0266



#### Saliva Ejector Tip Holder – A-dec

Non-autoclavable  
 Grey .....68070297 1@ \$5.90

# Evacuation Products

## HVE TIPS & ACCESSORIES (CONT.)

### Bullfrog E-VAC-U Vacuum Handpiece – Hager Worldwide



*Autoclavable*

Designed with a metal body, a thermoplastic lever for smooth lever action and Viton O-rings. Available in standard and extended configurations, the valve assembly accepts standard 7/16" O.D. tips and incorporates a swivel quick-disconnect: the connector adapts to any 1/2" or 5/8" I.D. central vacuum tubing.

**Handpiece – Standard with Swivel**  
Each .....68070614 1@ \$76.15

**Handpiece – Extended with Swivel**  
Each .....68070628 1@ \$86.80

#### Combination Kit

Contains: 3 Lever HVE valves, 1 hose adaptor swivel, 360° swivel – adapts to 1/2" or 5/8" I.D. tubing. Fits Universal holders, accepts standard vacuum tips.

**Combination Kit – Standard**  
Each .....68070632 1@ \$175.50

**Combination Kit – Extended**  
Each .....68070633 1@ \$180.45

### Bullfrog Ultra-Vac Handpiece (Medium) – Hager Worldwide



Aluminum Body and Swivel  
Each .....68070575 1@ \$88.40

### Bullfrog H.V.E. Handpiece – Hager Worldwide

*Aluminum Body and Swivel*

Designed with an aluminum body, an Ultem lever for smooth lever action and Viton O-rings. The quick-disconnect central vacuum valves can be autoclaved while fully assembled. They can be disassembled without tools for thorough cleaning and lubrication.

- One piece construction
- Clinically unbreakable lever w/smooth action
- 360° swivel quick-disconnect
- Adapts to 1/2" or 5/8" tubing
- Accepts 7/16" O.D. tips
- Fits all new holders
- Autoclavable

**Standard Handpiece without Swivel**  
Each .....68070547 1@ \$62.05

### Bullfrog H.V.E. Handpieces – Hager Worldwide



*Aluminum/Plastic*

Designed with a delrin body and a polysulfone lever, this quick aluminum disconnect central vacuum valve can be ultrasonically cleaned while fully assembled. The valve accepts standard 7/16" tips and the connector adapts to any 1/2" or 5/8" ID central vacuum tubing.

- Ultrasonically cleanable
- 360° swivel quick-disconnect
- Clinically unbreakable smooth lever action
- Easily disassembled for service and cleaning
- Cold sterilizable

**Extended Lever Handpiece with Swivel**  
Each .....68070754 1@ \$82.25

### Bullfrog Replacement Parts – Hager Worldwide

*Adapter Swivel*

**1/2" Tubing**  
#B11109.....68101688 1@ \$18.85



**Lever, 3/Pkg**  
#B11113.....68070556 1@ \$9.25



**O-Ring Standard, 6/Pkg**  
#B11114.....70050190 1@ \$15.10



**O-Ring Autoclavable, 6/Pkg**  
#B11115.....70050188 1@ \$19.55



**Valve Spool**  
#B11112.....68070746 1@ \$11.65

**Adapter With Swivel**  
#B11121.....68070693 1@ \$19.00



**Lever**  
#B11122.....68070554 1@ \$8.90



### Byrd Self-Cleaning Aspirator – Hu-Friedy



B1.....367-7747 1@ \$318.60  
B2.....367-7754  
B3.....367-7762

### #3 Cogswell Suction Tube – Hu-Friedy



Each .....367-7788 1@ \$125.50

### Coupland Aspirator – Hu-Friedy



Accepts 4 interchangeable tips.  
9 1/2" Handle Only.....368-4941 1@ \$179.20

#### Tips Only

#1.....368-5104 1@ \$32.45  
#2.....368-5138  
#3.....368-5161  
#4.....368-5195

#### Tips with Holes

#1.....368-4966 1@ \$40.45  
#2.....368-4974  
#3.....368-4982  
#4.....368-4990

### Disposable Aspirator Tip – A-dec

*Non-autoclavable*

Lightweight, polyethylene material.  
7/16" Outer Diameter x 5" Long  
100/Pkg.....103-1046 1@ \$77.00



### E-Vac Evacuation Tips – E-vac, Inc.



E-Vac Evacuation Tips prevent painful and traumatic tissue plugs and soft-tissue hematomas, while enabling the assistant to retrieve foreign objects with greater ease and safety. Prevents restorative pieces from being lost into the evacuation system. Fits all standard plastic and metal evacuation tubes.

100 Tips/Bag.....320-9939 1@ \$31.35

## Frazier Aspirator Tips – Hu-Friedy

#6 .....	367-7812	1@ \$126.85
#8 .....	367-7820	
#10 .....	367-7838	
#12 .....	367-7846	



## Grab-Free HVE High Volume Evacuation Tips Kit – Clive Craig

- Will not grab or pull tissue
- Curved tip allows easier access to posterior recesses
- Use as a cheek retractor

Contains: 2 nozzles w/100 tips.  
Kit .....228-1509 1@ \$79.25

Grab-Free Nozzle  
Each .....228-1517 1@ \$30.65

Grab-Free Replacement Tips  
100/Pack .....228-1525 1@ \$33.05



## HVE Tips – Heraeus Kulzer

Disposable  
100/Pkg.....449-7632 1@ \$75.95



## Max Vac HVE Tips – BDM

Evacuation tips  
100/pack .....134-4159 1@ \$5.10  
6@ \$4.90



## Maxvac™ Plus – Microbrush

Aspirator Tips  
Contains: 200 each of yellow, pink, purple, blue and green.  
1,000/Case  
Assorted Colours .....383-4447 1@ \$72.10

## Oral Evacuator Tips – Monoject

Disposable  
Soft opaque plastic. Can be used straight or easily shaped to desired angle. Vented tip.  
100/Bag  
Bag.....632-1632 1@ \$28.25  
3@ \$26.95



## Oratip® I – Pascal

For general-purpose removal of saliva and debris with high volume evacuation systems.  
Kit Contains: 1 adapter and 72 tips.  
Each .....555-8820 1@ \$46.20



## Oratip® II – Pascal

Disposable  
Designed especially for use in all oral surgery and endodontic procedures. 1/8" white opaque polyethylene.  
25/Pkg.....555-9315 1@ \$21.00  
50/Pkg.....555-8333 1@ \$42.45  
2@ \$41.30



## Oratip® III – Pascal

Fits All High Volume Evacuation Systems  
One piece construction. The soft tip design has a port near the attachment base for suction release. Sturdy construction allows cheek retraction.  
50/Pkg  
Blue.....555-8093 1@ \$41.55



## Roto-Vac – Young

Rotating Aspirator  
High-volume evacuator system with rotating control of the aspiration angle. Neoprene or plastic disposable sleeves prevent thermal shock and provide a soft, comfortable tip. Fits all high volume evacuation systems.  
Each .....473-3887 1@ \$120.40  
Neoprene Sleeves, Regular  
6/Pkg .....472-9380 1@ \$23.10  
Neoprene Sleeves, Vented  
6/Pkg .....472-9877 1@ \$24.50

## Vent-O-Vac Tips – Young

Disposable  
Side vented HVE tips.  
100/Pkg  
Reg. 5/8" .....473-8787 1@ \$43.70  
Short 3/8" .....473-8290 1@ \$33.10




## HVE Valve Brush – A-dec

6/Pkg  
3/4" Brush .....103-0147 1@ \$21.70  
1/2" Brush .....102-9081 1@ \$26.45



## SURGICAL TIPS & ACCESSORIES



## PATTERSON

### Surgical Aspirator Tips

These autoclavable, plastic surgical aspirator tips provide a secure fit in all evacuation systems. They're rigid enough to aspirate and retract cheek without bending, and the 30° permanent angle allows better visual contact of the site. Opaque colour masks fluids. Virgin-grade plastic has no impurities – we guarantee the highest quality.

25/Pkg  
Regular 1/8" Tip  
White.....089-4303 1@ \$9.45  
3@ \$8.25  
Large 1/4" Tip  
Green.....089-4329 1@ \$9.30  
3@ \$8.15

## Roeko Surgitip Micro – Coltène/Whaledent

Designed for use with microscopes and loupes. Ideal for use in endodontics. Extra fine tip ensures precise aspiration without blocking the field of vision.  
20/Pkg  
1.2mm .....766-0178 1@ \$55.65



# Evacuation Products

## SURGICAL TIPS & ACCESSORIES (CONT.)



### Surgical Aspirator Tip – A-dec Stainless Steel

Adapts to A.V.S. and Central Vacuum HVE valve, autoclavable.

3/16" Diameter ..... 103-2044 1@ \$51.45  
1/8" Diameter ..... 103-2069



### Surgical Aspirator Tips – Monoject

Disposable

Smooth opaque plastic. Individually packaged, sterile.

25/Box ..... 632-1145 1@ \$45.65  
3@ \$42.80



### Surgical Tip Brush – A-dec

3/16" Brush  
6/Pkg ..... 102-9586 1@ \$33.60

1/8" Brush  
6/Pkg ..... 102-9602 1@ \$46.25



### Surg-O-Vac Aspirator Tips – Young

Autoclavable – .125" 3mm Diameter Tip  
25/Pkg ..... 473-4869 1@ \$34.75

### Surg-O-Vac II Aspirator Tips – Young

Autoclavable – .25" 6mm Diameter Tip  
25/Pkg ..... 473-4935 1@ \$34.75

### Surg-O-Vac III Aspirator Tips – Young

Autoclavable  
0.375" (9mm) tip opening, blue.  
25/Pkg ..... 473-4976 1@ \$34.75  
0.375" (9mm) tip opening, vented, blue.  
25/Pkg ..... 473-4984 1@ \$34.75

## EVACUATOR TRAPS

### Dispos-A-Bowl – Apollo Dental Products, Inc.

Vacuum Filter



- Stainless steel screen is rolled at each end to prevent accidental puncture of protective gloves and skin
- Waste canister and lid protect against spills and leakage and are premarked with hazardous warning stickers
- Designed to ensure minimized exposure to infectious waste material and eliminate cross-contamination of the workplace

Contains waste canister and lid, filter bowl, filter screen, replacement rubber gasket and hazardous warning sticker

Each  
3/4" ..... 54020190 1@ \$23.50  
1" ..... 54060588 1@ \$29.60  
1-1/4" ..... 54020394 1@ \$37.00  
8/Pkg  
1-1/2" ..... 54060689 1@ \$242.90

### Dispose-A-Screens – Zirc

Colour	Size	Fits
Blue	21/8"	A-dec Model with plastic canister, Biotec and others
Yellow	21/4"	A-dec Model with metal canister
Green	17/8"	Forest, Belmont, Chayes, Den-tal-ez, Dentech, Pelton, Prisma and others
Coral	25/8"	Pelton and Crane Spirit unit
Beige	27/8"	Den-tal-ez with black filters
Grey	23/4"	Midwest 210 models

- Avoid contact with potentially hazardous material
- Reduce possibilities of infection and contamination
- Improve performance of suction system

25/Pkg 1@ \$17.05  
Blue ..... 771-1815 Yellow ..... 771-1823  
Green ..... 771-1831 Coral ..... 771-1849  
Beige ..... 771-1856 Grey ..... 771-1864

144/Pkg 1@ \$64.00  
Blue ..... 771-2669 Yellow ..... 771-2177  
Green ..... 771-2185 Coral ..... 771-2193  
Beige ..... 771-2201 Grey ..... 771-2219  
Maue ..... 771-2227



### Dispos-A-Trap – Pinnacle Disposable Vacuum System Traps

Model 5500  
Fits cascade/perform models. Diameter: 1 1/4"  
144/Pkg ..... 568-3487 1@ \$81.00

Model 5501  
Fits A-dec models with plastic canisters.  
Diameter: 2 1/2"  
144/Pkg ..... 568-3503 1@ \$81.00

Model 5502  
Fits A-dec models with metal canisters.  
Diameter: 2 1/4"  
144/Pkg ..... 568-3511 1@ \$81.00

Model 5504  
Fits Den-Tal-Ez models with wire screen filters.  
Diameter: 1"  
144/Pkg ..... 568-3537 1@ \$84.40

Model 5503  
Fits models by Belmont, Biotec, Chayes, Virginia, Dentech, Pelton & Crane, Prisma & other manufacturers. Diameter: 1 1/8"  
144/Pkg ..... 568-3529 1@ \$81.00

Model 5505  
Fits Pelton & Crane Spirit with white plastic canisters. Diameter: 3 3/4"  
144/Pkg ..... 568-3586 1@ \$84.40

Model 5506  
Fits Den-Tal-Ez models with black plastic filters. Diameter: 2 15/16"  
144/Pkg ..... 568-3545 1@ \$84.40

Model 5507  
Fits Midwest 210 models. Diameter: 2 1/4"  
144/Pkg ..... 568-3552 1@ \$84.40

Model 5508-S  
Fits all saliva ejectors.  
72/Pkg ..... 568-3560 1@ \$47.10

Model 5509  
Fits DCI/Belmont Reality Models, and System 6000 models.  
144/Pkg ..... 568-3594 1@ \$84.40

Model 5510  
Fits Sirona C8.  
144/Pkg ..... 568-4105 1@ \$84.40

Model 6100-C  
Fits most cuspidor bowls.  
144/Pkg ..... 568-3578 1@ \$86.00

Model 6200-C  
Fits European units, older style U.S. units with openings 1/2" and less.  
144/Pkg ..... 568-3461 1@ \$85.75

Model 6300-C  
Fits Planmeca units.  
144/Pkg ..... 568-3479 1@ \$84.40

Model 6400  
Fits A-dec Cascade and Performer cuspid bowls.  
144/Pkg ..... 568-3495 1@ \$85.75

**Evac-U-Trap – Pinnacle**

Disposable  
Central Suction  
System Traps



**Model 2200**

Fits Air Techniques (Vacstar models), Den-Tal-Ez (1 & 2 hp single models), Ohmeda/Matrx (all models), Turbine Industries (Whirl Vac models), and other manufacturers.  
12/Box.....568-3602 1@ \$141.70

**Model 2300**

Fits Den-Tal-Ez (dual models made after 1978).  
8/Box.....568-3628 1@ \$141.70

**Model 2350-FS**

Fits Apollo 2 hp, 3hp single and tandem models, Air Techniques-1997 and later dual models, Vacudent-Workhorse dual, new finer mesh screen.  
8/Box.....568-3636 1@ \$141.70

**Model 2400**

Fits Dentsply (MVS .75, 1.5, 2.0, 4.0 models made after 1978), Den-Tal-Ez (dual models made before 1978).  
8/Box.....568-3644 1@ \$141.70

**Model 2500**

Fits Dentsply-MVS Jr. models.  
8/Box.....568-3651 1@ \$141.70

**Model 2600**

Fits Ritter Vacuum Air models.  
8/Box.....568-3669 1@ \$141.70

**ORAL CUPS / LINERS**

**Molded Oral Cup – A-dec**

Plastic

Adapts to central vacuum HVE valve or A.V.S. handpiece. 4" diameter cup with snap-on disposable liners.  
Each .....6807008 1@ \$38.45



**Oral Cup Liners 6SR – Tidi**

Cut rolled rims, 6oz paper to fit A-dec units.  
250/Pkg.....699-3687 1@ \$23.10  
3@ \$21.70



**EVACUATION CLEANERS**

**PATTERSON**

**Enzymatic Tablets**

Enzyme cleaning tablets for ultrasonic cleaner, evacuation system and pre-soak. Two tablets make one gallon.



64/Box.....090-9796 1@ \$54.35  
3@ \$52.60

**Evacuation Cleaner Powder**

HVE powder sanitizes and deodorizes. Small particle size gives you better scrubbing action, and the easy-mix concentrate goes a long way.



2lb Canister (Makes 32 Gallons)....084-9091 1@ \$11.30  
3@ \$10.35

**Super Concentrate Evacuation Cleaner/Deodorizer**

This dental vacuum line cleaner dissolves all debris in the solid waste trap, evacuation tubing and vacuum pump. Add 1/2oz of solution to 1 quart of water.  
Quart .....090-9804 1@ \$29.80  
4@ \$28.80



**AlproCleaner – Alpro Medical**

Cleaning of aspiration systems

Alkaline liquid ready for use for the cleaning of spittoon bowls in all dental treatment units supplied with and without amalgam separators. Removes easily the residues of saliva, blood, etc.



500ml.....636-6652 1@ \$11.90  
5L Cubitainer .....636-6660 1@ \$89.50

**AlproJet D – Alpro Medical**

Cleaning of aspiration systems



Highly efficient liquid concentrate for the daily cleaning and maintenance of dental aspiration and separation systems. Deodorizes and dissolves the residues which have been aspirated, such as saliva, drilling dust, etc. FREE of phenols, aldehydes, phosphates and chlorine

1L Dosage bottle .....636-6520 1@ \$16.65  
3@ \$15.85  
10L Can.....636-6538 1@ \$85.75  
3@ \$83.00

**AlproJet W – Alpro Medical**

Cleaning of aspiration systems



Highly efficient liquid concentrate for weekly contra-cleaning and maintenance of dental aspiration and separation systems. Dissolves chalk (acid pH value), neutralizes hydrogen peroxide H<sub>2</sub>O<sub>2</sub> and sodium bicarbonate (from airflow). FREE of phenols, aldehydes, phosphates and chlorine

1L Dosage bottle .....636-6579 1@ \$16.65  
3@ \$15.85  
5L Cubitainer .....636-6587 1@ \$50.80  
3@ \$48.40



**Auto Walkabout™ – Sultan Healthcare**

Automatically cleans system without hassle, mess or work. Air and solution mixture creates a turbulence against tubing walls to detach and dissolve debris.

Each .....665-8678 1@ \$90.00

# Evacuation Products

## EVACUATION CLEANERS (CONT.)



### Citrzyme™ – Pascal

*Ultra Concentrated Enzyme Cleaner and Instrument Pre-Soak*

- Premium non-foaming formula
- Fast acting dual enzymes remove protein and biofilms quickly and leave no residue
- Powdered concentrate dissolves easily in warm tap water and has fresh citrus scent and colour
- Each 300g canister makes 35 gallons
- Dissolves instantly. No mixing required!
- Safe for environment, staff and equipment

300g/Canister .....555-2666 1@ \$49.70



### CleanStream™ – Air Techniques

*Evacuation System Cleaner*

- Fresh scent
- Cleans proteinaceous and synthetic material from all surfaces inside the evacuation system
- Non-foaming, biodegradable
- Convenient and easy-to-use dispenser system

#### Starter Kit

Contains 32 pouches of CleanStream Evacuation System Cleaner (sufficient for 64 applications) and a CleanStream Dispenser System with accessories

Each .....104-8248 1@ \$74.55  
2 Liter Bottle .....104-8263 1@ \$51.50  
3@ \$69.30



### Enzymax™ – Hu-Friedy

As a dual-enzyme formula, Enzymax is maximum cleaning and care for instruments and IMS cassettes and accessories. It can be used as an instrument and linen pre-soak, ultrasonic detergent and evacuation cleaner. Enzymax has a neutral pH. Does not contain caustic or hazardous chemicals, is phosphate-free and biodegradable, for hassle-free, earth friendly disposal.

#### Liquid

Quart .....371-0381 1@ \$53.95  
Gallon .....371-1710 1@ \$105.25  
**Pump Only for Quart**  
Bottle .....371-2643 1@ \$5.50

#### Powder

800gm .....371-2650 1@ \$66.10  
1/3oz Liquid Packets  
40/Box .....371-0373 1@ \$109.35

#### Enzymax PAX

32/Pkg .....371-2668 1@ \$70.20  
96/Pkg .....371-2676 1@ \$169.45

### E-VAC – L&R Manufacturing

Premium quality, nonfoaming, highly concentrated evacuation cleaner and sanitizer. 1oz provides 1qt of cleaning solution.

Gallon .....474-9594 1@ \$50.05



### FRESH-Vac – ECOLAB

*Daily Use Evacuation System Cleaner*

A daily evacuation system cleaner, containing powerful concentrated enzymes to dissolve organic material from lines and systems. Pleasant fresh scent, which eliminates unpleasant odours.

32oz .....376-4008 1@ \$30.95  
4@ \$29.95



### MERCONgel™ – Ross Dental

*Amalgam Trap & Evacuation Sanitizer*

Sustained Release Decontaminant/Cleanser eliminates and prevents the further build-up of mercury vapour in waste traps and evacuation systems.

Patented solution coats and protects lines from build-up of mercury-rich waste. Also contains Chlorane-T for added bacterial and viral control.

4L /Bottle .....317-3895 1@ \$144.95



### Orotol® Plus – Durr/SciCan

*Evacuation System Cleaner*

A mechanical and chemical vacuum system cleaner. Orotol Plus contains two chemical cleaners.

#### Orotol Plus Solution Only

2.5L Bottle .....638-7393 1@ \$160.00



### Purevac – Sultan

*Daily Evacuation System Cleaner*

Dissolves solid debris. Single application dissolves solids and keeps suction lines clean and running at peak efficiency. Dissolves calcium build-up in vacuum pumps. Non-foaming, lemon-scented and biodegradable.

2 Liter .....665-8652 1@ \$63.95  
6@ \$62.00

5 Liter .....665-8660 1@ \$141.00  
4@ \$138.95

#### Purevac Walkabout Unit

Each .....665-8686 1@ \$92.00





## PUREVAC® HG – Sultan Healthcare

### Dual Enzymatic Cleaner

- Patented formulation coagulates small particles and binds dissolved mercury, improving mercury capture in amalgam separators
- Fresh citrus lime scent
- Biodegradable formula
- Non-foaming

<b>Installation Kit</b>		
Each .....	665-9023	1@ \$100.00
<b>Quart Bottle</b>		
Each .....	665-9007	1@ \$48.75

## Sani-Treet Plus® – Enzyme Industries, Inc.

Natural enzymatic evacuation system cleaner specially formulated for use in ultrasonic cleaners and as an instrument presoak. Non-foaming, nontoxic, biodegradable. Makes 128 treatments of evacuation system cleaner solution or 512qt of instrument presoak or ultrasonic cleaning solution.



<b>Gallon</b>		
Peppermint Scent .....	557-6905	1@ \$83.95
Country Meadow Scent .....	299-0000	1@ \$72.80



## Turbo-Vac – KERR Total Care

### Daily Evacuation System Cleaner

Liquid concentrate dissolves debris, sanitizes and deodorizes. Is biodegradable and non-foaming, with a fresh scent. Makes 128 quarts. Contains: one gallon with dispensing pump. Model 5050 .....

568-4006 1@ \$67.50



## Turbo-Vac Dispenser – Kerr TotalCare

### Dispensing Unit for Turbo-Vac

Unique design mixes air and Turbo-Vac solution to create a powerful, non-foaming, turbulent mist that dislodges debris from the walls of dental tubings, greatly improving vacuum performance. This self-contained gallon container allows you to move freely between operatories and store any unused solution for next day use. The dispenser can be left unattended while in use.

<b>Dispensing Unit Only</b>		
Model 5000 .....	568-4014	1@ \$48.75

## Vac Attak™ – Premier

### Evacuation System Cleaner

- Unique organic bleach component is significantly more effective in removing debris and bioburden from evacuation lines
- Safe and environmentally friendly – biodegradable and noncorrosive to metals
- Nonfoaming formula is safe for use with all evacuation pumps
- Easy to use, fresh lemon scent



800gm Jar.....	585-8592	1@ \$73.50
----------------	----------	------------



## VacuKleen E2 – Heraeus Kulzer

Provides deep cleaning of blood, saliva and tissue. Convenient unit-dose, non-liquid package. Fresh, lemon scent. Contains alkaline protease enzymes, nonionic so it is easy on evacuation equipment. Biodegradable with neutral pH.

<b>Unit-Dose Pouch for Quarts</b>		
60/Pkg.....	448-9506	1@ \$87.95
<b>Unit-Dose Pouch for Quarts</b>		
140/Pkg.....	448-9514	1@ \$182.00
<b>Unit-Dose Pouch for Gallons</b>		
25/Pkg.....	448-9522	1@ \$82.45



## Vacusol Ultra – Biotrol

### Dental Vacuum Line Cleaner

Unique chelating formula that removes both synthetic and organic debris from the dental vacuum line tubing. Vacusol Ultra™ is non-foaming, anti-corrosive and non-acidic. 32oz bottle.....

123-8989	1@ \$52.50
	4@ \$51.15

## Cuspidor Screens – Zirc

### White Cuspidor

Fits: A-dec Cascade and Performer cuspidor bowls.

25/Pkg.....	771-0916	1@ \$17.05
144/Pkg.....	771-0957	1@ \$66.00

### Grey Cuspidor

Fits: A-dec, Ampxo, Belmont, Dabi Atlante, Dansereau, DCI, DentalEZ, Forest, Knight, Marus, MDT, Proma Ritter and Schein.

25/Pkg .....	771-0924	1@ \$17.05
144/Pkg .....	771-0965	1@ \$66.00

