

Implantology

01 PLACEMENT INSTRUMENT KITS

Omni-Tight™ Placement Instrument Kit	271
Midi™ Implant Placement Instrument Kit	271
Tool Kit Tray	271

02 OMNI-TIGHT™ IMPLANT SERIES

Fixed Restorative Implant Kits	272
Denture Retention Implant Kits	273
Add-on Packages	274

03 MIDI™ IMPLANT SERIES

Denture Retention	275
Fixed Restorative	275

04 ABUTMENT FIXED RESTORATIVES

Abutment Fixed Restoratives	275
-----------------------------	-----

05 INDIVIDUAL COMPONENTS

Individual Components	276
-----------------------	-----

06 INDIVIDUAL SURGICAL INSTRUMENTS

Individual Surgical Instruments	278
---------------------------------	-----

07 ANCILLARY SURGICAL PRODUCTS

Ancillary Surgical Products	280
-----------------------------	-----

08 LABORATORY PRODUCTS

Laboratory Products	280
---------------------	-----

09 MARKETING PRODUCTS

Marketing Products	281
--------------------	-----

10 EDUCATION

Education	282
-----------	-----

11 IMPLANTS MAINTENANCE INSTRUMENTS

Maintenance Instruments	282
-------------------------	-----

12 MISCELLANEOUS

Miscellaneous	283
---------------	-----

PLACEMENT INSTRUMENT KITS



Omni-Tight™ Placement Instrument Kit

Basic Dental Implants' Placement Instruments are logically arranged in the sequence used during placement. The Instrument Kit contains all instruments necessary to place and restore Fixed, Denture Retained, and Midi Implants. Placement instruments for Threaded Abutments and Bolt-Retained Abutments are available separately and are not included in this kit.

PRODUCT	PART #	DESCRIPTION	PRICE
	PLCMNTKIT <i>(Also contains Midi instruments)</i>	The Omni-Tight Instrument Kit contains: 2 Depth Gauges (13mm, 18mm) 1 Pilot Hole Soft Tissue Trephine 1 Pilot Hole Drill 4 Soft Tissue Cutters (3.5mm Midi, 4.5mm, 5.0mm, 5.5mm) 4 Osteotomy Drills (2.3mm Midi, 3.5mm, 4.0mm, 4.5mm, 5.7mm) 2 Healing Cap Tools (13mm, 18mm) 2 Omni-Tight Implant Driver/Carriers (13mm, 20mm) 1 Midi Abutment Driver 1 Midi Ball Driver 1 Ratchet Wrench 1 Drill Extender 1 Hard Bone Drill 1 Torque Wrench 1 Guide Pin (for Stent Fabrication) Coloured o-rings 1 Bone Mapping Instrument	(127-9215) \$2,100.00

Midi™ Implant Placement Instrument Kit

The Midi Instrument Kit contains the instruments necessary to place and restore BASIC Midi Dental Implants i.e. one-piece Midi Ball Implants or Midi Abutment Implants. It is recommended that these implants be placed laying a flap, but can be done flapless by experienced clinicians. In utilizing the Midi Abutment, the flapless approach is best accomplished with the three dimensional drilling guide stent as outlined for the standard Omni-Tight™ series implants.

PRODUCT	PART #	DESCRIPTION	PRICE
	MIDIPLCMNT	The Midi Instrument Kit contains: 1 Pilot Hole Soft Tissue Trephine 1 Pilot Hole Drill 1 Soft Tissue Cutter – 3.5mm 1 Midi Osteotomy Drill – 2.3mm 1 Midi Abutment Driver 1 Midi Ball Driver 1 Ratchet Wrench	(127-8241) \$964.65
	TRAY	Tool Kit Tray	(127-8274) \$220.80

OMNI-TIGHT™ IMPLANT SERIES

Fixed Restorative Implant Kits

BASIC Dental Implants' Fixed Restorative Implant Kits make ordering simple. Just specify the width and length of the desired implant.

Fixed Restorative Implant Kits come with every component necessary to restore most single or multi-unit cemented restorations. It is not necessary to select individual components for most of your cases.

Each Fixed Restorative Implant Kit includes:



PRODUCT	PART #	DESCRIPTION	PRICE
	35IMP11	3.5mm X 11mm Implant Kit: Fixed	(127-9017)
	35IMP13	3.5mm X 13mm Implant Kit: Fixed	(127-9025)
	35IMP15	3.5mm X 15mm Implant Kit: Fixed	(127-9033)
	40IMP11	4.0mm X 11mm Implant Kit: Fixed	(127-9041)
	40IMP13	4.0mm X 13mm Implant Kit: Fixed	(127-9058)
	40IMP15	4.0mm X 15mm Implant Kit: Fixed	(127-9066)
	45IMP11	4.5mm X 11mm Implant Kit: Fixed	(127-9074)
	45IMP13	4.5mm X 13mm Implant Kit: Fixed	(127-9082)
	45IMP15	4.5mm X 15mm Implant Kit: Fixed	(127-9090)
	50IMP9	5.0mm X 9mm Implant Kit: Fixed *	(128-0684)
	50IMP11	5.0mm X 11mm Implant Kit: Fixed *	(127-9108)
	50IMP13	5.0mm X 13mm Implant Kit: Fixed *	(127-9116)
	50IMP15	5.0mm X 15mm Implant Kit: Fixed *	(127-9124)
	60IMP8	6.0mm X 8mm Implant Kit: Fixed *	(128-0692)
	60IMP11	6.0mm X 11mm Implant Kit: Fixed *	(128-0700)
	60IMP13	6.0mm X 13mm Implant Kit: Fixed *	(128-0718)
60IMP15	6.0mm X 15mm Implant Kit: Fixed *	(128-0726)	

* Note: 5.0 mm and 6.0 mm diameter implants work best in soft bone.

Individual components are also available separately.

See Add-on Packages section for other restorative options including Bolt-retained restorations and Threaded abutments.

OMNI-TIGHT™ IMPLANT SERIES (CONT).

Denture Retention Implant Kits

BASIC Dental Implants' Denture Retention Implant Kits make ordering simple. Just specify the width and length of the desired implant.

Our Denture Retention Implant Kits come with every component necessary to allow the laboratory and clinician to fabricate an o-ring retained prosthesis over a standard implant. It is not necessary to select individual components for most of your cases.

Each Fixed Restorative Implant Kit includes:



PRODUCT	PART #	DESCRIPTION	PRICE
	35DIMP11	3.5mm X 11mm Implant Kit: Denture Retention	(127-9942) \$525.00
	35DIMP13	3.5mm X 13mm Implant Kit: Denture Retention	(127-9959)
	35DIMP15	3.5mm X 15mm Implant Kit: Denture Retention	(127-9967)
	40DIMP11	4.0mm X 11mm Implant Kit: Denture Retention	(127-9975)
	40DIMP13	4.0mm X 13mm Implant Kit: Denture Retention	(127-9983)
	40DIMP15	4.0mm X 15mm Implant Kit: Denture Retention	(127-9991)
	45DIMP11	4.5mm X 11mm Implant Kit: Denture Retention	(128-0007)
	45DIMP13	4.5mm X 13mm Implant Kit: Denture Retention	(128-0015)
	45DIMP15	4.5mm X 15mm Implant Kit: Denture Retention	(128-0023)
	50DIMP9	5.0mm X 9mm Implant Kit: Denture Retention *	(128-0734)
	50DIMP11	5.0mm X 11mm Implant Kit: Denture Retention *	(128-0031)
	50DIMP13	5.0mm X 13mm Implant Kit: Denture Retention *	(128-0049)
	50DIMP15	5.0mm X 15mm Implant Kit: Denture Retention *	(128-0056)
	60DIMP8	6.0mm X 8mm Implant Kit: Denture Retention *	(128-0742)
60DIMP11	6.0mm X 11mm Implant Kit: Denture Retention *	(128-0759)	
60DIMP13	6.0mm X 13mm Implant Kit: Denture Retention *	(128-0767)	
60DIMP15	6.0mm X 15mm Implant Kit: Denture Retention *	(128-0775)	

The following items are custom-packaged and require a minimum lead time of up to four weeks



30IMP11AW	3.0mm X 11mm Midi Implant, Analog & Waxing Sleeve	(127-9314) \$161.00
30IMP13AW	3.0mm X 13mm Midi Implant, Analog & Waxing Sleeve	(127-9322)
30IMP15AW	3.0mm X 15mm Midi Implant, Analog & Waxing Sleeve	(127-9330)

* Note: 5.0 mm and 6.0 mm diameter implants work best in soft bone.

Individual components are also available separately.

See individual components for other restorative options including Bolt-retained bar and Overdenture restorations.

OMNI-TIGHT™ IMPLANT SERIES (CONT).

The following items are custom-packaged and require a minimum lead time of up to four weeks.



30MAIMP11WSP	3.0mm X 11mm Midi Abutment Implant & Waxing Sleeve Profiler	(000-0000)
30MAIMP11WSP	3.0mm X 13mm Midi Abutment Implant & Waxing Sleeve Profiler	(000-0000)
30MAIMP11WSP	3.0mm X 15mm Midi Abutment Implant & Waxing Sleeve Profiler	(000-0000)

Price on request

Add-on Packages

Add-on Packages expand the Omni-Tight™ Implant Series by allowing you to customize your restorative options. Add-on Packages are ordered **in addition** to the Implant Package. Some components are available separately.

Implant Packages – Denture Retention

Our Denture Retention Implant Packages allow the laboratory and clinician to fabricate an O-ring retained prosthesis over a standard implant. The Denture Retention Package includes the following components: Ball Post, Delrin Ball Core, Titanium Ball Core, O-ring Cap, and O-rings. Denture Retention Kits are available as complete implant kits. (See Implant Kits Section - Pg. 3)

PRODUCT	PART #	DESCRIPTION	PRICE
	35DSK	3.5mm Denture Retention Kit (No Implant)	(127-9397)
	40DSK	4.0mm Denture Retention Kit (No Implant)	(127-9405)
	45DSK	4.5mm Denture Retention Kit (No Implant)	(127-9413)
	50DSK	5.0mm Denture Retention Kit (No Implant)	(127-9421)
	60DSK	6.0mm Denture Retention Kit (No Implant)	(127-9439)

Bolt Retained Post and Core Package

For dentists interested in the various overdenture designs, BASIC's Bolt Retained Post and Core Package is the answer. The larger diameter bolt retains a post and the restoration can be customized with a waxing sleeve. The Bolt Retained Post and Core Package includes: Bolt, Bolt Retained Post, and Waxing Sleeve.


PRODUCT	PART #	DESCRIPTION	PRICE
	35SP&C2	3.5mm Straight Post & Core Kit (No Implant)	(127-9470)
	40SP&C2	4.0mm Straight Post & Core Kit (No Implant)	(127-9488)
	45SP&C2	4.5mm Straight Post & Core Kit (No Implant)	(127-9496)
	50SP&C2	5.0mm Straight Post & Core Kit (No Implant)	(127-9504)
	60SP&C2	6.0mm Straight Post & Core Kit (No Implant)	(128-0817)

MIDI™ IMPLANT SERIES


Our Midi™ Implants were designed specifically for narrow bone ridges. Unlike our Omni-Tight™ Implants, our Midi™ Implants are a one-piece design which comes in two different types. Our Midi™ Ball Implant incorporates an O-ring retaining ball machined onto the top for denture stabilization, and our Midi™ Abutment Implant has an abutment formed on its top for fixed restorations of small teeth.

This minimizes the components and steps necessary. A 1.8mm pilot hole and a 2.3mm sizing Osteotomy Drill are used to drill a channel 1/2 to 2/3 of the length of the implant. The aggressive threads on the implant are 3mm in diameter providing good immediate fixation in the bone. The platform of the implant flares to 3.5mm.


Denture Retention

PRODUCT	PART #	DESCRIPTION	PRICE
	30DIMP11	3.0mm X 11mm Midi Ball Implant, O-ring cap and O-rings (127-9165)	\$128.25
	30DIMP13	3.0mm X 13mm Midi Ball Implant, O-ring cap and O-rings (127-9173)	
	30DIMP15	3.0mm X 15mm Midi Ball Implant, O-ring cap and O-rings (127-9181)	

Fixed Restorative

PRODUCT	PART #	DESCRIPTION	PRICE
	30ABUTIMP11	3.0mm X 11mm Midi Abutment Implant, Analog, Waxing Sleeve-Profilers (2 each), Guide Cylinders (128-0866)	\$149.45
	30ABUTIMP13	3.0mm X 13mm Midi Abutment Implant, Analog, Waxing Sleeve-Profilers (2 each), Guide Cylinders (128-0874)	
	30ABUTIMP15	3.0mm X 15mm Midi Abutment Implant, Analog, Waxing Sleeve-Profilers (2 each), Guide Cylinders (128-0882)	

ABUTMENT FIXED RESTORATIVES

PRODUCT	PART #	DESCRIPTION	PRICE
	38ABUTIMP11	3.8mm X 11mm Abutment Implant, Analog, Waxing Sleeve-Profilers (2 each), Guide Cylinders (128-0890)	\$149.45
	38ABUTIMP13	3.8mm X 13mm Abutment Implant, Analog, Waxing Sleeve-Profilers (2 each), Guide Cylinders (128-0908)	
	38ABUTIMP15	3.8mm X 15mm Abutment Implant, Analog, Waxing Sleeve-Profilers (2 each), Guide Cylinders (128-0916)	
	43ABUTIMP11	4.3mm X 11mm Abutment Implant, Analog, Waxing Sleeve-Profilers (2 each), Guide Cylinders (127-8316)	\$170.45
	43ABUTIMP13	4.3mm X 13mm Abutment Implant, Analog, Waxing Sleeve-Profilers (2 each), Guide Cylinders (127-8324)	
	43ABUTIMP15	4.3mm X 15mm Abutment Implant, Analog, Waxing Sleeve-Profilers (2 each), Guide Cylinders (127-8332)	

INDIVIDUAL COMPONENTS

PRODUCT	PART #	DESCRIPTION	PRICE
	CAP35	3.5mm Healing Cap	(128-0098) \$67.20
	CAP40+	Multi-size Healing Cap (Fits 4.0, 4.5, 5.0 and 6.0mm Diameter Omni-Tight Implants)	(128-0148)
	TRNFR35	3.5mm Transfer Device	(128-0171) \$28.00
	TRNFR40+	Multi-size Transfer Device (Fits 4.0, 4.5, 5.0 and 6.0mm Diameter Omni-Tight Implants)	(128-0213)
	30BALLEMPRO	Midi Ball Emergence Profiler	(128-0247) \$33.60
	35EMPRO	3.5mm Diameter Emergence Profiler	(128-0254)
	40EMPRO	4.0mm Diameter Emergence Profiler	(128-0262)
	45EMPRO	4.5mm Diameter Emergence Profiler	(128-0270)
	50EMPRO	5.0mm Diameter Emergence Profiler	(128-0288)
	60EMPRO	6.0mm Diameter Emergence Profiler	(128-0296)
	ANALG30BALL	3.0mm Midi Ball Lab Analog	(128-0379) \$46.20
	ANALG30ABUT	3.0mm Midi Abutment Lab Analog	(128-0957) \$46.20
	ANALG38ABUT	3.8mm Abutment Lab Analog	(128-0965) \$46.20
	ANALG43ABUT	4.3mm Abutment Lab Analog	(128-0973) \$46.20
	ANALG35	3.5mm Diameter Lab Analog	(128-0320) \$46.20
	ANALG40	4.0mm Diameter Lab Analog	(128-0338)
	ANALG45	4.5mm Diameter Lab Analog	(128-0346)
	ANALG50	5.0mm Diameter Lab Analog	(128-0353)
	ANALG60	6.0mm Diameter Lab Analog	(128-0361)
	PC35	3.5mm Cemented Post & Core Abutment	(128-0395) \$89.60
	PC40+	Multi-size Cemented Post & Core Abutment (Fits 4.0, 4.5, 5.0 and 6.0mm Diameter Omni-Tight Implants)	(128-0403)
	25PC35	3.5mm Cemented 25° Post & Core Abutment	(128-1005) \$142.10
	25PC40+	Multi-size Cemented 25° Post & Core Abutment (Fits 4.0, 4.5, 5.0 and 6.0mm Diameter Omni-Tight Implants)	(128-1013)
	BPOST35	3.5mm Ball Post Abutment	(127-8019) \$86.30
	BPOST40+	Multi-size Ball Post Abutment (Fits 4.0, 4.5, 5.0 and 6.0mm Diameter Omni-Tight Implants)	(127-8027)
	THRDBPOST35	3.5mm Threaded Ball Post Abutment	(127-8415) \$93.45
	THRDBPOST40+	Multi-size Threaded Ball Post Abutment (Fits 4.0, 4.5, 5.0 and 6.0mm Diameter Omni-Tight Implants)	(127-8423)







INDIVIDUAL COMPONENTS (CONT.)

PRODUCT	PART #	DESCRIPTION	PRICE
	CCORE	Castable Ball Core (Acrylic) (Fits 3.5, 4.0, 4.5, 5.0 and 6.0mm Diameter Omni-Tight Implants)	(127-8084) \$24.55
	WCORE	Weldable Ball Core (Titanium) (Fits 3.5, 4.0, 4.5, 5.0 and 6.0mm Diameter Omni-Tight Implants)	(127-8092) \$30.15
	GUIDE23	3.0mm Midi Guide Cylinders – Pilot and Osteotomy	(127-8118) \$21.55
	GUIDE35	3.5mm & 3.8mm Guide Cylinders – Pilot and Osteotomy	(127-8126)
	GUIDE40	4.0mm Guide Cylinders – Pilot and Osteotomy	(127-8134)
	GUIDE45	4.5mm/5.0mm Guide Cylinders – Pilot and Osteotomy (Fits 4.5 and 5.0 Diameter Omni-Tight Implants)	(127-8142)
	GUIDE57	6.0mm Guide Cylinders – Pilot and Osteotomy	(127-8159)
	ORINGCP	O-ring Cap with O-Rings	(127-9454) \$46.70
	ORINGS	O-rings Only: 1 Regular (Black), 1 Try-in (Red), 2 Blockout	(127-9512) \$18.65
	ORINGSR	O-rings Only: 5 Regular (Black), and 5 Try-in (Red)	(127-9520) \$32.70
	35SLEEVE	3.5mm Waxing Sleeve for Cemented Post & Core Abutment	(128-0619) \$22.40
	40SLEEVE	4.0mm Waxing Sleeve for Cemented Post & Core Abutment	(128-0627)
	45SLEEVE	4.5mm Waxing Sleeve for Cemented Post & Core Abutment	(128-0635)
	50SLEEVE	5.0mm Waxing Sleeve for Cemented Post & Core Abutment	(128-0643)
	60SLEEVE	6.0mm Waxing Sleeve for Cemented Post & Core Abutment	(128-0650)
	30BALLSLEEVE	3.0mm Midi Ball Waxing Sleeve	(128-1039) \$26.25
	30ABUTSLEEVE	3.0mm Midi Abutment Waxing Sleeve-Profiler	(128-1047) \$26.25
	38ABUTSLEEVE	3.8mm Abutment Waxing Sleeve-Profiler	(128-1054) \$26.25
	43ABUTSLEEVE	4.3mm Abutment Waxing Sleeve-Profiler	(128-1062) \$26.25
	35HEALABUT35	3.5mm Diameter Threaded Healing Abutment	(128-1112) \$105.00
	40HEALABUT35	4.0mm Diameter Threaded Healing Abutment	(128-1120)
	45HEALABUT35	4.5mm Diameter Threaded Healing Abutment	(128-1138)
	50HEALABUT41	5.0mm Diameter Threaded Healing Abutment	(128-1146)
	60HEALABUT39	6.0mm Diameter Threaded Healing Abutment	(128-1153)
	PCTRD35	3.5mm Threaded Post & Core Abutment	(128-0544) \$89.60
	PCTRD40+	Multi-size Threaded Post & Core Abutment (Fits 4.0, 4.5, 5.0 and 6.0mm Diameter Omni-Tight Implants)	(128-0593)
	35TTRANSFR	3.5mm Diameter Threaded Post Transfer Device	(128-1195) \$46.70
	40TTRANSFR	4.0mm Diameter Threaded Post Transfer Device	(128-1203)
	45TTRANSFR	4.5mm Diameter Threaded Post Transfer Device	(128-1211)
	50TTRANSFR	5.0mm Diameter Threaded Post Transfer Device	(128-1229)
	60TTRANSFR	6.0mm Diameter Threaded Post Transfer Device	(128-1237)




















INDIVIDUAL COMPONENTS (CONT.)

PRODUCT	PART #	DESCRIPTION	PRICE
	35TSLEEV	3.5mm Diameter Threaded Post Waxing Sleeve	(128-1278) \$26.25
	40TSLEEV	4.0mm Diameter Threaded Post Waxing Sleeve	(128-1286)
	45TSLEEV	4.5mm Diameter Threaded Post Waxing Sleeve	(128-1294)
	50TSLEEV	5.0mm Diameter Threaded Post Waxing Sleeve	(128-1302)
	60TSLEEV	6.0mm Diameter Threaded Post Waxing Sleeve	(128-1310)
	35ANALGT	3.5mm Diameter Threaded Post & Core Analog	(128-1369) \$46.70
	40ANALGT	4.0mm Diameter Threaded Post & Core Analog	(128-1377)
	45ANALGT	4.5mm Diameter Threaded Post & Core Analog	(128-1385)
	50ANALGT	5.0mm Diameter Threaded Post & Core Analog	(128-1393)
	60ANALGT	6.0mm Diameter Threaded Post & Core Analog	(128-1401)


INDIVIDUAL SURGICAL INSTRUMENTS

PRODUCT	PART #	DESCRIPTION	PRICE
	PPUNCH	Pilot Punch (1.8mm Diameter)	(127-9538) \$87.50
	TPUNCH35	3.5mm Tissue Punch (Use 1mm larger than Osteotomy Drill)	(127-9553) \$96.25
	TPUNCH40	4.0mm Tissue Punch (Use 1mm larger than Osteotomy Drill)	(127-9561)
	TPUNCH45	4.5mm Tissue Punch (Use 1mm larger than Osteotomy Drill)	(127-9579)
	TPUNCH50	5.0mm Tissue Punch (Use 1mm larger than Osteotomy Drill)	(127-9587)
	TPUNCH55	5.5mm Tissue Punch (Use 1mm larger than Osteotomy Drill)	(127-9595)
	TPUNCH65	6.5mm Tissue Punch	(128-1443) \$144.55
	HCITL13	13mm Length Healing Cap Insertion Tool – Short	(127-9611) \$72.80
	HCITL18	18mm Length Healing Cap Insertion Tool – Long	(127-9629)
	DRVBT13	13mm Length Drive Bit (Driver Carrier) for Implant – Short	(127-9660) \$99.15
	DRVBT18	18mm Length Drive Bit (Driver Carrier) for Implant – Long	(127-9678)
	MIDIBALLDRIVS	Midi Hex Driver for Midi Ball Implant – 18.5mm Short Length	(128-1484) \$102.05
	MIDIBALLDRIVL	Midi Hex Driver for Midi Ball Implant – 24.5mm Long Length	(128-1492) \$102.10
	MIDIABUTDRIVS	Midi Driver for Midi Abutment Implant – 18.5mm Short Length	(128-1534) \$107.40
	MIDIABUTDRIVL	Midi Driver for Midi Abutment Implant – 24.5mm Long Length	(128-1542) \$107.40




**INDIVIDUAL SURGICAL INSTRUMENTS
(CONT.)**

PRODUCT	PART #	DESCRIPTION	PRICE
	HEXDRIV-10	Hex Driver for Bolt Retained Post / Healing Abutment, 10mm length	(127-9728) \$70.00
	HEXDRIV-16	Hex Driver for Bolt Retained Post / Healing Abutment, 16mm length	(128-1575)
	PCDRIV	Threaded Post and Core Driver	(127-9736) \$105.00
	PDRILL-16	Pilot Drill - 1.8mm diameter X 16mm length	(127-9744) \$87.50
	PDRILL-21	Pilot Drill - 1.8mm diameter X 21mm length	(128-1583)
	DRIL23-16	Midi Osteotomy Drill - 2.3mm diameter X 16mm length (128-1625) (used exclusively for Midi Implants)	\$160.30
	DRIL35-16	Osteotomy Drill - 3.5mm diameter X 16mm length	(128-1633)
	DRIL40-16	Osteotomy Drill - 4.0mm diameter X 16mm length	(128-1641)
	DRIL45-16	Osteotomy Drill - 4.5mm/ 5.0mm diameter X 16mm length (same drill may be used for either 4.5 or 5.0mm)	(128-1658)
	DRIL57-16	Osteotomy Drill - 6.0mm diameter X 16mm length	(128-1666)
	DRIL23-21	Midi Osteotomy Drill - 2.3mm diameter X 21mm length (127-9769) (used exclusively for Midi Implants)	\$160.30
	DRIL35-21	Osteotomy Drill - 3.5mm diameter X 21mm length	(127-9777)
	DRIL40-21	Osteotomy Drill - 4.0mm diameter X 21mm length	(127-0795)
	DRIL45-21	Osteotomy Drill - 4.5mm/ 5.0mm diameter X 21mm length (same drill may be used for either 4.5 or 5.0mm)	(127-9793)
	DRIL57-21	Osteotomy Drill - 6.0mm diameter X 21mm length	(127-9801)
	HARDRIL	Hard Bone Drill	(127-9827) \$113.75
	DRLEXT	Drill Extender	(127-9835) \$72.95
	GAGE 1	Depth Gage Short	(127-9850) \$43.75
	GAGE 2	Depth Gage Long	(127-9868)
	CRPIN	Chrome Guide Pin (for 3D Guide Stent Fabrication)	(127-9918) \$11.69
	WRENCH	Wrench for Placement of Implants	(127-9900) \$291.90
	TORQUE	Torque Wrench	(127-9926) \$110.95







ANCILLARY SURGICAL PRODUCTS

PRODUCT	PART #	DESCRIPTION	PRICE
	BONEMAP	Bone Mapping Instrument: Includes 6 acupuncture needles, 8 stops and 1 needle holder	(127-8290) \$65.80

LABORATORY PRODUCTS

PRODUCT	PART #	DESCRIPTION	PRICE
	LABKIT	Laboratory Kit Includes One Each: Post Holder, Chrome Guide Pin, Midi Abutment Analog, Midi Ball Analog, 3.5 Analog, 4.0 Analog, 4.5 Analog, 5.0 Analog, 6.0 Analog	(127-9298) \$206.80
	CRPIN	Chrome Guide Pin For Guide Stent Fabrication	(127-9918) \$11.70
	PDRILL-21	Model Drill	(128-1583) \$87.50

MARKETING PRODUCTS

PRODUCT	PART #	DESCRIPTION	PRICE
	D2311	Dr. Gordon Christensen DVD Implant Replacement for Single Teeth – The BASIC System	(127-9249) \$73.10
	18691	Patient Education Flyers	(Contact Marketing)
	19204	Quick Reference Guide	(Contact Marketing)
	DENTDY05	Reprint: Bohunicky, Dr. Bruce. "Non-submerged Implant Surgery for the General Practitioner." Oral Health Aug. 2005	(Contact Marketing)
	DENTDY02	Reprint: Fernandes, Dr. Americo. "From Bridges to Implants." Dentistry Today, Aug. 2002	(Contact Marketing)
	DENTDY01	Reprint: Fernandes, Dr. Americo. "Combining the Single Implant With a CAD/CAM Restoration." Dentistry Today Dec. 2001	(Contact Marketing)

EDUCATION



Implant Prosthetic Model Kit

Three-model set includes one clear hadar bar model with a clear lower base with four-implant cast hadar bar and overdenture containing EDDS clips; one clear lower base model with implants and ball abutments, overdenture with O-rings; and one crown and bridge model with single unit anterior crown and three unit posterior bridge supported by implants and cementable abutments. Each set comes in a black protective plastic, hard-shelled case with padded interior. Models are made of high-impact plastic resins and denture teeth. Model: hadar bar: 2-1/2" W x 1-1/2" H x 2" D; clear 2 ball and o-ring: 2-1/2" W x 1-3/8" H x 2" D; case: 8-1/4" W x 2" H x 5" D. Each055-5292

1@ \$877.95

Exclusive Education Brochure

Easy-to-understand information and an eye-catching contemporary design come together in these cost-effective educational brochures to keep your patients up-to-date on the latest dental issues and help them make educated decisions about their oral care. Each six-panel brochure with space provided back for you to add practice information and special instructions. 3 5/8" W x 8 7/8" H, 50/Pkg Each Pkg050-6642

1@ \$21.95
2@ \$20.95
4@ \$19.95



ADA Flip Guide to Dental Implants

Covers implant types and benefits and guides patients through the procedure with large illustrations and easy-to-read text.

Each051-0602 1@ \$69.95



Single Pocket Literature Holder

Holding literature up to 4 1/8" W x 1 3/4" D, this single pocket unit stands on your countertop or mounts on the wall (mounting hardware not included) to make it stand out.

4 3/8" W x 7 3/4" H x 3 1/4" D
Each058-7634 1@ \$8.45
3@ \$7.45




Single Pocket Business Card Literature Display Rack - O®

Leaflet-size clear plastic pocket displays literature with a place for business cards in front. This display stands alone on countertop for easy accessibility. It's a great way to make sure clients take your card. Holds literature up to 4 1/8" W x 1 3/4" D.

4 3/8" W x 7 3/4" H x 4 1/8" D.
Each.....052-6152 1@ \$13.95
3@ \$12.95

IMPLANTS MAINTENANCE INSTRUMENTS



Implacare™ Implant Maintenance Instruments - Hu-Friedy

Used to remove plaque and calculus without damaging titanium abutments.

Implacare Tip Refill Packages

Package/10 pair
#4R/4L371-1629 1@ \$78.40
#204S371-1637 3@ \$75.80
#H6/H7371-1645

Implacare Implant Maintenance Kit #6 handle

Kit 1@ \$91.60
#16 371-5521
#26371-5539
#36371-5547
#A6371-5554

Implacare #6 Handle

Each371-5505 1@ \$20.65

Implant Cleaning Instrument - Brevet



144-8760




144-8778



144-8786



144-8794



144-8802

EachSee Item 1@ \$12.85



Implant Scalers – Premier Dental Products

Universal (Columbia 4L/4R)		
2/Pkg	575-8990	1@ \$86.75
Universal (Columbia 4L/4R)		
5/Pkg	575-8982	1@ \$184.75
Facial (Goldman-Fox)		
2/Pkg	585-3726	1@ \$88.15
Facial (Goldman-Fox)		
5/Pkg	585-3734	1@ \$177.75

Implant-Prophy+™ Plastic Implant and Composite Instrument System Kit – AIT

Contains one each Gracey 5/6, Gracey 11/12, Gracey 13/14, Columbia 13/14 Universal and Special Synthetic Sharpening Stone
Each..... 104-9808 1@ \$130.90

MISCELLANEOUS

PATTERSON



Resin Die Stone

Use for the construction of working models and dies for implant cases, precision attachment cases, and all crown and bridge cases. Provides smooth mixing and excellent stackability while promoting free material flow during vibration.

25lb.....085-1972 1@ \$58.35



Die Stone

Use for construction of working models and dies for implant cases, precision attachment cases, and all crown and bridge cases. Provides smooth mixing and ample setting time. Suitable for use with all impression materials.

25lb.....085-1956 1@ \$29.00

Implant Cement™ – Premier®

Designed for Implants; Per
Premier Implant Cement is resin-based cement for implant restorations and provisionals when retrievability is as important as retrievability.

- High retention eliminates patient recalls for re-cementation
- Easily retrievable when needed
- Low solubility with good marginal seal
- Initial gel-phase in 2.5 minutes for easy removal of excess cement



- Rigid, final set for excellent marginal seal and firm retention
- Automix syringe saves time
- Aesthetic gingival shading for natural appearance
- Virtually no taste or odour for increased patient comfort

Kit

Includes: 1 - 5ml Automix dual-barrel syringe, 10 mix tips, 1 mix pad.

Kit580-0834 1@ \$74.10

Mix Tips Refill

25 Automix dual-barrel syringe mix tips.

Kit580-0842 1@ \$36.75
6@ \$35.70



Cavitron® SofTip™ Ultrasonic Implant Insert – Dentsply Professional

- Benefit of ultrasonic lavage during implant maintenance
- Compatible with 30 kHz Cavitron® ultrasonic scalers

Contains Implant Insert, 15 SofTip disposable prophylaxis tips and 1 tip tool (wrench)

Each.....273-9852 1@ \$257.50

SofTip™ Disposable Prophylaxis Tips

100/Pkg494-8204 1@ \$133.95

Tip Tool (wrench)

Each.....494-8253 1@ \$27.15

Implant Organizer – Zirc

- Separate location for the Driver Tips and Torque Wrench; instruments have less chance of getting damaged
- Individual compartments for implants - quick and easy recognition
- The arch shape of the organizer matches the arch shape of the mouth for implant placement
- Cover prevents instruments or implants from falling out
- Available in blue only
- Not autoclavable

Each771-2888 1@ \$28.50



Miratray Implant – Hager Worldwide

Simplifying Impressions for Implants

When stability is important, Miratray Implant meets the demands of fit, rigidity and retention. When time is important, a dental practitioner can take a final implant impression in one session for dentulous and edentulous patients. Starter Sets include: 4 upper & 4 lower trays, calipers & measuring template

Dentulous Miratray Implant Set

Set.....777-0407 1@ \$135.80

Edentulous Miratray Implant Set

Set.....777-0415 1@ \$136.50

SensiTemp™ Resin – Sultan

Resin Desensitizing Temporal Cement

All-purpose desensitizing temporary cement for short and long-term cementation of temporary restorations and implant-retained crowns.

- Soft initial set for easy cleanup
- Final rigid set with high compressive strength
- Potassium nitrate formula to help reduce sensitivity
- Eugenol-free formula
- Hand-held automix delivery and automix cartridge

Hand Syringe

Contains: 1-4ml hand-held automix system and 20 mini tips.

Kit665-9759 1@ \$40.95

