

Instruments, Surgical

ELEVATORS

PATTERSON

					
1 Heidbrink 086-2011	1L Seldin 086-2029	1R Seldin 086-2037	2 Heidbrink 086-2045	2S Seldin 086-2052	3S Seldin 086-2060
					
4L Seldin 086-2078	4R Seldin 086-2086	8 Crane Pick 086-2110	12M 086-2128	13M 086-2136	15 Woodward 086-2144
					
16 Woodward 086-2151	25 Cryer 086-2169	26 Cryer 086-2177	31 Cryer 086-2185	32 Cryer 086-2193	34 Seldin 086-2201
					
34S Seldin 086-2219	40 Heidbrink 086-2227	41 Heidbrink 086-2235	44 Cryer 086-2243	45 Cryer 086-2250	46 086-2268
					
46R 086-2276	71 Miller Apexo 086-2284	72 Miller Apexo 086-2292	73 Miller Apexo 086-2300	74 Miller Apexo 086-2318	80 086-2326
					
81 086-2334	81A 086-2342	81E 086-2359	301 Apexo 086-2367	303 Apexo 086-2375	6 Potts T Handle 086-2094
					7 Potts T Handle 086-2102



Large Handle

1@ \$22.50
3@ \$21.05

Root Tip Elevators



1 Heidbrink
086-2532



2 Heidbrink
086-2540



3 Heidbrink
086-2557



4 Heidbrink
086-2565



5 Heidbrink
086-2573



1 Howard
086-2631



2 Howard
086-2649



3 Howard
086-2656



4 Howard
086-2664



5 Howard
086-2672

1@ \$22.90

3@ \$20.60

Hu-Friedy



Apexo Surgical Elevators



E301
369-6804



E302
369-6838



E303
369-6861

1@ \$82.50

Cogswell Elevators



EA
369-3108

1@ \$82.55



EB
369-3116

1@ \$82.55

Cryer Root Elevators



E25
369-5160



E26
369-5178



E44
369-7380



E45
369-7414

1@ \$82.55 3@ \$80.00

Heidbrink Elevators



E1
369-5764



E2
369-5798



E40
369-5806



E41
369-5814

1@ \$82.55 3@ \$80.00

Miller Apexo Elevators



E71
369-5319



E72
369-5327



E73
369-5335



E74
369-5343

1@ \$82.55 3@ \$80.00

Schmeckebier Elevators



E4
369-8271



E5
369-8289

1@ \$82.55 3@ \$80.00

Seldin Elevators



E34
369-8214



E34S
369-8222

1@ \$82.55

Serrated Elevators



E77
369-7562



E77R
369-7570

1@ \$82.55

Surgical Elevators



E190
369-3439



E191
369-3447

1@ \$82.55 3@ \$80.00

Instruments, Surgical

ELEVATORS (CONT.)

Hu-Friedy (continued)

Potts Surgical Elevators



E6X
369-3462




E7X
369-3470

1@ \$82.55
3@ \$80.00

Davis Root Tip Elevators




ED1 373-0389 ED2 373-0413 ED3 373-0447 ED4 373-0470



ED5 373-0504 ED8 373-0512 ED9 373-0520

1@ \$55.40 3@ \$54.50



ED10 373-0553 ED11 373-0561

1@ \$69.30

West Root Tip Elevators




EW1 373-0215 EW2 373-0223 EW3 373-0231



EW4 373-0249 EW5 373-0256 EW6 373-0264


1@ \$55.40 3@ \$53.35

Heidbrink Root Tip Picks



EHB13/14 373-0595 EHB1 373-0611 EHB2 373-0629 EHB3 373-0637

1@ \$75.15




EHB4 373-0645 EHB5 373-0652

1@ \$65.80 6@ \$64.50



Premier

Periosteal Elevators




9 577-1605 MOLT 2 577-1621

1@ \$50.00 3@ \$48.85

Hu-Friedy

Molt Periosteal Elevator




P9 372-1529

1@ \$69.65

Hu-Friedy

Periosteal Elevator, Prichard




PPR-3 372-1503

1@ \$77.65
3@ \$75.10

Hu-Friedy

Miller Bone Files




BF21 368-2903 BF45 368-2937 BF52 368-2960 BF64 368-2994


1@ \$86.45

FORCEPS

Extracting Forceps, Lower Molar



	15	16	17	23
	1st and 2nd Lower Molar Universal	1st and 2nd Lower Molar Cow Horn	1st and 2nd Lower Molar Universal	1st and 2nd Lower Molar Cow Horn
Patterson	086-4066	086-4074	086-4082	086-4132
Hu-Friedy		370-6439	370-5605	370-6488
Miltex	497-8797	497-8805	497-8813	497-8847



	217	222
	1st and 2nd Lower Molar Universal	3rd Lower Molar Universal
Patterson	086-4397	086-4405
Hu-Friedy		370-9029
Miltex	497-9084	497-9092

PATTERSON

1@ \$46.25

Hu-Friedy


1@ \$202.55 3@ \$196.10

Miltex


1@ \$178.20

FORCEPS (CONT.)


Extracting Forceps, Upper Molar



	10H	10S	18L	18R
	Upper Molar Hook Handle Universal	Upper Molar Universal	1st and 2nd Upper Molar Left	1st and 2nd Upper Molar Right
Patterson	086-4033	086-4041	086-4108	086-4090
Hu-Friedy			370-5662	370-5654
Miltex	497-8771	497-8789	497-8821	497-8839



	24	32	32A	53L
	Upper Molar Universal	Upper Cuspid, Bicuspid and Molar Universal	Upper Cuspid, Bicuspid, Molar and Root Universal	1st and 2nd Upper Molar Left
Patterson	086-4140	086-4157	086-4165	086-4199
Hu-Friedy				370-6629
Miltex	497-8854	497-8862	497-8870	497-8888



	53R	88L	88R	89
	1st and 2nd Upper Molar Right	1st and 2nd Upper Molar Left	1st and 2nd Upper Molar Right	Upper Molar Right
Patterson	086-4181	086-4264	086-4256	
Hu-Friedy	370-6611	370-6603	370-6595	370-6579*
Miltex	497-8896	497-8946	497-8953	497-8961



	90	210	210S
	Upper Molar Left	3rd Upper Molar Hook Handle Universal	3rd Upper Molar Universal
Patterson		086-4371	086-4389
Hu-Friedy	370-6587*		
Miltex	497-8979		497-9076

PATTERSON

1@ \$46.25

Hu-Friedy


1@ \$202.55 3@ \$196.10

*1@ \$335.40

Miltex


1@ \$178.20

Extracting Forceps, Lower Anterior




	85A	103	151	151A
	Lower Molar, Cuspid and Bicuspid Universal	Lower Bicuspid, Incisor and Root Universal	Lower Incisor, Cuspid, Bicuspid and Root Universal	Lower Incisor, Cuspid and Bicuspid Universal
Patterson	086-4249	086-4314	086-4348	086-4355
Hu-Friedy			370-9003	370-9011
Miltex	497-8938		497-9043	497-9050

Extracting Forceps, Upper Anterior



	1	62
	Upper Incisor and Cuspid	Upper and Lower Incisor, Cuspid, Bicuspid and Deciduous
Patterson	086-4017	086-4207
Hu-Friedy	370-6355	370-6512
Miltex	497-8763	497-8904



	203
	Lower Bicuspid, Cuspid, Incisor and Root Universal
Patterson	086-4363
Hu-Friedy	
Miltex	497-9068




	65	99A	99C
	Upper Incisor, Root and Fragments	Upper Incisor, Cuspid and Bicuspid	Upper Incisor, Cuspid and Bicuspid
Patterson	086-4215	086-4280	086-4298
Hu-Friedy	370-6645		
Miltex	497-8912		497-8987

PATTERSON
1@ \$46.25

Hu-Friedy
1@ \$202.55 3@ \$196.10
**1@ \$205.45


Miltex
1@ \$178.20



	101	150	150A	150AS
	Upper and Lower Bicuspid and Deciduous	Upper Incisor, Cuspid, Bicuspid and Root	Upper Incisor, Cuspid and Bicuspid	Upper Anterior Split Beak Serrated
Patterson	086-4306	086-4322	086-4330	
Hu-Friedy		370-8989	370-8997	370-8955**
Miltex	497-8995	497-9019	497-9027	497-9035

FORCEPS (CONT.)

Extracting Forceps, Mead



	MD1	MD2	MD3	MD4
	Upper Incisor, Cuspid, Bicuspid and Root	Upper 1st, 2nd and 3rd Molar Universal	Lower Incisor, Cuspid Bicuspid and Root	Lower 1st, 2nd and 3rd Molar Universal
Patterson	086-4710	086-4728	086-4736	086-4744
Hu-Friedy	370-6215	370-6223	370-6231	370-6249
Miltex	497-9217	497-9225	497-9233	497-9241

PATTERSON

1@ \$45.85


Hu-Friedy

1@ \$251.55


Miltex

1@ \$237.55

Extracting Forceps, Pado



	16S	62	101
	Lower Molar Cow Horn	Upper and Lower Incisor, Cuspid Bicuspid and Deciduous	Upper and Lower Bicuspid and Deciduous
Patterson	086-4918	086-4207*	086-4306*
Hu-Friedy	370-6454	370-6512	370-6710
Miltex	497-9399	497-8904**	497-8995**



	150S	151S
	Upper Incisor, Bicuspid and Root	Lower Incisor, Bicuspid and Root
Patterson	086-4926	086-4934
Hu-Friedy	370-6314	370-6322
Miltex	497-9407	497-9415

PATTERSON

1@ \$46.25

* 1@ \$46.25

Hu-Friedy

1@ \$202.55

3@ \$196.10


Miltex

1@ \$178.20

** 1@ \$178.20



Extracting Forceps, English Pattern



	FX13	FX22	FX33	FX74N
	Lower Bicuspid Universal	Lower Molar Universal	Lower Root	Lower Root Narrow Beak
Patterson	086-4058	086-4116	086-4173	086-4223
Hu-Friedy				
Miltex	497-9274	497-9282	497-9290	497-9308



	22A
	Lower Molar Universal
Patterson	086-4124
Hu-Friedy	
Miltex	


PATTERSON

1@ \$46.25

Miltex

1@ \$220.40

Extracting Forceps, Lower Molar Root



	69	286	46L	97
	Upper and Lower Fragment and Root	Upper Bicuspid, Incisor and Root	Lower Root	Upper Root
Patterson	086-5063	086-5071		
Hu-Friedy			370-6546*	370-6678*
Miltex	497-8920	497-9100		



	300	301
	Upper Root Universal	Lower Root Universal
Patterson		
Hu-Friedy	370-9169	370-9136
Miltex		

PATTERSON

1@ \$46.25

Hu-Friedy

1@ \$251.55







*1@ \$251.55

Miltex

1@ \$178.20

RONGEURS, SURGICAL CURETTES AND KNIVES

Apical Forceps

						
	FAF1	FAF150	FAF151	FAF222	FAFX67A	FAFX74N
	Upper Anterior	Upper Universal	Lower Universal	Lower Molar	Upper Molar	Lower Anterior
Hu-Friedy	370-5712**	370-5720**	370-5738**	370-5746**	370-5811***	370-5829***

Hu-Friedy

**1@ \$184.80

***1@ \$199.90

Rongeurs - Hu-Friedy



R1 Hu-Friedy
16.5cm/6½"
372-9274
1@ \$348.30



R4A Cleveland
17cm/6¾"
372-9852
1@ \$348.30



R1A Hu-Friedy
16.5cm/6½"
372-9373
1@ \$348.30



RBL Blumenthal Standard 30°
14.5cm/5¾"
373-1544
1@ \$348.30



R4 Cleveland
14cm/5½"
373-1056
1@ \$348.30



RF Friedman
14cm/5½"
373-2039
1@ \$325.10




PATTERSON
OFFICE SUPPLIES

PATTERSON
OFFICE SUPPLIES
offers everything
for your front office
including marketing
tools and patient
education
materials.

Contact your
Patterson representative
or local branch to place
an order or for
more information.

Rongeurs - Miltex




Blumenthal
6" 30°
499-1550
1@ \$314.90

Blumenthal
6" 45°
499-1576
1@ \$330.00

Curettes - Hu-Friedy

Lucas Surgical Curettes




CL84 368-5831	CL85 368-5807	CL86 368-5815	CL87 368-5823	CL88 368-5849
------------------	------------------	------------------	------------------	------------------

1@ \$65.70 3@ \$63.65

Curettes - Hu-Friedy

Miller Surgical Curettes




CM8 368-6698	CM9 368-6722	CM10 368-6755	CM11 368-6789	CM12 368-6813
-----------------	-----------------	------------------	------------------	------------------

1@ \$70.90 3@ \$68.70

Instruments, Surgical


RONGEURS, SURGICAL CURETTES AND KNIVES (CONT.)

Molt's Surgical Curettes



Knives - Hu-Friedy

Surgical Knives



Discover our full line
of Patterson Brand
Quality Products.
See our Patterson Private
Catalogue for more details.

We offer a full line of
Quality Equipment.
See our Equipment section
for more details.

SCISSORS

PATTERSON



Crown & Collar Scissors
Curved, Smooth
4 1/4"
089-3842
1@ \$62.55



Dean Scissors
6 3/4"
089-3800
1@ \$117.40



Goldman Fox Scissors
089-3859
1@ \$108.00



Iris Scissors, Curved
4 1/2"
089-3750
1@ \$32.40



Iris Scissors, Straight
4 1/2"
089-3743
1@ \$32.55



Kelly Gum Scissors
Angled, Serrated
089-3826
1@ \$107.65



Kelly Gum Scissors
Curved, Serrated
089-3818
1@ \$100.00



Kelly Scissors, Curved
6 1/4"
089-3792
1@ \$95.60



Lister Bandage Scissors
5 1/2"
089-3768
1@ \$40.70



Littauer Stitch Scissors
4 1/2"
089-3784
1@ \$36.65



Mayo Scissors, Curved
5 1/2"
089-3735
1@ \$31.10



Mayo Scissors, Straight
5 1/2"
089-3727
1@ \$32.55



Operating Scissors, Straight
5 1/2" - Sharp/Blunt
089-3719
1@ \$32.55



Operating Scissors, Straight
5 1/2" - Stainless Steel
089-3701
1@ \$33.55



Spencer Stitch Scissors
3 1/2"
089-3776
1@ \$32.25



Sullivan Scissors
Curved, Serrated
Stainless Steel
089-3834
1@ \$79.70

SCISSORS (CONT.)

Tissue Scissors - Hu-Friedy



S1
Kelly
Curved/Sharp
374-3556

1@ \$108.35



S2
Kelly
Straight/Sharp
374-4042

1@ \$108.35



S5
Wagner
Straight
374-0131

1@ \$82.50



S6
Wagner
Curved
374-0123

1@ \$91.60



S9
Dean
373-7616

1@ \$128.35



S11
Locklin
Curved Handles
374-0636

1@ \$127.70



S12
Locklin
374-0644

1@ \$127.70




S14
LaGrange
Double-Curved
374-2756

1@ \$127.70



S16
Goldman-Fox
374-3069

1@ \$121.25



S17
Iris
Straight/Delicate
373-8655


1@ \$85.15



S18
Iris
Curved/Delicate
373-9141

1@ \$85.15

Crown and Collar Scissors - Miltex




Smooth
Straight 4 1/4"

499-3168

1@ \$63.70 3@ \$61.70


Crown and Collar Scissors, Vantage - Miltex



Smooth
Curved 4 1/4"

499-2681

1@ \$9.50



Smooth
Straight 4 1/4"

499-2657

1@ \$11.10

Crown & Gold Scissors - Hu-Friedy



Crown & Gold
Straight
SCGS

373-8135



Crown & Gold
Curved
SCGS

373-8168

1@ \$80.65

Perma-Sharp™ Surgical Scissors - Hu-Friedy



Iris (18C)
Curved
S5083

373-8705

1@ \$277.35

Suture Scissors - Hu-Friedy



S13
374-1105


1@ \$108.35



S13S
374-1592

1@ \$100.00

Suture Scissors, Vantage - Miltex



Spencer
3 1/2"

499-8050

1@ \$15.60

Suture Scissors - Miltex



Littauer
4 1/2"

499-6773


1@ \$80.20



Littauer
5 1/2"

499-6807


1@ \$81.50



Long
6 1/4"

499-6906


1@ \$118.50



Spencer
3 1/2"

499-8043


1@ \$71.95



Littauer
4 1/2"

499-6781

1@ \$11.40




Littauer
5 1/2"

499-6815

1@ \$11.60

SCISSORS (CONT.)

Tissue Scissors, Carb-N-Sert – Miltex




Iris
Curved 4 1/2"

499-7060

1@ \$228.30


Tissue Scissors, Vantage – Miltex



Iris
Straight 4 1/4"

499-6963

1@ \$12.95




Iris
Curved 4 1/4"

499-7052

1@ \$11.75


Tissue Scissors – Miltex



Iris
Straight 4"

499-7003


1@ \$68.30
3@ \$65.00



Iris
Curved 4"

499-7029

1@ \$68.65
3@ \$66.00




Iris
Curved 4 1/2"

499-7045

1@ \$63.75
3@ \$61.70


Tissue Scissors – Miltex



Kelly
Straight 6 1/4"

499-5148


1@ \$107.25



Locklin
Curved 6 1/4"

499-6922


1@ \$130.05



Locklin
Straight 6 1/4"

499-6930

1@ \$121.10



LaGrange
Double-Curved 4 1/4"

499-5403


1@ \$126.40



Wagner
Curved 4 3/4"

499-8647

1@ \$86.80



Wagner
Straight 4 3/4"


499-8654

1@ \$86.80

NEEDLE HOLDERS, HAEMOSTATS AND TOWEL CLAMPS

PATTERSON

Needle Holder



Mayo-Hegar
6"
088-7208
1@ \$17.00 3@ \$15.25

Needle Holder




Crilewood
6" Surgical grade - Stainless Steel
088-7588
1@ \$16.95 3@ \$15.40

Castroviejo Needle Holder




Vantage	5½" Smooth Jaw 498-6899 1@ \$80.20
Hu-Friedy	(NHCV) 5½" Knurled Jaws 371-5794 1@ \$406.35
Miltex	5½" Smooth Jaw 498-6816 1@ \$449.45

Haemostats - Halsted - Mosquito



Patterson	5" Straight 086-3472 5" Curved 086-3480 1@ \$19.00 3@ \$17.30
Hu-Friedy	5" Straight 371-0787 5" Curved 371-0753 1@ \$71.50
Miltex	5" Straight 498-6808 1@ \$56.10 5" Curved 498-6790 1@ \$58.10

Haemostats - Kelly



Patterson	5½" Straight 086-3415 5½" Curved 086-3423 1@ \$19.20 3@ \$16.90
Hu-Friedy	5½" Straight 371-0530 5½" Curved 371-0548 1@ \$73.55
Miltex	5½" Straight 498-6758 1@ \$58.00 5½" Curved 498-6741 1@ \$60.55

Haemostats - Kelly-Rankin



Hu-Friedy
6¼" Straight 371-0605
6¼" Curved 371-0597
1@ \$69.70

NEEDLE HOLDERS, HAEMOSTATS AND TOWEL CLAMPS (CONT.)

Haemostats - H5

Needle Holders - Hu-Friedy



Hu-Friedy

1x2 teeth Straight 371-0803

1@ \$69.65



Crile-Wood
Grooved
NHCW

371-5828

1@ \$116.10



Hegar-Baumgartner
(Gardner 5")
NHHB

371-6313

1@ \$129.00




Mayo-Hegar
(Gardner 6")
NHMH

371-7295

1@ \$130.30

Perma Sharp™ Needle Holders - Hu-Friedy


With Tungsten Carbide Inserts



Crile-Wood
NH5038

371-5810

1@ \$192.20



Mayo-Hegar

371-5802

1@ \$192.20




Vascular 7"
NH5084

371-7360


1@ \$264.45

Needle Holders - Miltex




Crile-Wood
6"

Miltex	498-6824
	1@ \$94.40
Vantage	498-6832
	1@ \$19.30



Baumgartner
5 1/2"


Miltex	498-6857
	1@ \$97.60
Vantage	498-6865
	1@ \$18.90



Mayo-Hegar
6"


Miltex	498-6873
	1@ \$76.95

Needle Holders
Carb-N-Sert Jaws - Miltex



Crile-Wood
6"


498-6840
1@ \$171.30



Mayo-Hegar
6"


498-6881
1@ \$171.60

Backhaus Towel Clamps - Hu-Friedy



Backhaus
3/4"

375-3381
1@ \$77.35




Backhaus
5"

375-3415
1@ \$77.35


NEEDLE HOLDERS, HAEMOSTATS AND TOWEL CLAMPS (CONT.)

Molt Mouth Gag



Pediatric

Miltex 498-6014
1@ \$339.20




Child

Hu-Friedy 371-3229
1@ \$331.55

Miltex 498-6022
1@ \$339.20

Replacement Tips for Child Mouth Gag - Hu-Friedy

Silicone Rubber 1@ \$13.10
2/Pkg373-1031



Adult


Hu-Friedy 371-3211
1@ \$332.85

Miltex 498-5958
1@ \$336.60

Replacement Tips for Adult Mouth Gag - Hu-Friedy

Silicone Rubber 1@ \$13.10
2/Pkg373-1049

Seldin Retractor - Hu-Friedy



P23

369-0013
1@ \$96.10
3@ \$93.85

Mead Surgical Mallet - Hu-Friedy



1 Mead MAL1

371-2254
1@ \$199.95



Introducing Patterson's Technical Service Centre

Fast - handpieces repaired by the next business day

Reliable - fully trained and experienced in the products we service

Comprehensive - we offer the most repair services available

Convenient - as easy as a phone call

GUARANTEED!

Call Toll-Free 1-800-240-4604

SUTURES

PATTERSON



Each suture is individually wrapped in a sterile pack. The stainless steel needle is swagged directly to an 18" length of suture. 12 sterile packs to a box.



Suture Silk

CP-5.....089-9203 1@ \$31.90
100yd roll 3@ \$30.40

Suture Needles

Eyed, non-sterile suture needle. 1/2 circle, reverse cutting.

12/Package

Extra Small, CP-7.....089-9310 1@ \$31.20
Small, CP-6089-9328 3@ \$28.35

Dental Surgical Sutures

Actual Needle Size	Code	Needle Type	Size	Material	Length	Pkg Size	Item #	Price
 HS-16 1/2 Circle Reverse Cutting	CP62	HS-16 (X-8, C26)	3-0	Black Braided Silk	18"	12/Box	089-9609	1@ \$35.00 3@ \$33.55
 HRT-17 1/2 Circle Point Cutting	CP73	HRT-17 (V-5, Y5)	5-0	Green Braided Dacron*	18"	12/Box	089-9674	1@ \$46.05 3@ \$44.25
 DS-18 3/8 Circle Reverse Cutting	CP66	DS-18 (FS-2, C6)	4-0	Black Braided Silk	18"	12/Box	089-9641	1@ \$35.00 3@ \$33.55
	CP72	DS-18 (FS-2, C6)	3-0	Black Braided Silk	18"	12/Box	089-9666	1@ \$35.00 3@ \$33.55
	CP25	DS-18 (FS-2, C6)	3-0	Chromic Gut	18"	12/Box	089-9567	1@ \$38.80 3@ \$37.35
	CP58	DS-18 (FS-2, C6)	4-0	Chromic Gut	18"	12/Box	089-9559	1@ \$38.80 3@ \$37.35
	CP59	DS-18 (FS-2, C6)	4-0	Chromic Gut	27"	12/Box	089-9526	1@ \$33.45 3@ \$32.20
	CP60	DS-18 (FS-2, C6)	3-0	Chromic Gut	27"	12/Box	089-9500	1@ \$38.80 3@ \$37.35
	CP51	DS-18 (FS-2, C6)	4-0	Plain Gut	27"	12/Box	089-9542	1@ \$38.80 3@ \$37.35
	CP54	DS-18 (FS-2, C6)	3-0	Plain Gut	27"	12/Box	089-9518	1@ \$38.80 3@ \$37.35
	CP71	DS-18 (FS-2, C6)	4-0	Plain Gut	18"	12/Box	089-9658	1@ \$39.10 3@ \$37.65
 DS-24 3/8 Circle Reverse Cutting	CP53	DS-24 (FS-1, C7)	3-0	Plain Gut	27"	12/Box	089-9534	1@ \$38.80 3@ \$37.35
 HS-24 1/2 Circle Reverse Cutting	CP64	HS-24 (X-1, C9/C31)	3-0	Black Braided Silk	18"	12/Box	089-9625	1@ \$32.65 3@ \$29.75
	CP65	HS-24 (X-1, C9/C31)	3-0	Plain Gut	18"	12/Box	089-9633	1@ \$39.50 3@ \$37.70

*Trademark of E.I. DuPont



Silk



Synthetic












Absorbable Plain or Chromic Gut

Instruments, Surgical

SUTURES (CONT.)

Perma Sharp™ Sutures – Hu-Friedy


3/8 Circle	Code	Needle Type	Size	Material	Length	Pkg. Size	Item #	Price
 C-1 Premium Reverse Cutting	PSN1816C	C-1	6-0	Chromic Gut	18"	12/Box	371-3930	1@ \$160.95
 C-3 Premium Reverse Cutting	PSN640S	C-3	5-0	Black Braided Silk	18"	12/Box	371-4110	1@ \$158.05
	PSN641S	C-3	4-0	Black Braided Silk	18"	12/Box	371-4102	1@ \$155.80
	PSN8695P	C-3	6-0	Polypropylene	18"	12/Box	371-7501	1@ \$166.35
	PSN8698P	C-3	5-0	Polypropylene	18"	12/Box	371-7519	1@ \$181.85
	PSN686A	C-3	5-0	Plain Gut	18"	12/Box	371-4037	1@ \$164.50
	PSN687C	C-3	5-0	Chromic Gut	18"	12/Box	371-3971	1@ \$164.50
 C-6 Reverse Cutting	PSN1654C	C-3	4-0	Chromic Gut	18"	12/Box	371-3955	1@ \$163.20
	PSN682S	C-6	5-0	Black Braided Silk	18"	12/Box	371-4128	1@ \$65.20
	PSN683S	C-6	4-0	Black Braided Silk	18"	12/Box	371-4094	1@ \$65.20
	PSN7772S	C-6	3-0	Black Braided Silk	18"	12/Box	371-4078	1@ \$68.65
	PSN392V	C-6	4-0	PGA (Polyglycolic Acid)	18"	12/Box	371-7477	1@ \$80.95
	PSN393V	C-6	3-0	PGA (Polyglycolic Acid)	18"	12/Box	371-7493	1@ \$80.95
	PSN397V	C-6	4-0	PGA (Polyglycolic Acid)	27"	12/Box	371-7485	1@ \$80.95
	PSN8661P	C-6	5-0	Polypropylene	18"	12/Box	371-7527	1@ \$87.50
	PSN635C	C-6	4-0	Chromic Gut	27"	12/Box	371-3963	1@ \$73.20
	PSN636C	C-6	3-0	Chromic Gut	27"	12/Box	371-3948	1@ \$73.20
	PSN7771A	C-6	4-0	Plain Gut	18"	12/Box	371-4011	1@ \$73.20
 C-7 Reverse Cutting	PSNH821A	C-6	4-0	Plain Gut	27"	12/Box	371-4029	1@ \$73.20
	PSNH822A	C-6	3-0	Plain Gut	27"	12/Box	371-3997	1@ \$73.20
 C-22 Reverse Cutting	PSN684S	C-7	3-0	Black Braided Silk	18"	12/Box	371-4086	1@ \$65.20
	PSNH810A	C-7	3-0	Plain Gut	27"	12/Box	371-4003	1@ \$73.20
 C-22 Reverse Cutting	PSN681S	C-22	6-0	Black Braided Silk	18"	12/Box	371-4136	1@ \$56.90
1/2 Circle								
 C-9 Reverse Cutting	PSN632S	C-9	3-0	Black Braided Silk	18"	12/Box	371-4052	1@ \$65.15
	PSN612A	C-9	3-0	Plain Gut	18"	12/Box	371-3989	1@ \$73.20
 C-43 Reverse Cutting	PSN7762S	C-43	3-0	Black Braided Silk	18"	12/Box	371-4060	1@ \$60.45
 D-14 Cutting End Taper	PSN7773L	D-14	5-0	Green Braided Polyester	18"	12/Box	371-4045	1@ \$140.95

 Silk
  Synthetic
  Absorbable Plain or Chromic Gut

Suture needles – Hu-Friedy

Actual size

12/pkg 1@ \$40.65



371-4367 (SN1) 371-4391 (SN2) 371-4458 (SN3) 371-4516 (SN4) 371-4573 (SN5) 371-4631 (SN6) 371-4722 (SN10) 371-4755 (SN16) 371-4789 (SN17) 371-4813 (SN18) 371-4847 (SN20)

Ethicon Silk Sutures – Ethicon

Each Atraloc™ needle is permanently attached to the suture material. Every suture is packed in a sterile pack for convenience and reduced concern of cross infection.

Ethibond Excel Sutures – Ethicon A Coated, Braided Polyester Suture That Is Nonabsorbable

The coating eases the passage of the braided strands through the tissue and provides excellent pliability, handling qualities, and smooth tie-down with each throw of the knot. Compared to silk sutures, Ethibond elicits lower tissue reaction.

Monocryl Undyed Ethicon

Monofilament Suture With No Knot, No Needles, No Ties, No Gut Suture

A monofilament, synthetic absorbable suture. At 7 days post-implantation, approximately 55% of the initial tensile strength remains. All original strength will be lost after 21 days. Because of its monofilament structure, the smooth passage helps reduce trauma to tissue.



Vicryl Rapide Sutures – Ethicon

A Braided, Synthetic, Absorbable Suture

At 5 days post-implantation, only approximately 50% of the tensile strength remains. The suture begins to fall off in just 7 to 10 days as the wound heals. Lower tissue reaction than gut suture.

Vicryl Sutures – Ethicon

A Braided, Synthetic, Absorbable Suture

At 14 days post-implantation, approximately 65% of the tensile strength of coated vicryl remains. Absorption is essentially complete between days 56 and 70. Lower tissue reaction than gut suture.

Ethicon Sutures - Ethicon

3/8 Circle	Code	Needle Type	Size	Material	Length	Pkg. Size	Item #	Price
 C-3 Reverse Cutting	736G	C-3	4-0	Perma-Hand Silk Black Braided	18"	12/Box	408-8936	1@ \$73.15
	 FS-1 Reverse Cutting	684G	FS-1	3-0	Perma-Hand Silk Black Braided	18"	12/Box	408-8928
		H810H	FS-1	3-0	Plain Gut	27"	36/Box	408-9819
 FS-2 Reverse Cutting	682G	FS-2	5-0	Perma-Hand Silk Black Braided	18"	12/Box	408-8944	1@ \$53.50
	683G	FS-2	4-0	Perma-Hand Silk Black Braided	18"	12/Box	408-9801	1@ \$53.50
	661G	FS-2	5-0	Ethilon Nylon	18"	12/Box	408-8852	1@ \$54.40
	662G	FS-2	4-0	Ethilon Nylon	18"	12/Box	408-8860	1@ \$54.40
	635H	FS-2	4-0	Chromic Gut	27"	36/Box	409-0817	1@ \$225.00
	636H	FS-2	3-0	Chromic Gut	27"	36/Box	409-0791	1@ \$225.00
	H821H	FS-2	4-0	Plain Gut	27"	36/Box	409-0874	1@ \$225.00
	H822H	FS-2	3-0	Plain Gut	27"	36/Box	408-9835	1@ \$225.00
	J392H	FS-2	4-0	Vicryl Violet Braided	18"	36/Box	408-9108	1@ \$212.95
	J393H	FS-2	3-0	Vicryl Violet Braided	18"	36/Box	408-9090	1@ \$212.95
 P-3 Reverse Cutting	J422H	FS-2	4-0	Vicryl Undyed Braided	27"	36/Box	408-9124	1@ \$212.95
	J423H	FS-2	3-0	Vicryl Undyed Braided	27"	36/Box	409-3001	1@ \$212.95
	640G	P-3	5-0	Perma-Hand Silk Black Braided	18"	12/Box	409-3118	1@ \$119.95
	641G	P-3	4-0	Perma-Hand Silk Black Braided	18"	12/Box	408-8902	1@ \$117.50
	698G	P-3	5-0	Ethilon Nylon	18"	12/Box	409-3134	1@ \$116.50
	686G	P-3	5-0	Plain Gut	18"	12/Box	408-9231	1@ \$145.00
	687G	P-3	5-0	Chromic Gut	18"	12/Box	408-9207	1@ \$150.95
	1654G	P-3	4-0	Chromic Gut	18"	12/Box	408-9181	1@ \$138.10
	J463G	P-3	5-0	Vicryl Violet Braided	18"	12/Box	408-9140	1@ \$124.80
	J464G	P-3	4-0	Vicryl Violet Braided	18"	12/Box	408-9157	1@ \$124.80
 PC-3 Conventional Cutting	J492G	P-3	6-0	Vicryl Undyed Braided	18"	12/Box	409-3019	1@ \$133.75
	J493G	P-3	5-0	Vicryl Undyed Braided	18"	12/Box	408-9165	1@ \$124.80
	J494G	P-3	4-0	Vicryl Undyed Braided	18"	12/Box	408-9173	1@ \$124.80
	J845G	PC-3	4-0	Vicryl Undyed Braided	18"	12/Box	408-9132	1@ \$150.95



Silk



Synthetic



Absorbable Plain or Chromic Gut

(continued on next page)

SUTURES (CONT.)

Ethicon Sutures - Ethicon

3/8 Circle	Code	Needle Type	Size	Material	Length	Pkg. Size	Item #	Price
 PS-2 Reverse Cutting	1679H	PS-2	3-0	Perma-Hand Silk Black Braided	18"	36/Box	408-8977	1@ \$356.00
	J497G	PS-2	3-0	Vicryl Undyed Braided	18"	12/Box	408-9066	1@ \$124.80
	VR496	PS-2	4-0	Vicryl Rapide	18"	12/Box	409-3324	1@ \$136.50
	VR497	PS-2	3-0	Vicryl Rapide	18"	12/Box	409-3332	1@ \$136.50
1/2 Circle								
 J-1 Reverse Cutting	734G	J-1	4-0	Perma-Hand Silk Black Braided	18"	12/Box	408-9249	1@ \$73.15
 PS-4 Reverse Cutting	1632G	PS-4	5-0	Plain Gut	18"	12/Box	408-9868	1@ \$139.95
 V-5 Taper Cut	X951H	V-5	5-0	Ethibond Excel Polyester	30"	36/Box	408-8894	1@ \$307.95
 X-1 Reverse Cutting	632G	X-1	3-0	Perma-Hand Silk Black Braided	18"	12/Box	409-1765	1@ \$53.50
	612G	X-1	3-0	Plain Gut	18"	12/Box	409-2250	1@ \$82.50

Look Sutures - Surgical Specialties

3/8 Circle	Code	Needle Type	Size	Material	Length	Pkg. Size	Item #	Price
 C3 (P-3) Reverse Cutting	1248B	C3	5-0	Chromic Gut	18"	12/Box	667-6357	1@ \$85.15
	595B	C3	4-0	Chromic Gut	18"	12/Box	667-6308	1@ \$79.15
 C6 (FS-2) Reverse Cutting	500B	C6	5-0	Plain Gut	18"	12/Box	667-6142	1@ \$49.50
	551B	C6	4-0	Plain Gut	18"	12/Box	667-6126	1@ \$49.50
	558B	C6	4-0	Chromic Gut	18"	12/Box	667-6175	1@ \$49.50
	560B	C6	3-0	Chromic Gut	18"	12/Box	667-6167	1@ \$49.50
 C7 (FS-1) Reverse Cutting	786B	C7	3-0	Silk Black Braided	18"	12/Box	667-6043	1@ \$33.40
	553B	C7	4-0	Plain Gut	27"	12/Box	667-6118	1@ \$49.50
 C17 (C-3) Reverse Cutting	774B	C17	5-0	Silk Black Braided	18"	12/Box	667-6019	1@ \$33.40
	535B	C17	5-0	Plain Gut	18"	12/Box	667-6159	1@ \$49.50
 C22 (FS-3) Reverse Cutting	750B	C22	6-0	Silk Black Braided	18"	12/Box	667-6050	1@ \$33.40
 C19 (J-1) Reverse Cutting	783B	C19	3-0	Silk Black Braided	18"	12/Box	667-6035	1@ \$33.40
 C26 (X-8) Reverse Cutting	598B	C26	5-0	Chromic Gut	18"	12/Box	667-6134	1@ \$49.50
	777B	C26	4-0	Silk Black Braided	18"	12/Box	667-6027	1@ \$33.40



Silk



Synthetic



Absorbable Plain or Chromic Gut

Sutures – Tyco



Code	Needle Type	Size	Material	Length	Pkg. Size	Item #	Price
S1173	CE2	5-0	Silk	18"	Box/36	234-5015	1@ \$143.95
SG752	CE4	4-0	Plain Gut	18"	Box/36	234-5544	1@ \$161.25
SS694	CE4	3-0	Silk	18"	Box/36	234-5502	1@ \$136.95
LS638		2-0	Silk	Roll	Each	234-5494	1@ \$3.95
SG636G	C13	3-0	Chromic Gut	18"	Box/12	234-5684	1@ \$57.25
SG635G	C13	4-0	Chromic Gut	18"	Box/12	234-5676	1@ \$57.25
SG822G	C13	3-0	Plain Gut	18"	Box/12	234-5734	1@ \$57.25
SG822	C13	3-0	Plain Gut	18"	Box/36	234-5742	1@ \$174.95
SG821G	C13	4-0	Plain Gut	18"	Box/12	234-5726	1@ \$57.25
SG820G	C13	5-0	Plain Gut	18"	Box/12	234-5718	1@ \$57.25
SG637G	C14	2-0	Chromic Gut	18"	Box/12	234-5783	1@ \$57.25
SG823	C14	2-0	Plain Gut	18"	Box/36	234-5759	1@ \$174.95
SS684G	C14	3-0	Sofsilk	18"	Box/12	234-5809	1@ \$39.25
SS632G	C23	3-0	Sofsilk	18"	Box/12	234-5817	1@ \$39.25
SG5687G	P13	5-0	Chromic Gut	18"	Box/12	234-5700	1@ \$88.95
SG2686	P13	5-0	Plain Gut	18"	Box/12	234-5692	1@ \$98.75
SS641G	P13	4-0	Sofsilk	18"	Box/12	234-5791	1@ \$109.00
SS1640G	P13	5-0	Sofsilk	18"	Box/12	234-5825	1@ \$99.95



Silk



Synthetic



Absorbable Plain or Chromic Gut



Indermil – Tyco


Liquid suture/Tissue Adhesive

12 x 0.5 g/Pack.....234-5007

1@ \$53.50

SURGICAL BLADES

PATTERSON



Surgical Blades

These blades are manufactured to the highest quality control standards in the industry, which assures uniform cutting edges and a perfect fit on scalpel handles. Individually wrapped in foil and hermetically sealed. Sterilized using 2.5 Mrads gamma radiation.

Stainless Steel, 100/Box

#10.....	089-5268	1@ \$37.60
#11.....	089-5102	
#12.....	089-5110	
#15.....	089-5128	
#12B.....	089-5136	
#15C.....	089-5144	

Carbon Steel, 10/Box

#10.....	089-5276	1@ \$6.15
#11.....	089-5201	5@ \$4.65
#12.....	089-5219	
#15.....	089-5227	
#12B.....	089-5235	
#15C.....	089-5243	
#25.....	089-5250	

Disposable Scalpels

Each scalpel is individually packaged in medical grade Tyvek blister pack. Plastic handle includes millimeter scale and the stainless steel blade is covered with plastic blade protector. 10 individually wrapped scalpels per box. Sterilized using ETO gas.

10/Box

#10.....	089-5300	1@ \$12.95
#11.....	089-5318	3@ \$10.50
#12.....	089-5326	
#15.....	089-5334	
#12B.....	089-5342	
#15C.....	089-5292	

Scalpel Blades - Bard-Parker

Carbon Steel, Sterile (each)		
#12.....	129-0568	1@ \$2.60
		50@ \$2.40

Carbon Steel, Sterile (each)		
#10.....	129-1236	1@ \$0.85
#11.....	129-1269	50@ \$0.80
#15.....	129-1293	

Stainless Steel, Sterile (each)		
#10.....	128-9008	1@ \$0.85
		50@ \$0.80

#11.....	129-4289	129-0022	
#22.....	128-9032		
#20.....	128-9024	1@ \$0.95	

Stainless Steel, Sterile (each)		
#12.....	128-9016	1@ \$1.85
		50@ \$1.75

#23.....	128-9040	1@ \$2.05
----------	----------	-----------

Stainless Steel, Sterile (50/Box)		
#12.....	128-9081	1@ \$87.85

Stainless Steel, Sterile (each)		
#15.....	129-0022	1@ \$0.65
		50@ \$0.60

Stainless Steel, Sterile (each)		
#15C.....	129-4321	1@ \$2.60
		50@ \$2.45

Carbon Steel, Non-Sterile (6/Pkg)		
#12.....	129-3885	1@ \$6.35
#23.....	129-3893	
#24.....	129-3901	

Carbon Steel, Non-Sterile (6/Pkg)		
#10.....	129-3786	1@ \$4.50
#11.....	129-3794	25@ \$4.25
#15.....	129-3802	
#20.....	129-3810	
#21.....	129-3828	
#22.....	129-3836	
#25.....	129-3844	

Scalpel Blades - Hu-Friedy

Carbon Steel, Sterile (100/Box)		
#12.....	367-8414	1@ \$138.05

Carbon Steel, Sterile (100/Box) (Double-Edge, Push-Pull Cut)		
#12D.....	367-8422	1@ \$229.65

Carbon Steel, Sterile (100/Box)		
#15.....	367-8430	1@ \$138.05

Scalpel Blades Paragon - Medcom

Carbon Steel, Sterile (100/pkg)		
#12.....	114-8931	1@ \$40.55
		3@ \$37.55

Carbon Steel, Sterile (100/pkg)		
#10.....	114-8915	1@ \$40.55
#11.....	114-8923	3@ \$37.55
#15.....	114-8949	

Stainless Steel, Sterile (100/pkg)		
#12.....	114-8881	1@ \$39.95
		3@ \$37.55
#23.....	114-9079	1@ \$40.55
		3@ \$38.55

Stainless Steel, Sterile (100/pkg)		
#11.....	114-9020	1@ \$40.55
		3@ \$37.55
#15.....	114-8899	1@ \$39.95
		3@ \$37.00

ELECTROSURGERY



MC-4A Electrocauter Unit - Macan

A compact, powerful and economical unit, utilizing fully-rectified current for optimal hemostasis with 100 watts of cutting power. Comes complete with five electrodes and "The Electrosurgery Handbook" manual.

Each040-8500 1@ \$2,443.45

Radio Surge - Macan

The small size, 6¼" x 2¾" x 11½", fits in any drawer or under a cabinet. Weighs only five pounds. Handpiece, cord and electrodes are autoclavable. 100 watt power cuts smoothly. Comes complete with six electrodes and *The Electrosurgery Handbook* manual.

Each040-8518 1@ \$3,101.40

Electrosurgical Electrodes

	Loops										
	476-5053	476-5194	476-5152	476-5228	476-5236	476-5202	476-5145	476-5244	476-5251	476-5269	476-5277
Macan	M22	M23	M25	M16	M17	M18	M19	M20	M21	M26	M27
Coles	E-20-AX	E20AX90	E-20AXL	E-2-BHL	E-2-BHR	E-19-A	E-20-B	E-20-B	E-20-D	E-11-B	E-7-A
Ellman (old)	2B		2BL	2BJ		2A				9A	
Ellman (new)	108	108H	114	103	102	104	105	106	107	124	
Parkell	T-8		T-15	T-9		T-10	T-11				125

	Loops										
	476-5285	476-5293	476-5061	476-5079	476-5087	476-5343	476-5301	476-5319	476-5327	476-5095	476-5103
Macan	M30	M31	M32	M33	M34	M45	M35	M36	M37	M38	M39
Coles		E-12-B	E-13-A	E-13-B	E-13-C		E-6-A	E-6-B		E-9-A	E-10-A
Ellman (old)	4A	4B	3B	3B	3C				8C	8A	8B
Ellman (new)	119	120	126	127	128	129			123	121	122
Parkell		T-7	T-4	T-5	T-5					T-14	T-13

	Needles							Wire	
	476-5046	476-5178	476-5350	476-5038	476-5210	476-5137	476-5368	476-5129	476-5384
Macan	M15	M10	M12	M11	M13	M61	M62	M54	M55
Coles	E-AP1-1/2	E-1-A	E-2-C	E-2-B	E-3-A	E-2-D	E-21-B	E-16-B	
Ellman (old)		1A		1B	1C		7A	6A	6B
Ellman (new)	101	109	111	110	112	116,113F	117	133	134
Parkell	AP1-1/2	T-1		T-2	T-3	C1	C-2	C-6	

	Loops		Ball	
	476-5160	476-5335	476-5186	476-5111
Macan	M40	M42	M41	M51
Coles	E-10-B		E-10-C	E-18-A
Ellman (old)				5A
Ellman (new)		130	131	135
Parkell				C-3



Electrosurgical Electrodes

Fits all electrocauter units listed.

2/Pkg 1@ \$58.00
2@ \$56.45



ellman
Dento-Surg 90 FFP

Dento-Surg® 90 F.F.P. Radiosurgical Instrument
Uses a high frequency radio signal of 3.8 MHz that produces a pressureless, micro-smooth, bloodless incision.
Includes: Handpiece, Antenna Plate, Six Electrodes, Operation Manual & Built in Footswitch.
Each040-0143 1@ \$3,334.85

Disposable Electrodes 3/32" Shaft-Rigid Only

25/Package	
A2D Fine Wire	
3/8" Needle.....315-0612	1@ \$177.55
B1D Fine Wire	
1/4" Loop315-0638	

Retractable Electrode Protector Stand™

- Holds 12 electrodes
- Autoclavable
- Protects electrodes from damage during sterilization
- Will stand on any flat surface for easy identification and access to each electrode

Each314-9804	1@ \$61.55
--------------------	------------

Wipe-Its

For easy removal of tissue from radiosurgery tips.

4/Package315-5025	1@ \$26.35
-------------------------	------------

Radiosurgery Handpiece

For 1/16" - 3/32" shaft electrodes.
2/Package315-1404 1@ \$693.00

Dento-Surg® Antenna Plate-Coated

Each314-8533 1@ \$223.45

Perma Ground

Gold315-2451 1@ \$211.65
Silver315-2469

Dento-Surg® Operation Manual

Each314-9994 1@ \$43.85

Bipolar Forceps

8.6cm (3 7/16")
Each314-8525 1@ \$559.70

Bipolar Silicone Cord 10'

Each314-8517 1@ \$167.20

Electrode Block™

Made of durable lucite. Holds 18 electrodes upright for simple identification and selection.

1/16" Shaft315-0760	1@ \$56.25
3/32" Shaft315-0778	

Electrosurgical Electrodes - Ellman

Set of 6 Standard Electrodes

Contains: 1 each of T118, T108, T121, T124, T127 and T135.
Set314-9796 1@ \$295.35

Sherman Electrode Set

Contains: 1 each of T113F, T114, T117 and T128.
Set314-9846 1@ \$195.30

Flocken Electrode Set

Contains: 1 each of T118, T124, T126 and T127.
Set314-9812 1@ \$196.70

Elliot Feinberg Electrode Set


Contains: 1 each of EFL & EFR.
Set314-9838 1@ \$63.30

Ellman Micro-Tip Needles

2/Package314-9788 1@ \$211.15

Disposable Electrodes 1/16" Shaft-Rigid Only

25/Package	
A2D Fine Wire	
3/8" Needle.....315-0604	1@ \$177.55
A8D Vari-Tip.....315-0653	
B1D Fine Wire	
1/4" Loop315-0620	
C3D Fine Wire 3/16"	
Diamond Loop.....315-0679	
D8D 2.3mm Ball.....315-0646	
F1D Broad Needle -	
Bent Tip315-0687	
113ED315-0661	



Oral Radiosurgery - An Illustrated Clinical Guide (Second Edition) by Dr. Jeffrey A. Sherman

- Step by step instructional surgical techniques
- Over 350 colour photographs and line drawings
- Numerous charts summarizing waveforms, electrode tips and post-operative dressings
- Laboratory practice exercises included
- In depth discussion of radiosurgical theory
- Review of radiosurgical equipment and accessories available

Each315-2485	1@ \$331.65
--------------------	-------------

Instruments, Surgical

Electrosurgical Electrodes - Ellman

Frequency-matched to the Dento-Surg® 90 FFP for minimal tissue alteration and scar-free healing. These electrodes will increase the performance of all electrosurgical machines.



T101
Bendable
Gingival
Troughs
Electrode
1/16" Shaft
2/Package
315-0026



T102
Bendable
Gingival
Troughs
Electrode
1/16" Shaft
2/Package
315-0034



T103
Bendable
Gingival
Troughs
Electrode
1/16" Shaft
2/Package
315-0042



T104
Bendable
Gingival
Troughs
Electrode
1/16" Shaft
2/Package
315-0059



T105
Bendable
Gingival
Troughs
Electrode
1/16" Shaft
2/Package
315-0067



T106
Bendable
Gingival
Troughs
Electrode
1/16" Shaft
2/Package
315-0075



T107
Bendable
Gingival
Troughs
Electrode
1/16" Shaft
2/Package
315-0083



T108
Bendable
Gingival
Troughs
Electrode
1/16" Shaft
2/Package
314-9929



T108H
Bendable
Gingival
Troughs
Electrode
1/16" Shaft
2/Package
315-0091



T109
Bendable
Incising
Electrode
1/16" Shaft
2/Package
315-0109



T110
Bendable
Incising
Electrode
1/16" Shaft
2/Package
315-0117



T111
Bendable
Incising
Electrode
1/16" Shaft
2/Package
315-0125



T112
Bendable
Incising
Electrode
1/16" Shaft
2/Package
315-0133



T113
Bendable
Incising
Electrode
1/16" Shaft
2/Package
315-0141



T113F
Bendable
Incising
Electrode
1/16" Shaft
2/Package
315-0158



T114
Bendable
Periodontal
Electrode
1/16" Shaft
2/Package
315-0166



T115
Bendable
Periodontal
Electrode
1/16" Shaft
2/Package
315-0174



T116
Bendable
Fulgurating
Electrode
1/16" Shaft
2/Package
314-9895



T117
Bendable
Fulgurating
Electrode
1/16" Shaft
2/Package
315-0182



T118
Vari-Tip
Multi-
Purpose
Electrode
1/16" Shaft
2/Package
315-4812



T119
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
315-0190



T120
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
315-0208



T120A
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
314-9903



T121
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
315-0216



T122
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
315-0224



T123
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
315-0232



T124
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
315-0240



T125
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
314-9937



T126
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
315-0257



T127
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
315-0265



T128
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
315-0273



T129
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
315-0281




T130
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
315-0299



T131
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
315-0307



T132
Bendable
Tissue
Removal
Electrode
1/16" Shaft
2/Package
315-0315



T133
Bendable
Coagulating
Electrode
1/16" Shaft
2/Package
315-0323



T134
Bendable
Coagulating
Electrode
1/16" Shaft
2/Package
315-0331



T135
Bendable
Coagulating
Electrode
1/16" Shaft
2/Package
315-0349



T136
Bendable
Coagulating
Electrode
1/16" Shaft
2/Package
315-0356



T136D
Bendable
Coagulating
Electrode
1/16" Shaft
2/Package
314-9971



T137
Small Spoon
Electrode
1/16" Shaft
2/Package
315-0364



T138
Wide Spoon
Electrode
1/16" Shaft
2/Package
315-0372



T151
Lampe Loop
Electrode
1/16" Shaft
2/Package
314-9820

Electrosurgical Electrodes - Ellman

Fits all electrosurgical units listed.

#T118 1@ \$172.15
#T151 1@ \$70.13
Other numbers 1@ \$113.35

#T116
#T120A
#T136D
#T108
#T125

1@ \$156.45

DRY SOCKET DRESSING

Alvogyl – Septodont

Dry Socket Dressing

Incorporates itself with the clot that forms. No suturing required. Provides a moderate anaesthetic action that is effective for several hours after extraction. Contains a sedative and antiseptic.
11gm Jar631-1906 1@ \$87.25



Bone Wax – Johnson & Johnson Medical

Opaque, individually wrapped, sterilized package. A sterile mixture of beeswax and isopropyl palmitate (wax softening agent). Achieves local homeostasis of bone by acting as a mechanical barrier.
12/Box.....405-3807 1@ \$90.45

Dry Socket Paste – Sultan

(Gualocol, Balsam Peru, Eugenol, 1.6% Chlorobutanol)

For instant pain relief. A one-visit treatment with no packing to remove.
1oz Jar665-3836 1@ \$93.95



Hemcon® – Hemcon
Dental Dressing

- Physical barrier that provides pain relief
- Once in contact with blood, the dressing adheres to oral mucosal tissue and seals the site
- Designed to fit extraction sockets or the dressing can be cut to fit

10mm x 12mm
1/Pkg.....362-7791 1@ \$17.45
4/Pkg.....362-7809 1@ \$53.95



PACKING MATERIALS

Spongostan – Johnson & Johnson Medical

Absorbable Gelatin Sponges

To be used in oral surgery as a haemostatic agent and in the obliteration of “dead space” created by simple extraction, roots amputations and removal of cysts, tumors and impacted teeth.
Size: 1cm x 1cm x 1cm
24/Box407-5016 1@ \$56.25



Surgicel – Johnson & Johnson Medical

Absorbable Haemostat (Oxidized Regenerated Cellulose)

Sterile – patient ready. 1/2" x 2"
12/Box408-8837 1@ \$208.95
4@ \$202.50



Surgifoam Absorbable Gelatin Sponge, U.S.P. – Johnson & Johnson Medical

Oral Sponge

- Sterile
- Water-soluble
- Malleable
- Porcine gelatin absorbable sponges intended for haemostatic use by applying to a bleeding surface

2cm x 6cm x 7mm
Box of 12.....408-8704 1@ \$103.95



PERIODONTAL DRESSING

Barricaid® – Dentsply Caulk

Visible Light Cure Periodontal Surgical Dressing

4/Pkg Syringes.....169-1153 1@ \$65.50



Coe-Pak™ – GC America

Periodontal Dressing Standard Package

Contains: 90gm base, 90gm catalyst, and 1/4oz retarder.
Regular.....205-9194 1@ \$83.20
Hard & Fast Set205-9202 3@ \$80.40



Coe-Pak™ Automix NDS Cartridges – GC America

Periodontal Dressing in Cartridge Dispensing Form

2-Pak
Contains: Two double barrel cartridges (50ml each), each containing base and catalyst, 12 mix tips and instructions.
Pkg205-0938 1@ \$102.10



Periocard Periodontal Dressing – Pulpdent

Two-paste, highly elastic periodontal dressing which sets resiliently hard and does not chip or fall apart in the mouth. Assists in tissue placement after periodontal surgery and provides durable protection of tissue. It has a 7-minute working time and sets in 15 minutes. Contains no eugenol or asbestos.
Pkg.....590-1756 1@ \$46.90
3@ \$45.54

Zone Periodontal Pak Kit – Dux

Contains: 1 tube each of base and catalyst, lubricant, mix pad.
Kit.....162-9112 1@ \$36.10
3@ \$34.65

RETRACTORS AND MOUTH PROPS

PATTERSON



Mouth Props

Natural Rubber (non-latex)

Black (2/Pkg)		
Small.....085-4265	1@	\$19.50
Child.....085-4273	3@	\$17.55
Adult.....085-4281		

Silicone Rubber (non-latex)

Tan (2/Pkg)		
Small.....085-4315	1@	\$19.50
Child.....085-4323	3@	\$17.55
Adult.....085-4331		



372-6437 372-8888 372-8987

Cheek Retractors – Hu-Friedy

Useful for intraoral photography. Cold sterilization only.

University of Minnesota (CRM)
Each.....372-6437 1@ \$48.70

Plastic Adult (CRPA)
2/Pkg.....372-8888 1@ \$55.50

Plastic Child (CRPC)
2/Pkg.....372-8987 1@ \$55.50



ComfortView™ – Premier Dental Products

Lip and Cheek Retractor

- Retracts lips and cheeks for optimal buccal and gingival access, maximum visibility and patient comfort
- Suitable for use in diagnostic, preventive and therapeutic dental applications
- Autoclavable

2/pkg		
Universal.....575-7307	1@	\$55.00
Small.....575-7323	2@	\$53.75

Cushion		
10/Pkg.....575-7315	1@	\$14.80

EZ-Prop Mouth Prop – Dentsply/Rinn

The EZ-Prop is an easy way of maximizing access and visibility without sacrificing patient comfort. Made of non-elastic styrofoam, the EZ-Prop lets the teeth dig in and hold firmly.

100/Pkg.....614-0008	1@	\$36.25
	3@	\$34.25

Markel Mouth Props – Hu-Friedy

2/Package		
Small Child.....371-3419	1@	\$54.85
Child.....371-3427		
Adult.....371-3435		



OptraGate™ – Ivoclar Vivadent

- Complete patient comfort
- Easy access to a much larger work space
- Quick and easy placement
- High flexibility and elasticity
- Latex-Free

80/pkg		
Junior.....714-1666	1@	\$108.00
Small.....714-1633		
Regular.....714-1625		

Self-Span™ Cheek & Lip Retractors – Ellman

Self-retaining clear lip and cheek plastic retractors.

- Autoclavable (up to 250° F)
- For patients 7 years and older
- Slips into place easily and comfortably

3/Package.....315-4200	1@	\$92.35
------------------------	----	---------

Spandex Cheek Retractor – Hager Worldwide

Easy to use, more comfortable for the patient. Can be autoclaved or dry sterilized.

Spandex (pair)		
Child.....765-5327	1@	\$51.85
Adult.....765-5319		

Mirahold (pair)

Clear plastic cheek retractors with long handles. Can be cold sterilized.

Child.....765-5178	1@	\$32.00
Adult.....765-5160		

BONE GRAFTING MATERIALS

Bioplast® HTR® – Sybron Kerr

A synthetic bone substitute that promotes new bone growth.

HTR-24

Used for ridge preservation and ridge augmentation procedures

Patented, non-resorbable copolymer that promotes new bone growth. Its unique microporous design allows surrounding bone to grow into and integrate with the Bioplast HTR. Ideal bone substitute for ridge preservation cases in order to prevent the loss of jawbone following tooth extraction. HTR-24 Straight Syringe (0.25 gram each)



Curved syringes	.25g refill	
6/Pkg.....631-2144	1@	\$574.05
2/Pkg.....126-4159	1@	\$217.25

Curved syringes	.125g refill	
6/Pkg.....126-4142	1@	\$366.20

MISCELLANEOUS

Biosplint – Acteon

A Multiple-Use Dental Reinforcement Ribbon

- Stabilization and immobilization of loose teeth
- Immediate replacement of avulsed teeth
- Reinforcement of construction of both provisional and permanent bridges
- Maintenance of interdental spaces

Kit

Contains eight Biosplint splinting ribbons, one Biosplint thermal cutter, three Biosplint brushes and probes

Each.....102-7002	1@	\$435.60
-------------------	----	----------

Splinting Ribbons

Contains eight ribbons individually wrapped in aluminum bags, eight stickers for "Patient Traceability"

Each.....102-7036	1@	\$340.60
-------------------	----	----------

Thermal Cutter

Each.....102-7010	1@	\$110.90
-------------------	----	----------

Accessories

Contains eight brushes and three probes

Each.....102-7028	1@	\$103.00
-------------------	----	----------

Medicament Loop – Hu-Friedy

For subgingival application of medicaments.

Each.....371-2908	1@	\$88.40
-------------------	----	---------

Periowave – Ondine

A NEW Treatment for Periodontal Disease

Periowave uses a novel non-antibiotic approach for treating the underlying infection responsible for periodontal disease. *The procedure is expected to take 20 minutes to complete for the entire mouth.*



Starter Kit Complete

Contains: 1 PPD System, 1 Large Treatment Kit, 2 Small Treatment kit and 1 extra handpiece.
 Kit263-0085 1@ \$6,305.00

Periowave Photodynamic System

System263-0010 1@ \$5,450.00

Large Treatment Kits

Box of 10.....263-0036 1@ \$995.00



Small Treatment Kits

Box of 10.....263-0028 1@ \$495.00

Handpiece

Each263-0044 1@ \$430.00

Fibre Optic Cable

Each263-0051 1@ \$320.40

Fibre Optic Cleaning Tape

Each263-0069 1@ \$34.95

Safety Glasses

Each263-0077 1@ \$149.95



ViziLite™ – Zila

ViziLite™ is a non-toxic, diffuse chemiluminescent light source that when used in combination with conventional visual oral mucosal examination by health care providers, improves identification, evaluation and monitoring of oral mucosal abnormalities in a patient population at increased risk for oral cancer.

Each kit contains 1 ViziLite Device, 1 ViziLite Retractor, 30ml ViziLite Acetic Acid Solution (raspberry flavour), and instructions for use.
 40 Kits/Pkg.....770-6542 1@ \$1,053.40

WEBCOL® Alcohol Preps – Kendall

Premium pad material provides maximum absorbency for scrubbing and cleansin

- 70% isopropyl alcohol
- Sterile, nonwoven sponges designed for moderate cleansing
- Economical



Large, 2 ply

200/Box.....413-1025 1@ \$7.95

MEDICAL DEVICES

VELscope® – LED Dental

Mucosal Examination System

A revolutionary handheld device used for the early detection of abnormal tissue. When used as an adjunctive aid in combination with traditional oral cancer examination procedures, VELscope facilitates the early discovery and visualization of mucosal abnormalities that may lead to oral cancer.



System Package

Contains operation manual, user's guide, light source unit with base stand, handpiece with attached fibre-optic light guide, one pair deluxe safety glasses

Each458-1104 1@ \$4,748.85

Starter Kit

Contains 100 VELcaps, 250 VELsheaths and 100 VELtractors

Each458-1112 1@ \$475.00

VELcap

100/Pkg.....458-1120 1@ \$247.00

VELtractor

100/Pkg.....458-1153 1@ \$102.00

VELscope Vantage® System – LED Dental

An advance version of the revolutionary VELscope Oral Cancer Screening System.

Features a more powerful lamp technology that provides better fluorescence and improved camera imaging; a 2-year comprehensive warranty; enough disposable products (VELcaps and VELsheaths) to cover 200 screenings; a camera adapter kit for digital photography; and access to several continuing education courses related to use of the VELscope system.



System Pack458-1195 1@ \$6,775.00

LASER*



Odyssey® 2.4G Diode Laser – Ivoclar Vivadent

Its portable size, ease of use and affordable price have made the Odyssey® Diode Laser the #1 selling soft tissue laser in dentistry. The Odyssey® 2.4G now delivers the added convenience of wireless power control technology!

With the simple use of the Odyssey® 2.4G, perform Total Soft Tissue Management.

- Eliminate packing retraction cord
- Quick and easy aesthetic gingival sculpting
- Create and maintain haemostasis for a dry operating field
- Provide the most gentle, precise and least invasive dental experience for your patients

Odyssey® Navigator™ – Ivoclar Vivadent

Portable Soft Tissue Laser



- Ergonomic handpiece features an arc design to rest comfortably in your hand
- Pre-cleaved, procedure-specific unidose tips makes preparation fast and easy
- Battery operated for ultimate portability in the operatory
- Intuitive touch screen design for easy activation
- Over 30 pre-set clinical procedures installed
- Wireless foot control
- 810 nm wavelength - The Proven Diode Technology



SIROLaser – Sirona

- Less than one pound, fits in palm of hand
- Can easily be transported from operatory to operatory
- .05 - 7 watts of power
- Easy to read display indicates treatment modality and delivered power
- Pre-programmed and programmable treatment settings
- Foot switch
- Finger switch on handpiece
- Diode is completely maintenance-free

*Please call for pricing.