

Laboratory

MODEL STONES

PATTERSON



Die Stone

Use for construction of working models and dies for implant cases, precision attachment cases, and all crown and bridge cases.

Provides smooth mixing and ample setting time. Suitable for use with all impression materials.

25lb085-1956 1@ \$29.00



Resin Die Stone

Use for the construction of working models and dies for implant cases, precision attachment cases, and all crown and bridge cases. Provides smooth mixing and excellent stackability while promoting free material flow during vibration.

25lb085-1972 1@ \$58.35



Bitestone - Whip Mix

Extremely fast mounting cases.

White 1@ \$62.55
25lb Carton725-9708 4@ \$58.10

CoeCal™ - GC America, Inc.

The original hard, quick-setting, accurate stone. Unexcelled for its minimum processing change.

25lb Carton
Buff204-5342 1@ \$27.35
White204-5375 4@ \$26.00

50lb Drum
Buff204-6381 1@ \$38.50
White204-6415 4@ \$37.50



Denstone® - Heraeus Kulzer

Scientifically compounded to produce an extremely smooth, accurate cast with strength that gives an added factor on difficult cases.

25lb Carton
Golden513-8581 1@ \$46.05
White513-8599 4@ \$44.80

50lb Carton
Golden513-9258 1@ \$63.15
White513-9241 2@ \$61.65
Pink513-9266



Die-Keen® - Heraeus Kulzer Die Stone

Very hard but not brittle. For dies and crown and bridge materials.

Green
25lb Carton514-1551 1@ \$88.40
50lb Carton514-1783 1@ \$116.85

Blue
25lb Carton514-1577 1@ \$88.40
50lb Carton514-1809 1@ \$116.85
3@ \$113.85

Pink
25lb Carton514-1569 1@ \$88.40
50lb Carton514-1791 1@ \$116.85



Die-Stone™ - Heraeus Kulzer

Formulated especially as a master model material. For fixed and removable bridgework as well as individual dies.

Peach
25lb Carton514-2286 1@ \$62.85
4@ \$60.75
50lb Carton514-2807 1@ \$70.30
2@ \$68.65

FlowStone - Whip Mix

A fluid, accurate and dense base stone.

- Pours without vibration - saves time filling bases
- Low setting expansion ensures dimensional stability
- Dense surface resists abrasion from dowel pins
- Available in white or buff to provide excellent contrast to die stones

50lb Canister
White727-7049 1@ \$88.40
Buff727-7031
Black727-7056



Hydrock - KerrLab Model Stone

33lb/15kg Carton
Yellow426-5435 1@ \$50.85
White426-5419
Ivory426-5427

Rapid Stone
33lb
Yellow426-6102 1@ \$50.65
White426-6110
Green426-6128



Jade Stone – Whip Mix
Die Stone

Compatible with all types of impression material.



25lb Carton		
Green.....	728-3716	1@ \$71.20
Blue.....	728-3724	4@ \$69.90
50lb Carton		
Green.....	728-3906	1@ \$86.45
Blue.....	728-3914	2@ \$83.85

Labstone – Heraeus Kulzer

Blended for All-Around Laboratory Use



25lb Carton		
Buff.....	514-8580	1@ \$42.20
Blue.....	514-8598	
50lb Carton		
Buff.....	514-8853	1@ \$56.10
Blue.....	514-8887	2@ \$54.80
Green.....	514-8895	

Microstone – Whip Mix
Model Stone

Especially suited for processing acrylic dentures.



25lb Carton		
White.....	728-5992	1@ \$42.75
Golden.....	728-6008	4@ \$41.85
50lb Carton		
White.....	728-6214	1@ \$65.00
Golden.....	728-6222	4@ \$63.35

Mounting Stone – Heraeus Kulzer
Model Stone

White		
50lb Carton.....	514-4282	1@ \$59.70

Mounting Stone – Whip Mix
Model Stone

Very hard, fast setting for mounting casts to articulators.



25lb Carton		
White.....	728-7501	1@ \$48.30
Blue.....	728-7519	4@ \$46.75

Prima Rock – Whip Mix
Die Stone

Compatible with all types of impression material.



25lb Carton		
Violet.....	729-4101	1@ \$78.55
Yellow.....	729-4119	4@ \$75.25

Preweighed Envelopes		
70gm per envelope, 120/package		
Violet.....	729-4168	1@ \$136.65
Yellow.....	729-4176	

Quickstone – Whip Mix
Model Stone

A quality, economical dental stone.



25lb Carton		
Buff.....	729-9126	1@ \$39.95
Blue.....	729-9183	4@ \$38.60
50lb		
Buff.....	729-9910	1@ \$49.00

ResinRock – Whip Mix

Resin fortified, low expansion die material. Tough, non-brittle model resists chipping when trimming, fitting or finishing. Smooth, dense surface for superior detail reproduction. Highly abrasion resistant. We've added ivory to our selection. This natural colour die stone complements all ceramic restorations.



25kg (55lb)		
Blue.....	730-9024	1@ \$207.25
Grey.....	730-9032	5@ \$203.95
Peach.....	730-9016	
Ivory.....	730-9040	
White.....	730-9057	
Green.....	730-9065	
25 - 70gm packages (case of 4)		
Peach.....	727-7965	1@ \$152.45



ResinRock XL5 – Whip Mix

Has the same smoothness and high abrasion resistance you've come to appreciate in ResinRock. However, this ISO Type 5 die stone gives you increased expansion to compensate for shrinkage of vinyl polysiloxane and other rubber-based impression material.

- Ideal for use with rubber-based impression materials, which may contract over time
- Increased setting expansion of 0.2%
- Resin-fortified for advanced abrasion resistance
- Smooth, dense surface for superior detail reproduction
- Thixotropic, creamy mix stacks and vibrates easily into impression

25kg (55lb)		
Blue.....	727-8021	1@ \$207.25
Green.....	727-8013	5@ \$203.95
Ivory.....	727-8047	
Pink.....	727-8039	

MODEL STONES (CONT.)

Silky-Rock – Whip Mix

Die Stone

Works well with all types of impression material.

Preweighed Envelopes

70gm per envelope, 120/package.	1@ \$87.45
Violet.....	731-3208
Yellow.....	731-3232
White.....	731-3265

25lb Carton

Violet.....	731-4941	1@ \$60.40
Yellow.....	731-4974	
White.....	731-5005	

55lb Carton

Violet.....	731-5146	1@ \$80.45
White.....	731-5153	

Snap-Stone – Whip Mix

Model Stone

Exceptionally fast setting stone with an early high compressive strength.

Compatible with all impression materials.

25lb Carton

Pink.....	731-5799	1@ \$67.20
		5@ \$64.40

50lb Carton

Pink.....	731-5831	1@ \$97.65
-----------	----------	------------

Preweighed Envelopes

80 - 140gm.....	731-5823	1@ \$87.05
-----------------	----------	------------

StatStone – Heraeus Kulzer

- Fast and accurate
- Ready for removal from the impression in just 5 minutes
- Provides high-, early-compression strength and features low expansion
- Ideal for restorations when precision is needed for fabricating custom trays, mouthguards, vacuum-formed splints and denture repairs
- Extremely well suited for bleaching application

25lb Carton

Pink.....	514-9802	1@ \$70.15
-----------	----------	------------

Preweighed Packets – Pink

24/Pkg		
120gm.....	514-9810	1@ \$34.00

Super-Cal™ Plus – GC America

Die and Model Material

- Minimal shrinkage
- Pours easily
- Hard, smooth surface
- High strength to resist breakage and chipping
- Low setting expansion for greater accuracy
- Smooth consistency w/excellent flow and fit
- Superior performance characteristics

25lb Carton

Blue.....	204-4220	1@ \$51.00
Green.....	204-4212	

50lb Carton

Blue.....	204-4303	1@ \$80.35
Green.....	204-4295	

Super-Die – Whip Mix

Model and Die Material

Fast setting and requires special handling when pouring models. Vivid yellow.

25lb Carton.....	732-3082	1@ \$81.15
------------------	----------	------------

Suprastone – KerrLab

Die Stone

High strength, high expansion.

33lb Bag

Green.....	437-5572	1@ \$55.80
Blue.....	437-5580	

Tru-Stone™ – Heraeus Kulzer

High Production, Great Value - Type III

- Combines strength, accuracy, and low expansion
- Ideal for crown and bridge, dies and cast partials
- Excellent surface reproduction
- Quick working time for the most production-minded laboratory

Pink

25lb.....	516-8570	1@ \$41.80
-----------	----------	------------



Vanrock – Van R

- Provides the hardest of stone surfaces without brittleness
- Produces outstanding results when used with any alginate, polyvinylsiloxane, polyether or polysulfide impression material

4lb Carton

Green.....	705-3291	1@ \$26.60
------------	----------	------------

25lb Carton

Green.....	705-3325	1@ \$108.60
------------	----------	-------------

Preweighed Packets – 100gm

25/Pkg.....	705-3333	1@ \$47.45
-------------	----------	------------



Vel-Mix – KerrLab

Die Stone

33lb/15kg Carton

Pink.....	439-5133	1@ \$57.40
White.....	439-5141	

PLASTERS

Lab Plaster – Heraeus Kulzer

General Use Plaster

25lb Carton

Regular Set.....	514-6774	1@ \$35.50
Fast Set.....	514-6782	

50lb Carton

Regular Set.....	514-7251	1@ \$47.35
Fast Set.....	514-7269	





Lab Plaster – KerrLab

A high-quality, dental grade plaster. Most commonly used for articulation and general model work.

33lb		
Regular Set.....	431-1510	1@ \$47.35
Fast Set	431-1528	



Laboratory Plaster – Whip Mix

An economical white plaster formulated for general-purpose laboratory use. Improved to provide a creamier mix consistency.

20lb Carton		
Regular Set.....	729-1909	1@ \$27.65
Fast Set	729-1891	

50lb Carton		
Regular Set.....	729-1982	1@ \$42.20
Fast Set	729-1990	

Model Plaster – Heraeus Kulzer/Jelenko

White		
25lb.....	515-4091	1@ \$40.30
50lb.....	515-4158	1@ \$41.80



Mounting Plaster – Whip Mix

A specially formulated plaster which offers a low setting expansion and fast setting time. Its low compressive strength allows mounting plates to be easily cleaned.

20lb Carton	729-3970	1@ \$42.10
-------------------	----------	------------



Snow White No. 1 Plaster – KerrLab

Model Plaster – Natural Set

33lb Carton	431-2765	1@ \$66.85
		5@ \$65.55

Snow White No. 2 Plaster – KerrLab

Impression Plaster – Rapid Set

33lb Carton		
Plain.....	431-2773	1@ \$66.85
Mint	431-2781	5@ \$65.55

ORTHODONTIC PLASTERS/STONES



Orthodontic Plaster – Heraeus Kulzer

For orthodontia casts and study models.

Super-White		
25lb.....	516-0973	1@ \$39.35
50lb.....	516-0981	1@ \$71.30
		5@ \$70.05



Orthodontic Plaster – Whip Mix

Exceptional hardness and particular whiteness. Developed especially for orthodontic work.

Super-White		
20lb Carton	728-7964	1@ \$34.25
50lb Carton	728-8269	1@ \$52.40



Orthodontic Stone – Whip Mix

Hard, accurate, medium setting stone made with the orthodontist in mind.

Super-White		
25lb Carton	732-5046	1@ \$42.25
50lb Carton	732-4544	1@ \$58.90

O-67 Snow White Orthodontic Stone – Heraeus Kulzer

Specifically formulated for orthodontic casts where whiteness, accuracy, and durability are essential.

25lb.....	514-9554	1@ \$45.75
-----------	----------	------------



Pouring A Model

A thick, creamy mix of dental stone is made by using the proper water/powder proportion. Ideally, the mix should be mechanically spatulated under vacuum. The next best procedure is mechanical spatulation on a vibrator. Hand spatulation should be accomplished by “rubbing” the mix against the walls of the mixing bowl until completely smooth, and then vibrating thoroughly. The pour is made by placing a small amount of the mix in one corner of the impression and placing it on a vibrator, tipping it so the stone flows into the recesses within the impression.

INVESTING MATERIALS

Beauty-Cast – Whip Mix Inlay Investment



- For all precision techniques
- Produces extremely smooth castings

Envelopes

144-50gm	725-8080	1@ \$200.65
144-75gm	725-8635	1@ \$245.55
24-50gm	725-7108	1@ \$48.85

Cera-Fina Investment – Whip Mix Complete with Liquid



- Smooth casting surface
- Extended working time
- Very smooth texture
- Carbon-free

Envelopes

144-60gm	726-2355	1@ \$230.95
144-90gm	726-2405	1@ \$311.55

With Liquid

20lb Carton (2 liters)	726-2256	1@ \$159.65
------------------------------	----------	-------------

Ceramigold® – Whip Mix Phosphate Casting Investment With Carbon



Preweighed Envelopes with Liquid

60gm Per Envelope (1 litre and 1-340ml)		
144/Box.....	726-6471	1@ \$227.05
90gm Per Envelope (2 litres)		
144/Box.....	726-7453	1@ \$274.60

Cristobalite – Whip Mix

Gives high thermal expansion, even at comparatively low temperatures. It is possible to use the same creamy mix by simple standard vacuum or painting techniques for all types of inlays.



Preweighed Envelopes

50gm Per Envelope		
24/Box.....	726-7941	1@ \$41.15
144/Box.....	726-8436	1@ \$188.80

FastFire 15 Investment – Whip Mix

Universal phosphate investment for all crown and bridge alloys. Featuring freeze stable liquid and a 15-minute benchset. Produces smooth, accurate castings in less than an hour. Includes FastFire liquid which allows year-round problem-free shipping and controlled expansion. Reliable performance using many casting techniques and alloys.



Envelopes

150-60gm	727-6447	1@ \$271.95
100-90gm	727-6454	1@ \$271.95
25-60gm	727-6413	1@ \$68.35
25-90gm	727-6421	1@ \$89.10

20lb Carton

Each	727-6371	1@ \$165.00
------	----------	-------------

Fujivest II – GC America

GC Fujivest II is a carbon-free phosphate bonded investment for precision casting of all crown and bridge alloys. Fujivest II can be utilized in either a fast-fire or traditional step technique and with or without the use of a casting ring. Fujivest II is consistent, easy to use and has a long working time for easy investing while providing smooth, reliable castings.



Powder

100 x 60gm	326-4454	1@ \$118.75
67 x 90gm	326-4462	1@ \$107.20

Liquid

900ml.....	326-4470	1@ \$29.40
------------	----------	------------

Fujivest Platinum – GC America Inc.

All-Around Carbon-Free Phosphate-Bonded Investment

- For both quick heating and conventional step-heating process
- Devests easily
- Specifically developed for high-precision castings for crown and bridge work
- Extremely smooth casting surface

Promo Pack

(10x150gm)	203-9980	1@ \$130.00
------------------	----------	-------------

6kg

(100x60gm)	203-9949	1@ \$109.00
(40x150gm)	203-9956	5@ \$110.00
(67x90gm)	203-9964	

Liquid

900ml.....	203-9972	1@ \$33.05
------------	----------	------------

Fujivest Premium – GC America Inc.

All-Around Carbon-Free Phosphate-Bonded Investment

- Suitable for all types of dental alloys
- Accommodates ringless techniques in both quick- and slow-heating procedures
- Optimal expansion control by liquid dilution ensures a highly accurate fit
- Implant suprastructures with non precious alloys is now safe and straightforward

Promo Pack

(10x150gm)	204-0046	1@ \$130.25
------------------	----------	-------------

6kg

(100x60gm)	204-0004	1@ \$107.70
(40x150gm)	204-0012	5@ \$105.00
(67x90gm)	204-0020	1@ \$110.50
		5@ \$107.00

Liquid

900ml.....	204-0038	1@ \$30.55
------------	----------	------------

G-Cera Ceramic Pins – GC America, Inc.



- High-strength ceramic material, ideally used in refractory technique
 - Smooth, durable and can be used in an extended number of firings
 - Will not produce an oxide layer
- 10/Box.....326-4306 1@ \$75.85

G-Cera® Vest – GC America, Inc.

Refractory Investment

- Unrivalled strength, accuracy and fit
- Easier to mix and pour into impressions
- Its hard, smooth finish provides for an improved model surface for laminates and other all-ceramic restorations

200ml
Liquid Only.....326-4298 1@ \$23.95

Hi-Heat – Whip Mix

Soldering Investment

For high-fusing precious metals.
4½lb Can.....727-8161 1@ \$33.70



Hi-Temp – Whip Mix

Phosphate Casting Investment – Carbon Free

- Compensates for alloy shrinkage
- Withstands high-fusing temperatures
- Provides a non-contaminating atmosphere
- Allows casting of a lighter-density metal

25lb Carton
(1 Litre and 1-340ml)..728-0159 1@ \$157.05

Preweighed Envelopes with Liquid
60gm Per Envelope (1-340ml)
144/Box.....728-1116 1@ \$209.90

90gm Per Envelope (2 Litres)
144/Box.....728-0613 1@ \$282.50

Preweighed Envelopes without Liquid
90gm Per Envelope
144/Box.....728-1652 1@ \$220.45



Novocast – Whip Mix

For Low Fusing Crown and Bridge Alloys

An advanced inlay investment that provides increased accuracy for casting conventional gold alloys at a mold temperature of 650° C (1,200°F) without the use of a water bath.

Preweighed Envelopes
50gm Per Envelope
144/Box.....728-7683 1@ \$188.45



Polyvest – Whip Mix

Refractory Investment

- Ample working time to pour models
 - Fine grained with fluid consistency for sharp reproduction of detail
 - High strength, expansion controlled to match a wide range of porcelains
- 24/60gm729-4200 1@ \$128.70



Powercast Investment – Whip Mix

Fine-Grained, Carbon Free for Rapid, Ringless Casting

- Crack resistant and tolerant to rapid temperature rise
- Can be used both with ring or ringless technique
- Easy breakout for faster cleanup
- Smooth, accurate castings reduce finishing time

Preweighed Envelopes with Liquid
60gm Per Envelope
150/Box.....729-1537 1@ \$277.20
20lb Carton
Each729-1511 1@ \$179.20



Soldering Investment – Whip Mix

For all Soldering Purposes

Ideal for soldering all types of dental alloys, including chrome-cobalt. Its fine texture mixes readily with water to a smooth, stackable consistency that invests easily. Heat-shock resistance and high strength enable it to withstand being dried out or flushed with boiling water.

4½lb Jar.....731-5989 1@ \$34.00
25lb Jar.....731-6474 1@ \$71.65



Speed-E – Whip Mix

For Low-Fusing Precious Metals

Fast-setting investment with very high permeability and outstanding heat-shock resistance. This enables it to be dried significantly faster after setting. It shapes and trims easily, yet has adequate strength.

4½lb Jar.....731-6961 1@ \$44.65



V.H.T. Investment – Whip Mix

- Very high temperature
- Refractory material for fabricating porcelain veneers, jacket crowns and inlays

Envelopes w/Liquid
24/60gm/Liquid.....732-9121 1@ \$131.95

INVESTING ACCESSORIES



A B C

Crucibles – Dixon

Casting

Suitable for high temperature melting of platinum, chromium, gold alloy, etc.

B. Large "A"	284-3985	1@ \$45.95
For Large Machines		1@ \$47.25
C. Slotted "A"	284-3944	

Crucibles – Kerr

Casting

Clay		1@ \$36.30
Small – 1oz		
2/Pkg	417-7929	
Large – 2½oz		1@ \$24.40
1/Pkg	417-8414	
Quartz		1@ \$55.15
Slotted – 2½oz		
1/Pkg	417-8505	



Crucible Formers – Whip Mix

Casting Crucible Formers

Diameter Number	Product	Description Method	Casting
1¼" (31.6mm)	4081	Conical	Centrifugal
1½" (38mm)	4179	Conical	Centrifugal
1¾" (44.5mm)	4177	Conical	Centrifugal

1¼" Formers 1@ \$20.45
Conical-Large #4081...727-0457

1½" Formers 1@ \$22.90
Conical #4179727-1612

1¾" Formers 1@ \$22.80
Conical #4177727-1067

Casting Rings – Whip Mix

A Special Pyronel Alloy Metal Construction

Will remain precision-fitted and stable in size even after repeated heatings.

1¼" x 1¾"	730-0833	1@ \$18.70
1¼" x 1½"	730-1328	1@ \$20.35
1¼" x 1⅞"	730-1815	1@ \$23.10
1¾" x 2¾"	730-1823	1@ \$35.00
1½" x 1¾"	730-3282	1@ \$25.20
1¾" x 1¾"	730-3779	1@ \$28.50

Casting Liner – GC America

10m/Roll
Roll.....326-4314 1@ \$80.95

PLASTER BINS

Triplex – Nevin

Plaster Bin


- 21" x 18" x 8"
 - Three compartments with cover
 - 75lb capacity
 - Gravity feed
- French Grey.....536-5374 1@ \$376.20
Beige.....536-5366



81C – Handler

Plaster Bin

- 75lb capacity
- Heavy gauge sheet metal
- Grey enamel finish
- Weight: 20lb



* Please call for pricing.



Alginator Bowls – Dux

Blue

Medium.....161-8065 1@ \$28.00
Large.....161-8107 1@ \$32.35



Dispos-A-Bowl – Zirc

- Disposable, saves wear on mixing bowl
- Capable for use with electric alginate mixers
- Free-standing for use with or without mixing bowl

36/Pkg
White.....771-0312 1@ \$24.00

MIXING BOWLS

PATTERSON



Flex Bowls

Self-cleaning, non-stick mixing bowls are made of rigid yet flexible rubber. Available in 3 sizes.

Medium 350cc	093-6161	1@ \$8.30
		6@ \$8.15
Large 600cc	093-6179	1@ \$8.30
		6@ \$8.15
X-Large 850cc	093-6187	1@ \$10.90
		6@ \$10.60



Mini Mixing Bowl

Non-Stick, Silicone Bowls

For mixing acrylic. Designed with a suction cup bottom for convenient one-handed mixing.

Small, 3/Set.....	088-8917	1@ \$16.40
		3@ \$14.50
Medium, 2/Set	088-8909	1@ \$16.40
		3@ \$14.50
Large, 1/Set	088-8925	1@ \$13.75
		3@ \$12.25



FlexForm™ – BDM

Silicone Flex Mixing Bowls

Made of high quality silicone-based material, these colourful large 600ml FlexForm Mixing Bowls will not only enhance the style of your office and its environment, but make it fun!!! Clear and very flexible, making for fast, easy clean-up, yet tough and durable for mixing virtually every type of material.

Each	1@ \$7.35
Rose	145-6292 3@ \$7.00
Yellow	145-6300
Light Green	145-6318
Sea Blue	145-6326
Light Purple	145-6334



Flexibole Mixing Bowls – Hygenic

Flexible, Self-Cleaning Bowls

Plaster, stone, investment and alginates will not stick to these self-cleaning mixing bowls or to the nylon spatulas (not recommended for acrylic mixes). More flexible than rubber bowls for more efficient mixing.

Small, 150cc.....	377-8180	1@ \$8.85	3@ \$8.70
Medium, 350cc	377-8677	1@ \$11.95	2@ \$11.50
Large, 600cc	377-9162	1@ \$15.00	3@ \$14.60
Extra Large, 850cc	377-9659	1@ \$17.30	3@ \$16.90
Jumbo, 1,700cc	377-9667	1@ \$21.85	2@ \$21.30

SPATULAS AND KNIVES

PATTERSON



Magnabar

Magnetic Knife and Utensil Holder

Grasps objects firmly, yet releases without effort. Wood with magnetic strips. Installs easily. (Instruments shown not included.)

18"	093-6369	1@ \$48.20
-----------	----------	------------




Plaster Nipper

Stainless Steel

Cuts plaster or stone. May be used for trimming models.

Each	093-6781	1@ \$43.15
------------	----------	------------



Plaster Spatulas

Stainless Steel Blades

#2R Flexible.....	089-9690	1@ \$11.95
#3R Stiff.....	089-9724	3@ \$11.15
#4R Rigid	089-9732	
#8R Stiff.....	089-9708	1@ \$12.35
#10R Extra Stiff.....	089-9682	3@ \$11.65
#11R Flexible.....	089-9716	1@ \$12.60
		3@ \$11.95



Steel Spatulas

Stainless with octagonal handle for better control — less slippage when mixing.

#324	092-8184	1@ \$21.00	3@ \$20.30
#24	089-5706	1@ \$6.55	3@ \$6.25



Murphy Knives

Ground from 15-gauge chrome vanadium steel with brown oval hardwood handles.

Murphy 134		
1-3/4" Blade	093-6286	1@ \$11.50
Murphy 258		
2-5/8" Blade.....	093-6294	1@ \$11.20



Buffalo Style Lab Knives

Quality construction with wood handles. Carbon premium quality stainless steel razor-sharp knife blade.

Green Enamel Handle		
#3	088-7414	1@ \$8.00
#6	088-7422	1@ \$9.80
#7	088-7430	1@ \$14.70
#8	088-7448	1@ \$10.35
Rosewood Handle		
#6R.....	088-7455	1@ \$10.25
#7R.....	088-7463	1@ \$11.15
#12R.....	088-7661	1@ \$13.70

SPATULAS AND KNIVES (CONT.)



Form-Fit™ – BDM

Delrin Mixing Spatulas

Made of high quality, durable "Delrin" material, these ergonomically designed mixing spatulas will make any mixing procedure, less stressful and easier on the hands and arms!

Available in five pastel colours!

Each		1@ \$3.55
Cream	145-6573	3@ \$3.40
Rose	145-6581	
Sea Blue	145-6599	
Violet	145-6607	
White	145-6615	

Spatula 8X curved – Dentsply

Each184-1204 1@ \$54.75

Spatula – Miltex

Stainless Steel

All metal spatula for mixing alginate, stone and plaster.

No. 3 Flare 3 1/2"	528-1803	1@ \$19.95
		3@ \$19.15

Spatula Plastic Broad – Zirc

Alginate Spatula

White plastic spatula with broad, medium flex blade ideal for alginates, plaster and impression materials.

Each (50Z501)	771-8950	1@ \$9.00
---------------	----------	-----------



Lab Knife – Buffalo

Carbon Steel Blades

With Green Enameled Wood Handles

#3 – 2-1/2" Blade	154-0442	1@ \$12.50
#4 – 2-3/8" Blade	154-0939	1@ \$14.15
#6 – 1-5/8" Blade	154-1424	1@ \$13.95
#8 – 2" Blade	154-3388	1@ \$15.25

For Compound & Wax

#7 – 1-1/2" Blade	154-2406	1@ \$13.35
-------------------	----------	------------

With Rosewood Handles for Plaster and Miscellaneous

#6R – 1-5/8" Blade	154-1911	1@ \$18.50
#7R – 1-1/2" Blade	154-2893	1@ \$14.75
#12 – 2" x 3/4" Blade	154-3875	1@ \$22.45



Nievert Whittler – Temrex

Knife Sharpener

Each	385-1672	1@ \$77.30
------	----------	------------

CHISELS AND SCRAPERS


PATTERSON



Surgical Scalpel Handles

Stainless steel in a popular shape for comfortable fit. Provides secure blade fit with easy removal.

Metal #3	089-4501	1@ \$3.85
Metal #4	089-4519	3@ \$3.65



Vul-crylic Chisels – Dixon

Hardwood Handles

Can be resharpened.

A. #13 – 6"	292-4868	1@ \$19.25
#14 – 6"	292-4850	1@ \$19.25
(not pictured)		
B. #16 – 4"	292-3837	1@ \$17.00
C. #17 – 4"	292-4322	1@ \$17.70



Plaster Spatulas – Buffalo

Stainless Steel Blades With Rosewood Handles

#14R Tapered Rigid

Silicone and Rubber Base

4 1/4" Blade	157-3278	1@ \$17.35
--------------	----------	------------

#16R Extra-Stiff

Alginates, Stone and Plaster

3 1/2" Blade	157-3963	1@ \$16.75
--------------	----------	------------

#11R Flexible

Alginates, Stone and Plaster

3 1/2" Blade	157-2783	1@ \$17.80
--------------	----------	------------

#10R Extra-Stiff

Stone and Plaster

1 1/2" Blade	157-2296	1@ \$16.95
--------------	----------	------------

#15R Tapered Rigid

Large Mixes of Silicone and Rubber Base

6" Blade	157-3765	1@ \$21.55
----------	----------	------------

#8R Stiff

Plaster and Stone

4" Blade	157-0787	1@ \$17.25
----------	----------	------------

#4R Rigid

Plaster and Stone

4 1/2" Blade	157-0704	1@ \$16.15
--------------	----------	------------

#3R Stiff

Plaster and Stone

3 3/4" Blade	157-0670	1@ \$16.15
--------------	----------	------------


#2R Flexible

Impression Cream and Light Plaster Work

4" Blade	157-0605	1@ \$16.50
----------	----------	------------

Chisels & Scrapers – Buffalo


Carbon Steel Blades with Enameled Wood Handles



Pearson Chisels		
No. 16	150-7169	1@ \$7.50
No. 17	150-7177	
No. 18	150-7185	




Vulcanite Chisels		
No. 20	150-7672	1@ \$9.95
No. 21	150-7631	
No. 22	150-7664	



Kingsley Type Scrapers		
No. 1	156-5183	1@ \$9.25
No. 2	156-5175	1@ \$16.50
No. 4	156-5209	1@ \$13.65
No. 5	156-5217	

Blade Sizing Chart



Scalpel Blades – Bard-Parker

Carbon Steel, Non-Sterile, (6/Pkg)		
#12	129-3885	1@ \$6.35
#23	129-3893	
#24	129-3901	

Carbon Steel, Non-Sterile, (6/Pkg)		
#10	129-3786	1@ \$4.50
#11	129-3794	25@ \$4.25
#15	129-3802	
#20	129-3810	
#21	129-3828	
#22	129-3836	
#25	129-3844	

Disposable Scalpels – Bard-Parker

10/Box		
#10	128-1724	1@ \$15.75
#11	128-1757	
#15	128-1781	
#20	128-1815	
#21	128-1849	
#12B	128-2334	1@ \$36.00
#12	128-2854	1@ \$24.60
#23	128-2888	

Scalpel Handle – Bard-Parker

Metal		
#3	128-6038	1@ \$22.15
#4	128-6061	
Plastic		
#5	128-6582	1@ \$19.50
#6	128-6616	

Scalpel Handle – Hu-Friedy

Metal		
#3	373-7855	1@ \$34.20
Metal		
#5	373-7863	1@ \$85.80

Scalpel Blades – Hu-Friedy

Carbon Steel, Sterile-100/Box		
#12	367-8414	1@ \$138.05
Carbon Steel, Sterile-100/Box		
Double-Edge, Push-Pull Cut		
#12D	367-8422	1@ \$229.65
Carbon Steel, Sterile-100/Box		
#15	367-8430	1@ \$138.05

Scalpel Blade Safe Remover – Miltex

An easy-to-use instrument, for the safe removal of scalpel blades from the handle.

Each	496-3005	1@ \$24.90
		3@ \$23.70

PIN DRILLS

Pindex Accessories – Coltène/Whaledent

PX-114 - Long Pins and Sleeves		
1,000/Pkg	723-9916	1@ \$131.50 3@ \$128.45
PX-112 - Medium Pins and Sleeves		
1,000/Pkg	723-9890	1@ \$131.70 3@ \$126.40
PX-115 - Short Pins and Sleeves		
1,000/Pkg	724-0401	1@ \$129.15 3@ \$123.95
PX-116 - Long Pins and Sleeves		
100/Pkg	724-1383	1@ \$30.75 4@ \$29.85
PX-102 - Medium Pins and Sleeves		
100/Pkg	724-0427	1@ \$31.35
PX-117 - Short Pins and Sleeves		
100/Pkg	724-2365	1@ \$29.70
PX-118 - Articulating Pins		
100/Pkg	724-0898	1@ \$56.75
PX-119 - 3 Carbide Drills		
Pkg	724-0997	1@ \$116.15
PX-120 - 1 Reamer		
Pkg	724-1037	1@ \$46.40
PX-124 - 10 Bulbs		
Pkg	724-1953	1@ \$90.90
PX-153 - Dual Pins and Sleeves		
1,000/Pkg	724-3538	1@ \$381.50
PX-154 - Dual Pins and Sleeves		
100/Pkg	724-3546	1@ \$50.15
PX-144 - Self-Articulating Pins & Sleeves		
1,000/Pkg	724-3363	1@ \$410.50
PX-145 - Self-Articulating Pins & Sleeves		
100/Pkg	724-1052	1@ \$60.05
PX-145S - Self-Articulating Sleeves		
1,000/Pkg	724-1094	1@ \$198.00
PX-146 - 3 Flexible Molds		
Pkg	724-3421	1@ \$105.55
PX-147 - 1 Full Arch Mold		
Each	724-3405	1@ \$59.10
PX-148 - 2 Flexible Molds, Right/Left		
Pkg	724-3413	1@ \$75.90
P-48 - Pinmates		
100/Box	723-9429	1@ \$80.25

PIN DRILLS (CONT.)

Pindex Laser – Coltène/Whaledent

- Offers better visibility in daylight
 - Small laser beam makes drill hole positioning even more accurate
 - Laser has longer life span than regular bulbs
- Each048-8486 1@ \$1,397.50



DOWEL PINS AND ACCESSORIES

PATTERSON



Dowel Pins

Brass dowel pins.
1,000/Pkg.....086-9024 1@ \$38.85
Medium. #2 2@ \$36.10



Stabilizing Pins

Plastic
Micro-rib retention. Super accurate die seat. Helps eliminate problems of poor contact. No glue needed, just seat and pour.
1,000/Pkg.....087-7001 1@ \$28.65

Dowel Pins – Accubite

Ceramic Z
10/Pack101-0214 1@ \$52.20

Dowel Pins – Blitz Micro Turning

1,000/Pack132-6008 1@ \$75.40

Dental Positioning Pins – Blitz Micro Turning

1,000/Pack132-6057 1@ \$42.75



Dowel Pin Indexer – KerrLab

- Eliminates double pinning
- No rotation of dies
- Pin sleeves increase stability, prevent rotation and help identify the position of each model section by colour
- Dowel Pin Indexer kit contains: 1,600 sets of indexers and collars, sorting tray w/clear plastic cover, 1/2oz of dowel pin indexer adhesive, 1oz of indexer adhesive set, 250 indexer dowel pins and 50 adhesive wipes

Kit131-4236 1@ \$178.30

Brass Indexer™ Dowel Pins – KerrLab

Polished finish resists corrosion and oxidation.
1,000/Pkg.....131-4921 1@ \$96.70

ParaPost® Burnout Posts – Coltène/Whaledent

10/Package
P-51-3 (.036 Brown)...723-9015 1@ \$65.00
P-51-4 (.040 Yellow)..723-9023 2@ \$63.05
P-51-4 .5 (.045 Blue)..723-9064
P-51-5 (.050 Red)723-9031
P-51-5.5 (.055 Purple)723-9072
P-51-6 (.060 Black).....723-9049
P-51-7 (.070 Green) ...723-9056



Bi-Pin With Sleeves – Renfert

High luster nickel-plated pin with two guidance legs. Exact fitting metal sleeve with rectangular profile. Guide length for base heights of 17.5mm or 13.5 mm.

1,000/Pkg
Short267-1659 1@ \$296.35
Long.....267-1691 1@ \$267.95

Bi-V-Pin With Sleeves – Renfert

Narrow model pin with two guidance legs. Particularly suitable for narrow positioning in cramped space conditions.
1,000/Pkg.....267-1691 1@ \$267.95



WOW Dowel Pins – Premier Dental

300/Bag.....586-2057 1@ \$16.25



Accu-Trac Precision Die System – Coltène/Whaledent

Ideal for:

- Lower anterior cases
- No pins laminate veneer
- 1 pour technique longspan bridgework
- Excellent stability ceramic crowns
- Superior accuracy rush cases
- Simple fast technique

Introductory Kit – AT-6

Each719-2503 1@ \$191.40

Accessories

Articulation Keepers – ATK-25

Each719-2776 1@ \$248.15

Model Former – Columbia Dentoform

Model Base Formers of Rubber

For making neat, sturdy models easily and quickly.

2/Set
3" x 4" x 1"
#R1 Large213-8659 1@ \$61.40

2½" x 3½" x 1"
#R2 Medium213-8683 1@ \$61.40

3½" x 4½" x 1½"
#R1X Extra-Large.....213-9178 1@ \$64.70

2¼" x 3¼" x 3/4"
#R3 Small.....213-9665 1@ \$58.10

#R24.....213-9707 1@ \$85.80

#R22.....214-0150 1@ \$85.80



SAWS AND BLADES

PATTERSON



Saw Frame

Upper bow of frame is offset to allow a clear view of cutting motion.

Regular Pin 5"093-7102 1@ \$29.15

Accepts 5" pinned saw blades.

Offset Pin 5"093-7094 1@ \$29.60



Saw Blades

Pinned

5" length, .055" wide, .010" thick.

10/Pkg.....093-7052 1@ \$8.05

100/Pkg.....093-7060 1@ \$49.60

Saw Blades – Heraeus Kulzer

12/Pkg

Flat

.007"454-7105 1@ \$31.05

.010"454-7592 3@ \$30.75

Pinned

.007"454-8186 1@ \$42.05

.010"454-8087



Laboratory Saw Kit – Dentsply Neytech

For Hand Cutting of Stone & Acrylic

Contains: frame and 10 blades.

Kit541-2424 1@ \$28.70

Saw Blades – Dentsply Neytech

"Skip Tooth" Blades

10/Pkg.....541-1939 1@ \$13.45

100/Pkg.....541-1947 1@ \$130.60

Pindex Handsaw Kit – Coltène/Whaledent

Kit724-3603 1@ \$125.35

Replacement Blades

PX-201 .007" (Blue Colour) 12/Pkg722-4926 1@ \$30.75

PX-202 .010" (Gold Colour) 12/Pkg722-4926 1@ \$18.80

PX-301 .007" (Blue Colour) 100/Pkg722-4934 1@ \$123.45

PX-302 .010" (Gold Colour) 100/Pkg722-4942 1@ \$106.85



Ribbon Saw Blades – Pfingst

One-Sided Stainless Steel

12/Pkg

4mm.....566-9163 1@ \$22.60

6mm.....566-9130 3@ \$21.50

8mm.....566-9148



Swiss Cut Saw Blades and Frame – Keystone Industries

Made of the finest alloy steel using a unique tempering process for incredible flexibility and a long blade life.

Frame

.010" W x 5" L.....534-7778 1@ \$29.00

Blades

10/Pkg.....534-7794 1@ \$6.75

DIE PREPARATION MATERIALS

PATTERSON



Die Spacer

Consistent film thickness. One even coat provides 12 microns of controlled relief. Durable, tough coating is heat-resistant. Thinner is also available, to remove die spacer without distortion or damage.

1/Pkg

Kit contains 1oz die spacer and 1oz thinner.

Yellow.....093-5213 1@ \$20.30

Blue.....093-5221

1oz Die Spacer Refill Bottle

Yellow.....093-5296 1@ \$18.85

Blue.....093-5304

Die Spacer Thinner

1oz Bottle.....093-5379 1@ \$6.75



**Aron Alpha® – Taub
Cyanoacrylate Adhesive**

201 Type

Light viscosity for coating and sealing dies, fixing articulation systems, assembling crowns and bridgework for soldering.

5/2gm.....707-9460 1@ \$24.50

Die Hardener – Keystone Industries

A high performance plaster treatment and not comparable to die spacers. Extremely low viscosity liquid soaks into plaster and cures quickly.

20ml.....534-2514 1@ \$33.20

250ml.....534-2522 1@ \$222.45

DIE PREPARATION MATERIALS (CONT.)

Dry Teflon Spray – Handler

Each357-2203 1@ \$19.85



Fill-It Kit – American Dental Supply

A second generation cyanoacrylate system for repair, gap-filling, tacking, hardening surfaces, and bonding porous and dissimilar materials. Kit contains: Fill-it base and pinpoint applicators, 2oz Fill-it activator, 1oz Fill-it solvent. Economy kit contains: 61/2oz base and pinpoint applicators, 2-2oz activator spray, 1-8oz activator refill, and 1-1oz clean solvent.

Kit 109-7500 1@ \$28.60

Fill-It Base Refill – American Dental Supply

Activator
 2oz 109-7534 1@ \$9.65
 8oz 109-7567 1@ \$19.15
 32oz 109-7575 1@ \$69.90

Foilcoat – Whip Mix

Liquid Foil Substitute

When applied to plaster or stone, it seals the pores and forms a smooth surface that gives excellent separation of denture from the mold.

Gallon 727-6801 1@ \$68.00



Gator Die Lube – Whip Mix

- Provides a long-lasting and uniformly thin release film

Bottle w/Brush Cap
 1/2oz 727-7122 1@ \$12.55



Gi-Mask – Coltène/Whaledent

Silicone-Based

Gingival mask material for perfect aesthetics. Precise reproduction of gingival contours. For more restorations that fit the first time. Functional design of interdental spaces. Excellent dimensional stability. Quick, simple, reliable application. Kit contains 150ml each base and activator, 90ml separator, accessories.

Gi-Mask Starter Kit
 Kit 718-1571 1@ \$211.20

Gi-Mask Single Pack
 Contains base, activator and accessories.
 Each (168ml) 718-1563 1@ \$154.45

Gi-Mask Single Pack Refill
 Contains base and activator.
 Each (168ml) 718-1589 1@ \$142.25

Gi-Mask Separator
 1-90ml Bottle 718-1597 1@ \$41.60

Gi-Mask Automix Introductory Kit
 Contains: 1 - 75ml Gi-Mask cartridge, Gi-Mask separator, mixing tips, application tips, plunger, and jet dispenser.
 Kit 718-1282 1@ \$144.90

Gi-Mask Automix Single Pack
 Contains: 1 - 75ml Gi-Mask cartridge, mixing tips, and application tips.
 Each 718-1290 1@ \$57.75

Micro Mixing Tips
 40/Pkg 718-2371 1@ \$40.25
 3@ \$38.95

Gingifast – Zhermack

Gingival Reproduction Material

- Mixing ratio Base: Catalyst 1:1
- Total working time (73° F/23° C) 2'
- Setting time (73° F/23° C) 10'
- Hardness (Shore-A after 1h) 40

Contains 2 cartridges 50ml each, 1 bottle, separator 10ml, 12 mixing tips yellow and 12 intraoral tips yellow.

Gingifast Elastic
 2 - 50ml 770-3408 1@ \$68.85

Gingifast Rigid
 2 - 50ml 770-3416

GI-Mask New Formula – Coltène/Whaledent

- Precise reproduction of soft tissue contours for crown, bridge and combination case work
- Economical system that offers accurate mixing and dispensing and simple application
- Accurate colouring for a lifelike appearance
- Short, timesaving, setting times
- Delivers a positive appearance on the model
- Packaged in 50ml automix cartridges
- Self-leveling

GI-Mask Automix NF Starter Kit
 Contains 2 - 50ml cartridges, 30ml separator, 12 mixing tips (yellow), 12 oral tips (yellow)
 Kit 718-6471 1@ \$123.75

GI-Mask Automix NF Refill Kit
 Contains 2 - 50ml cartridges, 12 mixing tips (yellow), 12 oral tips (yellow)
 Kit 718-6489 1@ \$112.95
 3@ \$110.30

GI-Mask Separator
 30ml Bottle 718-6414 1@ \$32.35

GI-Mask Silicone Trimmer
 Each 718-6422 1@ \$25.45



Softtissue Moulage™ – KerrLab

The natural selection for replication of the soft gingival tissue on a stone cast in the laboratory. It is ideal for use in the fabrication of implants to determine emergence profiles, to check the fit of abutments on the implants and for all gingival simulation.

Contains 3 - 50ml cartridges, 24 mixing tips, separator medium, extruder and 24 intraoral tips
 Intro Kit 435-5004 1@ \$224.40

Softtissue Moulage Refill
 2/Pkg 435-5012 1@ \$68.70

Mixing Tips – Large
 48/Pkg 428-9039 1@ \$62.95

Intraoral Syringe Tips
 100/Pkg 428-8023 1@ \$61.30

Gypsum Hardener – Whip Mix

Replaces water as the gauging liquid for mixing with gypsum materials, where increased hardness and reduced brittleness are desired.



Gallon.....727-7189 1@ \$60.40
4@ \$58.40
Pint727-7676 1@ \$16.30

Lab Separator Extraoral – Coltène/Whaledent

20ml.....718-1985 1@ \$116.50

Lubritex No. 12 Die Lubricant – Whip Mix

- Water-soluble lubricant that is equally effective on all types of die materials
- Will not react w/die or pattern materials and will not affect investments



Pint728-5026 1@ \$47.20

Margin Liner – KerrLab

Refill 36/Pkg
Red.....131-4376 1@ \$20.20
Blue.....131-4384

Microfilm – KerrLab

Die Lubricant
For easy separation of wax patterns.
8oz Bottle.....428-5409 1@ \$33.85
32oz Bottle.....428-5417 1@ \$54.10

Model Gloss – Columbia Dentoform

Adds a smooth, lustrous, washable finish to stone or plaster models.
Quart214-1133 1@ \$35.00



Permabond 910 – BDM

General-purpose metal bonding adhesive with cyano-acrylate.
1oz Bottle.....145-7241 1@ \$32.95
6@ \$31.25
3 g Tube145-6763 1@ \$5.00



Potassium Sulfate – Dux Fixing Solution (K₂SO₄)

Improves stone hardness and keeps moisture from evaporating from the impression.
3/4lb Pkg704-6501 1@ \$22.65

Silicone Spray – Dentsply Trubyte All-Purpose Release Agent

For use on metals, many plastics, plasters, stone, amalgams and compounds. Contains no petroleum oil.
16oz Can276-5774 1@ \$28.30

Silicone Spray – T&S Dental Multi-Purpose Lubricant

16oz Can678-4151 1@ \$19.50



Stalite Special Solution – Buffalo Gypsum Hardener

Replaces water in gypsum mix for consistently stronger, harder, smoother, abrasion-resistant stone models and dies.
Gallon157-6891 1@ \$32.95

Stone Die & Plaster Hardener – Taub

- Seals and hardens stone
- Sets in 3-5 minutes
- A liquid plastic applied to any stone die, plaster, or stone model
- Protects margins and reduces abrasion
- Includes resin and thinner

Kits
2oz670-9570 1@ \$37.75



Strip-It – Keystone Industries

Acid Substitute
Gallon534-9881 1@ \$51.75



Super-Sep™ – KerrLab

An alcohol-based gypsum separating fluid that leaves no film and dries on contact. Plaster, stone, or investments can be poured against it immediately after application.
8oz437-6455 1@ \$42.20
32oz437-6463 1@ \$79.65



Tru-Fit Kit – Taub

- Die spacer and visual aid kit
- Controlled technique for relieving dies before wax up
- Use on stone, amalgam, epoxy and plated dies
- Contains 1/2oz gold and silver colours, plus 1oz thinner

Kit671-1048 1@ \$37.30
Thinner
2oz671-1527 1@ \$11.80

DIE PREPARATION MATERIALS (CONT.)



Zapit – Dental Ventures

Die Prep Problem Solver

Zapit consists of a special cyanoacrylate material of base and accelerator allowing a controlled setting time. Can be manipulated, carved and shaped into an appropriate design and then caused to harden by spraying with a slight amount of accelerator.


Introductory Kit

Contains: base with two tips, accelerator with pump, solvent.

Kit	279-7595	1@ \$44.10
Base	279-7579	1@ \$24.95
Accelerator w/Pump	279-7587	1@ \$16.90
Solvent.....	279-7561	1@ \$10.35
Accelerator, Refill	279-7553	1@ \$70.75
Zapit Accelerator Pens.....	279-7504	1@ \$27.70

ARTICULATORS AND ACCESSORIES

PATTERSON



Articulators

Crown and bridge articulator is popular Foster design, available with or without springs to allow protrusive and retrusive movements. Also available in full arch denture design. Chrome-plated construction at an economical price. Or choose the disposable clear, plastic model — eliminates the need to break off models.

A. Crown and Bridge	093-6005	1@ \$9.45
		6@ \$8.80
B. Denture High Arch.....	093-6039	1@ \$11.35
		6@ \$10.05
C. Spring.....	093-6096	1@ \$10.25
		6@ \$10.00
50/Pkg		
D. Disposable Plastic	083-8151	1@ \$40.95
		3@ \$37.45

Articulator Mounting Plates – Whip Mix

For Whip Mix Articulator Systems

Metal – Standard

Durable and reusable.

2/Set	725-2695	1@ \$15.50
-------------	----------	------------

Plastic

An economical alternative to metal mounting plates.

20/Pkg	725-2711	1@ \$60.40
100/Pkg	725-2729	1@ \$236.20



Articulator and Flask Cleaner


All-in-one, general-purpose cleaner removes plaster, stone and investments without ultrasonics. Heavy-duty concentrate penetrates joints, crevices and delicate areas. Rinses clean with water, no harsh chemicals.

Gallon	093-5007	1@ \$52.85
--------------	----------	------------

Blockout Compound

When adapting baseplates, use this compound to block out remaining teeth and undercuts on models.

Pkg.....	091-0893	1@ \$12.50
----------	----------	------------



Articulators – Foster

Apex® #2

Bronze

Full-Arch Articulator

A #50A.....	324-4522	1@ \$23.75
		3@ \$22.45

Foster® #8

Bronze

Denture Articulator

B #51A.....	324-5503	1@ \$25.90
-------------	----------	------------

Foster® #3

Bronze

Denture Articulator with Anterior Pin

C #54.....	324-5651	1@ \$32.70
------------	----------	------------

Ask your Patterson representative!

We carry many more laboratory products

Articulators – Shofu

Handy IIM
(without Occlusal Plate)
Each634-3826 1@ \$169.65



Handy IIM
(with Occlusal Plate)
Each634-3842 1@ \$194.05

Shofu Metal Mounting Plates #5072 for Handy IIM and Handy IIIM
Each634-6050 1@ \$26.35

Articulators WOW – Premier

Disposable
Intro kit: 20 anterior/30 posterior
Kit586-2008 1@ \$90.60

Kit anterior (50)
Kit586-2016 1@ \$81.10

Kit posterior (50)
Kit586-2032 1@ \$71.30

Deluxe Magnetic Denture Articulator – Keystone Industries

- Interchangeable
- Incisal pin included
- Lightweight and Teflon coated
- Magnetic split casting ensures easy, removable and exact setting
- Lateral excursions

Each534-0849 1@ \$66.35



Curved Discs 6/Pkg
Each534-3223 1@ \$1.50



Disposable Plastic Mounting Plates – Whip Mix

Designed to hold models securely in place. Metal insert minimizes cross threading and facilitates ease of use. Fits both HANAU and DENAR articulators.

10/Pack337-2091 1@ \$49.55

100/Pack337-2083 1@ \$252.00



Galetti Articulator – Keystone Industries

This quality articulator provides for the easy, speedy, firm grasp of models of any size, controlling the correct occlusion with maximum precision. Neither plaster work nor accessories required. Easy to clean.

Each131-4210 1@ \$580.50

Hanau® 96H2 Articulator – Whip Mix

- Semi-adjustable, non-arcon type articulator
- It has enclosed guidances; fixed intercondylar width of 110mm and lateral adjustments from 0° to 30° and protrusive adjustments from 0° to 75°
- It accepts a wide range of accessories from the Long Centric Incisal Guidance or Acrylic Guidance for those wishing greater freedom in centric or custom incisal fabrication to protrusive and lateral excursions
- Comes with black TFE finish and Orbitale indicator standard (additional modifications may be specified)



Standard Equipment Includes:

- 2 pair metal mounting plates (#015010-000)
 - 1 extension pin (#005243-000)
 - 1 standard incisal pin (#003330-000)
 - 1 standard adjustable incisal guide
 - 1 Orbitale indicator (#010046-000)
 - 1 technique manual
- Each354-6801 1@ \$1,411.40

Mounting Plates – Whip Mix

Metal
2/Set354-6843 1@ \$22.40



Trubyte® Fox Occlusal Plane Plate – Dentsply

Chairside Aid
Each276-1195 1@ \$39.45

Stephan Articulator – Dentsply Model “A”



Features include: ease of access to all surfaces, ease of which all adjustments can be made with fingers, and ease with which lateral motion is obtained. Can be taken apart and assembled quickly without any difficulty.
Each229-5764 1@ \$37.65
6@ \$36.40

W.O.W™ and W.O.W. 2 Articulators – Premier Dental Products

Self-Use Articulators

- Clear plastic view allows for accurate pin placement
- Pin placement design eliminates rotation
- Shallow retention grooves minimize dimensional change from stone expansion
- Fixed hinge allows for positive occlusal relationship between upper and lower models and appropriate lateral movement
- Remount bite, if necessary, with heat in less than a minute
- W.O.W. 2 features Flash Air Dam technology

Intro Kit

Contains 20 full arch, 30 posterior and 300 dowel pins
Each586-2008 1@ \$90.60

Kits

Contains 50 articulators and 300 dowel pins
Full Arch586-2016 1@ \$81.10
W.O.W.2 Full Arch586-2024 1@ \$81.10
Posterior586-2032 1@ \$71.30
W.O.W.2 Posterior586-2040 1@ \$71.30

Dowel Pins

300/Pkg586-2057 1@ \$16.90

DENTURE RESINS

High Impact-45 - Lang

1lb P + L		
Each	458-6558	1@ \$67.95
5lb P + L		
Each	458-6525	1@ \$315.70
25lb Powder only		
Each	458-6665	1@ \$1,149.45
5lb Powder Only		
Each	458-6590	1@ \$221.70
32oz Liquid Only		
Each	458-6707	1@ \$43.30



Lucitone® 199 - Dentsply

Denture Base Resin

120-Unit Package

Contains: 4-630gm powder, 3-430ml liquid, accessories.

Original	276-0890	1@ \$451.95
----------------	----------	-------------

30-Unit Powder Only

Contains: 630gm powder.

Original	276-0957	1@ \$132.50
Light	276-0924	2@ \$130.50
Lt. Reddish Pink	276-0932	
Dark Pink	276-0940	

540-Unit Powder Only

Contains: 25lb powder.

Original	276-0965	1@ \$824.00
Light	276-0999	
Lt. Reddish Pink	276-1005	

Lucitone® Liquid

430ml Bottle	275-9934	1@ \$43.25
		3@ \$42.35
Gallon	276-0189	1@ \$220.00



Lucitone® Clear Resin - Dentsply

30-Unit Powder Only

Contains: 630gm powder.

Each	275-9009	1@ \$74.40
------------	----------	------------

Contains: 4lb powder.

Each	275-9017	1@ \$143.75
------------	----------	-------------

Lucitone® FRS™ Flexible Dental Resin - Dentsply

Contains: 12 tubes, heat seal bags, instructions.

Original	276-1575	1@ \$207.25
Lt. Pink	276-1583	
Lt. Reddish Pink	276-1591	
Dark Pink	276-1609	
Clear	276-1617	

Lucitone® FRS™ Success® Starter Kit - Dentsply

Kit contains: Lucitone® FRS™ 12-pack original shade, 12-pack clear, cartridge furnace, safety gloves, flask heat lamp, injection insert, cartridge sleeve, knockout rod, knockout base, timer, video, technique and manual.

Each	276-9701	1@ \$2,520.95
------------	----------	---------------



Lucitone® Characterized - Dentsply

Denture Base Resin

10-Unit Package

Contains: 10-21gm powders, 120ml liquid, 10gm gingival tone powder and accessories.

Light	276-0015	1@ \$143.95
Lt. Reddish Pink	276-0023	
Reddish Pink	276-0031	

50-Unit Package

Contains: 50-21gm powders, 430ml liquid, 120ml liquid, 2-10gm gingival tone powders and accessories.

Light	276-0064	1@ \$493.00
Lt. Reddish Pink	276-0072	

250-Unit Package

Contains: 250-21gm powders, 10-10gm gingival tone powders and accessories.

Light	276-0122	1@ \$2,600.00
Lt. Reddish Pink	276-0130	
Reddish Pink	276-0148	
Bluish Pink	276-0155	

Characterized Lucitone Gingival Tone Powder - 10gm

Self-Curing	276-0205	1@ \$14.80
Heat-Curing	276-0213	1@ \$18.10

Lucitone® Hy Pro Denture Base Resin - Dentsply

11-Unit Package

Contains: 231gm powder, 120ml liquid and accessories.

Pink	276-0361	1@ \$55.40
Fibered	276-0379	
Fibered Light	276-0387	
Fibered Dark	276-0395	

120-Unit Package

Contains: 4-630gm powders, 3-430ml liquid and accessories.

Pink	276-0429	1@ \$403.00
Fibered	276-0437	
Fibered Light	276-0445	
Fibered Dark	276-0452	

30-Unit Powder Only

Contains: 630gm powder.

Pink	276-0528	1@ \$86.95
Fibered	276-0536	
Fibered Light	276-0544	
Fibered Dark	276-0551	

540-Unit Powder Only

Contains: 25lb powder.

Pink	276-0627	1@ \$1,030.00
Fibered	276-0635	
Fibered Light	276-0643	
Fibered Dark	276-0650	



Microlon Denture Resin - Hygenic

A premium denture material that responds beautifully to characterization.

Powder & Liquid

1lb Pkg

Veined	379-3593	1@ \$73.25
Special Veined	379-3601	
Heavy Veined	379-3619	
Clear	379-3635	

Powder & Liquid

5lb Pkg

Veined	379-3759	1@ \$274.55
Special Veined	379-3783	
Heavy Veined	379-3817	
Clear	379-3874	

Powder Only

1lb Pkg

Veined	379-3932	1@ \$63.05
Special Veined	379-3940	
Heavy Veined	379-3957	
Clear	379-3973	

Orthodontic Resin – Dentsply/Caulk

Orthodontic Appliances

For fabrication of maintainers or retainers. Material never runs or slumps.

Regular Package

Contains: 605gm powder, 430ml liquid and liquid dispenser.

Pink.....	179-8750	1@ \$160.00
Clear	179-8784	2@ \$157.00

Powder Refill Only 1@ \$108.00

Contains: 16oz powder. 3@ \$106.00
Pkg..... 179-9824

Liquid Refill Only 1@ \$52.25

12oz Bottle 3@ \$50.95
Pink..... 180-0341
Clear

Orthoresin – DETREY/Dentsply

Powder
3 Kg 180-0408 1@ \$336.00

Liquid Clear
1 Litre 180-0424 1@ \$109.00

Liquid Pink
1 Litre 180-0432 1@ \$109.00

Paragon Denture Resin – Hygenic

Economically priced resin featuring uniform shade and working properties.

Powder & Liquid
1lb Pkg

Veined	379-5432	1@ \$57.55
Special Veined	379-5465	
Heavy Veined.....	379-5499	
Clear	379-5523	

Powder & Liquid
5lb Pkg

Veined	379-6042	1@ \$221.10
Special Veined	379-6059	
Heavy Veined.....	379-6067	
Clear	379-6075	

Liquid Only
Quart 379-6505 1@ \$42.25

Perma-Cryl CF – GC America

Cadmium Free

- Economical cross-linked resin
- Cure 20 minutes at high temp & 8 hours at lower temperatures
- Easily finished & polished

Powder Only - 5lb

Light Fibered Plus.....	208-3574	1@ \$156.00
Clear	208-3616	
Ethnic Mild.....	208-3624	
Ethnic Heavy	208-3640	

Powder Only - 25lb

Light Fibered Plus.....	208-3863	1@ \$568.50
Clear	208-3905	

Pour-N-Cure™ – GC America, Inc.

Self-Cure Denture Base and Repair

Pour type resin for complete dentures. May be used with alginate or hydrocolloid. Introductory pack contains: 1lb powder, 12oz liquid, 2oz COE-Sep, accessories.

Laboratory Package
Contains: 5lbs powder, 2qts liquid

Pink Fibered.....	208-8078	1@ \$263.50
Light Fibered +	208-8045	

25lb

Pink Fibered.....	208-8474	1@ \$599.50
Light	208-8417	
Light Fibered +	208-8425	
Clear	208-8441	

POURABLE RESINS

Lucitone® Fas-Por™ E Pourable Denture Resin – Dentsply

32-Unit Powder Only
Contains: 640gm Powder
Light 276-1229 1@ \$60.95

Lucitone® Fas-Por™ E Liquid/Repair Material
175ml..... 276-5154 1@ \$30.55
946ml..... 276-5105 1@ \$78.00

Lucitone® Fas-Por™ + Premium Pourable Denture Resin – Dentsply

32-Unit Powder Only
Contains: 640gm powder.
Original 276-1427 1@ \$96.50
Light 276-1435
Lt. Reddish Pink 276-1443

125-Unit Powder Only
Contains: 5.5lbs powder.
Original 276-1518 1@ \$278.00
Light 276-1526
Lt. Reddish Pink 276-1534
Dark Pink..... 276-1542

Lucitone® Fas-Por™ + Liquid
Contains: 1qt bottle.
Each 276-1187 1@ \$95.95

www.pattersondental.ca

PATTERSON DENTAL/DENTAIRE

479

REPAIR MATERIALS



Jet Repair Acrylic – Lang

Self Cure Denture Repair Acrylic

Crosslinked and colourfast for easy, reliable denture repairs.

Powder & Liquid, 1lb Package	1@ \$91.10
Pink	459-2630
Fibered Pink	459-2614
Clear	459-2648

Powder Only, 4oz Bottle	1@ \$33.55
Pink	459-2382
Fibered Pink	459-2390
Translux Pink	459-2408
Clear	459-2416

Powder Only, 1lb Bottle	1@ \$68.65
Pink	459-2903
Fibered Pink	459-2994
Translux Pink	459-3000
Clear	459-2879

Liquid Only, 120ml Bottle	
120ml Bottle	1@ \$26.40
240ml Bottle	1@ \$39.95
Quart	1@ \$80.65



Lucitone 199® Repair Material – Dentsply

Contains: 1lb matching material in all four aesthetic Lucitone 199® acrylic shades.

Original	276-1021	1@ \$86.75
Light	276-1070	
Lt. Reddish Pink	276-1088	
Dark Pink	276-1096	

Powder Only – 4lb		
Pink	276-5055	1@ \$223.25
Pink Free-Flow	276-5063	
Pink Fibered	276-5071	
Light Fibered	276-5089	
Original	276-0809	1@ \$288.00



Repair Material – Dentsply

Self-Cure Denture Repair Acrylic

Provides an accurate, stable colour match and assures a strong, long-lasting repair.

Starter Kit

Contains: 21gm each of pink free-flow, pink fibered and light fibered; 58ml liquid, 65ml Al-Cote separator and 1 brush.
Kit276-5162 **1@ \$84.60**

Labpak Kit

Contains: 454gm powder, 21gm pink free-flow, 2-180ml liquid and accessories.
Pink276-4934 **1@ \$157.50**
Pink Free-Flow276-4942
Pink Fibered276-4959
Light Fibered276-4967

Labpak Refill Kit

Contains: 454gm powder, 2-180ml liquid.
Pink276-4876 **1@ \$128.00**
Pink Free-Flow276-4884
Pink Fibered276-4892
Light Fibered276-4900

Lab Powder Only

Contains: 454gm powder.
Pink276-4991 **1@ \$84.35**
Pink Free-Flow276-5006
Pink Fibered276-5014
Light Fibered276-5022

Repair Material Liquid Only

180ml Bottle276-5154 **1@ \$30.55**
946ml Bottle276-5105 **1@ \$78.00**



Perm Reline & Repair Resin – Hygenic

A superfine, non-slumping, multipurpose resin formulated for exacting repairs.

1lb Pkg	
Lt Fibered	379-8840 1@ \$83.85
Lt Reddish Pink	379-8857
Clear	379-8824
Pink	379-8642
Veined	379-8709
Special Veined	379-8733
Heavy Veined	379-8766

5lb Pkg	
Lt Fibered	379-8956 1@ \$311.55
Veined	379-8923
Special Veined	379-8949
Heavy Veined	379-8915

Perm Liquid Only	
1qt Pkg	380-2394 1@ \$54.00



TruRepair™ – Bosworth

Fast Bonding Repair Material

A cross-linked, colour stable, high strength bonding repair material. A special internally pigmented spherical globe of powder will blend with denture base materials.

Standard Kit

Includes: 6oz pink powder, 6oz veined powder, 8oz liquid and liquid dispenser.
Each142-9885 **1@ \$66.00**

RELINE MATERIALS (CONT.)



Molloplast-B® – Buffalo Dental Denture Reliner

Lab processed silicone soft reline material. Microwaveable consistency offers updated fast curing and added plaque inhibitors. Convenient lab procedure for use with existing dentures or during construction of new dentures.

Regular Size		
45gm Jar.....	155-0045	1@ \$148.95
Lab Size		
170gm Jar.....	155-0052	1@ \$341.00

Mollosil - BDM

Non-hardening, self-curing, silicone based material for direct rebasing in the mouth.

Pkg.....	133-1164	1@ \$165.95
----------	----------	-------------

Mollosil® Plus – Buffalo Dental Soft Reline Material

Long-term soft lining material for chairside use, silicone based, cold curing.

Complete Kit (Automix I)		
Each.....	133-1198	1@ \$245.00
Primer 5ml		
Each.....	133-1123	1@ \$31.50
50ml Cartridge		
Each.....	133-1115	1@ \$223.95
Complete Kit (Automix II)		
Each.....	133-2212	1@ \$319.00
Refill Cartridge		
Each.....	133-2220	1@ \$207.95



PermaSoft® Denture Liner – DENTSPLY

Quick Cure Chairside and Laboratory Denture Liner

Non-methyl methacrylate, self-bonding to acrylic surfaces, 15-minute cure in hot water, resists staining, cracking and hardening; remains soft for 18 months and longer.

Single Kit (60gm)

Contains: 60gm powder, 60ml liquid, 30ml sealer, and accessories.
Pink.....545-3006 1@ \$214.00

Double Kit (120gm)

Contains: 120gm powder, 120ml liquid, 30ml sealer, and accessories.
Pink.....545-3105 1@ \$346.00

Primo – Buffalo Dental

Contains: 15ml bottle and brush
Each.....156-0390 1@ \$40.95

Rebase II – Tokuyama

Chairside Hard Denture Reline Material

- Lower chemical irritation
- Lower heat generation
- Durability
- Easier polishing
- One-appointment procedure
- Accurate adaptation of the denture Kit

Contains 1 - 80g powder, 1 - 50ml liquid, 1 - 15ml adhesive, 1 measuring cup, 1 dropper, 1 rubber cup, 1 spatula, 1 plastic cup, 1 brush, 1 Tokuso Resin Hardener II and 1 spoon
Each.....674-9816 1@ \$133.35

Powder

80g Bottle.....674-9832 1@ \$76.30

Liquid

50ml bottle.....674-9824 1@ \$72.35

Resin Hardener

48/g Bottle.....675-0087 1@ \$23.85

Adhesive

15ml Bottle.....674-9840 1@ \$28.40



SOFRELINER-S – Tokuyama

SOFRELINER technology improves upon TSR to provide more patient comfort. The new delivery system reduces wasted material.

Starter Kit

Contains: same as SOFRELINER-S Combination Kit, minus TD II Mixing Dispenser.
12 Full Lower Relines...675-0210 1@ \$223.10

SOFRELINER-S Refill Cartridge

54gm x 1675-0228 1@ \$201.95

SOFRELINER-S Twin-Pack

54gm x 2675-0236 1@ \$296.95

SOFRELINER Primer

10ml x 1.....675-0244 1@ \$53.80

SOFRELINER-MS – Tokuyama

Combination Kit

Contains: TD II Mixing Dispenser 1, SOFRELINER-MS Cartridge (54gm) 1, SOFRELINER Primer (10ml), TD II Mixing Tip (size S) 5, TD II Mixing Tip (size XS) 5, Plastic Cup 2, Disposable Brush 2, Spatula 1, Steel Bur (Coarse) 1, Finishing Wheel 1.

Starter Kit

Contains: same as SOFRELINER-MS Combination Kit, minus TD II Mixing Dispenser.
12 Full Lower Relines.675-0319 1@ \$223.10

SOFRELINER-MS Refill Cartridge

54gm x 1675-0327 1@ \$203.25

SOFRELINER-MS Twin-Pack

54gm x 2675-0335 1@ \$297.70

Sofreliner Tough – Tokuyama

Soft Denture Reline Material

An addition-cured silicone chairside soft lining material for dentures. Has no odour, taste and chemical irritations. It offers a high tear resistance, high stain resistance and outstanding adhesion to the denture base. Designed to be used when durability is the primary concern. The paste cures in five minutes intraorally, and the kit comes with two different finishing points.

Kit

Contains paste/base paste (27g), catalyst paste (27g), 1 - 10ml primer, 10 TD II mixing tips XS, finishing point (brown), shape adjustment point (white), diamond round bur, spatula, drip-cap for primer, 2 plastic cup, 2 brush no. 5 and marking pen
Each.....675-0152 1@ \$241.50

Softline – Kerr

- Offers ease of use
- Cures in the mouth in under three minutes
- Versatile and easy to colour match
- Automix cartridge delivers an accurate and pre-measured mixture each time

Intro Kit

Contains 1 - 50ml cartridge relining material, 6 mix tips, 1 - 5ml bottle bonding agent, 1 - 5ml bottle varnish/base, 1 - 5ml bottle varnish/catalyst, 25 four-cavity mixing wells and 38 applicators

Light Pink435-5038 1@ \$247.70
 Dark Pink.....435-5046

Refill

Contains 1 - 50ml cartridge relining material, 6 mix tips, 1 - 5ml bottle bonding agent, 25 four-cavity mixing wells and 38 applicators

Light Pink425-5212 1@ \$201.75
 Dark Pink.....425-5220
 Clear Pink.....425-5238



Softone – Bosworth

Tissue Conditioner

Especially developed as a treatment for toning of the tissue and generating a functional impression.

Standard Package

Contains: 3oz powder, 4oz liquid and plastic mixing cup.

Pkg.....140-7725 1@ \$77.50

Special Bulk Package

Contains: 1lb powder and liquid. (Powder packaged in 2-8oz bottles.)

Pkg.....141-0174 1@ \$206.50

Powder Refill Only

3oz Bottle.....140-8210 1@ \$44.90

Powder Refill Only

16oz Pkg.....140-9192 1@ \$163.05

Liquid Refill Only

4oz Bottle.....140-9689 1@ \$40.65

8oz Bottle.....141-0661 1@ \$71.30

Super-Soft™ – GC America, Inc.

Denture Reliner

A heat-processed permanent liner.

Professional Package

Contains: 6oz powder and 6oz liquid.

Each206-6033 1@ \$110.45



Tempo - Lang

Tissue conditioner and reline material

Prof. Pkg Each 1@ \$79.90

Clear460-1100

Pink.....460-1076

White.....460-1043

1 lb Pkg Each 1@ \$150.75

Clear460-1183

Pink.....460-1175

White.....460-1167



Truliner – Bosworth

Denture Reliner

Original Formula Truliner

Has no sting, colour stable, and remains aesthetic indefinitely.

Standard Package

Contains: 8oz powder and 6oz liquid.

Pink.....143-2749 1@ \$116.20

Fibered143-2772

Clear143-2806



Truliner New Formula – Bosworth

Completely bland, no burn, no sting formula for complete patient acceptance.

Standard Package

Contains: 8oz powder and 6oz liquid.

Pink.....143-6534 1@ \$112.00

Fibered143-6567

Clear143-6591

Trusoft – Bosworth

Intermediary Denture Reliner

An intermediary liner for prolonged adjustment periods. Will remain soft for a minimum of 12 months. Its elasticity will ensure against pressure pain, tender ridges, sensitive areas, and irritated tuberosities.

Standard Package

Contains: 3oz pink powder, 4oz liquid, 1/3oz bonding liquid, powder measure, liquid measure and mixing cups.

Each143-9041 1@ \$73.95

3@ \$71.95

Bulk Package

Contains: 2-8oz pink powder, 8oz liquid, 1/3oz bonding liquid, powder measure, liquid measure and mixing cups.

Each143-9538 1@ \$199.35

Powder Refill Only

3oz Bottle.....144-0023 1@ \$44.30

1lb Bottle144-1989 1@ \$163.15

Liquid Refill Only

4oz Bottle.....144-0510 1@ \$33.00

8oz Bottle.....144-1005 1@ \$65.35



Ufi Gel Hard (Hand Mix) – VOCO

Powder / Liquid - Hand Mix

Ufi Gel Hard Kit

Includes 1-60gm powder, 40ml liquid, 20ml conditioner and accessories.

Intro Kit.....706-1328 1@ \$156.50

Ufi Gel Hard Refill

60ml Powder Refill....711-9027 1@ \$74.70

3@ \$73.40

40ml Liquid Refill711-9019 1@ \$52.75

20ml Conditioner711-9001 1@ \$28.80



RELINE MATERIALS (CONT.)

Ufi Gel P - VOCO

Hand Mix - Paste to Paste

Ufi Gel P Kit

Contains: 50ml base, 50ml catalyst, 10ml ea. glazing base/catalyst, 10ml adhesive and accessories.

Kit706-1302 1@ \$160.95

Ufi Gel P Refills

50ml Each Base

and Catalyst706-1310 1@ \$125.50

2@ \$120.00

10ml Each Glazing

Base/Catalyst712-0017 1@ \$34.75



Ufi Gel Hard C Intro Kit with Dispenser - Voco

Includes: 80gm SC, 10ml adhesive, 10 mixing tips (type 13), dispenser (type 4).

Each706-0767 1@ \$251.00

Without Dispenser

Each706-1906 1@ \$167.00

Refills

80gm SC and

10 Tips (Type 13).....706-1914 1@ \$146.95

10ml Adhesive706-1922 1@ \$37.75

Dispenser

(Type 4).....709-3487 1@ \$124.95



Ufi Gel SC - VOCO

A silicone-based permanently soft denture relining material.

- Complete chairside reline in minutes, in one appointment
- High adhesion between denture and relining
- Aesthetic, stable colour
- Available in 2 automix styles (SC/safety cartridge and standard) and a hand mix (paste to paste)

Versa-Soft Denture Liner - Sultan Chemists

Soft, Resilient
Denture Material

- Adhesive, longest lasting automix chairside permanent denture relining material
- Easy-to-use, quick setting, durable and remains soft indefinitely
- Adheres tenaciously to acrylic resin
- Methyl methacrylate-free and very comfortable for denture wearers
- Perfect for accurate indirect relines, exhibiting negligible curing shrinkage
- Long-lasting, with a life-span of 2 years or more!

1-25ml cartridge, applicator gun, 8 tips, primer, varnish, brushes, trimmer
Introductory Kit.....666-3702 1@ \$215.00



Visco-gel - Dentsply

Tissue conditioner and functional impression material

Kit119-6849 1@ \$139.95

3@ \$136.00

Automix - SC (Safety Cartridge)

Ufi Gel SC (Safety Cartridge) Kit

Contains: 50ml SC cartridge, 10ml adhesive, 10ml ea. of glazing base/catalyst, 10 mixing tips (type 8) and accessories.

Kit706-1955 1@ \$232.75

Ufi SC Refills

50ml SC Cartridge

and 10 Tips.....706-1963 1@ \$142.50

50 Mix Tips, Type 8.....706-0551 1@ \$62.95

3@ \$60.95

Ufi Gel SC and C Refills

10ml Adhesive706-1286 1@ \$45.20

20ml Glazing706-1237 1@ \$44.75

VACUUM FORMERS

PATTERSON

Blockout Compound

Blocks out remaining teeth and undercuts on models when adapting baseplates.

Pkg.....091-0893 1@ \$12.50

Air Vac QXQ - Keystone Industries

Results rival those of expensive European pressure forming machines. Produces lab quality in-office mouthguards and bleaching trays in 5 minutes!

The Air Vac QXQ plugs into your standard air line in seconds. Simply seat your model in the patented Bead Well, heat your EVA material, and process. This compact unit is ideal for your operatory and lab area with its virtually silent operation.

- Processes all sizes of round and square EVA and baseplate/tray material
 - Simple to operate
 - Solid metal construction
 - Free installation kit included
 - 5 year warranty
 - UL/cUL/CE Rated
- Each510-2736 1@ \$659.35



Econo-Vac - BDM

A sturdy vacuum unit for the production of bleaching trays, mouthguards, nitelguards, splints, etc. It takes standard 5" x 5" sheets of thermoplastic and is easily operated.

Each039-6580 1@ \$525.00



Sta-Vac II - BDM

This deluxe vacuum former is engineered for precision adaptation of bleaching trays, mouthguards, custom trays, base plates, temporary splints, bruxism splints, copings, etc, quickly and easily.

Each039-6598 1@ \$900.00



VACUUM FORMING

PATTERSON



Tray-Vac - BDM

The perfect unit for a small lab or office needing a solid, reliable vacuum former for constructing bleaching trays, custom impression trays, mouthguard and splints
 Each039-6341 1@ \$495.00



The Machine III - Keystone Industries

With the most powerful 2-stage electric motor, finished appliances are produced with ease. Other quality features include:
 • All aluminum construction
 • Processes all sizes of square EVA and baseplate/tray material
 • Heavy-duty silicone gaskets
 • 5 year warranty
 • UL/cUL/CE Rated
 Each510-2728 1@ \$455.40



RazerTip - BDM

Electric Trimming Device. For vacuum formed material. Variable heat control. Handpiece selector for 2 handpieces. 7 different cutting handpieces available.

Unit039-6424 1@ \$289.95

Perfecta® Tray Magic Trimmer - American Dental Hygienics



Heated, curved, double-edged blade makes trimming trays fast and easy. Adjustable heat setting allows for use with denture relines, mouthguards, nightguards and other applications.
 Each 104-0799 1@ \$489.95



Tray Material

Ideal for custom impression trays. Blue, high-impact styrene, available in two thicknesses.

100/Pkg		
.125091-0927	1@ \$112.95 3@ \$110.05
.150091-0968	1@ \$128.60 3@ \$124.15

Baseplate Material

Ideal for wax biteplates and denture bases. Pink high-impact styrene, available in three thicknesses.

100/Pkg		
.060091-1008	1@ \$64.95 3@ \$62.75
.080091-1040	1@ \$83.40 3@ \$81.20
.100091-1081	1@ \$89.55

Mouthguard Material

Ideal for mouthguards, nightguards and model duplication. Soft EVA material for energy absorption, available in three thicknesses.

25/Pkg		
.040089-5003	1@ \$27.50 3@ \$26.05
.080091-1180	1@ \$29.75 3@ \$29.00
.150091-1214	1@ \$46.05

Coping Material

Ideal for temporary crown and bridge splints. Clear material, available in two thicknesses.

50/Pkg		
.020091-1123	1@ \$23.35 3@ \$21.75
.030091-1149	1@ \$23.15 3@ \$21.70

Temporary Splint Material

Ideal for hand splints, ortho retainers and surgical trays. Stiff, clear material available in three thicknesses.

50/Pkg		
.020091-1263	1@ \$21.40 3@ \$19.95
.060091-1305	1@ \$46.20
.080091-1347	1@ \$56.50 3@ \$52.95



1. Heat setting 7 or 8. Remove excess material from model. Trim tray while still on model.



2. Start at most posterior molar. Allow Tray Magic blade to rest on free gingival margin. Start interproximally.



3. Keep blade 90 degrees to facial surface of tooth. Leave 1/2 to 1mm short of gingiva.



4. Remove excess. Tray Magic blade makes smooth cut, so no additional finishing is required.

Vacuum Forming Material – BDM

Clear material .060 25/Pack	133-0513	1@ \$27.60
Clear material .080 25/Pack	133-0505	1@ \$33.85
Dual laminate 12/Pack	133-0802	1@ \$113.40
Mouthguard strap 10/Pack	155-0466	1@ \$13.40

Vacuum Forming Material – Dental Resources

Mouthguard Material

5" x 5"	12/Pkg	1@ \$82.20
Clear	279-5219	Red
Orange	279-5227	White
Blue	279-5235	Yellow
Green	279-5243	Assorted
Black	279-5250	Blue/White
Purple	279-5276	Black/Red

Laminate Mouthguard Material

5" x 5" – .150	25/Pkg	1@ \$44.70
		6@ \$42.05
Yellow	279-6142	Black
Blue	279-6159	Green
Orange	279-6191	

Materials	Characteristics	Applications	Techniques
PROFORM Mouthguard Material 12/Package 279-5219 Clear 279-5227 Orange 279-5235 Blue 279-5243 Green 279-5250 Black 279-5276 Purple 279-5284 Red 279-5292 White 279-5300 Yellow 279-5318 Assorted 279-4287 Blue/White 279-4311 Black/Red	<ul style="list-style-type: none"> • Soft EVA • 5" x 5" laminated sheets • Variety of colours • Anterior brace mechanism • Strap available 	<ul style="list-style-type: none"> • Professionally made mouthguards for athletics 	<ul style="list-style-type: none"> • Heat until 1" droop develops • Vacuum form with brace behind anteriors and trim with scissors
PROFORM Niteguard Material 12/Package 279-5508	<ul style="list-style-type: none"> • Two layers of EVA – one dense and one softer • 5" x 5" laminated sheets • 30% denser than mouthguards 	<ul style="list-style-type: none"> • Soft bruxing appliance • Custom fluoride trays • Ortho repositioner • Surgical applications 	<ul style="list-style-type: none"> • Place niteguard name up • Heat until 1" droop • Vacuum form and trim with scissors
PROFORM Dual Laminate Material 12/Package 279-5052	<ul style="list-style-type: none"> • Hard upper level • Soft retentive lower layer • 5" x 5" laminated sheets • Acrylic bonds to both layers 	<ul style="list-style-type: none"> • Hard bruxing appliance • Ortho repositioning appliance • TMJ splint 	<ul style="list-style-type: none"> • Heat soft side first until material is rubbery • Flip and heat hard side • Vacuum form
PROFORM Bleaching Laminate Material 50/Package 279-5003	<ul style="list-style-type: none"> • Stiff upper thin material for margin control • Foam liner for bleach absorption • 5" x 5" laminated sheets 	<ul style="list-style-type: none"> • Foam lined bleaching tray • Fluoride tray 	<ul style="list-style-type: none"> • Heat material until foam is softened (does not droop) • Vacuum form and trim with scissors
PROFORM .020 Temporary Crown and Bridge Material 50/Package 279-6274	<ul style="list-style-type: none"> • Thin opaque material • 5" x 5" sheets • Acrylic does not stick to material 	<ul style="list-style-type: none"> • For production of temporary crown and bridge forms 	
PROFORM .030 Temporary Matrix Retainer Material 50/Package 279-5938	<ul style="list-style-type: none"> • Stiff, clear material • 5" x 5" sheets of PolyCarb • Acrylic sticks to material • High melt temperature 	<ul style="list-style-type: none"> • Invisible orthodontic retainers 	<ul style="list-style-type: none"> • Heat until bubbles appear • Drop onto model and vacuum form and trim with scissors
PROFORM .150 Mouthguard Material 25/Package 279-6142 Yellow 279-6159 Blue 279-6167 Red 279-6175 Black 279-6183 Green 279-6191 Orange	<ul style="list-style-type: none"> • EVA • Soft for energy absorption • 5" x 5" sheets 	<ul style="list-style-type: none"> • Mouthguards • Niteguards • Model duplication 	<ul style="list-style-type: none"> • Heat until 1" droop develops • Vacuum form and trim with scissors • Bevel edges with micro torch
PROFORM Foam Liner 50/Package 279-5029	<ul style="list-style-type: none"> • 5" x 5" thin packing foam • Air bubbles trapped in material 	<ul style="list-style-type: none"> • Bleaching liner 	<ul style="list-style-type: none"> • Lay on top of model • Vacuum form .040 or .080 soft EVA down on top of liner

LIGHT CURING



Eclipse Prosthetic Resin System – Dentsply Prosthetics

Handles like a dental wax, yet cures to the standards of a conventional denture base.

- Monomer-free
- New indirect buildup technique providing a processed baseplate for every patient. Now “Fit at try-in equals fit at delivery”
- Final denture cure time of less than 13 minutes, increasing capacity and turnaround time potential
- Eliminates flasking and investing
- Meets or exceeds ADA and ISO specifications

Starter Kit 1

Contains Eclipse processing unit with lamps and filters, Dentsply conditioning oven with heat resistant gloves, Dentsply Wax Pencil Pro, Dentsply hot air gun, Dentsply melting pot, tooth cage, tooth storage boxes, 10 shade guides, 12/pkg boxes of resins (18) and miscellaneous accessories
Kit274-0983 1@ \$20,698.50

Starter Kit 2

Contains Eclipse processing unit with lamps and filters, Dentsply conditioning oven with heat resistant gloves, Dentsply Wax Pencil Pro (2), Dentsply hot air gun (2), Dentsply melting pot (2), 12/pkg boxes of resins (11) and miscellaneous accessories
Kit274-0991 1@ \$20,598.50

Eclipse Processing Unit II (120V)

Each274-1106 1@ \$9,999.00

Conditioning Oven (120V)

Each274-1023 1@ \$3,840.00

Melting Pot

Each274-1031 1@ \$1,216.00

Hot Air Gun

Each274-1049 1@ \$321.50

Wax Pencil Pro

Each274-1056 1@ \$767.00



Eclipse Resin Material – Dentsply Prosthetics

Upper Baseplate

12/pkg		
Clear	274-1353	1@ \$110.25
Dark Pink	274-1346	
Light Pink	274-1320	
Light Reddish Pink	274-1338	
Original	274-1312	

60/pkg

Clear	274-1759	1@ \$521.50
Dark Pink	274-1742	
Light Pink	274-1726	
Light Reddish Pink	274-1734	
Original	274-1718	

Lower Baseplate

12/pkg		
Clear	274-1452	1@ \$86.20
Dark Pink	274-1445	
Light Pink	274-1429	
Light Reddish Pink	274-1437	
Original	274-1411	

60/pkg

Clear	274-1858	1@ \$407.00
Dark Pink	274-1841	
Light Pink	274-1825	
Light Reddish Pink	274-1833	
Original	274-1817	

Set-Up Resin

12/pkg		
Dark Pink	274-1544	1@ \$67.25
Light Pink	274-1528	
Light Reddish Pink	274-1536	
Original	274-1510	

60/pkg

Dark Pink	274-1940	1@ \$321.00
Light Pink	274-1924	
Light Reddish Pink	274-1932	
Original	274-1916	

Contour Resin

12/pkg		
Dark Pink	274-1643	1@ \$67.35
Light Pink	274-1627	
Light Reddish Pink	274-1635	
Original	274-1619	

60/pkg

Dark Pink	274-2047	1@ \$321.00
Light Pink	274-2021	
Light Reddish Pink	274-2039	
Original	274-2013	



iCure Colour Curing Unit – BDM

The iCure is the perfect addition to any dental office or laboratory. It rapidly cures all light cured Custom Tray Material, Temporary Material, Baseplate Material, and Block-Out Materials in mere minutes without distortion, heat or unhealthy vapours. CSA approved.

Each133-2238 1@ \$479.95

Triad® 2000™ – Dentsply Light Curing Unit

The Triad 2000 unit features a rotating, adjustable turntable, roomy chamber, automatic timer, and a smaller “footprint.” Its powerful halogen tungsten lamp and an advanced cure-management system augment the 2000's performance. Triad sheet, rope, and clear gel materials allow quick and easy fabrication of appliances for use in removable, maxillofacial, and C & B applications, such as:

- Provisional restorations
 - Partial denture saddles and repairs
 - Flippers
 - Record bases
 - Denture repairs and relines
 - Occlusal splints
 - Implant stents and trays
 - Identification procedures
 - Custom impression trays
 - Orthodontic appliances
 - Emergency dentures
- 115 Volt276-6368 1@ \$1,325.00



Triad® Accessories – Dentsply

For Use with Triad VLC Materials

Contains all the accessories that are used with Triad VLC materials. Contains: 1-5oz air barrier coating, 1-4oz model release agent, 1-20ml VLC bonding agent, 1 modeling tool, 2 brushes (#4 flat, #6 round).

Each276-7515 1@ \$82.90

Air Barrier Coating (ABC)

5oz Bottle.....276-6459 1@ \$19.25

Model Release Agent

4oz Bottle.....276-6855 1@ \$17.20

VLC Bonding Agent

20oz Bottle.....276-6632 1@ \$17.20

Triad® VLC Modeling Instrument – Dentsply

Used with Triad® VLC materials.

Each276-6954 1@ \$18.10

Triad® DuaLine™ – Dentsply

A dual-cure full or partial reline material which combines the accuracy of self-cure, the strength and dimensional stability of light cure, and an easy time-saving technique. Extrudable through a dual paste cartridge.

DuaLine is Methyl and Butyl Methacrylate Monomer-free.

Intro Package

Contains: 2-50gm DuaLine cartridges, 8 mixing tips, 30ml bottle bonding agent, 15ml bottle air barrier coating, 25 brushes, instructions, MSDS.

Pkg.....276-8786 1@ \$170.00

Refill Package

Contains: 2-50gm DuaLine cartridges, 8 mixing tips, instructions, MSDS.

Pkg.....276-8794 1@ \$129.70

Triad® Gel – Dentsply

Ready-To-Use, Flowable VLC Material

For orthodontic, implant, and other removable appliances.

Standard Package

Contains: 4-22gm tubes.

Clear Colourless276-7176 1@ \$179.00

Clear Pink276-7184 3@ \$176.50

Assortment.....276-7143

Assortment Package

Contains: 1-22gm tube each of clear colourless, clear pink, blue and red gels.

Pkg.....276-7143 1@ \$179.00

3@ \$177.00

Triad® Materials – Dentsply

Visible Light Cure Materials

For denture base fabrication, relines, repairs, as well as orthodontic appliances, custom trays, provisional crown and bridge and many other applications.

VLC Denture Base & Orthodontic Material

Sheets 6-Unit Package

Light Pink Fibered276-6822 1@ \$48.25

Reg. Pink Fibered276-6830

Reg. Pink Unfibered..276-6848

Dark Pink Fibered.....276-6798

30-Unit Package

Light Pink

Fibered276-6665 1@ \$191.95

Reg. Pink Fibered276-6673 3@ \$188.95

Reg. Pink Unfibered..276-6681

Ropes 6-Unit Package

Light Pink Fibered276-7069 1@ \$48.25

Reg. Pink Fibered276-7077

Reg. Pink Unfibered..276-7085

Dark Pink Fibered.....276-7093

Triad® Provisional Material – Dentsply

Tooth Colour

Ropes, 6-Unit Package

Ivory Extra Light276-7861 1@ \$47.70

Ivory Light276-6566

Ivory Medium276-6574

Ivory Medium

Yellow.....276-7879

Ivory Dark.....276-6582

Ivory Dark Brown276-7887

Enamel.....276-6590

30-Unit Package

Ivory Extra Light276-7952 1@ \$144.00

Ivory Light276-8000

Ivory Medium276-8018

Ropes Assortment Package

Contains: 1 rope of each ivory shade plus an enamel.

Pkg.....276-7531 1@ \$48.95

Triad® Reline Material – Dentsply

Original Flow

Sheets 6-Unit Package

Light Pink

Fibered276-6871 1@ \$63.75

Regular Pink

Fibered276-6889

30-Unit Package

Light Pink

Fibered276-6905 1@ \$193.50

Regular Pink

Fibered276-6913 3@ \$188.95

VLC Reline Material

Ropes, 6-Unit Package

Light Pink

Fibered276-6939 1@ \$48.20

Regular Pink

Fibered276-6947

Hi-Flow VLC Reline Material

6-Unit Package

Light Pink

Fibered276-7317 1@ \$125.50

Regular Pink

Fibered276-7325 3@ \$123.50

30-Unit Package

Light Pink

Fibered276-7366 1@ \$362.95

Regular Pink

Fibered276-7374

VLC Custom Tray Material

Blue Sheets

6-Unit Pkg276-6756 1@ \$37.00

3@ \$36.50

30-Unit Pkg276-6640 1@ \$122.55

3@ \$119.65

Triad® Resiline™ VLC Reline Material – Dentsply Trubyte

Resilient Reline Material

Brings convenience, control, and quality to extended-term full and partial denture relines, implant prostheses, and numerous prosthetic treatment options. Benefits include chairside delivery and two-minute curing.

Intro Kit

Contains: 3-46gm tubes Resiline, 1 application key, 1-20ml bonding agent, 1 modeling tool, MSDS and instruction sheet.

Kit276-8059 1@ \$185.50

Refill Kit

Contains: 2-46gm tubes Resiline, 1 application key, MSDS and instruction sheet.

Kit276-8075 1@ \$103.50



No. 6NC Bunsen Burner – Buffalo Dental

Chrome finish. For natural gas.
Each 148-3445 1@ \$29.95



Kisag Burner – Pflingst

Complete With Refill
Each 562-0208 1@ \$224.40

Refill for Kisag Burner
Each 562-0190 1@ \$26.40

Laboratory Burner – Wall Lenk

Small, compact, portable.
Complete Package
Contains: 1 burner head, gas refill.
Pkg 709-4600 1@ \$54.05



Burner Head Only
Up to 1/2" to 6" adjustable flame.
Each 709-4618 1@ \$61.20

Gas Refill
5.5oz Can 709-4691 1@ \$15.85

MiraTorch – Hager Worldwide

Each 765-7976 1@ \$81.85



Touch-O-Matic® Burner – Waterpik

Natural Gas
Complete with Pilot Shield.
Each 354-9292 1@ \$315.50



Waxing Cup
Each 354-9367 1@ \$28.40

Loop Heat Conductor
Each 354-9375 1@ \$26.75

Waxing Cup – Buffalo Dental

Designed to fit 5N and 5U burners, has a 3" OD and is constructed of one-piece case zinc aluminum with a polished finish
Each 151-8059 1@ \$17.50



TORCHES AND ACCESSORIES

Alcohol Lamp – Buffalo Dental

Chrome Plated Brass
Base diameter – 3 3/4". Holds 5 fluid ounces.
Each 154-4857 1@ \$39.25



Alcohol Torch – Buffalo Dental

Produces a needle-point flame with gentle pressure on the weighted plastic bottle. May also be used as an alcohol lamp; 5" high.
Each 158-1875 1@ \$8.85



Alcohol Lamp Wicks – Buffalo Dental

Each 154-5367 1@ \$0.80
12/Pkg 154-5342 1@ \$3.75

Alcohol Torch – Waterpik

Alcohol Torch
Each 354-3907 1@ \$298.25



#3 Nozzle Tip
Each 354-4889 1@ \$16.30

Nozzle Holder Assembly, complete
Each 354-3956 1@ \$60.65

Alcohol Wicks
12/Pkg 354-5860 1@ \$2.25

Glass Alcohol Lamp – Dixon

Nickel-plated brass burner and ratchet. Safety cap prevents fuel evaporation. Includes 3" wick.
Each 286-9972 1@ \$12.50



Hand Torch – Blazer

Multi-purpose, refillable, self-igniting.
Each 132-7709 1@ \$80.45



Gas Refill – Blazer
150gm Can 132-7717 1@ \$5.90
4@ \$5.30

Hot Shot MT 3000 – Blazer

Automatic, self-igniting, refillable, Piezo electronic ignition.
Each 132-7691 1@ \$81.50



Metal Alcohol Lamp – Dixon

Polished chrome-plated brass with ratchet. Safety cap on chain prevents evaporation. Supplied with 3/16" wick. 3oz capacity.
Each 286-9980 1@ \$19.95



Replacement Wicks – Dixon

For alcohol torches. 3/16" diameter, supplied.
Each 287-1507 1@ \$1.30

BRUSHES

PATTERSON



Disposable Brush Applicators

These bendable brushes have a wire core for easier bending and retain shape once bent. Ideal for use with etchants, bonding agents or other resins. Available in both regular bristle and fine bristle in an assortment of colours.

Gross Pack: 144 Brushes (Black-Regular Bristle)

144/Box.....083-8961 1@ \$20.25

Gross Pack: 144 Brushes (White-Fine Bristle)

144/Box.....083-8953 1@ \$20.25

Assorted (72 standard and 72 fine)

144/Box.....083-8136 1@ \$22.25

Economy Pack: 432 Brushes (Black-Regular Bristle)

Each084-2443 1@ \$55.35

Economy Pack: 432 Brushes (White-Fine Bristle)

Each084-2450 1@ \$55.35

Assorted Economy Pack: 432 Brushes (White and Black)

Each084-2468 1@ \$55.35



Composite Brush Handle and Tips

Applicator Brush Handle

White

2/Pkg.....084-3383 1@ \$6.70
6@ \$6.10

Brush Tips

100/Pkg.....084-3391 1@ \$14.65
3@ \$13.20

Sable Brushes

A. #0 Right Angle.....	093-8209	1@ \$4.10
B. #00 Right Angle....	093-8217	1@ \$4.10
C. #000 Right Angle..	093-8225	1@ \$4.10
D. #0 Straight.....	093-8258	1@ \$3.40
E. #00 Straight	093-8266	1@ \$3.40
F. #000 Straight.....	093-8274	1@ \$3.40
G. #1 Straight.....	093-8282	1@ \$3.70
H. #2 Straight.....	093-8290	1@ \$4.75
I. #3 Straight	093-8308	1@ \$4.95
J. #4 Straight	093-8316	1@ \$5.95
K. #5 Straight.....	093-8324	1@ \$6.60
L. #6 Straight	093-8332	1@ \$6.30
M. #7 Straight.....	093-8340	1@ \$7.20



Porcelain Brushes

Golden Takalon bristles hold their shape long after sable brushes have been discarded. Straight brushes in assorted sizes.

N. #2 Flat.....	093-8001	1@ \$4.65
O. #4 Flat.....	093-8019	1@ \$5.00
P. #6 Flat	093-8027	1@ \$7.60
Q. #0 Round.....	093-8050	1@ \$4.15
R. #00 Round	093-8068	1@ \$4.10
S. #000 Round.....	093-8076	1@ \$4.10
T. #1 Round	093-8084	1@ \$4.10
U. #2 Round	093-8092	1@ \$4.95
V. #3 Round.....	093-8100	1@ \$5.25
W. #4 Round	093-8118	1@ \$5.80
X. #5 Round	093-8126	1@ \$6.30
Y. #6 Round.....	093-8134	1@ \$7.55
Z. #7 Round.....	093-8142	1@ \$8.00



Composite Brush Handle – Hu-Friedy

#HCB1.....368-4701 1@ \$22.25

BRISTLE BRUSHES

PATTERSON



Bristle Brushes

12/Box

Soft #9083-1628
#11083-1792
#12083-1826

Stiff #9083-1636
#11083-1800
#12083-1834

X-Stiff #9083-1644
#11083-1818
#12083-1842

144/Box

Soft #9083-1693
#11083-1883
#12083-1917

Stiff #9083-1701
#11083-1891
#12083-1925

X-Stiff #9083-1719
#11083-1909
#12083-1933

	12/Box	144/Box
#9	\$9.10	\$87.05
#11 & #12	\$9.10	\$85.15



Abbott-Robinson Bristle Brushes - Buffalo Dental

Handpiece Shank

144/Package

#9 Extra Stiff		1@ \$113.00
Black147-3446	3@ \$107.50
#9 Standard Stiff		
Silver Grey147-3479	
#11 Extra Stiff		
Black147-3537	
#11 Standard Stiff		
Silver Grey147-3560	
#12 Extra Stiff		
Black147-3651	
#12 Standard Stiff		
Silver Grey147-3685	

Abbott-Robinson™ Bristle Brushes - Buffalo Dental

144/package

Small Cup Standard Stiff		1@ \$113.00
Silver Grey147-3867	3@ \$107.50

Large Cup Standard Stiff		
Silver Grey147-3776	

Large Cup Extra Stiff		
Black147-3743	

#9 Soft		
Grey147-2711	1@ \$106.00

#11 Soft		
Grey147-2778	

#12 Soft		
Grey147-2836	

Large Cup Soft		
Grey147-2893	1@ \$106.00

Abbott-Robinson™ Bristle Brushes - Buffalo Dental

12/Package

#9 Standard Stiff		1@ \$10.75
Silver Grey147-1861	

#9 Extra Stiff		
Black147-1929	

#11 Standard Stiff		
Silver Grey147-1960	

#11 Extra Stiff		
Black147-2018	

#12 Standard Stiff		
Silver Grey147-2034	

#12 Extra Stiff		
Black147-2059	

Small Cup Standard Stiff		
Silver Grey147-2190	

Large Cup Standard Stiff		
Silver Grey147-2109	

Large Cup Extra Stiff		
Black147-2075	

#9 Soft		
Grey147-4386	1@ \$12.00

#11 Soft		
Grey147-4451	

#12 Soft		
Grey147-4501	

Large Cup Soft		
Grey147-4568	1@ \$12.00

Small Cup Right Angle Soft		
Grey147-4816	1@ \$20.35

BRUSH WHEELS

PATTERSON



Brush Wheels

Wooden Center

2-Row Converging		1@ \$4.70
B12 - 1 1/8"091-0109	12@ \$3.80

2-Row Converging		1@ \$2.50
B14 - 2"091-0224	12@ \$2.25

2-Row Converging		1@ \$2.50
B20 - 2 1/2"091-0125	12@ \$2.25

3-Row Converging		1@ \$3.85
B23 - 2 3/4"091-0141	12@ \$3.30

3-Row Straight		1@ \$3.85
B25 - 2 3/4"091-0166	12@ \$3.30

4-Row Straight		1@ \$3.85
B27 - 3"091-0182	

4-Row Straight		1@ \$4.85
B29 - 3"091-0208	12@ \$4.25

Brush Wheels

Plastic Center

2-Row Converging.		
B12P - 1 1/8"		
Each091-0216	1@ \$2.50

B20P - 2 1/2"		1@ \$3.20
Each091-0232	12@ \$2.65

BUFFS, FELT

Bufs, Felt - Buffalo Dental

Soft miniature felt wheels are 1/8" thick. Pinhole centers for mounting on tapered lathe chucks and HP mandrels as appropriate.




1/2" Soft Square Edge Wheels		
24/Pkg152-2580	1@ \$5.50

5/8" Soft Square Edge Wheels		
24/Pkg152-2606	1@ \$6.40

3/4" Soft Square Edge Wheels		
24/Pkg152-2622	1@ \$6.25

BUFFS, FELT (CONT.)



Pointed Blunt

Felt Cones - Buffalo

White Wool Felt

Compressed and hard-finished, providing durable service, even when wet.

1/2" Diameter x 1" Long

#1C Pointed	152-0519	1@ \$3.20
#4C Blunt.....	152-0543	6@ \$3.10

3/4" Diameter x 1-1/2" Long

#2C Pointed	152-1061	1@ \$4.60
#5C Blunt.....	152-1095	6@ \$4.15

1" Diameter x 2" Long

#3C Pointed	152-1616	1@ \$6.10
#6C Blunt.....	152-1640	6@ \$5.85

BUFF WHEELS

PATTERSON

Buff Wheels

Muslin buff wheels available in a variety of multi-ply sizes.

Available in high-quality muslin with uniform stitching for long-lasting durability. 1 per package.



2" x 24-Ply.....	091-0406	1@ \$2.05	10@ \$1.65
3" x 24-Ply.....	091-0422	1@ \$1.65	10@ \$1.45
3" x 30-Ply.....	091-0448	1@ \$1.75	10@ \$1.55
3" x 35-Ply.....	091-0463	1@ \$2.00	10@ \$1.75
4" x 24-Ply.....	091-0489	1@ \$1.90	10@ \$1.75
4" x 30-Ply.....	091-0505	1@ \$2.10	10@ \$1.90
4" x 42-Ply.....	091-0521	1@ \$2.95	10@ \$2.60
5" x 45-Ply.....	091-0547	1@ \$4.25	10@ \$3.90
6" x 30-Ply.....	091-0562	1@ \$3.80	10@ \$3.45
1".....	091-0687	1@ \$1.20	10@ \$1.00
1 1/2".....	091-0703	1@ \$1.40	10@ \$1.10

Buff Wheel - Dixon

Lead Center

4" x 36-Ply	
Each	283-8555 1@ \$4.25

BASE PLATES

Baseplates - Miltex Moyco

Regular/Pink

45/Pkg		
Lower.....	524-2037	1@ \$24.75
Upper.....	524-2060	

400/Pkg

Lower.....	524-2722	1@ \$192.75
Upper.....	524-2714	2@ \$188.75

Double Thick/Pink

25/Pkg		
Lower.....	524-2581	1@ \$26.10
Upper.....	524-2615	

200/Pkg

Lower.....	524-2664	1@ \$153.15
Upper.....	524-2656	

Trubase® Baseplates - Dentsply

Single Thickness

50/Package		
Upper Pink	269-3901	1@ \$46.10
Lower Pink	269-3935	

Double Thickness

25/Package		
Upper Pink	269-5583	1@ \$42.85
Lower Pink	269-5617	

WAXES

PATTERSON



Wax sheets, Canadian

5 lb/Box		1@ \$33.50
#400 Pink thick	164-6710	
#401 Red thick	164-6744	
#402 Pink thin.....	164-6728	
#403 Red thin	164-6736	



Baseplate Wax

All-Season #2 Pink		1@ \$73.80
5lb Box	091-1396	3@ \$68.25

Set-Up Wax

Pink #2		
5lb Box	091-1412	1@ \$83.05
		3@ \$80.95

PATTERSON



Boxing Wax

11" Long x 1.5" Wide

1lb Box

Regular.....	091-1446	1@ \$20.40
X-Thin.....	091-1453	3@ \$19.30



Occlusal Bite Blocks

Recessed, 3/8"

Pink wax for occlusal rims. U-shaped; pliable for custom reshaping at room temperature to fit the individual arch.

100/Box.....	091-0000	1@ \$28.35
		6@ \$27.55



Utility Wax Strips

Large - 11" Long x 3/16" Diameter

75 Strips/Box

White.....	091-1560	1@ \$14.00
Red.....	091-1578	3@ \$12.80

Small - 11" Long x 1/8" Diameter

114 Strips/Box

White.....	091-1586	1@ \$14.00
Red.....	091-1594	3@ \$12.80



Wax Square Ropes

11" Long x 3/16" Square

55 Strips/Box

White.....	091-1503	1@ \$14.15
Red.....	091-1511	3@ \$12.60



Wax sheets, Canadian – Carmel Industries

5 lb/Box		
#400 Pink thick	164-6710	1@ \$33.50
#401 Red thick	164-6744	3@ \$32.50
#402 Pink thin	164-6728	
#403 Red thin	164-6736	



Baseplate Wax – Heraeus Kulzer
Sheets

Modern Pink No. 3		
Soft, with high melting range.		1@ \$28.80
1lb Box	517-4917	5@ \$27.85
5lb Box	517-5401	1@ \$118.50
		3@ \$114.85

Red Base Plate		
All purpose, extra tough.		
1lb Box	517-5898	1@ \$28.05
		5@ \$27.20
5lb Box	517-6383	1@ \$123.30

Baseplate Wax – Coltène Hygenic
Sheets – .055" Thick

Extra Tough Pink		
1lb Box	382-5502	6@ \$16.05
5lb Box	382-5999	1@ \$71.40
Medium Soft No. 3 Pink		
1lb Box	382-8100	6@ \$16.05
5lb Box	382-8209	1@ \$71.40
Red Wax		
1lb Box	382-7615	6@ \$16.05
5lb Box	382-6484	1@ \$71.40

Beauty Pink Dental Wax – Miltex Moyco
Sheets

For extremely accurate impression and bites.		
5lb Box		
Extra Hard	529-7726	1@ \$97.35
Hard	529-7734	

Beeswax – Coltène Hygenic

<i>Sheets – Pure Yellow</i>		
2-7/8" x 5-3/4" x .055" Thick		
1/2lb Box	382-9587	1@ \$20.60

Bite & Impression Wax – Aluwax
Denture Wax

Sheets – 2-3/4" x 5-3/4" x 1/8"		1@ \$33.60
Thick Box	105-9500	3@ \$31.05

Scored Wax		
Sheets – 2-3/4" x 5-3/4" x 5/16". Thick identical formulation to denture material, but scored to allow easy use in small quantities.		
Box	106-0979	1@ \$25.55
		3@ \$23.55

Combination Pak		
Contains an assortment of denture and scored waxes in one package.		
Box	106-0482	1@ \$25.55
		3@ \$23.55

Denture Forms Wax <i>Cut-To-Shape Forms</i>		
Made of the standard denture material.		
Box	105-9997	1@ \$34.45
		3@ \$32.30

Waxed Cloth Forms <i>Cut-To-Shape</i>		
Permits an impression of the entire teeth. Contains coarse cloth at the centre of the material.		
Box	106-1951	1@ \$34.45
		3@ \$32.30

Waxed Cloth Sheets <i>2-3/4" x 5-3/4" x 3/16" Thick</i>		
Box	106-1464	1@ \$26.55
		3@ \$24.60

Boxing Wax – Coltène Hygenic

For boxing impressions, withstands extreme heat.		
1lb Regular	382-6229	1@ \$30.90
1lb Extra Thin	377-5129	1@ \$33.65
5lb Regular	382-6302	1@ \$127.00

Boxing Wax – KerrLab

<i>Melting Temperature: 176°F/71°C</i>		
For boxing impressions prior to pouring the mold.		
Sticks		
6.2oz Pkg	440-1113	1@ \$40.60
Strips		
1lb Pkg	440-1600	1@ \$70.20

Byte Ryte Wax – Mizzy

<i>Sheets – Bite Registration Wax</i>		
5lb Box	512-9390	1@ \$84.40



Check Bite Wafers – Heraeus Kulzer

<i>Pre-Formed Wafers – 1/8" Thick</i>		
White – for checking bites.		
Yellow – for taking bites.		
Trial Package		
Contains: 56 wafers.		
Yellow	517-3380	1@ \$31.10
Economy Package		
Contains: 192 wafers.		
Yellow	517-3935	1@ \$101.60

Coprwx Bite Wafers – Heraeus Kulzer

<i>Bite Registration Wax</i>		
Contains a thin sheet of foil to help prevent the teeth from cutting through and reduce distortion.		
84/Box	457-4240	1@ \$108.20
		3@ \$105.60

Inlay Casting Wax – KerrLab

<i>Melting Temperature: Regular 160°F/71°C, Hard 163°F/73°C</i>		
Standard Package		
12 Sticks/Package		
Regular Blue	440-8571	1@ \$22.25
Hard Blue	440-8605	3@ \$21.60
Regular Ivory	441-0650	1@ \$20.00
		3@ \$19.35
Bulk Package		
120 Sticks/Package		
Regular Blue	440-8738	1@ \$52.00
Hard Blue	440-8720	
Hard Green	441-0130	1@ \$66.75

Model Cement – Dentsply

12 sticks/Box	277-9502	1@ \$24.95
120 sticks/Box	277-9510	1@ \$189.50

Neowax® – Dentsply

<i>Flexible Baseplate Wax</i>		
.050" Thick		
5lb Pkg	278-1623	1@ \$47.30
.070" Thick		
5lb Pkg	278-1656	1@ \$47.40

WAXES (CONT.)

Orthodontic Tray Wax – Coltene Hygenic

Sticks – 6" x 1/4" Dia.

For peripheral lining of trays, increasing tray height, post-damming, etc. Soft and tacky, they may be adapted without heating. Pure white colour.

48/Box383-1450 1@ \$16.85

Orthodontic Tray Wax – Heraeus Kulzer

Mint Flavoured

60 Sticks/Box517-7506 1@ \$28.20

Orthodontic Tray Wax – Kerr

Melting Temperature: 164.3°F/73°C

White1@ \$30.50

48 Sticks/Box441-2276

Economy Pack1@ \$227.00

480 Sticks/Box441-2763

Periphery Wax – Heraeus Kulzer

Sticks

For use with all types of impression materials. For making wax stops, post-damming and lining the periphery of trays. Soft, tacky and easily shaped at room temperature. High melting point for use with hydrocolloid material. Minimum tissue displacement.

60 Sticks/Box457-5221 1@ \$30.85

3@ \$29.85

Plastic Wax Sticks – Hygenic

Scented

Sticks readily formed without heating for post damming, peripheral lining and making stops in impression trays. Adheres to hydrocolloid impression materials.

48/Box383-2086 1@ \$17.25

Pure Beeswax – Moyco

Sheets – Pure Yellow

For bite registration.

1lb Box529-8229 1@ \$30.55

Ready Made Wax Shapes – Kerr

Sprue Wax

1/2lb Spools, Round

6ga441-5675 1@ \$34.75

8ga441-5683

Regular U-Shaped Bite Blocks – Hygenic

Lab Pack

100/Box376-7753 1@ \$80.85

Rijform Bite Rims – Yates & Bird

- Unique pink ridge-formed bite rim saves time, wax and labour
- Highest-quality ingredients ensure the utmost accuracy in bite registration
- Sticks are easily shaped to conform to the dental arch

4lb Carton

Each767-3981 1@ \$53.80

Set-Up Wax – Dentsply

Specially formulated for arranging and articulating artificial teeth.

Sheets

5lb Pkg278-8511 1@ \$114.75



Shur Wax – Heraeus Kulzer

Baseplate Wax Sheets – Pink

Regular

1lb Box517-8835 1@ \$29.25

5lb Box517-9320 1@ \$105.55

Extra Hard

1lb Box517-9726 1@ \$31.10

5@ \$30.10

5lb Box517-9817 1@ \$113.85

Slaycris Wax – Slaycris

Casting Wax

Can

Regular636-9391 1@ \$15.30

Hard636-9409 3@ \$14.45

Special Formula Wax – Keystone

Sheets – 5-7/8" x 2-7/8"

Articulating wax made for Dr. Clyde Schuyler's technique. Dark red.

1lb Box512-9887 1@ \$37.30



Sticky-Wax – KerrLab

Melting Temperature: 136°F/58°C

Standard Package

12 Sticks/Box441-8232 1@ \$17.35

3@ \$16.85

Bulk Package

120 Sticks/Box441-9214 1@ \$57.10



TruWax® – Dentsply Trubyte

Baseplate Wax

1lb Box

Regular278-2639 1@ \$34.45

Extra Hard278-4106

5lb Box

Extra Hard278-4593 1@ \$89.50

Utility Wax – Coltene Hygenic

Shaped in Round Ropes

80/Box White383-3183 1@ \$15.30

80/Box Red383-3159

Utility Wax Ropes – Heraeus Kulzer

11" Long x 3/16" Square

55/Package

Red518-3488 1@ \$22.10

White518-3512 6@ \$21.40

Utility Wax Strips – Heraeus Kulzer

Large Diameter – 11" Long x 3/16" Diameter

75 Strips/Package

Large Red518-2878 1@ \$22.10

Large White518-2902 6@ \$21.40

Small Diameter – 11" Long x 3/16" Diameter

114 Strips/Package

Small Red518-2936 1@ \$22.10

Small White518-2969 6@ \$21.40

Wax Bite Blocks – Heraeus Kulzer

Pre-Formed – 7/16" Thick
Pink, Laboratory Package
120/Package

Hard.....	517-0352	1@ \$101.25
Soft.....	517-0386	3@ \$98.10

Wax Bite Wafers – Hygenic

U-Shaped Bite Wax
Double Thick-Foil Laminated
Light Blue

24/Box.....	383-0551	1@ \$21.95
		5@ \$21.05

Single Thickness-Without Foil Lamination

Lemon Yellow		
50/Box.....	383-0064	1@ \$21.95
		2@ \$21.05

Wax Occlusal Wafers – Miltex Moyco

45/box

Pink		
Soft.....	530-0751	1@ \$42.00
Medium.....	530-0785	
Hard.....	530-0819	
X-Hard.....	530-0843	

Yellow Bite Wax – Coltene Hygenic

Economical alternative to beeswax with many of the same characteristics.

1lb Box.....	382-9082	1@ \$16.70
		5@ \$16.25
5lb Box.....	382-9579	1@ \$72.05

Yellow Bite Wax – Heraeus Kulzer

Sheets – Rich In Beeswax

Softens quickly and flows easily, sets hard when chilled.

1lb Box.....	517-1368	1@ \$36.10
		3@ \$34.90
5lb Box.....	517-1855	1@ \$135.95

10X Light Pink Wax – Miltex Moyco

Thin, Pliable and Soft

For taking accurate bites and waxing setup.

1lb Box.....	530-1403	1@ \$31.70
5lb Box.....	530-1411	1@ \$137.30

CROWN & BRIDGE WAXES

PATTERSON



Dipping Wax

- Bright colour
- Excellent margin

Powdered 8oz

Red.....	091-1677	1@ \$23.90
Yellow.....	091-1685	1@ \$24.35
Green.....	091-1669	1@ \$23.90

Carving Wax

This high-quality medium-hard wax is used for modeling in connection with the crown and bridge technique. Easy to mill, but also to scrape, it does not smear and is dimensionally stable. An ideal wax for all milling work.

2oz

Green.....	091-0034	1@ \$9.40
Blue.....	091-0042	



Margin Wax

Special elastic, tension-free, soft and shrink-free wax for modeling cervical edges of crowns, precision parts and inlays. The edges do not break off!

2oz

Red.....	091-0018	1@ \$9.40
Ivory.....	091-1701	1@ \$10.65



Xtra Hard Carving Wax

High-quality hard wax for modeling in connection with crown and bridge technique. Ideal for scraping, quick hardening after application, easy to shape.

2oz

Green.....	091-0059	1@ \$8.90
Blue.....	091-0067	1@ \$9.45



BelleWax Dipping Wax Cubits™ – KerrLab

18oz Bottle

Green.....	131-5696	1@ \$53.95
Red.....	131-5712	
Yellow.....	131-5704	

9oz Bottle

Green.....	131-5613	1@ \$38.95
Red.....	131-5647	
Yellow.....	131-5639	
Blue.....	131-5621	

Inlay Casting Wax – KerrLab

Melting Temperature: Regular 160°F/71°C, Hard 163°F/73°C

Standard Package

12 Sticks/Package

Regular Blue.....	440-8571	1@ \$22.25
Hard Blue.....	440-8605	3@ \$21.60
Regular Ivory.....	441-0650	1@ \$20.00
		3@ \$19.35

Bulk Package

120 Sticks/Package		
Regular Blue.....	440-8738	1@ \$48.40
Hard Blue.....	440-8720	
Hard Green.....	441-0130	1@ \$66.75

Slaycris Wax – Slaycris

Casting Wax

Can

Regular.....	636-9391	1@ \$15.30
Hard.....	636-9409	3@ \$14.45

Special Formula Wax – Mizzy

Sheets – 5-7/8" x 2-7/8"

Articulating wax made for Dr. Clyde Schuyler's technique. Dark red.

1lb Box.....	512-9887	1@ \$37.30
--------------	----------	------------

Sticky-Wax – KerrLab

Melting Temperature: 136°F/58°C

Standard Package

12 Sticks/Box.....	441-8232	3@ \$16.85
--------------------	----------	------------

Bulk Package

120 Sticks/Box.....	441-9214	1@ \$57.10
---------------------	----------	------------

SPRUES AND RODS

PATTERSON



Wax Sprues

Direct Sprues

100/Pkg

Button Sprues091-1875 1@ \$14.40

MISCELLANEOUS

PATTERSON



Boil Out Solution

Concentrate dissolves and disperses wax completely with no residue. Nontoxic, non-flammable, safe for all materials. Rinses completely clean.

Gallon093-5056 1@ \$48.15



Debubblizer and Wax Pattern Cleaner

Formulated for investment to flow evenly over the wax pattern. Eliminates bubbles and reduces the surface tension of waxes. Produces castings with a satin smooth finish. Can be painted on, sprayed on, or used as a dip.

8oz Bottle.....093-5106 1@ \$12.30

Quart093-5122 1@ \$28.60

Gallon.....093-5148 1@ \$67.90

PATTERSON



Liquid Foil

Superior for handling all types of acrylic restorations. Non-stick formula cleanly separates acrylic and stone. Requires only one coat.

Clear, Gallon093-5502 1@ \$51.65
4@ \$50.00



Place-It 2

Just a dab holds your assembly during soldering, welding or firing. Applies easily to board or coping interior. Works instantly, no mixing or degassing – so you save time and costly repairs are done in seconds! Place-It 2 makes it easy to change shades or add on to porcelain restorations. Safe, asbestos-free formula comes in a large, economy syringe.

Each087-7019 1@ \$35.45



Wax Solvent

Smoothly dissolves wax in denture molds. Removes oils and die lubricants. Cleans away grinding residues, polishing compounds and greasy films on instruments and work pieces. Safe for use on acrylics.

Quart093-5601 1@ \$25.55

4@ \$24.10

Gallon093-5627 1@ \$67.30

4@ \$60.20



Reline Jig

Chrome

Each093-6898 1@ \$55.40



Alcote® Separator – Dentsply

Tinfoil Substitute – Pink

65ml Bottle269-1608 1@ \$19.95

3@ \$19.10

Gallon269-1657 1@ \$95.75



Coe-Sep™ – GC Lab Technologies

Tinfoil Substitute

Separating Medium

Quart206-2115 1@ \$36.85

Gallon206-3097 1@ \$74.75



Debubblizer – KerrLab

Surface tension reducing agent for hand-mixed investments.

Pump Spray Bottle

8oz Bottle.....420-2800 1@ \$38.65

Refill Bottle

32oz Bottle.....420-2818 1@ \$56.70



MULTI-SEP – GC America, Inc.
Separating Medium

Universal separating agent for multiple applications in today's dental laboratories. MULTI-SEP features a very low viscosity and leaves no residue on the model. Can be used to separate wax, stone and both composite and acrylic resins.

2oz326-3027 1@ \$25.50
12oz326-3035 1@ \$79.50



Silicone Emulsion – Mizzy

Debubbler

4oz Spray Bottle512-4862 1@ \$18.95
16oz Bottle.....512-5356 1@ \$54.15



Smoothex – Whip Mix

Debubbler and Wax Pattern Cleaner

A surface tension reducing liquid for treating patterns to permit better adherence and flow of investment. Also for treating synthetic impressions before pouring casts.

Pump

8oz 731-5518 1@ \$37.90

Refills

16oz731-5526 1@ \$47.55
Gallon731-5542 1@ \$271.85



Wax Remover – Buffalo Dental

Ideal to dress finished wax-ups, to completely eliminate wax from investments, to remove wax film from teeth, and to assure the bonding of plastic teeth. Contains only natural ingredients – nonchlorinated hydrocarbons, phosphates or petroleum distillates. It is water soluble, has no toxic fumes and has a pleasant citrus scent.

Gallon 159-0918 1@ \$97.50



Wax Solvent #20 – Justi
For Wax Removal

An efficient, safe and pleasant smelling solvent used for wax elimination and shining after waxing a prosthetic case.

Gallon412-5878 1@ \$72.00

Atlas Flask - Handler

Regular #44 1@ \$130.60

Lower357-2054

Upper357-2062

Giant #44G

Each357-2047 1@ \$465.95

Disk Ejector #44D

Each357-2021 1@ \$23.10



Flask Compresses – Whip Mix

As dentures cure, the acrylic materials have a tendency to expand. The HANAU Flask Compress exerts evenly distributed pressure up to 400 pounds on the flask during the curing cycle. This yields dense, compact, and accurate dentures.

Compress355-2726 1@ \$411.80



Plastipac – Yates & Bird

Plastic Separating Sheets

1,000/Box

4" x 4"767-2512 1@ \$10.30

5" x 5"767-3007 1@ \$11.80



Tin Foil – Buffalo Dental

Includes 6" wide, 1-lb rolls

.001 Thickness 158-1321 1@ \$70.25

.002 Thickness 158-1354

.003 Thickness 158-1388



Varsity Ejector Type Flask – Whip Mix

Upper355-3245 1@ \$215.10


Lower355-3278

Disc Only for Ejector Type Flasks


Each355-3765 1@ \$21.90

WAX CARVERS

Carvers – Premier



574-9445 574-9502 574-9569



574-9593 574-9627
1@ \$32.95 3@ \$30.75

PATTERSON

Akers Wax Carver

Double-ended instrument is specially designed for precision wax carving and "setup" work. Blade and spoon screw into handle and are easily replaceable — saves you money! This is the most widely used design in dental labs.












Complete Carver

(Includes Handle, Blade and Spoon)

Pkg.....083-1008 1@ \$10.90

Carvers – KerrLab

Titanium Nitride coated instrument is indicated by "TN" after the name. The second listing is the Stainless Steel instrument indicated by "SS". Each instrument half comes with a threaded connector.

TN	1@ \$32.20	SS	1@ \$28.15
			
Hollenback, Regular, TN	417-6210		Large Discoid, TN
Hollenback, Regular, SS	417-5980		Large Discoid, SS
			
Hollenback, Lateral, TN	417-6228		Small Discoid, TN
Hollenback, Lateral, SS	417-5998		Small Discoid, SS
			
1/2 Hollenback, Regular, TN	417-6236		Mini Discoid, TN
1/2 Hollenback, Regular, SS	417-6004		Mini Discoid, SS
			
1/2 Hollenback, Lateral, TN	417-6244		Large Scoop, TN
1/2 Hollenback, Lateral, SS	417-6012		Large Scoop, SS
			
Large SHOA, TN	417-6251		Small Scoop, TN
Large SHOA, SS	417-6020		Small Scoop, SS
			
Small SHOA, TN	417-6269		Large Cleoid, TN
Small SHOA, SS	417-6038		Large Cleoid, SS
			
2.0 mm Chisel, TN	417-6350		Small Cleoid, TN
2.0 mm Chisel, SS	417-6111		Small Cleoid, SS
			
1.5 mm Chisel, TN	417-6368		Le Cron, TN
1.5 mm Chisel, SS	417-6129		Le Cron, SS

WAX INSTRUMENTS

PATTERSON



Boley Gauge

With Lever Lock

Accurately calibrated and fitted.

Each084-3300 1@ \$28.55
 3@ \$27.55



Iwanson Spring Caliper

Precision Caliper

Blunt End084-3318 1@ \$27.40
 Sharp End.....084-3326 3@ \$26.45

Wax Spatula

Double-end #7 wax spatula is crafted of sturdy stainless steel.

Each092-8234 1@ \$21.05
 6@ \$20.30



Washout Brush

Stiff black bristles with a strong, wooden handle. This sturdy brush is available in flat or tufted design, each with four rows of bristles.

Flat090-9978 1@ \$6.15
 Tufted090-9986



Margin Carvers – KerrLab

This instrument has been specifically designated for carving the margins of dipped wax copings and includes: large scoop paired with 1/2 Hollenback, regular and threaded connector.

Margin Carver, TN, Titanium Nitride Coated Blades

Each414-2550 1@ \$55.15



PK Thomas – Premier

Wax Instruments


#1 - SE.....585-9301 1@ \$17.10
 #2 - DE.....585-9319 1@ \$28.80
 #3 - DE.....585-9335 1@ \$28.80
 #4 - DE.....585-9343 1@ \$28.80
 #5 - DE.....585-9327 1@ \$28.80



Wax Spatulas – Hu-Friedy

Stainless Steel

#7 Wax – Double End
 Each375-1591 1@ \$33.75



Wax Spatulas – Premier

Stainless Steel

Double End

#7A is identical to #7 except that one side of the small end is flat.

#7585-0235 1@ \$30.35
 #7A585-0268 3@ \$28.70

CROWN HOLDERS

Pic-n-Stic™ – Pulpdent


Pic-n-Stic is a 2 1/2" long plastic stick, 2mm in diameter, with an adhesive tip at one end. It easily picks up and places those small objects that fingers simply were not made to handle. Use Pic-n-Stic as a handle for crowns, inlays, veneers, bands, brackets, files, reamers, clasps, precision attachments, implant screws and numerous other objects.



Pic-n-Stic
 60/Box.....590-1053 1@ \$23.10
 3@ \$22.45

HAEMOSTATS

PATTERSON



Haemostats – Kelly


5 1/2" Straight086-3415 1@ \$19.20
 5 1/2" Curved086-3423 3@ \$16.90

Hu-Friedy

5 1/2" Straight371-0530 1@ \$73.55
 5 1/2" Curved371-0548

Miltex

5 1/2" Straight498-6758 1@ \$58.00
 5 1/2" Curved498-6741 1@ \$60.55



Halsted Mosquito Haemostats

5" Straight086-3472 1@ \$18.95
 5" Curved086-3480 3@ \$17.30

Hu-Friedy

5" Straight371-0787 1@ \$71.45
 5" Curved371-0753

Miltex

5" Straight498-6808 1@ \$56.10
 5" Curved498-6790 1@ \$58.10

TWEEZERS

PATTERSON



Cotton Pliers

Stainless steel with serrated tips to provide a secure grip. Guide pin properly aligns tips.

Standard 317

Each088-7687 1@ \$5.80
6@ \$4.75

SS Deluxe, Self-Locking 652

Each088-7711 1@ \$30.80
3@ \$29.95

Economy Self-Locking

Each088-7703 1@ \$8.40
3@ \$7.90

MODEL TRIMMERS

PATTERSON



Model Trimmer Wheel

10" x 1/4" - Coarse

Reversible. For Handler, Torit, Buffalo, Foster and Toothmaster units. Also fits old style (single-hole) Patterson model trimmers.

Each534-3595 1@ \$40.95

Diamondback™ - Premier

Diamond Model Trimming Wheel

The Diamondback™ trims models faster, runs quieter and is substantially longer lasting than traditional wheels. It has an extra-coarse diamond cutting surface, which is patterned for aggressive, yet cool cutting. It is light-weight, durable and perfectly balanced for concentric operation.

Thickness — 6.4mm (1/4")

Center Hole Diameter — 25.4mm (1")

12" Extra Coarse

Each580-6799 1@ \$747.15



DualTrim™ Dry Model Trimmer - BDM

New High Performance Rotary Dry Model Trimming Machine That Is Convertible For Wet Trimming Operation.

- Fast and Accurate Trimming
- Convertible Dry or Wet Operation
- New Exhaust System Design
- New Quick-Stick (TM) Abrasive System
- Improved Water Spray Bar Design
- New Noise Dampening Feet

Each039-6150 1@ \$765.00



Model Trimmer #10 - BDM

Buffalo Model Trimmers are engineered for greater visibility, plus efficient trimming of models without splashing or plaster build-up. Quiet and smooth running, they are ideal for any requirement - from occasional use to high production laboratory work.

Each039-6366 1@ \$720.00



Model Trimmer Wheels - Foster

12" x 1/4" - Coarse

Reversible. For Foster model trimmers.

Each324-8580 1@ \$63.05

10" x 1/4" - Coarse

Fits new style Patterson model trimmers.

Each324-8549 1@ \$51.80

Model Trimmer Wheels - Handler

For all makes of model trimmers.

12" Coarse.....357-7293 1@ \$50.50

12" Medium.....357-7301

12" Fine.....357-7319

12" X-coarse357-7327



Router - BDM

Premium trimmer provides high-torque cutting power and improved dust control! Fully equipped with built-in dust control, safety shield, and assortment of cutting tools

Each035-0959 1@ \$625.00

PLASTER TRAPS

PATTERSON



Plaster Trap

Eliminates odours, cleaning problems and mess. Built-in drain tube attached to disposable vinyl liner (included). Inlet and outlet are standard 1-1/2" FPT.

Each043-4464 1@ \$151.60



Disposable Plaster Trap

Complete kits include two hoses plus end caps.

Kits

3.5 gallon 083-1206 1@ \$61.95

5 gallon 083-1222 1@ \$66.50

Refill Components

3.5 gallon 083-1214 1@ \$30.50

5 gallon 083-1230 1@ \$36.95



Trap-EZE BT Bottle Trap – BDM

The new Trap-EZE BT provides a convenient and easy to use bottle-style trap for permanent plumbing into the existing sanitary drain beneath the lab sink. The replaceable bottles have excellent capacity for smaller labs and dental offices and are fast and easy to replace. Including 64oz bottle and 1 1/2 piping installation kit

System 133-6304 1@ \$81.50



General Purpose Vibrator – Whip Mix

The 3-speed vibrator provides quiet, long lasting performance. The vibrator features cast metal housing construction, shock absorbing base, and a removable 4-1/2" (11.5cm) diameter flanged rubber top. The flexible, easy to clean top has a raised perimeter lip that helps confine the objects being vibrated and prevents material spillage onto vibrator. An additional raised area provides delicate vibration.
White.....732-8677 1@ \$225.65



Combination Unit – Whip Mix

Two mixing speeds: vacuum invest at high speed or power mix with or without vacuum at low speed.

- Vacuvestor/Power Mixer, equipped with convenient vibrator arm
- For vacuum investing use the Vac-U-Spat, 200 or 300ml size, for inlays and crowns
- Use the Die Mixer, 25ml size (15gr cap.), for small amounts of die materials

Note: If the Combination unit, motor only, without accessories is ordered, the #5700 plastic tubing and connections will be included.

Combination Unit With Accessories Including 4100 Vac-U-Spat

White.....048-9559 1@ \$1,844.50

VIBRATORS



No. 1A Vibrator – Buffalo

Improved electronics and strong, dependable, performance make this our most popular small vibrator. Durable 4" diameter rubber cover and suction cup feet improve stability and ensure quiet operation.

4-1/2"L x 4-1/2"W x 4-1/2"H

Each 133-0752 1@ \$131.00



No. 2 Heavy Duty Vibrator – Buffalo

Ribbed 6" diameter rubber-covered platform to hold larger flasks. Three-speed control switch. Beige enameled housing.

11"W x 4-1/2"D x 6-1/2"H

Each 133-1891 1@ \$286.65

VACUUM MIXERS



Accessories – Whip Mix

Vac-U-Spat

#4120
200ml.....732-8214 1@ \$178.55

#6380
Bench Stand
White.....732-8461 1@ \$159.75

#5700
Tubing With
Connections729-1495 1@ \$81.80

O-Rings

#1 4008 5/16"
6/Pkg.....730-5246 1@ \$15.60

#2 5702 7/16"
6/Pkg.....730-7200 1@ \$18.95

#3 4005 1 1/2"
6/Pkg.....730-4264 1@ \$21.15

#4 4006 3 3/8"
1/Pkg.....730-5907 1@ \$12.10

#5 4110 5/8"
6/Pkg.....730-5733 1@ \$18.35

#6 6505 4"
1/Pkg.....730-7697 1@ \$12.75

Set 4020

Contains three of #1, one of #2, one of #3, and one of #5.

Plastic Trap #5705.....729-1453 1@ \$16.50

Metal Trap Cap #5701 ..728-5711 1@ \$36.45

Vac-U-Mixers – Whip Mix

For slow speed (425rpm) spatulation of gypsums and alginates under vacuum. Vacuum mixing of plasters, stones and die materials is the proven method to produce smoother, denser, "bubble-free" casts and dies.

#6400:
200ml, 25-100gm Gypsum, 1-2 Scoops
Alginate
Friction732-6515 1@ \$170.60
Slotted.....732-6523

#4450:
300ml, 25-150gm Gypsum, 1-3 scoops
Alginate
Friction732-6531 1@ \$170.60
Slotted.....732-6549

#6600:
875ml, 50-800gm Gypsum
Slotted.....732-7547 1@ \$192.75

#7600:
1,200ml, 100-1,000gm Gypsum
Slotted.....732-7562 1@ \$254.70

VACUUM MIXERS (CONT.)



Vacuum Power Mixer Plus - Whip Mix

Includes 300ml and 500ml Vac-U-Mixers.

- Low speed (425 rpm)
- Convenient vibrator arm
- Light tan in colour

Each048-9583 1@ \$1,677.75



Vac-U-Timer - Whip Mix

- Accurate mixing time for stone & investment
- Programmable up to 99 minutes
- Easy-to-read LCD display
- Touch tone buttons include countdown, alarm shutoff, memory activation
- Velcro strip for easy positioning

Each048-9682 1@ \$46.10

VV Lubricant for Oiler - Whip Mix

For Combination Unit, Power Mixer, Vac-U-Vestor


4oz Bottle.....728-4441 1@ \$15.70

WAXING UNITS

Mini-Waxer - Almore

Includes three PKT tips (#2 small, #3 medium and #7 spoon tip).

Each105-4436 1@ \$290.35



Mastertouch Tips - KerrLab

- | | | |
|----------------|----------|------------|
| 1. B20..... | 428-5011 | 1@ \$42.05 |
| 2. B22..... | 428-5029 | |
| 3. B24..... | 428-5037 | |
| 4. B28..... | 428-5045 | |
| 5. LB20..... | 428-5052 | |
| 6. LB22..... | 428-5060 | |
| 7. FS 20..... | 428-5078 | |
| 8. FS 22..... | 428-5086 | |
| 9. RN20..... | 428-5094 | |
| 10. RN22..... | 428-5102 | |
| 11. RN24..... | 428-5110 | |
| 12. S20..... | 428-5128 | |
| 13. RFS20..... | 425-5136 | |
| 14. RFS22..... | 428-5144 | |
| 15. P22..... | 428-5151 | |



Micro Slip-On Tips - KerrLab

- | | | |
|------------------------|----------|------------|
| 1. Needle Tip..... | 428-4568 | 1@ \$34.60 |
| 2. Small PKT..... | 428-4998 | 1@ \$34.60 |
| 3. Small #7..... | 428-4725 | 1@ \$41.05 |
| 4. Large PKT..... | 428-4733 | 1@ \$41.05 |
| 5. Large #7..... | 428-4576 | 1@ \$32.20 |
| 6. 1/2 Hollenback..... | 428-4741 | 1@ \$41.05 |
| 7. Large Denture..... | 428-4758 | 1@ \$41.05 |
| 8. Small Denture..... | 428-4766 | 1@ \$41.05 |
| 9. Beavertail..... | 131-5217 | 1@ \$34.60 |



System Waxing Station - KerrLab

Includes:

- 27" x 15.5" rigid aluminum base with composite surface
- Ultra-waxer w/three Ultra-spatulas™ heating tips
- Digital dipping unit
- Wet-Waxer with four belleWax samples
- Shortened power cords for equipment
- 16" goose neck lamp with 60-watt light capsule
- 4" x 2", 2X magnifier with 4X spot magnification
- Kleen Lube and Stick Lube die lubricants
- Four-outlet power strip for one-plug electrical supply

C & B Station439-9200 1@ \$1,735.05
Denture Station416-0818 1@ \$1,041.70



Ultra Waxter™ – KerrLab

State-of-the-art electronic wax pen. Total wax control. Exclusive feature offers 4 programmable °C/°F pre-sets.

- Four programmable temperature pre-sets
- “Quick Heat” cycle for rapid heating of tips
- Temperature range 60°C/100°F to 370°C/700°F
- Continuous heating or foot control mode
- Selectable °F or °C temperature display
- All aluminum lightweight handle
- Includes one Ultra-Waxter™ Tip #27599

W/ FC Sm PKT439-1017 1@ \$673.40



Ultra Spatula Heating Tips – KerrLab

Short

1. Sm. PKT, Blue.....	439-0894	1@ \$89.35
2. Lg. PKT, Teal	439-0902	
3. Beavertail, Purple .	439-0910	
4. Sm. #7, Grey	439-0928	
5. Lg. #7, Green	439-0936	
6. Needle, Red	439-0944	
7. 1/2 Hollenback, Black	439-0951	



Waxlectric Modeling Tips – Renfert

- 9 different tips
- Heating element is situated within the tip
- Quick and easy tip replacement

1. Probe, SM.....	267-4877	1@ \$53.20
2. Probe, MD	267-4885	
3. Probe, LG	267-4893	
4. Blade, Narrow	267-4901	
5. Blade, Wide	267-4919	
6. Blade, Chamfer.....	267-4927	
7. Hollenback, Horizontal.....	267-4935	
8. Hollenback, Vertical	267-4943	
9. Probe, Mini Angled.....	267-4950	



Hotty E – Renfert

Wax die pot with adjustable working temperature from 65-105°C. The temperature is set via a potentiometer on the top of the unit. As with all Hotty models the unit boasts minimal temperature fluctuations (+/-0.5°C).

Each267-3200 1@ \$248.65

DIPPING POTS



Digital Wax Pot – Yates & Bird

Precision controls take the guesswork out of dipping wax. Easy-to-use, accurate and efficient machine quickly reaches and maintains desired temperature, so it makes perfect copings every time. Temperature range 50°-115°C (122°-239°F). Includes green dipping wax.

Each766-7892 1@ \$283.75



Hotty LED – Renfert

Hotty is equipped with modern low voltage technology. This involves use of an external transformer.

- The unit itself is ergonomically compact
- The LED display shows the current temperature
- Heats up the post continuously until it reaches the desired temperature
- Touch control panel for easy adjusting and cleaning

110V.....267-1733 1@ \$450.15



Hotty – Renfert

The cost-effective wax die pot with preset temperature of 90°C. The temperature is set via a potentiometer on the back of the unit. The temperature range is 65-105°C. As with all Hotty models the unit boasts minimal temperature fluctuations (+/-0.5°C).

Each267-3218 1@ \$201.70

DIPPING POTS (CONT.)

Wax Dipping Unit – KerrLab

Digital Micro Dip

This electronic wax dipping unit precisely controls the thickness of dipping wax copings. LED digital readout displays the desired temperature setting. Its superior accuracy will maintain the dipping wax to within 1.0°F of the displayed setting.

(3.75" W x 3.75" D x 3" H)

Micro Dip Digital Unit Complete (with wax and instructions)

Each428-4501 1@ \$439.15

Standard Micro Dip

Precise temperature control prevents wax overheating. Since the wax working properties are never destroyed by overheating, you never have to discard dipping wax. Numbered dial increments allow temperature settings in 10°F increments.

Micro Dip Standard Unit Complete (with wax and instructions)

Each428-4519 1@ \$204.85

Dip Well Insert

Each420-0101 1@ \$14.35



Wet-Waxer™ – KerrLab

- Uniquely designed receptacles maintain various lab waxes at working temperatures including a warming bath for softening sprues and wax forms
- Does not overheat wax, thus ensuring consistent working properties
- Programmable temperatures
- Compact design
- Comes with four square tins with wax samples and one empty rectangular tin

Each439-9218 1@ \$364.80



STEAM CLEANERS

Super Steam Cleaner – Vident

- Self-contained 1/2 gallon water tank
- Three separate safety devices
- 70-75psi
- Full day of steam for most labs

Each048-8338 1@ \$1,485.00



WAX POTS

Waxing Unit – Buffalo

- Power cord switch
- Three compartment
- Adjustable temperature, 140°-210°F
- Front panel indicator light

Each039-6325 1@ \$98.75



AIR DRIVEN HANDPIECE

No. 220 Lab Air Turbine Handpiece – Buffalo

Produces 350,000 rpm for faster cutting. Complete 220 includes handpiece, foot control, regulator filter with bracket, automatic lubricator, tubing, fittings, replacement chuck, chuck tool and turbine oil.

Each039-6606 1@ \$525.00

220 Replacement Cartridge
Each61050080 1@ \$285.25

220 Metal Chucks
Each61050092 1@ \$28.75



ELECTRIC HANDPIECES

PATTERSON



Electric Lab Engine NC-350II

This new system is released a new with higher controllability, functionality and durability.

- Rotation speed: 1,000 - 35,000rpm
- Maximum torque: A@410gf - cm
- Quick bur change by twisting handpiece sleeve
- Lightweight – 6.95oz (without cable)
- Handpiece housing is thermally insulated and virtually indestructible
- Available on/off or variable foot pedals

Each043-4233 1@ \$1,256.75



M-35II Handpiece System – BDM

Provides excellent cutting power and smooth control in a compact system. Improved controls provide fully variable speed via console or foot control. System excels at providing a full range of usable operating speeds from ultra-low to 35,000 RPM.

M-35II Handpiece System, 120 V AC, including Control Console, Handpiece, Foot Control)
Each035-1148 1@ \$710.00



X35 Micro Engine System - 110v - Buffalo

- High torque, cool running handpiece
 - Built-in sensor prevents overloads, and provides more power
 - Convenient twist-lock collet
 - Solid state all metal non-slip foot control
- Each039-6358 1@ \$795.00

X35 Micro Engine Replacement Carbon Brushes
 Each61010021 1@ \$53.10

X35 Micro Engine Handpiece Only
 Each61050062 1@ \$462.30

SANDBLASTERS



Opti™ Workstation - Buffalo

- Cabinet for Optiblast or other blasters
- Available with or without built-in suction
- Remote outlet in rear for electric handpiece
- Compact size 12" x 11" x 8-1/2", portable

Opti Workstation w/o Suction & Light 120VAC
 Each039-6168 1@ \$425.00



Programat CS Furnace - Ivoclar Vivadent

Ivoclar Vivadent has developed a ceramic and crystallization furnace for the dentist. The Programat CS is an easy-to-operate, small furnace for the dental office. It is particularly suitable for glaze and crystallization firing of IPS e.max CAD LT blocks. The furnace features 20 programs and is very easy to use. Because of its small dimensions, the furnace occupies very little space in the dental office. Thanks to the vacuum built-in function, more complex firing cycles with high-quality results can be conducted.

- Easy operation
 - Automatic temperature calibration
 - State-of-the-art design
 - 3 pre-programmed standard programs
 - 17 individually created programs
 - QTK (Quartz Tube Kanthal) heating technology
 - Fast, turbo opening system (5 seconds)
 - Compact ergonomic design
 - 2-year warranty
- Each048-8858 1@ \$5,495.00

FURNACES



Centurion Q50 - Dentsply/Neytech

- 50-program memory
 - Quartz sleeve muffle
 - Continuous step technology
 - Enhanced temperature stability
 - Closed muffle drying
 - Programmable muffle movement
 - Fast Cooling
 - Bright, backlit LCD readout
 - Icon graphic controls
 - Auto calibration, auto purge, auto night mode
 - Program override
 - Large 10 cm (4") work platform
 - Internal surge/spike protection
 - Program advance key
 - °C to °F conversion
 - Full two-year warranty
- Each042-4010 1@ \$4,871.95

LATHES AND ACCESSORIES

PATTERSON



Arbor Band Chuck
For Standard Lathes with 3/8" Shaft
 Left Hand083-9316 1@ \$15.45
 Right Hand.....083-9324

Arbor Band Mandrel
Handpiece Shank
 For 3/8" band, complete with rubber.
 Each088-6929 1@ \$4.60

Stone Chuck
 Steel; for standard lathes with 3/8" shaft.
 Right Hand.....083-9357 1@ \$16.35

LATHES AND ACCESSORIES (CONT.)

Bench Polishing Lathes – BDM

High Boy	039-6382	1@ \$400.00
Low Boy	039-6390	1@ \$670.00

With Quick Chuck.



Bench Polishing Lathes – Handler

Handler lathes offer two-speed operation (1725/3450 rpm) with a full 1/4 hp motor. Manufactured to provide years of durable service. Two-year warranty against defects or workmanship.

Model 26A		
4 1/2" Shaft Height	039-8255	1@ \$352.75
Model 26		
6" Shaft Height	039-8263	1@ \$375.60

Lathe Splash Pans – Handler

<i>Aluminum</i>		
Splash Pan – Without Light		
12" high x 6" wide x 10" deep with insert pan.		
Each	357-3367	1@ \$87.20
		3@ \$82.80

Splash Pan – With Light & Sta-Klean Shield		
Standard #83WL splash pan with an #86B Sta-Klean Shield.		
Right	357-5362	1@ \$102.50
Left	357-5370	

Lathe Splash Pan – Clive Craig

<i>All Plastic, One Piece, Seamless White Colour</i>		
7-3/4" high x 6-7/8" wide x 10" deep.		
Each	228-1830	1@ \$56.90

Lathe Splasher – Buffalo Dental

Lathe Splasher only		
Each	157-6727	1@ \$92.30
Lathe Splasher with Tray and Light		
Each	157-6735	1@ \$106.50

Arbor Band Chucks – Buffalo

For Standard Lathes
 Emery Cloth Arbor – #100
 Has a rubber head that is adjusted with a turn of the screw in the end of the spindle.
 Right

145-8744	1@ \$14.25
Left	145-8777



Jacobs Type Chuck – Buffalo

Keyless, True-Running Chuck
 Completely sealed. Made for right side only.
 Each

150-8159	1@ \$99.00
----------	------------

Taper Chucks – Buffalo

For Brush Wheels and Buffs
 3 1/2" Length
Steel
 Right

150-8639	1@ \$9.40
----------	-----------

Left	150-8621
------------	----------

Brass	1@ \$11.45
--------------	------------

Right	150-8704
-------------	----------

Left	150-8670
------------	----------

Jacobs Type #18 Chuck – Handler

Atlas Precision
 Fits right shaft only.
 Each #18

357-2377	1@ \$210.70
----------	-------------

AIR ABRASION

PATTERSON



Aluminum Oxide

Moisture protected packaging. 50 micron.
White Abrasive
 50lb.....

083-1602	1@ \$143.05
----------	-------------



Aluminum Oxide 50 Micron – Danville

1lb (White).....	236-4099	1@ \$13.70
------------------	----------	------------

Aluminum Oxide – Deldent/Mant

50 Micron		
2lb	252-5566	1@ \$27.50
90 Micron		
2lb	252-5558	1@ \$31.25

Emery Powder – Protech Professional

25lb.....	587-9085	1@ \$60.65
50lb	587-9093	1@ \$101.80

KCP Evacuation Filter – ADT

Each	58060151	1@ \$59.10
------------	----------	------------

KCP Preventive Maintenance Kits – ADT

For KCP1000		1@ \$60.35
Each	58060232	

For KCP2000		1@ \$60.85
Each	58060210	

KCP Replacement Nozzles – ADT

Nozzle, .015 Standard	1@ \$155.90
2/Package	58010138

Nozzle, .015 - 20 degree	1@ \$155.85
2/Package	58010182

Nozzle, .015 - 100 degree	1@ \$159.00
2/Package	58010181

Air Scribe® – Chicago Pneumatic

Air chisel for perfect finishing of stone models, virtually dust free.

Complete Kit CP-9361-1
 Contains: Air Scribe, 5 points, air hose, air filter assembly and carrying case
 Kit Assembly

199-0001	1@ \$880.00
----------	-------------

Air Scribe only	
Each	199-0019
	1@ \$448.00





Microetcher II – Danville

The Microetcher II Intraoral Sandblaster dramatically enhances dental bonding by micro-roughening.

Precision air abrasive stream with finger button control. It is the most advanced air abrasive handpiece in the marketplace. Minimum pressure of 50 psi (3.5 bar) air supply is required.

Complete Kit

Contains: Microetcher II, 60° nozzle, 9 foot air line, two abrasive jars, sample abrasive and instructions.

Kit039-6614 1@ \$305.70



Optiblast – Buffalo

Microblaster & Workstation with Suction

- Low-cost, non-recycling blaster for use in cleaning alloy in the lab
- FDA approved for intraoral use in the operatory
- Precise fingertip controls for instant on/off
- Can use aluminum oxide abrasives and glass beads
- Built-in quiet suction + bright, even lighting
- Use as a cabinet for blasting or as a workstation to grind appliances
- Remote plug in rear for electric handpieces
- Compact 12"D x 11"W x 8-1/2"H

Optiblast.....039-6283 1@ \$605.00

Optiblast Handpiece

Each155-7354 1@ \$227.00

Pattern Resin™ – GC Lab Technologies

Acrylic Die Modifier

1:1 Package

Contains: 100gm powder, 100ml liquid, mixing cups, brush, measure.

Pkg326-3654 1@ \$103.00

DUST COLLECTORS

60U Dust Collector – Handler

- 1/2hp capacitor motor
- 2-3" suction outlets, 200 CFM each - 400 CFM when using only one outlet
- Shakeable filter bags (14) with special poly bags for aluminum oxide
- 24" clearance with silencer



60U-1 Dust Collector – Handler

- Special poly bags for aluminum oxide

66 Porta-Vac Dust Collector – Handler

- Bench top unit with splash hood
- 150 CFM

66A Porta-Vac Dust Collector – Handler

- Under bench unit with splash hood
- 150 CFM

66D Redi-Vac With Magnetic Flange – Handler

- 150 CFM



101 Dust Collector – Handler

A 1 hp unit designed to provide suction up to four 3" (7.6cm) work positions, simultaneously. The 101 is provided with four 3" (7.6cm) inlets. 26"W x 24"D x 41"H. Suction - 800 CFM.



62 Super Sucker – Handler

- Charcoal filtration reduces acrylic monomer odours
- Suction from two 3" inlets, ultra quiet operation
- Secondary filtration entraps fine dust particles to 0.1 microns
- Compact size of 14" x 18" x 8"

* Please call for pricing

POLISHING/GLAZING

PATTERSON



Arbor Bands

Coarse Grit

100/Box

3/8"083-8094 1@ \$18.50

3/4"083-8102 10@ \$18.20



Pumice

White Pumice Stone for Polishing

3lb (1.36 kg) Container

Flour (Prophy).....088-7554 1@ \$16.00

4@ \$15.00

Fine088-7521 1@ \$20.75

Medium.....088-7539 4@ \$19.05

25lb Container 1@ \$53.90

Fine088-7364

Medium.....088-7372

Acrilustre – Buffalo

Greaseless Polishing Compound

For use on acrylic denture materials.

Used dry on cloth buffs.

1lb Pkg145-6268 1@ \$21.25



Bands arbor – BDM

100/Box

1/4" Box 1@ \$9.50

Fine.....145-9551

Medium.....145-9569

Coarse.....145-9502

3/8" Box 1@ \$9.50

Fine.....145-9593

Medium.....145-9601

Coarse.....132-9820 1@ \$9.85

1/2" Box 1@ \$9.50

Fine.....145-9536

Medium.....145-9528

Coarse.....145-9544

3/4" Box 1@ \$9.50

Fine.....145-9577

Medium.....145-9585

Coarse.....132-9812 1@ \$9.85

POLISHING/GLAZING (CONT.)

Cerama Seal/Sep – KerrLab

Chemically pure, noncontaminating two-liquid system with brush cap applicators. Allows building of porcelain directly to model without porcelain discoloration or sticking. May also be used for porcelain margins or lapping.

#003000 Kit		
Each	416-6005	1@ \$35.05
#003001 1/2oz Cerama Seal		
Each	131-4327	1@ \$17.70
#003002 1/2oz Cerama Sep		
Each	131-5670	1@ \$22.50

Composite Accessories – Kerr

Luster Paste

Squeeze Bottle

10gm Bottle	419-2217	1@ \$80.40
		3@ \$78.00

Diaglaze – Antraco

Porcelain Polish

2gm Syringe	114-6208	1@ \$24.40
-------------------	----------	------------

Diaglaze Refill – Deldent/Mant

Porcelain Polish

Each		1@ \$25.15
X-Fine	252-5624	
Fine	252-5632	
Regular	252-5640	

Econo-Cutters – E.C. Moore

Cut-Off and Trimming Wheels

For plastic or metal. Fast-cutting, long-lasting; resists breakage.

2-5/16" diameter with 1/4" hole for use with Equa Flange mandrel.

Each	518-7513	1@ \$3.25
------------	----------	-----------

#2 1-1/2" diameter. Designed for use with standard handpiece mandrels.

Each	518-8008	1@ \$2.10
------------	----------	-----------

Faskut® Wheel – Dentsply

White Aluminum Oxide

Ideal for trimming denture base material, grinding stone casts and shaping baseplate material. 3" x 3/8" – 1/4" hole.

Each	275-8738	1@ \$14.20
------------	----------	------------



Insta-Glaze Kit – Taub

Porcelain Polishing Paste

Complete Package

Contains: 2gm dispensing syringe, felt wheel, cone, knife edge, and rubber wheel.

Kit	670-5404	1@ \$47.25
-----------	----------	------------



Metal Etch-It Kit – American Dental Supply

For Non-Precious or Semi-Precious Metal

Contains: 1/2oz Etch-it, 3 pkgs Dip-it, plastic spatula and Etch-it labels.

Each	109-6916	1@ \$45.00
------------	----------	------------



Lathe Wheels – Dedeco

3" x 3/8"

Soft

The traditional favourite for fastest grinding of porcelain teeth.

Green	255-9060	1@ \$20.85
-------------	----------	------------

Medium

For fast grinding of porcelain teeth.

Black	255-9078	1@ \$20.85
-------------	----------	------------

Medium

Rubber-bonded wheel for rapidly removing scratches from porcelain and metals.

Green	255-9110	1@ \$20.85
-------------	----------	------------

Coarse

Rubber-bonded wheel for fast, clog-free trimming of baseplates and acrylic dentures.

Acrylic Base Plate	255-9102	1@ \$20.85
--------------------------	----------	------------

Coarse

For fast, clog-free grinding of acrylics.

Red Acrylic	255-9128	1@ \$20.85
-------------------	----------	------------

Coarse

Coarsest for trimming plaster and stone models. Also suitable for acrylics.

Lead center.

White Knock Down	255-9136	1@ \$20.85
------------------------	----------	------------

Milling & Polishing Cream – Justi

All-Purpose Cream – Pink

For milling, finishing and intraoral polishing. Pleasant tasting.

3oz Tube	410-6282	1@ \$65.20
----------------	----------	------------

Moldent – Motloid

Buffing Agent – Red

1lb Bar	522-7145	1@ \$18.00
---------------	----------	------------

Porcelain Etch-It Kit – American Dental Supply

For all-porcelain crowns and veneers, inlays and onlays. Contains: 1/2oz

Etch-it, 3 pkgs Dip-it, plastic spatula and Etch-it labels.

Kit	109-6924	1@ \$46.20
-----------	----------	------------



Pumice – Moyco/Union Broach

For Polishing White Pumice Stone

4lb Bottle 1@ \$12.15

Coarse

527-9724

5lb Bottle 1@ \$12.05

Fine

527-8650

Flour

527-8684

Medium

527-8759

Rapid Glaze – Yates & Bird

Diamond Polishing Paste

Eliminates need for reglazing after adjusting porcelain.

- Eliminates costly refiring
- Saves chair & lab time
- Instant adjustments & insertions
- Safe, clean & easy to use
- Convenient syringe applicator

2gm	767-3940	1@ \$27.45
		6@ \$26.20

Red Rouge Stick – Buffalo

Final Polish On All Gold Work

Used on dry cloth buffs.

1/2lb Stick	156-2842	1@ \$7.00
-------------------	----------	-----------

The Porcelain Wash – American Dental Supply

- Cleans and decontaminates soiled porcelain restorations
- Will not affect shades
- Leaves porcelain microscopically clean

- Prevents fracturing and gassing of porcelain cases previously in the mouth

Quart	109-7724	1@ \$25.80
-------------	----------	------------



Tripoli – Buffalo

Secondary Polishing

For denture materials and metals. Will remove scratches and prepare the work for final polishing.

1lb Bar	158-7260	1@ \$7.55
---------------	----------	-----------



MISCELLANEOUS

PATTERSON



Denture Boxes

Plastic

12/Pkg

Beige.....084-3714	1@ \$12.45
Blue.....084-3722	6@ \$11.45
Mauve.....084-3730	



Lab Work Pans

Poly, Medium Back w/Clip

Red.....279-4543	1@ \$4.10
Light Blue.....279-4527	6@ \$3.60
Green.....279-4535	
Yellow.....279-4576	
White.....279-4568	
Tan.....279-4550	

Denture Bags – Buffalo

Clear – 5" x 6"

500/Box.....146-4627	1@ \$28.50
----------------------	------------



Denture Baths – BDM

For safe convenient cleaning and storage of dentures.

White Bath133-2543	12@ \$11.15
Beige Bath.....133-2576	
Hot Pink Bath133-2550	
Light Blue Bath.....133-2504	
Yellow Bath.....133-2568	
Dusty Rose Bath133-2592	

Den-Shur Cups® – Modalities Research

Poly

Assorted colours per package: white, blue and pink.

12/Pkg.....513-3400	1@ \$17.35
	3@ \$16.45



Denture Cups – Dental Resources

Each	1@ \$1.55
Dark Blue.....279-6613	36@ \$1.30
Dusty Rose.....279-6621	144@ \$1.15
French Vanilla.....279-6639	
Light Blue.....279-6647	

Denture Cups – Union Broach

Assorted	
12/Pack680-3423	1@ \$31.25



Denture Boxes – Zirc

3-3/4" long x 1-3/4" deep

12/Pkg	
Blue.....771-0999	1@ \$15.40
Mauve.....771-1005	
Beige.....771-1013	
Yellow.....771-1021	
Assorted771-1039	



Smile Again Denture Cleaner – BDM

A revolutionary denture cleaner and disinfectant all in one!

Pre-measured in convenient foil packets with a fresh mint flavour. Removes tough coffee, tea, iron and nicotine stains. Each packet lasts seven to ten days. Used as directed, safe for all dentures including soft relines, partial dentures and orthodontic appliances, orthodontic night guards and retainers.

* Not Recommended for use on partials fabricated from Ticonium® Alloy.

Intro pack of 12133-2451	1@ \$27.50
Economy Pack of 48...133-2527	1@ \$101.55

Individuo Lux – VOCO

Light-Cured Custom Tray Material

- Very stable and nonsticky trays
- Peppermint flavour and pleasant smell
- Short curing time under halogen light devices or UVA light

Blue Opaque, 50/Pkg	
Maxillae.....711-2020	1@ \$112.25
Mandibles711-3038	
Transparent Blue, 50/Pkg	
Maxillae.....706-1336	1@ \$112.25
Mandibles711-3061	



Lab Work Pan – T&S Dental

Poly

Red.....678-4011	1@ \$4.15
Blue.....678-4029	6@ \$3.30
Green.....678-4037	
Yellow.....678-4045	
Pink.....678-4052	
White.....678-4060	
Beige.....678-4078	

Lab Work Pans – Clive Craig

Poly

7-3/8" x 4-7/8" x 2-1/2" deep, 7-3/8" back

Blue.....228-0717	1@ \$7.65
White.....228-0725	6@ \$7.40
Yellow.....228-0733	
Green.....228-0741	
Red.....228-0758	
Cocoa.....228-0766	

Plastic Lab Job Pans – Zirc

- Organizes and transports impressions, models, dentures, and crown and bridge appliances from operator to lab to operator on a daily basis
- Colour-coded as per day usage
- Stainless-steel spring card holder
- Stacking guides – no sticking
- Available in six colours

Lab Job Pan Plastic Counter 1@ \$5.50

White771-3514	Green771-3548
Blue771-3522	Coral771-3563
Yellow.....771-3530	Beige771-3571



Wall Hanging Lab Job Pans – Zirc

- Same as standard version but with strip on top back for hanging on any vertical surface
- Saves on installing expensive shelving

Lab Job Pan Plastic Wall 1@ \$6.60

White771-3613	Green771-3647
Blue.....771-3621	Coral.....771-3662
Yellow.....771-3639	Beige771-3670

Wall Strips – Zirc

- Drilled and countersunk with screws and toggle bolts

Wall Strip, 5" Holds 1 Pan
Each771-3704 1@ \$4.50

Wall Strip, 18" Holds 3 Pans
Each771-3712 1@ \$10.00