

Orthodontics

PATTERSON

DIRECT BOND METAL BRACKETS



A Roth System



Vertical Slot



An Alexander System



Standard Edgewise

Patterson Brand Direct Bond Metal Brackets are available in four prescriptions: A Roth System, Vertical Slot, An Alexander System and Standard Edgewise. These systems offer:



- Quality brackets
- Unsurpassed bond strength
- Pre-adjusted systems feature hooks on cuspids and bicuspids
- Mini bracket systems
- Pre-adjusted systems feature Parallel Torque and Parallel Angulation Metal Injection Molded (MIM) stainless steel brackets

Patterson Brand Metal Brackets use parallel torque to achieve superior results.

Parallel Torque Brackets

Slot walls and base walls are parallel.

Facial surface and base are parallel.




Benefits of Parallel Torque

Slots are always centered over base for easier, more accurate bracket placement during bonding.

Archwire force is transmitted directly to center of tooth's crown.

Why is parallel torque superior to torque in slot and torque in base?


Either slot walls and base walls or facial surface and base are parallel, but not both.



Torque in Slot



Torque in Base




Problems with Torque in Slot

In order to line up slots, bracket bases must vary in height, increasing difficulty and decreasing accuracy during bonding.

Problems with Torque in Base

Torque in Base brackets stick out further facially increasing the chance for occlusal interference and patient discomfort.




Black: Parallel Torque
Red: Torque in Base

Patterson Brand Metal Brackets use parallel angulation to achieve superior results.


Parallel Angulation Brackets

Provide visual references which are parallel with tooth anatomy for easier, more accurate placement during bonding.



Why is parallel angulation superior to brackets with angulation in the face?

Conflicting visual references make it harder to accurately place brackets.



Vertical Slot System

With Hooks on Cuspids and Bicuspids

Five Patient Kit contains: 100 brackets total

Upper and lower, 2nd bicuspid to 2nd bicuspid, 5 sets

.018 Five Patient Kit093-9025

1@ \$429.15

.022 Five Patient Kit093-9033

TORQUE. ANGULATION, ROTATION - REFERENCE CHART

MAXILLARY

	A Roth System	V-Slot	An Alexander System	Standard Edgewise
Central				
Torque	+12	+14	+12	0
Angulation	+5	+5	+5	0
Lateral				
Torque	+8	+8	+8	0
Angulation	+9	+9	+9	0
Cuspid				
Torque	-2	-2	-3	0
Angulation	+10	+10	+10	0
Rotation	4	0	0	0
1st Bicuspid				
Torque	-7	-7	-7	0
Angulation	0	0	0	0
Rotation	2	0	0	0
2nd Bicuspid				
Torque	-7	-7	-7	0
Angulation	0	0	0	0
Rotation	2	0	0	0

MANDIBULAR

	A Roth System	V-Slot	An Alexander System	Standard Edgewise
Central				
Torque	0	0	-5	0
Angulation	0	0	0	0
Lateral				
Torque	0	0	-5	0
Angulation	0	0	0	0
Cuspid				
Torque	-11	-7	-7	0
Angulation	+7	+7	+6	0
Rotation	2	0	0	0
1st Bicuspid				
Torque	-17	-11	-11	0
Angulation	0	0	0	0
Rotation	4	0	0	0
2nd Bicuspid				
Torque	-22	-17	-17	0
Angulation	0	0	0	0
Rotation	4	0	0	0

Individual Brackets

Maxillary Arch, 5/pkg

1@ \$22.90


10@ \$22.35

Tooth		Torque	Angulation	.018	.022
Central	UR	+14	+5	093-9561	093-9710
	UL			093-9579	093-9728
Lateral	UR	+8	+9	093-9587	093-9736
	UL			093-9595	093-9744
Cuspid	UR	-2	+10	093-9603	093-9751
	UL			093-9611	093-9769
1 st Bicuspid	UR	-7	0	093-9629	093-9777
	UL			093-9637	093-9785
2 nd Bicuspid	UR	-7	0	093-9629	093-9777
	UL			093-9637	093-9785

Mandibular Arch, 5/pkg

Tooth		Torque	Angulation	.018	.022
Central	LR/LL	0	0	093-9645	093-9793
	LR/LL	0	0	093-9645	093-9793
Lateral	LR	-7	+7	093-9652	093-9801
	LL			093-9660	093-9819
1 st Bicuspid	LR	-11	0	093-9678	093-9827
	LL			093-9686	093-9835
2 nd Bicuspid	LR	-17	0	093-9694	093-9843
	LL			093-9702	093-9850

Bracket Identification



Patterson's direct bond metal brackets are colour-coded for easy identification. Green designates patient's upper right, Blue upper left, Yellow lower right, and Red lower left. Universal right/left is designated with no colour code.

DIRECT BOND METAL BRACKETS (CONT.)



A Roth System

With or Without Hooks on Cuspids and Bicuspids
 Five Patient Kit contains: 100 brackets total
 Upper and lower, 2nd bicuspid to 2nd bicuspid, 5 sets

With Hooks

.018 Five Patient Kit	093-9009	1@ \$390.05
.022 Five Patient Kit	093-9017	10@ \$369.05

Without Hooks

.018 Five Patient Kit	093-9082	1@ \$390.05
.022 Five Patient Kit	093-9090	10@ \$369.05

Individual Brackets


Maxillary Arch, 5/pkg	1@ \$20.50
	10@ \$20.15

Tooth		Torque	Angulation	Rotation	With Hooks		Without Hooks	
					.018	.022	.018	.022
Central	UR	+12	+5	0			093-9132	093-9280
	UL						093-9140	093-9298
Lateral	UR	+8	+9	0			093-9157	093-9306
	UL						093-9165	093-9314
Cuspid	UR	-2	+10	4	093-9173	093-9322	093-8563	093-8662
	UL				093-9181	093-9330	093-8571	093-8670
1 st Bicuspid	UR	-7	0	2	093-9199	093-9348	093-8589	093-8688
	UL				093-9207	093-9355	093-8597	093-8696
2 nd Bicuspid	UR	-7	0	2	093-9199	093-9348	093-8589	093-8688
	UL				093-9207	093-9355	093-8597	093-8696

Mandibular Arch, 5/pkg

Tooth		Torque	Angulation	Rotation	With Hooks		Without Hooks	
					.018	.022	.018	.022
Central	LR/LL	0	0	0			093-9215	093-9363
Lateral	LR/LL	0	0	0			093-9215	093-9363
Cuspid	LR	-11	+7	2	093-9223	093-9371	093-8605	093-8704
	LL				093-9231	093-9389	093-8613	093-8712
1 st Bicuspid	LR	-17	0	4	093-9249	093-9397	093-8621	093-8720
	LL				093-9256	093-9405	093-8639	093-8738
2 nd Bicuspid	LR	-22	0	4	093-9264	093-9413	093-8647	093-8746
	LL				093-9272	093-9421	093-8654	093-8753

Bracket Identification



Patterson's direct bond metal brackets are colour-coded for easy identification. Green designates patient's upper right, Blue upper left, Yellow lower right, and Red lower left. Universal right/left is designated with no colour code.

Bracket Organizer

Used chairside to organize brackets prior to bonding.

50/Pkg

Package	093-8530	1@ \$52.95
		3@ \$50.80



If you are treating patients for orthodontic reasons, consider updating or adding a panoramic/cephalometric X-ray unit to your practice.

Please contact your Patterson representative or local branch for more information.



An Alexander System

With Hooks on Cuspids and Bicuspids

Five Patient Kit contains: 100 brackets total

Upper and lower, 2nd bicuspid to 2nd bicuspid, 5 sets

.018 Five Patient Kit093-9041 **1@ \$390.05**

.022 Five Patient Kit093-9058 **10@ \$369.05**


Individual Brackets
Maxillary Arch, 5/pkg **1@ \$20.50**
10@ \$20.15

Tooth		Torque	Angulation	.018	.022
Central	UR	+12	+5	094-0015	094-0130
	UL			094-0023	094-0148
Lateral	UR	+8	+9	094-0031	094-0155
	UL			094-0049	094-0163
Cuspid	UR	-3	+10	094-0056	094-0171
	UL			094-0064	094-0189
1 st Bicuspid	UR/UL	-7	0	094-0072	094-0197
2 nd Bicuspid	UR/UL	-7	0	094-0072	094-0197

Mandibular Arch, 5/pkg

Tooth		Torque	Angulation	.018	.022
Central	LR/LL	-5	0	094-0080	094-0205
Lateral	LR/LL	-5	0	094-0080	094-0205
Cuspid	LR	-7	+6	094-0098	094-0213
	LL			094-0106	094-0221
1 st Bicuspid	LR/LL	-11	0	094-0114	094-0239
2 nd Bicuspid	LR/LL	-17	0	094-0122	094-0247

Bracket Identification
Brackets are colour coded by quadrant.



Patterson's direct bond metal brackets are colour-coded for easy identification. Green designates patient's upper right, Blue upper left, Yellow lower right, and Red lower left. Universal right/left is designated with no colour code.



Standard Edgewise System

Five Patient Kit contains: 100 brackets total

Upper and lower, 2nd bicuspid to 2nd bicuspid, 5 sets

.018 Five Patient Kit.....093-9066 **1@ \$313.15**

.022 Five Patient Kit.....093-9074 **10@ \$303.05**


Individual Brackets
Maxillary Arch, 5/pkg **1@ \$15.95**
10@ \$15.65

Tooth		Torque	Angulation	.018	.022
Central	UR/UL	0	0	094-0411	094-0460
Lateral	UR/UL	0	0	094-0429	094-0478
Cuspid	UR/UL	0	0	094-0437	094-0486
1 st Bicuspid	UR/UL	0	0	094-0445	094-0494
2 nd Bicuspid	UR/UL	0	0	094-0445	094-0494

Mandibular Arch, 5/pkg

Tooth		Torque	Angulation	.018	.022
Central	LR/LL	0	0	094-0452	094-0502
Lateral	LR/LL	0	0	094-0452	094-0502
Cuspid	LR/LL	0	0	094-0437	094-0486
1 st Bicuspid	LR/LL	0	0	094-0445	094-0494
2 nd Bicuspid	LR/LL	0	0	094-0445	094-0494

Bracket Identification
Brackets are not colour coded.



High Output - Fast Cure!

Curing Light Halogen

- Efficient - fast ten second cure
- Ergonomic - clean compact design
- Economical - high power, not price
- Extra long (2m) flexible cord
- Audible trigger start and stop
- Quick and simple bulb replacement
- Ten second beeper with auto-shut-off
- Thermal protection against over-heating
- Autoclavable light guides in various sizes
- Voltage stabilizer assures constant light output
- High efficiency fan for continuous cool operation
- Improved optical system provides high power from a halogen source
- 11mm light guide included.



Curing Light Halogen 085-1121

\$537.70

DIRECT BOND METAL BRACKETS (CONT.)



Ceramic System

With Hooks on Cuspids and Bicuspids

- A Roth prescription
- Translucent ceramic material is cosmetic yet very strong
- Mechanical bonding base for use with standard adhesives without the need for silane

Five Patient Kit contains: 100 brackets total

Upper and lower, 2nd bicuspid to 2nd bicuspid, 5 sets

.018 Five Patient Kit094-0700

1@ \$1,203.50

.022 Five Patient Kit094-0718

Individual Brackets

Maxillary Arch, 5/Pkg

1@ \$61.10
10@ \$60.05

Tooth	Torque	Angulation	Rotation	.018	.022
Central	UR	+12	+5	0	094-0759 094-0908
	UL				094-0767 094-0916
Lateral	UR	+8	+9	0	094-0775 094-0924
	UL				094-0783 094-0932
Cuspid	UR	0	+11	4	094-0791 094-0940
	UL				094-0809 094-0957
1 st Bicuspid	UR	-7	0	2	094-0817 094-0965
	UL				094-0825 094-0973
2 nd Bicuspid	UR	-7	0	2	094-0817 094-0965
	UL				094-0825 094-0973

Mandibular Arch, 5/Pkg

Tooth	Torque	Angulation	Rotation	.018	.022
Central	LR/LL	0	0	094-0833 094-0981	
Lateral	LR/LL	0	0	094-0833 094-0981	
Cuspid	LR	-11	+5	2	094-0841 094-0999
	LL				094-0858 094-1005
1 st Bicuspid	LR	-17	0	0	094-0866 094-1013
	LL				094-0874 094-1021
2 nd Bicuspid	LR	-22	0	0	094-0882 094-1039
	LL				094-0890 094-1047

Reinforced Plastic, Metal Slot Plastic System

With Hooks on Cuspids and Bicuspids

- A Roth prescription
- Translucent, fortified urethane material for strength and aesthetics
- Slot reinforced with metal for low friction and accurate control of torque
- Mechanical bonding base for use with standard adhesives without the need for plastic bracket primer

Five Patient Kit contains: 100 brackets total

Upper and lower, 2nd bicuspid to 2nd bicuspid, 5 sets

.018 Five Patient Kit094-1195

1@ \$1,203.50

.022 Five Patient Kit094-1203

Individual Brackets

Maxillary Arch, 5/Pkg

1@ \$61.10
10@ \$60.05


Tooth	Torque	Angulation	Rotation	.018	.022
Central	UR	+12	+5	0	094-1260 094-1419
	UL				094-1278 094-1427
Lateral	UR	+8	+9	0	094-1286 094-1435
	UL				094-1294 094-1443
Cuspid	UR	0	+11	4	094-1302 094-1450
	UL				094-1310 094-1468
1 st Bicuspid	UR	-7	0	2	094-1328 094-1476
	UL				094-1336 094-1484
2 nd Bicuspid	UR	-7	0	2	094-1328 094-1476
	UL				094-1336 094-1484

Mandibular Arch, 5/Pkg

Tooth	Torque	Angulation	Rotation	.018	.022
Central	LR/LL	0	0	094-1344 094-1492	
Lateral	LR/LL	0	0	094-1344 094-1492	
Cuspid	LR	-11	+5	2	094-1351 094-1500
	LL				094-1369 094-1518
1 st Bicuspid	LR	-17	0	0	094-1377 094-1526
	LL				094-1385 094-1534
2 nd Bicuspid	LR	-22	0	0	094-1393 094-1542
	LL				094-1401 094-1559

Bracket Identification

Brackets are colour-coded by quadrant.



Patterson's direct bond cosmetic brackets are colour-coded for easy identification. Green designates patient's upper right, Blue upper left, Yellow lower right, and Red lower left. Universal right/left is designated with Black colour code.

PLAIN MOLAR BANDS

Patterson Brand molar bands are available in either plain or prewelded styles with a buccal tube and lingual cleat.

Patterson's first molar bands feature:

- Right and left design, with anatomically precise contouring for exceptional fit
- 24 sizes
- Medium temper for ease of placement with no distortion
- Rolled occlusal edge for cement retention and strength during placement
- Reduced interproximal thickness and height for ease of banding, less separation between teeth, and less tissue impingement
- Permanent laser size marking

Plain bands are available in kits of 100 or individual sizes in bags of three per size.



Plain Band Kit, Upper, 100pc (1/2 right / 1/2 left)094-5006 1@ \$242.15

Size	d	c	b	a	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Qty	2	2	2	2	2	4	4	4	4	4	4	4	4	6	6	6	4	4	4	4	4	4	4	4	2	2	2	2

Plain Band Kit, Lower, 100pc (1/2 right / 1/2 left)094-5014 1@ \$242.15

Size	d	c	b	a	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Qty	2	2	2	2	2	4	4	4	4	4	4	4	4	6	6	6	4	4	4	4	4	4	4	4	2	2	2	2

Plain Band Individual Sizes

Size	Upper Right	Upper Left
1	095-1517	094-5030
2	095-1525	094-5048
3	095-1533	094-5055
4	095-1541	094-5063
5	095-1558	094-5071
6	095-1566	094-5089
7	095-1574	094-5097
8	095-1582	094-5105
9	095-1590	094-5113
10	095-1608	094-5121
11	095-1616	094-5139
12	095-1624	094-5147
13	095-1632	094-5154
14	095-1640	094-5162
15	095-1657	094-5170
16	095-1665	094-5188
17	095-1673	094-5196
18	095-1681	094-5204
19	095-1699	094-5212
20	095-1707	094-5220
21	095-1715	094-5238
22	095-1723	094-5246
23	095-1731	094-5253
24	095-1749	094-5261
Marking	Size+ 6UR	Size+ 6UL

Size	Lower Right	Lower Left
1	095-1756	094-5279
2	095-1764	094-5287
3	095-1772	094-5295
4	095-1780	094-5303
5	095-1798	094-5311
6	095-1806	094-5329
7	095-1814	094-5337
8	095-1822	094-5345
9	095-1830	094-5352
10	095-1848	094-5360
11	095-1855	094-5378
12	095-1863	094-5386
13	095-1871	094-5394
14	095-1889	094-5402
15	095-1897	094-5410
16	095-1905	094-5428
17	095-1913	094-5436
18	095-1921	094-5444
19	095-1939	094-5451
20	095-1947	094-5469
21	095-1954	094-5477
22	095-1962	094-5485
23	095-1970	094-5493
24	095-1988	094-5501
Marking	Size+ 6LR	Size+ 6LL

Sizes smaller than 1; d, c, b and a, with d being the smallest.

Size	Upper Right	Upper Left
d	095-2028	094-5543
c	095-2010	094-5535
b	095-2002	094-5527
a	095-1996	094-5519
Size	Lower Right	Lower Left
d	095-2069	094-5584
c	095-2051	094-5576
b	095-2044	094-5568
a	095-2036	094-5550
3/Bag:	1@ \$7.75	3@ \$7.25

WELDABLE TUBES

- Metal Injection Molded (MIM) for smooth, accurate, reduced size tubes
- Convertible tubes allow easier insertion of wires when 2nd molar tubes are also present
- Mesial hook cast as one piece with tube for strength and reduction of bulk

1st Molar Convertible Tubes

3/pkg 1@ \$23.40

Maxillary

Cast Convertible Triple

	Torque	Rotation	.018	.022
Right	-10	8	094-2094	094-2110
Left			094-2102	094-2128



Cast Convertible Double

	Torque	Rotation	.018	.022
Right	-10	8	094-2052	094-2078
Left			094-2060	094-2086



Cast Convertible Single

	Torque	Rotation	.018	.022
Right	-10	6	094-2011	094-2037
Left			094-2029	094-2045



Mandibular

Cast Convertible Double

	Torque	Rotation	.018	.022
Right	-25	0	094-2177	094-2193
Left			094-2185	094-2201




Cast Convertible Single

	Torque	Rotation	.018	.022
Right	-25	0	094-2136	094-2151
Left			094-2144	094-2169



Weldable Linguals

- Cleat or button for intraoral elastics or chain
- Sheath for use with palatal bars
- 10/pkg



Molar Cleat	094-2813	1@ \$8.50
Bicuspid Cleat.....	094-2805	1@ \$8.50
Button-flat.....	094-2821	1@ \$17.60
Button-curved.....	094-2839	1@ \$17.60
Sheath.....	094-2847	1@ \$21.95

1st Molar Non-Convertible Tubes

3/pkg 1@ \$19.95

Maxillary

Cast Non-Convertible Triple

	Torque	Rotation	.018	.022
Right	-10	0	094-2524	094-2557
Left			094-2540	094-2565



Cast Non-Convertible Double

	Torque	Rotation	.018	.022
Right	-10	0	094-2482	094-2508
Left			094-2490	094-2516



Cast Non-Convertible Single

	Torque	Rotation	.018	.022
Right	-10	0	094-2441	094-2466
Left			094-2458	094-2474



Mandibular

Cast Non-Convertible Double

	Torque	Rotation	.018	.022
Right	-25	0	094-2623	094-2649
Left			094-2631	094-2656



Cast Non-Convertible Single

	Torque	Rotation	.018	.022
Right	-25	0	094-2573	094-2599
Left			094-2581	094-2607



2nd Molar Non-Convertible Tubes

3/pkg 1@ \$12.25

Maxillary

	Torque	Rotation	.018	.022
Right	-10	6	094-2367	094-2383
Left			094-2375	094-2391



Mandibular

	Torque	Rotation	.018	.022
Right	-25	0	094-2409	094-2425
Left			094-2417	094-2433



DIRECT BOND TUBES AND PREWELDED MOLAR BANDS

1st Molar Convertible Tubes

5/pkg 1@ \$54.10
10@ \$53.20

Maxillary

Cast Convertible

	Torque	Rotation	.018	.022
Right	-10	8	094-2961	094-2987
Left			094-2979	094-2995



Mandibular

Cast Convertible

	Torque	Rotation	.018	.022
Right	-25	0	094-3001	094-3027
Left			094-3019	094-3035



Direct Bond Linguals

- Cleat or button for intraoral elastics or chain
- 10/pkg

Cleat.....	094-3456	1@ \$16.95 10@ \$16.70
Button-flat.....	094-3464	1@ \$30.40 10@ \$29.90
Button-curved.....	094-3472	1@ \$30.25 10@ \$29.65



1st Molar Non-Convertible Tubes

5/pkg 1@ \$48.50
10@ \$47.70

Maxillary

Cast Non-Convertible

	Torque	Rotation	.018	.022
Right	-10	0	094-3183	094-3209
Left			094-3191	094-3217



Mandibular

Cast Non-Convertible

	Torque	Rotation	.018	.022
Right	-25	0	094-3266	094-3282
Left			094-3274	094-3290



2nd Molar Non-Convertible Tubes

5/pkg 1@ \$48.50
10@ \$47.70

Maxillary

Cast Non-Convertible

	Torque	Rotation	.018	.022
Right	-10	0	094-3225	094-3241
Left			094-3233	094-3258



Mandibular

Cast Non-Convertible

	Torque	Rotation	.018	.022
Right	-25	0	094-3308	094-3324
Left			094-3316	094-3332



Patterson Brand Kit Box

- Features:
- Interlocking for convenient storage
 - Pull-out drawer for easy access
 - Plastic lid on drawers for improved hygiene
 - Can be used with included tray or with Patterson Band Kits
- First Molar, Upper Right094-5816 1@ \$14.65
First Molar, Upper Left094-5824
First Molar, Lower Right094-5832
First Molar, Lower Left094-5840



Patterson Brand Prewelded 1st Molar Bands

- Cast buccal tubes with ball hook
- Buccal tubes are convertible, allows easy archwire insertion when 2nd molars are present
- Lingual cleat for additional treatment options
- Prewelded bands save time
- Accurate tube placement provides more consistent treatment results

Prewelded bands are available in kits of 100 or individual sizes in bags of three per size.

Prewelded Band Kits, 100pc (1/2 right/1/2 left) 1@ \$1,086.05

Tooth	Tube	Torque	.018	.022
Upper 1 st Molar	Single	-10	094-6012	094-6020
Upper 1 st Molar	Double	-10	094-6038	094-6046
Upper 1 st Molar	Triple	-10	094-6053	094-6061
Lower 1 st Molar	Single	-25	094-6079	094-6087
Lower 1 st Molar	Double	-25	094-6095	094-6103

Kit Distribution

Size	1	2	3	4	5	6	7	8	9	10	11	12
Qty	2	4	4	4	4	4	4	4	4	6	6	6

Size	13	14	15	16	17	18	19	20	21	22	23	24
Qty	6	4	4	4	4	4	4	4	4	4	4	2

Orthodontics

PREWELDED MOLAR BANDS

.018 Prewelded Band Individual Sizes

Upper Single -10° Torque

Size	Right	Left
1	094-6319	094-6558
2	094-6327	094-6574
3	094-6335	094-6582
4	094-6343	094-6590
5	094-6350	094-6608
6	094-6368	094-6616
7	094-6376	094-6624
8	094-6384	094-6632
9	094-6392	094-6640
10	094-6400	094-6657
11	094-6418	094-6665
12	094-6426	094-6673
13	094-6434	094-6681
14	094-6442	094-6699
15	094-6459	094-6707
16	094-6467	094-6715
17	094-6475	094-6723
18	094-6483	094-6731
19	094-6491	094-6749
20	094-6509	094-6756
21	094-6517	094-6764
22	094-6525	094-6772
23	094-6533	094-6780
24	094-6541	094-6798



3/Bag: 1@ \$46.20
10@ \$44.80

Upper Double -10° Torque

Size	Right	Left
1	094-7291	094-7531
2	094-7309	094-7549
3	094-7317	094-7556
4	094-7325	094-7564
5	094-7333	094-7572
6	094-7341	094-7580
7	094-7358	094-7598
8	094-7366	094-7606
9	094-7374	094-7614
10	094-7382	094-7622
11	094-7390	094-7630
12	094-7408	094-7648
13	094-7416	094-7655
14	094-7424	094-7663
15	094-7432	094-7671
16	094-7440	094-7689
17	094-7457	094-7697
18	094-7465	094-7705
19	094-7473	094-7713
20	094-7481	094-7721
21	094-7499	094-7739
22	094-7507	094-7747
23	094-7515	094-7754
24	094-7523	094-7762

Upper Triple -10° Torque

Size	Right	Left
1	094-8257	094-8497
2	094-8265	094-8505
3	094-8273	094-8513
4	094-8281	094-8521
5	094-8299	094-8539
6	094-8307	094-8547
7	094-8315	094-8554
8	094-8323	094-8562
9	094-8331	094-8570
10	094-8349	094-8588
11	094-8356	094-8596
12	094-8364	094-8604
13	094-8372	094-8612
14	094-8380	094-8620
15	094-8398	094-8638
16	094-8406	094-8646
17	094-8414	094-8653
18	094-8422	094-8661
19	094-8430	094-8679
20	094-8448	094-8687
21	094-8455	094-8695
22	094-8463	094-8703
23	094-8471	094-8711
24	094-8489	094-8729

Lower Single -25° Torque

Size	Right	Left
1	094-9214	094-9453
2	094-9222	094-9461
3	094-9230	094-9479
4	094-9248	094-9487
5	094-9255	094-9495
6	094-9263	094-9503
7	094-9271	094-9511
8	094-9289	094-9529
9	094-9297	094-9537
10	094-9305	094-9545
11	094-9313	094-9552
12	094-9321	094-9560
13	094-9339	094-9578
14	094-9347	094-9586
15	094-9354	094-9594
16	094-9362	094-9602
17	094-9370	094-9610
18	094-9388	094-9628
19	094-9396	094-9636
20	094-9404	094-9644
21	094-9412	094-9651
22	094-9420	094-9669
23	094-9438	094-9677
24	094-9446	094-9685

Lower Double -25° Torque

Size	Right	Left
1	095-0170	095-0410
2	095-0188	095-0428
3	095-0196	095-0436
4	095-0204	095-0444
5	095-0212	095-0451
6	095-0220	095-0469
7	095-0238	095-0477
8	095-0246	095-0485
9	095-0253	095-0493
10	095-0261	095-0501
11	095-0279	095-0519
12	095-0287	095-0527
13	095-0295	095-0535
14	095-0303	095-0543
15	095-0311	095-0550
16	095-0329	095-0568
17	095-0337	095-0576
18	095-0345	095-0584
19	095-0352	095-0592
20	095-0360	095-0600
21	095-0378	095-0618
22	095-0386	095-0626
23	095-0394	095-0634
24	095-0402	095-0642

.022 Prewelded Band Individual Sizes

Upper Single -10° Torque

Size	Right	Left
1	094-6806	094-7044
2	094-6814	094-7069
3	094-6822	094-7077
4	094-6830	094-7085
5	094-6848	094-7093
6	094-6855	094-7101
7	094-6863	094-7119
8	094-6871	094-7127
9	094-6889	094-7135
10	094-6897	094-7143
11	094-6905	094-7150
12	094-6913	094-7168
13	094-6921	094-7176
14	094-6939	094-7184
15	094-6947	094-7192
16	094-6954	094-7200
17	094-6962	094-7218
18	094-6970	094-7226
19	094-6988	094-7234
20	094-6996	094-7242
21	094-7002	094-7259
22	094-7010	094-7267
23	094-7028	094-7275
24	094-7036	094-7283

Upper Double -10° Torque

Size	Right	Left
1	094-7770	094-8018
2	094-7788	094-8026
3	094-7796	094-8034
4	094-7804	094-8042
5	094-7812	094-8059
6	094-7820	094-8067
7	094-7838	094-8075
8	094-7846	094-8083
9	094-7853	094-8091
10	094-7861	094-8109
11	094-7879	094-8117
12	094-7887	094-8125
13	094-7895	094-8133
14	094-7903	094-8141
15	094-7911	094-8158
16	094-7929	094-8166
17	094-7937	094-8174
18	094-7945	094-8182
19	094-7952	094-8190
20	094-7960	094-8208
21	094-7978	094-8216
22	094-7986	094-8224
23	094-7994	094-8232
24	094-8000	094-8240



3/Bag: 1@ \$46.20
10@ \$44.80

Upper Triple -10° Torque

Size	Right	Left
1	094-8737	094-8976
2	094-8745	094-8984
3	094-8752	094-8992
4	094-8760	094-9008
5	094-8778	094-9016
6	094-8786	094-9024
7	094-8794	094-9032
8	094-8802	094-9040
9	094-8810	094-9057
10	094-8828	094-9065
11	094-8836	094-9073
12	094-8844	094-9081
13	094-8851	094-9099
14	094-8869	094-9107
15	094-8877	094-9115
16	094-8885	094-9123
17	094-8893	094-9131
18	094-8901	094-9149
19	094-8919	094-9156
20	094-8927	094-9164
21	094-8935	094-9172
22	094-8943	094-9180
23	094-8950	094-9198
24	094-8968	094-9206

Lower Single -25° Torque

Size	Right	Left
1	094-9693	094-9933
2	094-9701	094-9941
3	094-9719	094-9958
4	094-9727	094-9966
5	094-9735	094-9974
6	094-9743	094-9982
7	094-9750	094-9990
8	094-9768	095-0006
9	094-9776	095-0014
10	094-9784	095-0022
11	094-9792	095-0030
12	094-9800	095-0048
13	094-9818	095-0055
14	094-9826	095-0063
15	094-9834	095-0071
16	094-9842	095-0089
17	094-9859	095-0097
18	094-9867	095-0105
19	094-9875	095-0113
20	094-9883	095-0121
21	094-9891	095-0139
22	094-9909	095-0147
23	094-9917	095-0154
24	094-9925	095-0162

Lower Double -25° Torque

Size	Right	Left
1	095-0659	095-0899
2	095-0667	095-0907
3	095-0675	095-0915
4	095-0683	095-0923
5	095-0691	095-0931
6	095-0709	095-0949
7	095-0717	095-0956
8	095-0725	095-0964
9	095-0733	095-0972
10	095-0741	095-0980
11	095-0758	095-0998
12	095-0766	095-1004
13	095-0774	095-1012
14	095-0782	095-1020
15	095-0790	095-1038
16	095-0808	095-1046
17	095-0816	095-1053
18	095-0824	095-1061
19	095-0832	095-1079
20	095-0840	095-1087
21	095-0857	095-1095
22	095-0865	095-1103
23	095-0873	095-1111
24	095-0881	095-1129

ARCHWIRE

Exact Arch Form

Ideal for today's pre-adjusted bracket systems. Available in a broad range of materials for continuity of arch shape throughout all treatment phases. Use Super Elastic Nickel Titanium for leveling and early treatment, Beta Titanium for medium force during mid treatment and Stainless Steel for superior control during final detailing. Permanent midline marks for accurate alignment.



Super Elastic Nickel Titanium

- Exact arch form
- Leveling/Initial treatment stages
- Low, continuous force
- 10 pieces/pkg

Round Upper		Round Lower	
.012.....095-4016	1@ \$19.30	.012.....095-4024	1@ \$19.30
.014.....095-4032	10@ \$18.50	.014.....095-4040	10@ \$18.50
.016.....095-4057		.016.....095-4065	
.018.....095-4073		.018.....095-4081	
.020.....095-4099		.020.....095-4107	

Rectangular Upper		Rectangular Lower	
.016 x .016...095-4214	1@ \$30.15	.016 x .016...095-4222	1@ \$30.15
.016 x .022...095-4230	10@ \$29.60	.016 x .022...095-4248	10@ \$29.60
.017 x .025...095-4255		.017 x .025...095-4263	
.018 x .018...095-4271		.018 x .018...095-4289	
.018 x .025...095-4297		.018 x .025...095-4305	
.019 x .025...095-4313		.019 x .025...095-4321	

Thermal Sensitive Nickel Titanium

- Exact arch form
- Soft and flexible at room temperature for easier insertion
- Transforms at 35°C to provide low, continuous force
- 10 pieces/pkg

Round Upper		Round Lower	
.014.....095-4438	1@ \$30.80	.014.....095-4446	1@ \$30.80
.016.....095-4453	10@ \$30.30	.016.....095-4461	10@ \$30.30
.018.....095-4479		.018.....095-4487	

Rectangular Upper		Rectangular Lower	
.016 x .016...095-4560	1@ \$40.60	.016 x .016...095-4578	1@ \$40.60
.016 x .022...095-4586	10@ \$40.25	.016 x .022...095-4594	10@ \$40.25
.017 x .025...095-4602		.017 x .025...095-4610	
.019 x .025...095-4628		.019 x .025...095-4636	

Beta Titanium Wire

- Exact arch form
- Medium force
- Extra length allows loops and bends
- 10 pieces/pkg

Round Upper		Round Lower	
.016.....095-4719	1@ \$84.90	.016.....095-4727	1@ \$84.90
.018.....095-4735	10@ \$84.55	.018.....095-4743	10@ \$84.55

Rectangular Upper		Rectangular Lower	
.016 x .022...095-4818	1@ \$88.90	.016 x .022...095-4826	1@ \$88.90
.017 x .025...095-4834	10@ \$86.95	.017 x .025...095-4842	10@ \$86.95
.019 x .025...095-4859		.019 x .025...095-4867	

Six Stranded Wire

- Exact arch form
- Stainless steel
- Very flexible, gentle force
- Will not unravel when cut
- 10 pieces/pkg



Upper/Lower		
.015.....095-5567		1@ \$25.25
.0175.....095-5575		10@ \$24.80

Three Stranded Wire

- Exact arch form
- Stainless steel
- More flexible than solid wire
- Will not unravel when cut
- 10 pieces/pkg



Upper/Lower		
.015.....095-5484		1@ \$18.50
.0175.....095-5492		10@ \$17.90

Reverse Curve of Spee

- Super Elastic Nickel Titanium material
- Open/close a bite
- Make corrections to curve of Spee
- Close space after extraction
- Stabilize teeth while using plastic chain
- 10 pieces/pkg

Round Upper		Round Lower	
.016.....095-5260	1@ \$57.15	.016.....095-5278	1@ \$57.15
.018.....095-5286	10@ \$56.10	.018.....095-5294	10@ \$56.10

Rectangular Upper		Rectangular Lower	
.016 x .022...095-5377	1@ \$58.80	.016 x .022...095-5385	1@ \$58.80
.017 x .025...095-5393	10@ \$57.85	.017 x .025...095-5401	10@ \$57.85
.019 x .025...095-5419		.019 x .025...095-5427	

Stainless Steel Wire

- Exact arch form
- High force for complete control
- 10 pieces/pkg

Round Upper		Round Lower	
.014.....095-4933	1@ \$4.85	.014.....095-4941	1@ \$4.85
.016.....095-4958	10@ \$4.55	.016.....095-4966	10@ \$4.55
.018.....095-4974		.018.....095-4982	
.020.....095-4990		.020.....095-5005	

Rectangular Upper		Rectangular Lower	
.016 x .016...095-5062	1@ \$8.00	.016 x .016...095-5070	1@ \$8.00
.016 x .022...095-5088	10@ \$7.85	.016 x .022...095-5096	10@ \$7.85
.017 x .025...095-5104		.017 x .025...095-5112	
.018 x .025...095-5120		.018 x .025...095-5138	
.019 x .025...095-5146		.019 x .025...095-5153	
.021 x .025...095-5161		.021 x .025...095-5179	

Straight Wire

- 14" straight lengths
- Most popular sizes for forming arches, sectional wires and auxiliaries
- 10 pieces/pkg

Round		
.016.....096-0419		1@ \$5.60
.018.....096-0427		10@ \$5.30
.028.....096-0435		
.030.....096-0443		
.032.....096-0450		
.036.....096-0468		
.040.....096-0476		
.045.....096-0484		
.051.....096-0492		

WIRE PRODUCTS, STOPS AND HOOKS

Retainer Wire


- 14" length
- Slightly softer than regular stainless steel
- 20 pieces/pkg

.025	095-5633	1@ \$10.75
.030	095-5641	10@ \$10.50
.032	095-5658	
.036	095-5666	

Ball Retainer Clasps

- Increases retention of removable appliances
- Ball tip is twice the diameter of the wire
- 25 pieces/pkg

.028	095-5716	1@ \$7.75
.032	095-5724	10@ \$7.60



Open Coil Springs

- Nickel Titanium
- Light, continuous force
- Used for opening spaces
- 15" spool

.030 ID	095-5757	1@ \$35.40
.036 ID	095-5765	10@ \$35.10



Closed Coil Springs

- Nickel Titanium
- Light, continuous force
- Used for closing spaces
- Attach to brackets, tubes or posted wires
- 10 pieces/pkg

9mm, .030 ID.....	095-5773	1@ \$92.80
12mm, .030 ID.....	095-5799	10@ \$91.35



Preformed Ligature Ties

- Dead soft for easy tying
- Long ties for lacing and for use with Steiner or Coon style ligature tying instruments
- 1,000 pieces/pkg

.009	095-5823	1@ \$27.35
.010	095-5831	10@ \$26.85



Pre-Twisted Ligature Ties

- Convenient and efficient
- For use with haemostats
- 1,000 pieces/pkg

.010	095-5989	1@ \$53.50
		10@ \$52.65



Kobayashi Hook

- Adds a hook for intraoral elastics
- 100 pieces/pkg

.012	095-6003	1@ \$24.15
.014	095-6011	10@ \$23.70



Pre-Twisted Kobayashi Hook

- Adds a hook in just seconds
- 100 pieces/pkg

.012	095-6052	1@ \$35.50
.014	095-6060	10@ \$34.80



Separating Springs

- For separating teeth prior to banding
- Stainless steel with comfortable, efficient low force
- 25 pieces/pkg

Medium.....	095-6102	1@ \$41.75
Long.....	095-6110	10@ \$41.00



Split Rectangular Stop

- Easy to attach without removing arch wire
- 100 pieces/pkg

.018 x .025	095-6409	1@ \$53.50
.021 x .028	095-6417	10@ \$52.65

- Allows easy adjustments
- Hook adds mechanical options

10 Locks without Hooks.....	095-6458	1@ \$94.25
		10@ \$92.80



10 Locks with Hooks	095-6466	1@ \$115.50
		10@ \$113.55

One Hex Handle	095-6474	1@ \$8.55
		10@ \$8.40

Surgical Hooks

- Secure attachment point for elastomers
- Post surgical intermaxillary fixation
- 10 pieces/pkg

Slide-on hook

.018 x .025, 2mm Cast with Ball Hook.....	095-6516	1@ \$33.20
.021 x .025, 2mm Cast with Ball Hook.....	095-6524	10@ \$32.65



Split hook

.018 x .025, 2mm Cast with Ball Hook.....	095-6532	1@ \$33.20
.021 x .025, 2mm Cast with Ball Hook.....	095-6540	10@ \$32.65

Bondable Retainers for Lower 3 to 3

- Made from .032 wire
- Mesh pads for strong retention
- 5 pieces/pkg

22mm.....	095-6615	1@ \$22.90
24mm.....	095-6623	10@ \$22.50
26mm.....	095-6631	
28mm.....	095-6649	
30mm.....	095-6656	



Direct Bond Space Maintainers

- Made from .032 wire
- Mesh pads for strong retention
- Sizes measured from center of pad to center of pad
- Used to stabilize mixed or adult dentition
- Five pieces/pkg

17mm Right	096-0518
17mm Left.....	096-0559
19mm Right	096-0526
19mm Left.....	096-0567
21mm Right	096-0534
21mm Left.....	096-0575
23mm Right	096-0542
23mm Left.....	096-0583



	1@ \$22.90
	10@ \$22.50

FACE BOWS AND EXPANSION SCREWS

High Pull Safety Release Straps

- For use with safety high pull head caps
- Adjustable force
- One pair/pkg

Navy Blue095-6904 1@ \$16.15
 10@ \$15.65



Safety High Pull Head Caps

- Flat seams for patient comfort
- One piece/pkg

Navy Blue095-6912 1@ \$15.25
 10@ \$14.95

Cervical Safety Release Straps

- For use with safety cervical neck pad
- Adjustable force
- Medium, 20oz
- Five pairs/pkg

Navy Blue095-6920 1@ \$64.70
 10@ \$63.65

Cervical Safety Neck Pads

- Flat seams for patient comfort
- Five pieces/pkg

Navy Blue095-6938 1@ \$31.10
 10@ \$30.60



Face Bows

- Standard bow allows full adaptation to each patient
- Short outer bow for use with high pull system
- Five sizes of looped face bows for convenience
- Measure along upper arch from mesial to mesial of headgear tubes

Standard Face Bow - .045 inner
 Standard Bow095-7019


1@ \$9.50
 10@ \$9.30

Short Outer Bow - .045 inner
 Each095-7027

1@ \$9.50
 10@ \$9.30

Loop Style Face Bow - .045 inner

Size 2, 90mm.....095-7035 1@ \$9.50
 Size 3, 97mm.....095-7043 10@ \$9.30
 Size 4, 104mm.....095-7050
 Size 5, 111mm.....095-7068
 Size 6, 118mm.....095-7076




Warning: Safety device must be used to avoid injury.

Expansion Screws

- Standard for transverse expansion of the upper arch
- Medium for expansion in upper or lower arch
- Micro for narrow palates or for anterior expansion
- Five piece/pkg (screw keys included)


Standard

5mm Maximum Expansion095-7183 1@ \$47.80
 10@ \$47.00




Medium

5mm Maximum Expansion095-7191



Micro

4mm Maximum Expansion095-7209




7mm Maximum Expansion095-7217

Rapid Palatal Expansion Screw

- Weldable legs for soldering to bicuspid and molar bands
- Most widely used screw for palatal expansion
- One piece/pkg (screw key included)

8mm Expansion095-7316 1@ \$37.30
 9mm Expansion095-7324 10@ \$36.70
 11mm Expansion095-7332
 13mm Expansion095-7340



Expansion Screw Key

- Includes safety straps
- 25 pieces/pkg

Each095-7407 1@ \$113.30
 10@ \$111.35



Mesial Loop Palatal Bar

- Rotate, expand, anchor, torque, upright upper molars
- Avoids gag reflex
- Improved tongue activation for faster results
- 5 pieces/pkg

37mm095-6763 1@ \$11.30
 40mm095-6771 10@ \$10.75
 43mm095-6789
 46mm095-6797
 49mm095-6805
 52mm095-6813



Warning: Palatal bars must be securely tied into the lingual sheath with steel or plastic ligatures. Improper use may result in injury or death by aspiration or ingestion.

ELASTOMERICS

Adventure Series Elastics

Adventure Series reflects today's active lifestyles, popular with both teens and adults.



Latex Elastics

- Made from surgical grade latex tubing
- 100 elastics/bag, 30 bags/pkg

	2½oz Light	4½oz Medium	6½oz Heavy	
1/8"	095-7514 Motorcycling	095-7563 Scuba Diving	095-7613 Inline Skating	1@ \$27.40 10@ \$27.00
3/16"	095-7522 Skateboarding	095-7571 Rock Climbing	095-7621 Hang Gliding	
1/4"	095-7530 Bicycling	095-7589 Parachuting	095-7639 Auto Racing	
5/16"	095-7548 Kayaking	095-7597 Sailing	095-7647 Water Skiing	
3/8"	095-7555 Skiing	095-7605 Wind Surfing	095-7654 Surfing	

X-ray Opaque Plastic Separators

- Latex free material shows up on X-ray in case separator gets lodged below the gum line
- For use on posterior teeth prior to banding
- 1,000 separators/pkg

Blue.....095-8801	1@ \$18.35 10@ \$17.75
-------------------	---------------------------



Plastic Tie Tubing

- Latex free, strong material
- For retracting along the arch wire, consolidation of space, extrusion of impacted teeth

Clear	095-8900	1@ \$10.60 10@ \$10.30
-------------	----------	---------------------------



Arch Wire Sleeve

- Latex free material
- Slip over arch wire sections to isolate wires from sensitive cheek tissue
- .030 ID, .065 OD
- 10ft/spool

Grey	095-8827	1@ \$10.60
Tooth Colour.....	095-8835	10@ \$10.30



Ortho CleanNcare™

- Fits retainers, pacifiers, mouthgards, etc.
- Dishwasher safe
- Flexible
- Sanitizes
- Easy to use
- Convenient
- Could be colour coded for each family member
- FDA approved silicone

Ortho CleanNcare	093-7854	1@ 30.80\$ 3@ 29.30\$
------------------------	----------	--------------------------



Plastic Ligatures and Chain

- 12 eye-catching, patient pleasing colours
- Ligature and chain colours match
- Latex free medical grade polyurethane material

Plastic Ligatures

- Molded ligatures are strong, smooth and comfortable
- 20 ligatures/stick, 25 sticks/pkg



Black.....	095-7985	Light Pink....	095-8033	1@ \$11.25
Clear	095-8025	Purple.....	095-7936	10@ \$10.90
Green.....	095-7993	Red	095-7977	
Light Blue....	095-7951	Silver.....	095-7910	
Navy.....	095-7928	White.....	095-8017	
Orange	095-7969	Yellow	095-8009	
Pink.....	095-7944			

Plastic Chain

- Available in closed, medium and long filament gaps
- 15ft/spool



	Closed	Short	Long	
Black	095-8181	095-8389	095-8587	1@ \$23.20
Clear	095-8223	095-8421	095-8629	10@ \$22.50
Green	095-8199	095-8397	095-8595	
Light Blue	095-8157	095-8355	095-8553	
Navy	095-8124	095-8322	095-8520	
Orange	095-8165	095-8363	095-8561	
Pink	095-8140	095-8348	095-8546	
Light Pink	095-8231	095-8439	095-8637	
Purple	095-8132	095-8330	095-8538	
Red	095-8173	095-8371	095-8579	
Silver	095-8116	095-8314	095-8512	
White	095-8215	095-8413	095-8611	
Yellow	095-8207	095-8405	095-8603	

Rotation Wedges

- Clear, latex free material
- Attach to bracket tie wings for increased rotation control
- 100 wedges/pkg

Clear	095-8819	1@ \$31.45 10@ \$30.90
-------------	----------	---------------------------



Retainer Box

- Sturdy plastic
- Available in eight colours

12/Pkg		
Beige.....	084-3797	1@ \$11.25
Green.....	084-3805	6@ \$10.45
White.....	084-3839	
Blue.....	084-3847	
Yellow.....	084-3854	
Red.....	084-3862	
Pink.....	084-3821	



Band Comparison Chart

Patterson	Amer. Ortho	EOP	Sullivan Schein/ Masel	3/M Unitek	Rocky Mtn.	Ortho Org.	Ormco	TP Ortho	Dentaurum
1	1.5	31.5	31.5	31.5	9	30.5	7	29.5	5
2	2	32	32	31	10	31	8	30	6
3	2.5	32.5	32.5	32.5	11	31.5	9	30.5	7
4	3	33	33	33	12	32	10	31	8
5	3.5	33.5	33.5	33.5	13	32.5	11	31.5	9
6	4	34	34	34	14	33	12	32	10
7	4.5	34.5	34.5	34.5	15	33.5	13	32.5	11
8	5	35	35	35	16	34	14	33	12
9	5.5	35.5	35.5	35.5	17	34.5	15	33.5	13
10	6	36	36	36	18	35	16	34	14
11	6.5	36.5	36.5	36.5	19	35.5	17	34.5	15
12	7	37	37	37	20	36	18	35	16
13	7.5	37.5	37.5	37.5	21	36.5	19	35.5	17
14	8	38	38	38	22	37	20	36	18
15	8.5	38.5	38.5	38.5	23	37.5	21	36.5	19
16	9	39	39	39	24	38	22	37	20
17	9.5	39.5	39.5	39.5	25	38.5	23	37.5	21
18	10	40	40	40	26	39	24	38	22
19	10.5	40.5	40.5	40.5	27	39.5	24	38.5	23
20	11	41	41	41	28	40	26	39	24
21	11.5	41.5	41.5	41.5	29	40.5	27	39.5	25
22	12	42	42	42	30	41	28	40	26
23	12.5	42.5			31	41.5	29	40.5	27
24	13	43			32	42	30	41	28

Patterson Brand bands and the corresponding competitor bands are a close approximation. Due to manufacturing tolerances concerning band contour, stretching techniques and material, sizes may differ slightly.

HAND INSTRUMENTS – PATTERSON

Cutting Instruments



Hard Wire Cutter
096-0062
1@ \$140.30



Mini Pin and Ligature Cutter
096-0013
1@ \$140.30



Angled 15° Hard Wire Cutter
096-0070
1@ \$140.30



45° Ligature Cutter
096-0021
1@ \$140.30



Safety Shear and Hold Distal End Cutter
096-0039
1@ \$140.30



Safety Shear and Hold Distal End Cutter for Niti Wire
096-0047
1@ \$140.30

Bending/Utility Instruments



Nance Loop Pliers
095-9155
1@ \$129.05



Band Contouring Pliers
095-9163
1@ \$129.05



Bird Beak Pliers
095-9049
1@ \$129.05



Light Wire Pliers
095-9056
1@ \$129.05



Super Looper Light Wire Pliers
095-9064
1@ \$129.05



Three Jaw Straight Pliers
095-9080
1@ \$129.05



Tweed Style Rectangular Arch Forming Pliers
095-9122
1@ \$129.05



How Style Straight Utility Pliers
095-9015
1@ \$129.05



How Style Offset Utility Pliers
095-9023
1@ \$129.05



Weingart Pliers
095-9072
1@ \$129.05



Lingual Arch Forming Pliers
095-9114
1@ \$129.05



Tweed Style Loop Forming Pliers
095-9098
1@ \$129.05



Tweed Style Omega Loop Forming Pliers
095-9106
1@ \$129.05



V-Bending Pliers
095-9148
1@ \$129.05



Jarabak Pliers
095-9130
1@ \$129.05



Crimping Pliers for Surgical Hooks
095-9304
1@ \$140.30



Heavy Duty Wire Cutter
096-0088
1@ \$140.30



Convertible Cap Remover
095-9858
1@ \$41.05



Haemostat
095-9908
1@ \$48.60



Fine Tip Haemostat
095-9916
1@ \$48.60



Mathieu
095-9874
1@ \$72.95



Narrow Tip Mathieu
095-9882
1@ \$72.95



Perry Dressing Tweezers
095-9809
1@ \$43.30



Stainless Steel Mouth Mirror
096-0203
1@ \$19.65

HAND INSTRUMENTS (CONT.) – PATTERSON

Direct Bond and Band Removing Instruments



Direct Bond
Bracket Holder
095-9791
1@ \$37.35



Direct Bond Bracket
Holder, Curved
095-9775
1@ \$44.75



Direct Bond Bracket
Holder with
Alignment Tip
095-9817
1@ \$37.45



Direct Bond Bracket
Holder with Fine Tip
095-9825
1@ \$37.45



Positive Positioner
095-9759
1@ \$129.05



Band Seating Pliers
095-9031
1@ \$129.05



High Heat
Bite Stick, 3 pieces/pkg
095-9700
1@ \$72.90



Force Module
Separator
095-9890
1@ \$129.05



Boone Gauge
096-0229
1@ \$76.20



Long Beak Posterior
Band Remover
095-9536*
1@ \$138.25



Band Slitting Pliers
096-0054
1@ \$140.30



Adjustable Enamel
Stripper Kit
096-0245
1@ \$207.35



Direct Bond Bracket
Remover, Narrow
095-9502
1@ \$135.95



Direct Bond Bracket
Remover, Wide
095-9510
1@ \$138.25



Padded Direct Bond
Bracket Remover,
Narrow
095-9544**
1@ \$138.25



Padded Direct Bond
Bracket Remover
Wide
095-9551**
1@ \$138.25



Short Beak Posterior
Bracket Remover
095-9528*
1@ \$138.25

Adjustable Enamel
Stripper
Replacement Blades
Fine
(.005" single sided) -
096-0310
Medium
(.006" single sided) -
096-0328
Coarse
(.007" single sided) -
096-0336
1@ \$19.95

*Replacement Pads for Posterior Band Removers
3/16" Ultem
Replacement Pad095-9676 1@ \$3.20
1/4" Ultem
Replacement Pad095-9684

**Replacement Tips for Padded Direct Bond
Bracket Removers
Two Replacement
Tips and Wrench - .125"095-9601 1@ \$21.40
Two Replacement
Tips and Wrench - .170"095-9619



Pliers Rack
096-0237
1@ \$70.95

Instrument Racks

Archwire Holder 4"W x 11"L x 2"H.....	088-8628	1@ \$42.80
Archwire Holder with Instrument Holding Tube 6½"W x 4¼"L x 4¼"H.....	088-8636	1@ \$44.25
Straight Wire Holder 8"W x 3½"L x 4½"H.....	088-8644	1@ \$39.15
Upright Plier Rack 5"W x 9"L x 5½"H.....	088-8651	1@ \$23.90



Sample Instrument Set-up Trays

Bonding Appointments

Explorer/Scaler, Mirror, Direct Bond
Bracket Holder or Positive
Positioner, Boone Gauge

Banding Appointments
























Explorer/Scaler, Mirror, Band
Seating Pliers, Band Contouring
Pliers, Band Removing Pliers

Wire Change Appointments

Explorer/Scaler, Mirror, Weingart or
How Instrument, Ligature Cutter,
Distal End Cutter, Haemostat,
Ligature Director/Tucker

Debonding/Debanding Appointments

Explorer/Scaler, Mirror,
Appropriate Debond/Deband
Instrument and/or Band Slitting
Pliers, Weingart or How
Instrument, Ligature Cutter

					
Omega Loop Forming Pliers 375-8505	Pin & Ligature Cutter 375-8125	Micro-Mini Pin & Ligature Cutter 375-8158	Mini Pin & Ligature Cutter 375-8133	Micro-Mini Pin & Ligature Cutter Long Handle 375-8166	O'Brien Pliers 375-8588
1@ \$202.55 3@ \$194.80	1@ \$231.55 3@ \$221.90	1@ \$231.55 3@ \$221.90	1@ \$231.55 3@ \$221.90	1@ \$231.55 3@ \$221.90	1@ \$202.55 3@ \$194.80
					
Posterior Band Removing Pliers Long 375-8448	Posterior Band Removing Pliers Short 375-8455	Rectangular Arch Bending Pliers 375-8554	Three Jaw Pliers 375-8497	Mini Three Jaw Pliers 375-8596	Tweed Loop Forming Pliers 375-8539
1@ \$202.55 3@ \$194.80	1@ \$202.55 3@ \$194.80	1@ \$202.55 3@ \$194.80	1@ \$202.55 3@ \$194.80	1@ \$202.55 3@ \$194.80	1@ \$202.55 3@ \$194.80
					
Universal Cut & Hold Distal End Cutter 375-8000	Universal Cut & Hold Long Handle 375-8018	Weingart Pliers 375-8372	Weingart Pliers Long Handle 375-8380	Slim Weingart Pliers 375-8398	Adhesive Removing Pliers 375-8430
1@ \$242.35 3@ \$221.50	1@ \$231.55 3@ \$221.90	1@ \$211.95 3@ \$198.80	1@ \$202.55 3@ \$194.80	1@ \$202.55 3@ \$194.80	1@ \$202.55 3@ \$194.80
					
Arch Bending Pliers 375-8547	Bird Beak Pliers 375-8513	Bracket Removing 375-8463	Micro-Mini Ligature Cutter 15 Degree 375-8182	Crown and Band Contouring Pliers 375-8471	
1@ \$202.55 3@ \$194.80	1@ \$202.55 3@ \$194.80	1@ \$202.55 3@ \$194.80	1@ \$231.55 3@ \$221.90	1@ \$202.55 3@ \$194.80	

HAND INSTRUMENTS (CONT.)



**Coon Style
Ligature Tying
Pliers**

375-8331
1@ \$202.55



**Distal End Cutter
Flush Cut**

375-8026
1@ \$231.55
3@ \$221.90



**Force Module
Separating Pliers**

375-8307
1@ \$202.55
3@ \$194.80



**Hard Wire Cutter
Straight**

375-8117
1@ \$231.55
3@ \$221.90



**Hard Wire Cutter
15 Degree**

375-8109
1@ \$231.55
3@ \$221.90



**Hollow Chop
Contouring Pliers**

375-8570
1@ \$202.55
3@ \$194.80



How Pliers Offset

375-8414
1@ \$202.55
3@ \$194.80



How Pliers Straight

375-8406
1@ \$202.55
3@ \$194.80



Jarabak Pliers

375-8489
1@ \$202.55
3@ \$194.80



**Light Wire Bird Beak
Pliers with Groove**

375-8521
1@ \$202.55
3@ \$194.80



**Lingual Arch
Forming Pliers**

375-8562
1@ \$202.55
3@ \$194.80




Wire Cutting Scissors

374-6104
1@ \$170.95



Perma Sharp Mathieu Pliers

371-7436
1@ \$261.90



Twirl-on Ligature Applicator

371-1793
1@ \$48.95



Small Tip Band Pusher Scaler

367-3514
1@ \$71.60



Ligature Director

371-1819
1@ \$73.50



Single-ended Ligature Director

371-1991
1@ \$67.40



Double-ended Ligature Director

371-1801
1@ \$73.85



Band Pusher (Bp300)

367-3506
1@ \$71.55



Ligature Director Tucker

371-1975
1@ \$73.50



Mathieu Pliers

371-6800
1@ \$190.60



Band Pusher Scaler

367-3498
1@ \$42.50

Orthodontic Instruments - LM Instruments



Ligature Tucker

306-8566
1@ \$28.85 3@ \$26.25



Ligature Tucker - Scaler U15

306-8558
1@ \$28.85 3@ \$26.25



Ring Tie Placer

306-8905
1@ \$28.85 3@ \$26.25

DENTAL TEACHING AIDS



Fredd E. Flossisaurus
058-7477
021-5475, Braces



Al E. Gator
058-7444
021-5483, Braces



Brusher B. Boa
058-7485
021-5459, Braces



Molar Z. Monster
013-5459
021-5434, Braces



Mojo Monkey
013-5509
021-5442, Braces



Finn Z. Fish
013-5475
021-5509, Braces



Ollie Z. Mutt
058-7501
021-5491, Braces



Magi Z. Dragon
013-5467
021-5558, Braces



With Braces



Demonstration Mirror
021-5418



Orthodontic Consultation Model Kit

Kit is comprised of nine sets of orthodontic and dental malocclusions. Includes the following sets: Class I ideal, Class I anterior open bite, Class II division I blocked cuspids, Class II division I open bite narrow arches, Class II division I excessive overbite excessive overjet, Class II division II mixed dentition, Class II division II permanent dentition, Class II mixed dentition posterior crossbite, Class III mixed dentition. Each kit comes in a black, hard-shelled protective plastic case with padded interior. Models are made of high-impact plastic resins and denture teeth. Models approximately 2-7/8" W x 2" H x 2" D; case: 13-1/2" W x 4" H x 10" D.

Each.....055-5094 1@ \$938.95

Dental Puppet Teaching Aids

Brushing is fun! Animated puppets make it fun for children to learn about good brushing habits. Each puppet comes with an oversized, bristled toothbrush. Made of fire-retardant, poly-filled plush with thermoplastic resin mouth. Easy to clean and stain resistant – just spray and wipe clean. Puppets: 15"-17" H; Snake Puppet: 8" W x 40" L; Toothbrush: 11-1/2" H, Mirror 12". Please allow 3-4 weeks for delivery of Orthodontic Teaching Aids.

Puppet Teaching Aid		
Each	See item	1@ \$163.90
		3@ \$158.90
Additional Toothbrush		
Each	058-7519	1@ \$24.95
		3@ \$23.95
Orthodontic Puppet Teaching Aid		
Each	See item	1@ \$254.90
		3@ \$250.90
Demonstration Mirror		
Each	021-5418	1@ \$24.95
		3@ \$23.95

PATIENT EDUCATION BROCHURES



#050-6329 Orthodontic Care



#050-6600 Orthodontics

Educational Dental Brochures

Dental brochures are designed to address questions your patients may have on a variety of timely topics. Informative, colourful, and easy-to-read, these brochures can help clarify your recommendations and provide patients with the clear-cut information they need to make informed decisions about their care.

Perfect for use as handouts, billing and newsletter inserts, or for display in your office. Space to personalize with your practice information using a stamp or label. An inexpensive, economical way to educate your patients.


Dental brochures
50/pkg.
PkgSee item 1@ \$21.95
 2@ \$20.95
 4@ \$19.95




#055-9179

#055-9195

055-9203



#051-0172



#051-0206



#051-0149



#051-0156

Dental Education Booklets

- Space provided to stamp your practice information on back
- Larger size (7" W x 10" H folded)

25/pkg.
PkgSee item 1@ \$52.50
 2@ \$50.65
 4@ \$48.50

Kids' Mini Brochures

Colourful illustrations with clear and simple instructions that help teach children good oral hygiene habits at a young age. The small size is perfect to use as an office or school handout and as a statement insert. There's space available to stamp your practice information on the back of the brochure for personalization.

- 3-3/4" x 8-1/2"; one sided

Kid's Mini brochures
50/pkg.
PkgSee item 1@ \$16.45
 2@ \$15.45
 4@ \$14.45

Patient Education Brochure

This informative, illustrated brochure provides a fast and effective way to discuss treatment procedures, purposes and benefits with your patients. This great educational tool supports your recommendations and helps patients make treatment decisions.

- Cost-effective and easy-to-use
- 4" x 9", 6-panel brochure
- Space provided to stamp your practice information on back
- 50/pkg.

Pkg051-0388 1@ \$35.95
 2@ \$32.25
 4@ \$30.75



#051-0388

CEMENTS

GC Fuji ORTHO BAND Paste – GC America, Inc.



GC Fuji ORTHO BAND is a blue-coloured resin modified glass ionomer for bonding orthodontic bands. It can also be used for the cementation of appliances retained with stainless steel crowns (e.g., Herbst or expansion appliances). The new Paste Pak Dispenser saves time and reduces waste.

Paste Pak Starter Kit

Contains: Two Fuji ORTHO BAND Paste Pak Cartridges, one mixing pad, one Paste Pak Dispenser, and one mixing spatula.
Kit203-9436 1@ \$213.50

Paste Pak Refill

Contains: Two 13.3gm (7.2ml) Paste Pak Cartridges of Fuji ORTHO BAND.
Refill203-9444 1@ \$135.00
3@ \$129.95

Paste Pak Dispenser

Each203-8800 1@ \$81.35

GC Fuji Ortho Band LC for Automix – GC America, Inc.



GC Fuji ORTHO BAND is dual cured reinforced glass ionomer for bonding orthodontic bands. It can also be used for the cementation of appliances retained with stainless steel crowns (e.g., Herbst or expansion appliances).

Paste Pak Starter Set

Contains: 2 Paste pak Cartridge (13.3g ea.), 1 Paste Pak Dispenser, 20 Mixing Tips and 1 Mixing Pad
Set327-0816 1@ \$246.00

Paste Pak Refill

Contains: 2 Paste pak Cartridge (13.3g ea.), 1 Mixing Pad
Refill327-0824 1@ \$184.50
3@ \$178.50

Mixing Tips Refill

Contains: 20 Mixing Tips for LC
Pkg327-0832 1@ \$44.35

GC Fuji Ortho Self-Cure – GC America, Inc.



Self-Cured Resin Reinforced Orthodontic Cement

Excellent for bonding brackets, bands and acrylic appliances.

Self-Cure Powder/Liquid Delivery:

GC Fuji Ortho Self-Cure Standard Package
Contains: 1 - 40gm powder, 2 - 6.7ml liquid, scoop and mixing pad. Standard Pkg (439503)
Each203-9303 1@ \$317.70

GC Fuji Ortho Self-Cure Powder

Refill, 1 - 40gm (439521).....203-9329 1@ \$183.45

GC Fuji Ortho Self-Cure Liquid

Refill, 1 - 6.7ml (439591).....203-9311 1@ \$74.95



GC Fuji Ortho LC – America, Inc.

Light-Cured Resin Reinforced Orthodontic Cement

Fuji Ortho LC is a light-cured orthodontic adhesive that bonds metal, porcelain and polycarbonate brackets and metal bands to enamel. Fuji Ortho LC features low sensitivity to moisture and can be used in the presence of saliva. Available in both capsules and powder-liquid form.

Light-Cured Capsules Delivery

Fuji Ortho LC Refill
50 Capsules/Pkg
Box (439450)326-1732 1@ \$246.00

Fuji Ortho LC Capsule Starter Pack

Contains: 50 capsules, applicator and GC Ortho conditioner.
Applier Kit (439451) ..326-2243 1@ \$352.75

Capsule Applier III

Each203-9345 1@ \$92.25

Light-Cured Powder/Liquid Delivery

Fuji Ortho LC Standard Kit

Contains: 1 - 40gm bottle of powder, 2 - 6.8ml bottle of liquid, scoop and mixing pad.
Standard Pkg (439403).....326-1708 1@ \$311.75

Fuji Ortho LC Introductory Kit

Contains: 1 - 15gm bottle of powder, 1 - 6.8ml bottle of liquid, scoop and mixing pad.
Intro. Pkg (439401).....326-1690 1@ \$142.35

Fuji Ortho LC Powder Refill

40gm Refill (439421) ..326-1724 1@ \$179.85

Fuji Ortho LC Liquid Refill

6.8ml/8gm (439409) ..326-1716 1@ \$73.35

Fuji Ortho Conditioner

25gm (439425)208-0505 1@ \$50.10

Mixing Pad

3" x 3", #22
5 Pads326-3084 1@ \$55.10

GC Fuji Ortho LC (brackets) – GC America, Inc.



GC Fuji ORTHO BAND is dual cured reinforced glass ionomer orthodontic bracket cement.

It is recommended for bonding metal brackets and attachments, ceramic brackets and acrylic appliances (e.g. Rapid Palatal expander)

Paste Pak Automix Starter Set

Contains: 1 Paste pak Cartridge (13.5g), 1 Paste Pak Dispenser, 2 Ortho Gel Conditioner Syringes (1.4g ea.) 20 Mixing Tips and 15 Brush Tips
Set327-0774 1@ \$312.60

Paste Pak Automix Set

Contains: 1 Paste pak Cartridge (13.5g), 2 Ortho Gel Conditioner Syringes (1.4g ea.) 20 Mixing Tips and 15 Brush Tips
Refill327-0782 1@ \$286.95

Paste Pak Refill

Contains: 2 Paste pak Cartridge (13.5g ea.), Refill327-0790 1@ \$270.00

Fuji Ortho Gel Conditioner 20%

Contains: 3 Ortho Gel Conditioners (1.4g ea.) and 25 Brush Tips
Pkg.....208-0808 1@ \$9.65

Heliosit®-Orthodontic™ – Ivoclar Vivadent



Light Cure Composite

A light cured composite for the direct bonding of metal and plastic orthodontic brackets. Heliosit®-Orthodontic™ is a single component material in easy to use syringes. The material allows for virtually unlimited working time and is extremely fast curing.

Syringe Refill Package

Contains: 3 - 2.5gm syringes.
Pkg708-6721 1@ \$171.95

Ortho-Choice™ Band-Rite™ – Pulpdent



- Two-paste resin-ionomer band cement
- Dual cure – fluoride releasing

Band Kit contains 24gm base, 24gm catalyst, mixing pad, spatulas.
Kit589-5594 1@ \$125.35

Band-R Refill contains 12gm base, 12gm catalyst.
Refill589-5602 1@ \$62.00

Ortho-Choice™ Glass Ionomer Band Cement – Pulpdent



Fluoride releasing. Self-cure glass ionomer cement. Starter kit contains 30gm powder, 15ml liquid, mixing pad and scoop.

Kit589-8424 1@ \$68.75

Ortho Band Bulk Pack contains 240gm powder, 120ml liquid, scoop, mixing pad.

Kit589-8432 1@ \$304.90

Ortho-Choice™ OBA – Pulpdent



- Orthodontic bracket adhesive
- Light cure – fluoride releasing

OBA Syringe Kit

Contains 2 - 4gm syringes adhesive, 5ml bonding resin, 6ml etchant, plus accessories.

Kit589-8358 1@ \$78.60

Orthodontic Plaster Whip Mix



Exceptional hardness and particular whiteness. Developed especially for orthodontic work.

Super-White
20lb Carton728-7964 1@ \$35.30
3@ \$34.25

50lb Carton728-8269 1@ \$55.10
5@ \$52.35

Orthodontic Stone – Whip Mix



Hard, accurate, medium setting stone made with the orthodontist in mind.

Super-White
25lb Carton732-5046 1@ \$43.55
4@ \$42.25

50lb Carton732-4544 1@ \$61.80
5@ \$58.90

Orthodontic Resin – Dentsply/Caulk

Orthodontic Appliances

For fabrication of maintainers or retainers. Material never runs or slumps.

Regular Package

Contains: 605gm powder, 430ml liquid and liquid dispenser

Pink179-8750 1@ \$160.00
Clear179-8784

Powder Refill Only

Contains: 16oz powder

Pkg179-9824 1@ \$108.00
3@ \$106.00

Liquid Refill Only

12oz Bottle

Pink180-0341 1@ \$52.25
Clear180-0374



Orthodontic Resin – Hygenic

An auto-polymerizing resin used to construct orthodontic appliances and retainers. Cures in 6-8 minutes at room temperature.

Lab Pack (includes powder and liquid)

1lb

Pink379-4880 1@ \$82.90
Clear379-4898

PLASTERS / STONES

Orthodontic Plaster – Heraeus Kulzer



For orthodontia casts and study models.

Super-White

25lb516-0973 1@ \$39.35

50lb516-0981 1@ \$71.30
5@ \$70.05

Snow White Orthodontic Stone – Heraeus Kulzer



Specifically formulated for orthodontic casts where whiteness, accuracy, and durability are essential.

25lb514-9554 1@ \$46.90

ACRYLICS

NEOCRIL™ Splint Material – Bosworth



Bosworth's Neocril Splint material is an extra strong, high density acrylic resin suited for the fabrication of orthodontic splints. It will not chip, break or wear easily and is crystal clear in the patient's mouth.

Splint Kit

Contains: 1lb (454gm) powder, and 8oz (236ml) liquid.

Each138-5202 1@ \$88.30

Orthoresin – DETREY/Dentsply

Powder

3 Kg180-0408 1@ \$336.00

Liquid Clear

1 Liter180-0424 1@ \$109.00

Liquid Pink

1 Liter180-0432 1@ \$109.00

Ortho-Jet – Lang

A self-curing resin for fabrication of orthodontic appliances.



4oz1@ \$37.80

Red459-4180

Green459-4198

Blue459-4206

Orange459-4214

Yellow459-4222

Black459-4230

White459-4248

16oz1@ \$89.10

Snow White459-4081

White459-4099

240ml Liquid1@ \$37.80

Pink459-4420

Clear459-4412

Quart Liquid1@ \$81.50

Pink459-4321

Clear459-4313

TRAY WAXES

Orthodontic Tray Wax – Coltène/Hygenic

Sticks – 6" x 1/4"
Diameter

For peripheral lining of trays, increasing tray height, post-damming, etc. Soft and tacky, they may be adapted without heating. Pure white colour.

48/Box.....383-1450 1@ \$16.85



Orthodontic Tray Wax – KerrLab

Melting Temperature:
164.3°F/73°C

White

48/Box.....441-2276 1@ \$30.50

Economy Pack

480/Box.....441-2763 1@ \$227.00



Orthodontic Tray Wax – Heraeus Kulzer

Mint Flavoured

60/Box.....517-7506 1@ \$28.20



IMPRESSION MATERIAL

Alginmax – Bolton Dental

Chromatic Impression Material

Alginate designed for orthodontics, crowns and bridges and dentures, with these unique features:

- Complete impression in just 2 minutes!
- Dimensional Stability - 120+ hours!
- Pleasant vanilla scent with neutral taste.

1lb Pkg145-7365 1@ \$17.50



Orthoprint – Zhermack

- Highly elastic alginate
- Extra-rapid processing and setting times
- Thixotropic
- Vanilla flavour
- Yellow colour

500gm/Bag
Each

342-6277 1@ \$15.75
6@ \$15.40



MISCELLANEOUS

Ball Retainer Clasps – Keystone

Wire ball clasps are stainless steel to resist corrosion. Available in two sizes.

10/Pkg1@ \$6.35

Small .028"533-9015

Large .032"533-9023

100/Pkg1@ \$62.10

Small .028"533-9049 3@ \$59.75

Large .032"533-9056



Interra – Dentsply

The one visit, custom fit, value-packed Nightguard for dreamy results.

- Easy-to-use
- Less expensive than conventional nightguards
- Improved fit over laboratory nightguards
- Requires fewer adjustments
- Improved patient compliance
- Highly durable and wear resistant

Starter Package

1 Small Arch, 2 Medium Arches, 1 Large Arch, 1 Regisil Rigid 50ml Cartridge, 6 Mix Tips (yellow hub), Directions and Technique.

Pkg173-4201 1@ \$168.05



Niteguard Material – Dental Resources

Ortho repositioner/surgical applications.

12/Pkg.....279-5508 1@ \$97.60

Ortho-Coat™ – Pulpdent

Coats bracket and tooth. Keeps teeth clean under brackets.

Prevents decalcification, staining and discolouration under orthodontic brackets. Protects soft tissue from sharp edges of brackets, bands and wires.

Contains: 2 x 5ml syringes + applicator tips

Pkg.....590-1806 1@ \$65.95



Orthodontic Set LS-905 – Axis Dental

For fabricating and refining orthodontic appliances in three easy steps.

Set120-6952 1@ \$141.90



Orthodontic Travel – Oral-B

Compact convenience for traveling.

12/Pkg223-2189 1@ \$28.35
12@ \$26.95



Platypus™ Ortho Flosser – PDT

The easiest and most effective orthodontic flossing device ever?

- Increase patient compliance
- Better oral hygiene

Patient Sampler case (75 bags of 3 flossers)

Case561-1017 1@ \$113.75

Refills (Bag of 25 flossers)


Bag.....561-1009 1@ \$15.35



Intensiv Ortho-Strips


Thin, semiflexible strips for initial interproximal stripping and finishing in orthodontic procedures. Easily access areas that burs and discs cannot. Minimal heat generated, reducing patient discomfort and chance of damaging tooth structure or soft tissue. Use with the KaVo 61LRG LUX prep control head.

Double-Sided



MFG. NO	OS-15	OS-25	OS-40	OS-60	OS-90
Thickness (mm)	.10	.12	.13	.30	.35
Order No.	115-1257	115-1240	115-1232	115-1224	115-1216
Price	1@ \$174.10				
Assortment	(115-1265): 1 of each: \$174.10				
Package	4 Per Pack			4 Per Pack	

Single-Sided



MFG. NO	OS-15L-SS	OS-25L-SS	OS-40L-SS	OS-15R-SS	OS-25R-SS	OS-40R-SS
Thickness (mm)	.10	.11	.12	.10	.11	.12
Order No.	115-2768	115-2750	115-2743	115-2735	115-2727	115-2719
Type	Left Side			Right Side		
Price	1@ \$53.20			1@ \$53.20		
Assortment	(115-2693): 1 of each: \$316.05					
Package	Individual					

Ortho Triumph® with SmartGuide™ - Oral-B



The only brushing system with a wireless display that visibly guides patients to promote better brushing habits. The unique stand-alone SmartGuide provides real-time feedback and encourages longer brushing time.

Each221-2173 1@ \$124.00



Orthodontic Brush Refill
With special tufting pattern thoroughly cleans and gently removes plaque from teeth with braces.
1-pack223-1694 1@ \$9.50



Power Tip Brush Refill
Clean deep between teeth as well as around orthodontic appliances.
1-pack223-2171 1@ \$9.50



Oral-B® BracesCare Essentials
Ideal for Orthodontic Patients. Each kit contains: 2 Power Tip® brushheads and 2 Ortho brushheads
1-pack224-1081 1@ \$18.95