

Pins & Posts

PINS

Bondent® – Coltène/Whaledent Self-Threading Pin System

- Titanium Alloy or Gold Plated Stainless Steel
- 2 head sizes: Long = 2.5mm (yellow) and Short = 2.0mm (blue)
- .017"/0.425mm

Bondent Complete Kits

Contains: 10 long (yellow), 10 short (blue), 2 drills, 1 universal hand driver and FG carbide bur.

Gold-Plated Stainless Steel Complete Kit (DB-20SL)

Each 719-0507 1@ \$237.60

Titanium Alloy Kit

Complete Kit (DB-20SLT)
Each 719-0002 1@ \$237.60

Bondent Bulk Kits

Contains: 50 long (yellow) pins.

Long – Gold Plated Stainless Steel Kit

Bulk (DB-50L) 719-0531 1@ \$427.70

Long – Titanium Alloy Kit

Bulk (DB-50LT) 719-0093 1@ \$427.70

Contains: 50 short (blue) pins.

Short – Gold Plated Stainless Steel Kit

Bulk (DB-50S) 719-0549 1@ \$427.70

Short – Titanium Alloy Kit

Bulk (DB-50ST) 719-0101 1@ \$427.70



Bondent® Twist Drills – Coltène/Whaledent

Depth Limiting

6/Pkg
Bronze 719-0515 1@ \$158.35



Filpin™ System – Filhol

Pure Titanium

The Filpin System includes pure titanium self-shearing retention pins for both automatic and manual placement. The special 47° thread form crushes and displaces less dentin and causes less crazing. Contains: 38 pins and 1 twist drill.



Universal
.60mm/.021 708-5384 1@ \$132.70

Large
.76mm/.027 708-5392 1@ \$132.70

R-Series

The Filpin R-Series is used for manual placement.

Universal
.60mm/.021 708-5418 1@ \$132.70

Large
.76mm/.027 708-5426 1@ \$132.70

Filpin™ Twist Drill – Filhol

Universal – .60mm
5/Pkg 708-5459 1@ \$56.95

Large – .76mm
5/Pkg 708-5467 1@ \$56.45



Stabilok Pins – Stabilok

Reinforcement pins

Pkg of 100
Box 1@ \$198.00

Titanium
Small Yellow 703-0182
Medium Orange 703-0174

Pkg of 20
Titanium 1@ \$45.90

Small Yellow 703-0083
Medium Orange 703-0091



Max® Restorative Pins – Coltène/Whaledent



Titanium Alloy, Self-Threading Pin System

- Countersink shoulder stop reduces stress
- Two sizes and colour coded (.017"/blue and .021"/purple)
- Inverted thread design maximizes retention while minimizing the risk of dentin crazing
- Fewer threads displace less dentin

.021" Max Pin Bulk Kit
Contains: 50 pins.
Bulk Kit (M-23) 718-6570 1@ \$278.55

.021" Max Pin Complete Kit
Contains: 20 pins, 2 each .021" drills, 1 universal hand instrument.
Complete Kit (M-21) .. 718-6547 1@ \$188.10

.017" Max Pin Bulk Kit
Contains: 50 pins.
Bulk Kit (M-13) 718-6539 1@ \$278.55

.017" Max Pin Kit
Contains: 20 pins, 2 each .017" drills, 1 universal hand instrument.
Complete Kit (M-11) .. 718-6505 1@ \$188.10

Max® Restorative Drills – Coltène/Whaledent

- Depth limiting
- Blister package
- 2 sizes (.017" & .021")

.017"/0.425mm (blue)
5/Pkg 718-6588 1@ \$178.20

.021"/0.525mm (purple)
5/Pkg 718-6604



Stabilok Twist Drills – Stabilok

Each 1@ \$15.50

Small Blue 660-4938
Medium Green 660-4912
Small Yellow 660-4946
Medium Orange 660-4920

TMS® - Thread Mate System® - Coltène/Whaledent



Self-Threading Pin System

- The original TMS self-threading pins
- Gold plated stainless steel
- Available in 4 sizes
- Manual placement with a hand wrench or mechanical placement with a reduction gear contra angle

TMS - Manual Placement

.027 Regular Complete Kit (gold colour code)
 Contains: 20 pins, 1 Kodex® drill (K-96), and 1 universal hand drive.
 Complete Kit (T-01) ...719-4905 1@ \$151.80

.027 Regular Refill Kit
 Contains: 50 pins.
 Refill (T-02)719-5456 1@ \$242.25

.027 Regular Bulk Kit
 Contains: 100 pins.
 Refill (T-03)719-6371 1@ \$421.10

.021 Minim® Complete Kit (silver colour code)
 Contains: 20 pins, 1 Kodex® drill (K-93) and 1 universal hand drive.
 Complete Kit (T-04) ...719-6868 1@ \$151.80

.021 Minim® Refill Kit
 Contains: 50 pins.
 Refill (T-05)719-7403 1@ \$242.25

.021 Minim® Bulk Kit
 Contains: 100 pins.
 Bulk (T-06)720-7210 1@ \$421.10

TMS Two in One

.027 Regular Complete Kit (gold colour code)
 Contains: 20 double shear pins, 1 Kodex® drill (K-96) and 1 universal hand drive.
 Complete Kit (T-07) ...720-7707 1@ \$202.00

.027 Regular Refill Kit
 Contains: 50 double shear pins.
 Refill (T-08)720-8291 1@ \$347.20

.027 Regular Bulk Kit
 Contains: 100 double shear pins.
 Bulk (T-09)720-8689 1@ \$636.25

TMS Two in One

.021 Minim® Complete Kit (silver colour code)
 Contains: 20 double shear pins, 1 Kodex® drill (K-92), 1 universal hand drive.
 Complete Kit (T-10) ...720-9174 1@ \$202.00

.021 Minim® Refill Kit
 Contains: 50 double shear pins.
 Refill (T-11)720-9679 1@ \$347.20

.021 Minim® Bulk Kit
 Contains: 100 double shear pins.
 Bulk (T-12)721-0156 1@ \$634.25

TMS Self-Shearing

.027 Regular Complete Kit (gold colour code)
 Contains: 20 pins, 1 Kodex® drill (K-97), 1 universal hand drive.
 Complete Kit (T-13) ...719-7841 1@ \$202.00

.027 Regular Refill Kit
 Contains: 50 pins.
 Refill (T-14)720-5305 1@ \$347.20

.027 Regular Bulk Kit
 Contains: 100 pins.
 Bulk (T-15)720-6725 1@ \$636.25

.017 Minikin® Complete Kit (red colour code)
 Contains: 20 pins, 1 Kodex® drill (K-91) and 1 universal hand drive.
 Complete Kit (T-107) .724-8735 1@ \$202.00

.017 Minikin® Refill Kit
 Contains: 50 pins.
 Refill (T-108)724-9253 1@ \$347.20

.017 Minikin® Bulk Kit
 Contains: 100 pins.
 Bulk (T-109)724-9733 1@ \$636.25

.0135 Minuta® Complete Kit (pink colour code)
 Contains: 20 pins, 1 Kodex® drill (K-89), 1 universal hand drive.
 Complete Kit (T-111) .724-9741 1@ \$202.00

.0135 Minuta® Refill Kit
 Contains: 50 pins.
 Refill (T-112)724-9808 1@ \$347.20

.0135 Minuta® Bulk Kit
 Contains: 100 pins.
 Bulk (T-113)724-9824 1@ \$636.25



TMS Link Plus® - Coltène/Whaledent

Self-Threading Pin System

- Built-in shoulder stop guarantees shearing and proper pin placement
- Titanium alloy or gold plated stainless steel
- 4 sizes, colour coded
- Mechanical or manual placement

Stainless Steel - TMS Link Plus

.021 Minim Single Shear Complete Kit
 Contains: 20 pins, 1 Kodex® drill (K-92) and 1 hand driver.
 Kit (L-721)722-6806 1@ \$159.40

.021 Minim Single Shear Bulk
 Contains: 50 pins.
 Bulk (L-722)722-6814 1@ \$271.95

.021 Minim Double Shear Complete Kit
 Contains: 20 double shear pins, 1 Kodex® drill (K-92) and 1 hand driver.
 Kit (L-731)722-6830 1@ \$208.60

.021 Minim Double Shear Bulk
 Contains: 50 double shear pins.
 Bulk (L-732)722-6848 1@ \$370.95

.027 Regular Single Shear Complete Kit
 Contains: 20 pins, 1 Kodex® drill (K-96) and 1 hand driver.
 Kit (L-741)722-6855 1@ \$159.40

.027 Regular Single Shear Bulk
 Contains: 50 pins.
 Bulk (L-742)722-6822 1@ \$271.95

.027 Regular Double Shear Complete Kit
 Contains: 20 double shear pins, 1 Kodex® drill (K-96) and 1 hand driver.
 Kit (L-751)722-6897 1@ \$208.60

.027 Regular Double Shear Bulk
 Contains: 50 double shear pins.
 Bulk (L-752)722-6863 1@ \$370.95

Titanium Alloy - TMS Link Plus

.021 Minim Single Shear Complete Kit
 Contains: 20 pins, 1 Kodex® drill (K-92) and 1 hand driver.
 Kit (L-821)722-6889 1@ \$159.40

.021 Minim Single Shear Bulk
 Contains: 50 pins.
 Bulk (L-822)722-6913 1@ \$271.95

.021 Minim Double Shear Complete Kit
 Contains: 20 double shear pins, 1 Kodex® drill (K-92) and 1 hand driver.
 Kit (L-841)722-6939 1@ \$208.60

.021 Minim Double Shear Bulk
 Contains: 50 double shear pins.
 Bulk (L-842)722-6947 1@ \$370.95

TMS® Accessories

Kodex® Drills - Coltène/Whaledent



- Colour coded for use with Minuta®, Minikin®, Minim®, and Regular TMS pins
- Depth limiting

6/Pkg 1@ \$102.00
 3@ \$98.70

- K-89 Pink.....721-3895 K-94 Black ...721-4034
 - K-90 Red.....721-3911 K-95 Black ...721-4067 *
 - K-91 Red.....721-3945 K-96 Gold....721-4091
 - K-92 Silver...721-3978 K-97 Gold....721-4125 *
 - K-93 Silver...721-4000* K-98 Green...721-4158 *
- *not depth limiting

TMS Hand Wrench - Coltène/Whaledent

For TMS Regular, Regular 2-in-1, Minikin, Self-Shearing.
 #B-59-3721-1865 1@ \$19.70

For TMS Minim, Minim 2-in-1
 #B-59-5721-1899 1@ \$19.70

For TMS Minuta
 #B-59-10721-1931 1@ \$19.70

TMS Bending Tool - Coltène/Whaledent
 Regular Gold.....721-0651 1@ \$28.85
 Minim Silver721-0669

PINS (CONT.)



TMS Link Series® - Coltène/Whaledent Self-Threading Pin System

TMS (Thread Mate System) offers:

- Titanium alloy or gold-plated stainless steel
- 4 sizes, colour coded
- Ball joint link to shank provides smooth alignment during placement
- Mechanical or manual placement

Stainless Steel - TMS Link Series

.0135 Minuta Single Shear Complete Kit

Contains: 20 pins, 1 Kodex® drill (K-89) and 1 hand driver.

Kit (L-501).....722-6509 1@ \$180.20

.0135 Minuta Single Shear Bulk Kit

Contains: 50 pins.

Bulk (L-502).....722-6517 1@ \$367.00

.017 Minikin Single Shear Complete Kit

Contains: 20 pins, 1 Kodex® drill (K-91) and 1 hand driver.

Kit (L-511).....722-6525 1@ \$180.20

.017 Minikin Single Shear Bulk

Contains: 50 pins.

Bulk (L-512).....722-6533 1@ \$367.00

.021 Minim Single Shear Complete Kit

Contains: 20 pins, 1 Kodex® drill (K-92) and 1 hand driver.

Kit (L-521).....722-6541 1@ \$180.20

.021 Minim Single Shear Bulk

Contains: 50 pins.

Bulk (L-522).....722-6558 1@ \$277.20

.021 Minim Double Shear Complete Kit

Contains: 20 double shear pins, 1 Kodex® drill (K-92) and 1 hand driver.

Kit (L-531).....722-6566 1@ \$214.50

.021 Minim Double Shear Bulk

Contains: 50 double shear pins.

Bulk (L-532).....722-6574 1@ \$367.00

.027 Regular Single Shear Complete Kit

Contains: 20 pins, 1 Kodex® drill (K-96) and 1 hand driver.

Kit (L-541).....722-6582 1@ \$180.20

.027 Regular Single Shear Bulk

Contains: 50 pins.

Bulk (L-542).....722-6590 1@ \$277.20

.027 Regular Double Shear Complete Kit

Contains: 20 double shear pins, 1 Kodex® drill (K-96) and 1 hand driver.

Kit (L-551).....722-6608 1@ \$214.50

.027 Regular Double Shear Bulk

Contains: 50 double shear pins.

Bulk (L-552).....722-6616 1@ \$367.00

Titanium Alloy - TMS Link Series

.017 Minikin Single Shear Complete Kit

Contains: 20 pins, 1 Kodex® drill (K-91) and 1 hand driver.

Kit (L-811).....722-6921 1@ \$180.20

.017 Minikin Single Shear Bulk

Contains: 50 pins.

Bulk (L-812).....722-6905 1@ \$367.00



Introducing Patterson's Technical Service Centre

Fast - handpieces repaired by the next business day

Reliable - fully trained and experienced in the products we service

Comprehensive - we offer the most repair services available

Convenient - as easy as a phone call

GUARANTEED!

Call Toll-Free 1-800-240-4604

Post Chart

This chart is designed as a quick reference to popular products available in this section. Please check with your local Patterson representative for additional information.

Post	Pg #	Name	Manufacturer	Metal Free		Metal			Shape		Design	
				Fiber/Resin	Plastic	Titanium	Gold Plate	Stainless	Parallel Side	Tapered	Passive	Active
	-	Achromat Post	Axis	✓					✓		✓	
	-	Cytco - K	Dentsply/Maillefer			✓			✓	✓		✓
	553	Flexi-Flange	EDS			✓		✓	✓			✓
	553	Flexi-Overdenture	EDS			✓		✓	✓			✓
	553	Flexi-Post	EDS			✓		✓	✓			✓
	553	Filpost System	Filhol Dental USA			✓				✓	✓	
	552	FRC Postec	Vivadent	✓						✓	✓	
	552	Glass Fibre Post	Ellman International	✓						✓	✓	
	554	IntegraPost	Premier			✓			✓		✓	
	554	KG-Endowel	Star Dental		✓					✓	✓	
	554	Luminex Light Transmitting	Dentatus USA			✓		✓	✓	tip	✓	
	552	Luscent Anchors	Dentatus USA	✓						✓	✓	
	554	NuBond Post	Ellman International		✓				✓		✓	
	555	ParaPost Burnout Post	Coltène/Whaledent		✓				✓		✓	
	552	ParaPost Fiber Lux	Coltène/Whaledent	✓					✓		✓	
	552	ParaPost Fiber White	Coltène/Whaledent	✓					✓		✓	
	555	ParaPost Impression	Coltène/Whaledent		✓				✓		✓	
	554	ParaPost Plus	Coltène/Whaledent			✓		✓	✓		✓	
	555	ParaPost System	Coltène/Whaledent			✓		✓	✓		✓	
	-	ParaPost Temporary	Coltène/Whaledent			✓			✓		✓	
	556	ParaPost XH	Coltène/Whaledent			✓			✓		✓	
	555	ParaPost XP	Coltène/Whaledent			✓		✓	✓		✓	
	555	ParaPost XT	Coltène/Whaledent			✓			✓			✓
	556	Radix-Anchor	Dentsply/Maillefer			✓			✓			✓
	553	Snowpost	Danville	✓					✓	end	✓	
	556	Surtex	Dentatus USA			✓	✓	✓	✓	end	✓	
	-	Vented Posts	Roydent			✓		✓	✓		✓	

POSTS – METAL FREE



FRC Postec® – Ivoclar Vivadent

Glass Fibre Reinforced Composite Posts

FRC Postec provides a biocompatible, aesthetic, metal-free alternative to traditional metal posts.

- Outstanding aesthetics
- Elastic behavior similar to that of dentin
- Excellent fracture resistance
- Light transmission
- 2 sizes: size 1 (0.8mm), size 3 (1.0mm)

Refills

FRC Postec, 5 each		
Size 1	710-6222	1@ \$91.70
Size 3	710-6230	
FRC Postec Pilot Reamer		
Size 1	710-6321	1@ \$34.40
Size 3	710-6339	
FRC Postec Impression Posts, 5 each		
Size 1	710-6271	1@ \$34.80
Size 3	710-6289	

Luscent Anchors™ – Dentatus USA



The resin reinforced fibre glass Luscent Anchors fulfill dentistry's preference for very strong, non-metallic posts that are consistent with adhesive practices and procedures. Radiolucent.

Introductory Kit

Contains: 12 anchors (4 of each – sizes 1.4mm, 1.6mm, 1.8mm), 3 depth-marked reamers, 1 Probos II instrument and convenient core forms-to-fit.
Introductory Kit.....237-3405 1@ \$343.20

Luscent Anchors Refills

Contains: 12 Luscent Anchors, 1 matching reamer, 1 Probos II router and 12 core forms-to-fit.
Large 1.8mm.....237-4718 1@ \$250.75
Medium 1.6mm.....237-4726
Small 1.4mm.....237-4734

Glass Fibre Post – Ellman International



Kit contains 6 posts each (#1-4), 4 reamers (#1-4) and 1 screwdriver.
Kit315-2659 1@ \$442.20

Glass Fibre Post Refills

6/Pkg		
Size 1	315-2675	1@ \$80.85
Size 2	315-2683	
Size 3	315-2691	
Size 4	315-2709	

ParaPost® Fiber White – Coltène/Whaledent



ParaPost Fiber White is for all of your aesthetic needs. Features include: metal free aesthetics, uni-directional glass fibre/resin formulation, anti-rotational head and passive parallel design. ParaPost Fiber White Posts use standard Parapost drills located on page 555.

Introductory Kit

Contains: 10 assorted posts, 4 drills and 1 hand-driver.
Kit721-5817 1@ \$242.10

ParaPost Fiber White Refills

5/Pkg		
PF-161-4.5 (.045).....	721-5833	1@ \$71.25
PF-161-5 (.050).....	721-5841	3@ \$68.85
PF-161-5.5 (.055).....	721-5858	
PF-161-6 (.060).....	721-5866	

ParaPost Fiber Lux – Coltène/Whaledent



- Translucent
- Superior head design – double undercuts lock in core material
- Enhanced radiopacity
- Colour-coded – easy identification

Intro Kit

Contains: 15 posts (2 of sizes 3, 5.5, 6) (3 of sizes 4, 4.5, 5), 6 drills (1 for each size), 1 universal hand driver.
Introductory Kit.....722-3308 1@ \$290.80

Refills

6/Pkg		
Size 3	722-3324	1@ \$68.10
Size 4	722-3332	
Size 4.5	722-3340	
Size 5	722-3357	
Size 5.5	722-3365	
Size 6	722-3373	

ParaPost® Taper Lux – Coltène/Whaledent

Intro Kit

Contains 15 assorted posts 4 Taper Lux drills and 1 universal hand driver
Each722-3035 1@ \$288.50

Refills

5/Pkg		
4.5 (Blue).....	722-3050	1@ \$68.10
5.0 (Red).....	722-3068	3@ \$66.50
5.5 (Purple).....	722-3076	
6.0 (Black).....	722-3084	

Taper Lux Drills

3/Pkg		
.045" (Blue/Black).....	722-3126	1@ \$67.25
.050" (Red/Black).....	722-3134	1@ \$65.20
.055" (Purple/Black).....	722-3142	
.060" (Black/Black).....	722-3159	

P82A Assorted

Contains 4 drills (1 of each size)
Each722-3209 1@ \$89.70



RelyX™ Fiber Post – 3M ESPE

RelyX Fiber Post has elasticity more similar to dentin than to metal and ceramic posts – so there's a lower risk of root fracture due to the "wedge effect".

Intro Kit

Contains 5 of each: size 1 post (1.3mm dia.), size 2 post (1.6mm dia.), 4 drills: size 1, size 2, size 3 and a universal drill, 20 elongation tips, instructions and technique card
Each503-4210 1@ \$317.00

RelyX™ Fiber Post/RelyX™ Unicem Cement Combination Kit

Contains 5 of each size posts: 1, 2, 3; 4 of each size drills: 1, 2, 3 and a universal drill; 20 elongation tips, instructions, technique card and 1 RelyX Unicem cement trial kit including 20 assorted capsules, 1 aplicap activator/applier set, instructions and a technique card
Each503-4228 1@ \$410.00

Fiber Post Refill Packs

Contains 10 posts and 10 Aplicap elongation tips
Size 1

Size 1	503-4244	1@ \$134.00
Size 2	503-4251	1@ \$128.50
Size 3	503-4269	

Drills

Size 1	503-4319	1@ \$27.25
Size 2	503-4327	1@ \$26.50
Size 3	503-4335	

Drill Reamer

Each503-4376 1@ \$27.75

RelyX Unicem Aplicap Elongation Tips

10/Pkg.....503-3899 1@ \$22.50
1@ \$22.15

Snowpost – Danville Materials

Aesthetics Post

Snowposts are white, transparent, composite fibre posts that offer excellent aesthetics.

- Composite/silica-zirconia fibre post system
- 4 sizes: 1.0, 1.2, 1.4 and 1.6mm
- 60% filled for high strength
- Flexural modulus closely matches dentin to distribute stress

Introductory Kit

Contains: 20 posts (4 sizes) and matching reamers.

Kit236-2101 1@ \$245.55

Refills

Contains: 10 posts of one size and 1 reamer.

1.0mm (90082).....236-2127 1@ \$128.00
 1.2mm (90083)236-2135
 1.4mm (90088).....236-2143
 1.6mm (90089).....236-2150



Flexi-Flange® – Essential Dental Systems

Flexi-Flange is a post designed for situations where there is inadequate coronal dentin and excessive occlusal force.

Introductory Kit

Includes 12 posts (4 each sizes 0,1,2)
 Kit309-3507 1@ \$288.70

Refills

Contains: 10 posts, reamer and countersink drill. Available in Stainless Steel or Titanium.

Stainless Steel

Size 0 Yellow.....309-3549 1@ \$198.00
 Size 1 Red.....309-3556
 Size 2 Blue.....309-3564
 Size 3 Green.....309-3572

Titanium

Size 0 Yellow.....309-3713 1@ \$207.00
 Size 1 Red.....309-3721
 Size 2 Blue.....309-3739
 Size 3 Green.....309-3747



Flexi-Post® – Essential Dental Systems

Prefabricated Split Shank Post

- Split shank minimizes stress to the root
- Parallel-sided and threaded for retention
- Maximum retention with minimal stress

Introductory Kits

Contains: 12 serrated posts (4 of each size), 3 primary reamers, 3 secondary drills, 1 extender, 2 wrenches and 1 gauge.


Stainless Steel Kit, size #0, 1, 2
 Each309-5197 1@ \$288.70

Stainless Steel Kit, size #1, 2, 3
 Each309-5221 1@ \$288.70

Titanium Kit, size #0, 1, 2
 Each309-5205 1@ \$299.25



POSTS – METAL



EZ-Change™ – Essential Dental Systems

Flexi-Overdenture Accessories

Patented keeper and cap insert system allows for “quick and easy” nylon cap replacement. When worn, the cap insert may be easily removed from the permanent metal keeper, requiring only seconds for replacement.

EZ-Change Intro Kit

Contains: 2 keepers, 2 cap inserts, 2 rubber bands & EZ-Change Wrench
 Kit309-4547 1@ \$113.40

Combo Refill

Contains: 2 keepers, 2 caps, 2 rubber bands
 Pkg.....309-4554 1@ \$84.65

Cap Inserts

6/Pkg.....309-4562 1@ \$72.60

Filpost™ System – Filhol

Titanium Alloy

An anatomically designed, passive-fitting, prefabricated post that offers preservation of tooth structure and enhanced retention of build-up materials. Contains: 20 titanium posts, 1 FilLock reamer.

Standard

1.3mm x 16mm.....708-5335 1@ \$188.70

Large

1.6mm x 20mm.....708-5343 1@ \$188.45



Flexi-Overdenture® Attachment Kits – Essential Dental Systems

The unique split shank of the Flexi-Post, with a ball and socket attachment for the highest retention and ideal distribution of stress. For full or partial dentures.

Flexi-Overdenture Introductory Kit – Stainless Steel

Contains: 12 posts (6 each, denture post #1 and #2), and 12 nylon attachment caps, 2 primary reamers, 2 countersink drills, 1 wrench, 1 depth gauge.
 Stainless Steel309-5387 1@ \$459.40

Flexi-Overdenture Introductory Kit – Titanium

Contains: 4 posts (2 each, denture post #1 and #2), and 4 nylon attachment caps, 2 primary reamers, 2 countersink drills, 1 wrench, 1 depth gauge.
 Stainless Steel309-5379 1@ \$252.40

Flexi-Overdenture Refills – Stainless Steel

Contains: 6 denture posts, 6 attachment caps, 1 primary reamer, 1 countersink drill.
 Size 1 Red309-5403 1@ \$225.20
 Size 2 Blue309-5411
 Size 3 Green.....309-5429

Overdenture Attachment Caps® – Essential Dental Systems

Contains 6 nylon female attachment caps
 6/Pkg.....309-5742 1@ \$66.55



Refills – Flexi-Post

Contains: 10 posts, reamer and drill.

	Titanium	Stainless Steel
Size 00 White ..	309-4414	309-5312
Size 0 Yellow ..	309-4422	309-5304
Size 1 Red	309-4430	309-5270
Size 2 Blue	309-4448	309-5288
Size 3 Green	309-4455	309-5296
	1@ \$207.00	1@ \$197.80

Bulk Refills – Flexi-Post

Contains: 30 posts, reamer and drill.

	Titanium	Stainless Steel
Size 00 White ..	309-4497	309-5601
Size 0 Yellow ..	309-4505	309-5619
Size 1 Red	309-4513	309-5627
Size 2 Blue	309-4521	309-5635
Size 3 Green	309-4539	309-5643
	1@ \$525.90	1@ \$504.80

Flexi-Post® Primary Reamers – Essential Dental Systems

For Flexi-Post System

Size 00309-5767 1@ \$25.75
 Size 0309-5775
 Size 1309-5783
 Size 2309-5791
 Size 3309-5809

Flexi-Post® Secondary Drills – Essential Dental Systems

For Second Tier Preparation

Size 0, Yellow.....309-5874 1@ \$25.75
 Size 1, Red.....309-5882
 Size 2, Blue.....309-5890
 Size 3, Green.....309-5908



Pins & Posts

POSTS – METAL (CONT.)

Flexi-Post® Wrench – Essential Dental Systems

For Serrated Head Posts
For Sizes 00 through 3.

External	309-5494	1@ \$25.75
Internal.....	309-5247	1@ \$25.75



IntegraPost™ – Premier



The ultra-passive, ultra-secure integrating core retainer with the patented Flow-through head that ensures full integration of the post head and core buildup.

Intro Kit (in nonspillable, autoclavable Hot Box)

Contains: 18 posts (three each #3, #5; five each #4, #4.5; two each #6), 1 IP carrier, 1 IP probe, 1 assorted package of 5 drills.
Kit 579-3073 1@ \$409.20

Mini Intro Kit (in lucite block)

Contains: 6 posts (two each #3, #4, #4.5), 3 drills (one each #3, #4, #4.5), 1 IP carrier.
Kit 579-3005 1@ \$124.40

IntegraPost™ Refills – Premier

10/Pkg		
Size 3	579-3047	1@ \$121.45
Size 4	579-3054	
Size 4.5	579-3062	
Size 4.5L	579-3070	
Size 5	579-3088	
Size 6	579-3096	

IntegraPost™ Integradrills – Premier

3/Pkg		
Size 3	579-3260	1@ \$63.35
Size 4	579-3278	3@ \$62.05
Size 4.5	579-3286	
Size 5	579-3294	
Size 6	579-3302	

5/Pkg		
Assorted	579-3369	1@ \$93.55

Integra Instruments – Premier

IP Carrier	580-0206	1@ \$21.85
IP Probe	580-0214	1@ \$32.95



KG Endowel Posts – Star Dental

Plastic, 50/Pkg		
#80	645-7659	1@ \$33.80
#90	645-7683	
#100	645-7717	
#120	645-7741	
#140	645-7774	
Assorted	645-7808	



Luminex® Light Transmitting Post System – Dentatus USA

For restoring compromised thin-walled roots with strong adhesive material.

Luminex Introductory Kit

Stainless Steel	765-6507	1@ \$342.55
Titanium	765-6515	

Luminex Starter Kit

Stainless Steel	765-6523	1@ \$463.35
-----------------------	----------	-------------



NuBOND™ Post – Eilman International



- Cement-in system reduces root stress
- Self-limiting reamers prevent perforations
- Tapered posts and reamers prevent root fractures

NuBOND Post Kit

Contains 36 posts (assorted sizes 1-6) and 1 vial of Cyanodent - Slow Set. Kit (EPO165) 315-3236 1@ \$343.15

NuBOND Refills 6/Pkg

Size 1	315-3251	1@ \$50.95
Size 2	315-3269	
Size 3	315-3277	
Size 4	315-3285	
Size 5	315-3293	
Size 6	315-3301	

NuBOND Reamers

Tapered, self limiting.

Size 1	315-3400	1@ \$38.60
Size 2	315-3418	
Size 3	315-3426	
Size 4	315-3434	
Size 5	315-3442	
Size 6	315-3459	



ParaPost® Plus – Coltène/Whaledent

Cylindrical, passively cemented endodontic post system for direct post and core build-up technique. Spiral flutes provide excellent venting of excess cement and resistance to rotation. Available in Stainless Steel and Titanium Alloy.

ParaPost Plus, Stainless Steel

Refill (Bulk), 25/Pkg

P-244-3B	.036 Brown ..	722-2219	1@ \$312.85
P-244-4B	.040 Yellow ..	722-2227	
P-244-4.5B	.045 Blue	722-2318	
P-244-5B	.050 Red	722-2235	
P-244-5.5B	.055 Purple ..	722-2326	
P-244-6B	.060 Black	722-2243	
P-244-7B	.070 Green ..	722-2250	

Refill, 10/Pkg

P-244-3	.036 Brown ..	722-2136	1@ \$134.65
P-244-4	.040 Yellow ..	722-2144	3@ \$132.00
P-244-4.5	.045 Blue	722-2185	
P-244-5	.050 Red	722-2151	
P-244-5.5	.055 Purple ..	722-2193	
P-244-6	.060 Black	722-2169	
P-244-7	.070 Green ..	722-2177	

ParaPost Plus, Titanium

Refill (Bulk), 25/Pkg

P-284-3B	.036 Brown ..	722-2706	1@ \$303.60
P-284-4B	.040 Yellow ..	722-2714	
P-284-4.5B	.045 Blue	722-2755	
P-284-5B	.050 Red	722-2722	
P-284-5.5B	.055 Purple ..	722-2763	
P-284-6B	.060 Black	722-2730	
P-284-7B	.070 Green ..	722-2748	

Refill, 10/Pkg

P-284-3	.036 Brown ..	722-2417	1@ \$134.65
P-284-4	.040 Yellow ..	722-2425	3@ \$132.00
P-284-4.5	.045 Blue	722-2482	
P-284-5	.050 Red	722-2433	
P-284-5.5	.055 Purple ..	722-2490	
P-284-6	.060 Black	722-2441	
P-284-7	.070 Green ..	722-2458	

ParaPost® System – Coltène/Whaledent

Universal endodontic post system for direct post and core buildups and all casting techniques. Available in Stainless Steel and Titanium.

ParaPost® System, Stainless Steel



10/Pkg			
P-44-3	.036 Brown ..723-5781	1@	\$131.35
P-44-4	.040 Yellow ..723-5815		
P-44-5	.050 Red723-5849		
P-44-5.5	.055 Purple ..723-5914		
P-44-6	.060 Black723-5872		
P-44-7	.070 Green ..723-5906		

25/Pkg			
P-44-3B	.036 Brown ..723-6094	1@	\$310.20
P-44-4B	.040 Yellow ..723-6102		
P-44-4.5B	.045 Blue723-6144		
P-44-5B	.050 Red723-6110		
P-44-5.5B	.055 Purple ..723-6151		
P-44-6B	.060 Black723-6128		
P-44-7B	.070 Green ..723-6136		

ParaPost® System, Titanium, 10/Pkg			
P-84-3	.036 Brown ..723-5930	1@	\$131.35
P-84-4	.040 Yellow ..723-5948		
P-84-4.5	.045 Blue723-5989		
P-84-5	.050 Red723-5955		
P-84-5.5	.055 Purple ..723-5997		
P-84-6	.060 Black723-5963		
P-84-7	.070 Green ..723-5971		



ParaPost® Plastic Burnout Posts – Coltène/Whaledent

10/Pkg			
P-51-3	.036 Brown ..723-9015	1@	\$65.00
P-51-4	.040 Yellow ..723-9023	2@	\$63.05
P-51-4.5	.045 Blue723-9064		
P-51-5	.050 Red723-9031		
P-51-5.5	.055 Purple ..723-9072		
P-51-6	.060 Black723-9049		
P-51-7	.070 Green ..723-9056		

ParaPost® Plastic Impression Posts – Coltène/Whaledent

20/Pkg			
P-43-3	.036 Brown ..723-2614	1@	\$54.50
P-43-4	.040 Yellow ..723-2622		
P-43-4.5	.045 Blue723-2671		
P-43-5	.050 Red723-2630		
P-43-5.5	.055 Purple ..723-2663		
P-43-6	.060 Black723-2648		
P-43-7	.070 Green ..723-2655		

ParaPost® Titanium Alloy Temporary Posts – Coltène/Whaledent

For Temporary Use

20/Pkg			
P-46-3	.036 Brown ..721-8654	1@	\$92.35
P-46-4	.040 Yellow ..721-8688		
P-46-5	.050 Red721-8712		
P-46-6	.060 Black721-8746		
P-46-7	.070 Green ..721-8779		



ParaPost®, ParaPost® Plus, and ParaPost® X™ Drills – Coltène/Whaledent

Side-Cutting/Two-Fluted, 3/Pkg			
P-42-3	.036 Brown ..722-7499	1@	\$76.50
P-42-4	.040 Yellow ..722-7523	3@	\$72.60
P-42-4.5	.045 Blue722-7622		
P-42-5	.050 Red722-7556		
P-42-5.5	.055 Purple ..722-7630		
P-42-6	.060 Black722-7580		
P-42-7	.070 Green ..722-7614		

P-42-A, 6 sizes 3-6.....	722-6954	1@	\$92.35
--------------------------	----------	----	---------



ParaPost® XT™ – Coltène/Whaledent

One-Office Visit, Titanium Alloy Post



Introductory Kit
Contains: 25 assorted posts, 6 XT drills, 1 hand driver, metal wrench A & B.
Kit..... 721-8795 1@ \$434.75

Refills			
10/Pkg			
P-683-0	.036 Brown ..722-4611	1@	\$140.90
P-684-0	.040 Yellow ..722-4629		
P-684-5	.045 Blue722-4637		
P-685-0	.050 Red722-4645		
P-685-5	.055 Purple ..722-4652		
P-686-0	.060 Black722-4660		

30/Pkg			
P-683-0B	.036 Brown ..722-4686	1@	\$356.40
P-684-0B	.040 Yellow ..722-4694		
P-684-5B	.045 Blue722-4702		
P-685-0B	.050 Red722-4710		
P-685-5B	.055 Purple ..722-4728		
P-686-0B	.060 Black722-4736		



ParaPost® XT™ Drill – Coltène/Whaledent

Smooth Cutting and Laser Marked

3/Package			
Size 3 Brown.....	721-4216	1@	\$92.25
Size 4 Yellow.....	721-4224	2@	\$89.30
Size 4.5 Blue.....	721-4232		
Size 5 Red.....	721-4240		
Size 5.5 Purple.....	721-4257		
Size 6 Black.....	721-4265		

ParaPost® XT™ Metal Wrench – Coltène/Whaledent

Size A for 3,4 4.5.....	718-2058	1@	\$43.70
Size B for 5,5.5,6.....	718-2066		
Assorted A&B.....	718-2074	1@	\$86.80



ParaPost® XP™ – Coltène/Whaledent

- 47% more resistance to rotation
- 50% more tensile retention
- 68% more resistance to oblique forces than ParaPost or ParaPost Plus due to the unique patented raised diamond pattern on the post shaft



ParaPost XP Introductory Kits

One-Office-Visit Kit
25 assorted posts, 7 X drills and 1 hand driver.
Titanium Alloy Kit 721-4620 1@ \$377.50
Stainless Steel Kit 721-4612 1@ \$377.50

Refills – Stainless Steel

10/Pkg			
P-744-3	.036 Brown ..721-4885	1@	\$104.60
P-744-4	.040 Yellow ..721-4893		
P-744-4.5	.045 Blue721-4901		
P-744-5	.050 Red721-4919		
P-744-5.5	.055 Purple ..721-4927		
P-744-6	.060 Black721-4935		
P-744-7	.070 Green ..721-4943		

25/Pkg			
P-744-3B	.036 Brown ..721-4992	1@	\$254.75
P-744-4B	.040 Yellow ..721-5007		
P-744-4.5B	.045 Blue721-5015		
P-744-5B	.050 Red721-5023		
P-744-5.5B	.055 Purple ..721-5031		
P-744-6B	.060 Black721-5049		
P-744-7B	.070 Green ..721-5056		

Refills – Titanium

10/Pkg			
P-784-3	.036 Brown ..721-5114	1@	\$104.60
P-784-4	.040 Yellow ..721-5122		
P-784-4.5	.045 Blue721-5130		
P-784-5	.050 Red721-5148		
P-784-5.5	.055 Purple ..721-5155		
P-784-6	.060 Black721-5163		
P-784-7	.070 Green ..721-5171		

25/Pkg			
P-784-3B	.036 Brown ..721-5312	1@	\$254.75
P-784-4B	.040 Yellow ..721-5320		
P-784-4.5B	.045 Blue721-5338		
P-784-5B	.050 Red721-5346		
P-784-5.5B	.055 Purple ..721-5353		
P-784-6B	.060 Black721-5361		
P-784-7B	.070 Green ..721-5379		

ParaPost® Hand Driver Universal – Coltène/Whaledent

Each 722-6632 1@ \$7.00

POSTS – METAL (CONT.)

ParaPost XP™ Casting Technique – Coltène/Whaledent

Contains 25 plastic impression posts, 25 titanium alloy temporary posts, 25 plastic burnout posts, 7 ParaPost X Drills (1 of each size), universal hand driver.
 P-781.....721-4646 **1@ \$384.20**

Refills ParaPost XP™ Plastic Burnout Posts 10/Pkg
 P-751-3.....721-4661 **1@ \$56.65**
 P-751-4.....721-4679
 P-751-4.5.....721-4687
 P-751-5.....721-4695
 P-751-5.5.....721-4703
 P-751-6.....721-4711
 P-751-7.....721-4729



ParaPost® XH™ – Coltène/Whaledent

Titanium Alloy
 A post head with rounded edges to minimize stress points when using composite or glass ionomer cores. Multiple channels for improved cement venting.

ParaPost XH Introductory Kits

One-Office-Visit Kit
 25 assorted posts, 7 X drills and 1 hand driver.
 Kit.....721-4638 **1@ \$421.15**

Refills

10/Pkg
 P-88-3 .036 Brown ..721-5460 **1@ \$113.20**
 P-88-4 .040 Yellow ..721-5478
 P-88-4.5 .045 Blue721-5486
 P-88-5 .050 Red721-5494
 P-88-5.5 .055 Purple ..721-5528
 P-88-6 .060 Black721-5502
 P-88-7 .070 Green ..721-5510

25/Pkg
 P-88-3B .036 Brown ..721-5569 **1@ \$266.65**
 P-88-4B .040 Yellow ..721-5577
 P-88-4.5B .045 Blue721-5585
 P-88-5B .050 Red721-5593
 P-88-5.5B .055 Purple ..721-5601
 P-88-6B .060 Black721-5619
 P-88-7B .070 Green ..721-5627



SURTEX™ – Dentatus USA

Surface Treated, Classic Post
 The Surtex surface treatment provides greater surface area for micro-adhesion and an anti-rotational and a non-reflective surface to enhance aesthetic results.

Surtex Classic Bulk Posts Kit

Bulk kit include 120 assorted posts, 6 assorted reamers, 1 hollow key, 1 cross key and 20 stoppers. Available in gold plate, titanium and stainless steel.

Stainless Steel Bulk Kit (STS-B).....237-3314 1@ \$1,419.00

SURTEX™ Surface Treated Stainless Steel Post Refills – Dentatus USA

Size/ Diameter	1	2	3	4	5	6
Short	237-3868	237-3876	237-3884	237-3892	237-3900	237-3918
Medium	237-3926	237-3934	237-3942	237-3959	237-3967	237-3975
Long	237-3983	237-3991	237-4007	237-4015	237-4023	237-4031
X-Long	N/A	N/A	237-4049	237-4056	237-4064	237-4072
S-Long	N/A	N/A	237-4080	237-4098	237-4106	237-4114

Length: Short - 7.8mm, Medium - 9.3mm, Long - 11.8mm, Extra Long - 14.2mm, Super Long - 17mm
12/Package Each: 1@ \$159.75



Radix-Anchor® System – Dentsply Maillefer

Titanium Alloy
 A one-piece, fully vented parallel post means less chance of excess pressure at the apical portion. Self-tapping helical blades make these posts ideal for reconstruction of non-vital teeth without the excess stress of fully threaded posts. The Radix-Anchor core section has 20 undercuts that secure any core material to the new surface.

Introductory Kit

Contains: 12 posts of 1 size and accessories.
 Kit #1182-1412 **1@ \$234.50**
 Kit #2182-1420
 Kit #3182-1438

Radix Post Refills, 12/Pkg


Kit #1182-1495 **1@ \$193.95**
 Kit #2182-1503 **3@ \$188.95**
 Kit #3182-1511

Unity Post Drills – Coltène/Whaledent

Side-cutting, 3/package.
 .030 Brown.....724-7117 **1@ \$116.70**
 .040 Yellow724-7125 **3@ \$113.05**
 .050 Red724-7133
 .060 Black724-7141
 .070 Green.....724-7158

Unity Drills Assorted Package

Side-cutting. Contains: 1 drill of each size.
 5/Pkg.....724-7091 **1@ \$133.60**



Discover the variety of cements we offer in the Cements Section
 Starting at page 80




PATTERSON
DENTAL/DENTAIRE

CREATING THE PERFECT SMILE

Your smile...

Faster than ever, CEREC creates value, raising your practice to the next level. The rapid return on investment will bring a smile to your face.

Your patients' smile...

With great quality, precision colour matching and same day restorations, your patients will rave about CEREC dentistry!