



**EXCLUSIVE OFFER**

**BUY \$3,999**  
in C/W products,  
**GET A FREE 16G iPad!\***

## ParaPost®

Fiber Lux & Taper Lux Translucent  
Fiber Post Systems

- Outstanding strength
- Low-stress head design
- Light-transmitting for fast cementation
- Parallel or tapered


ParaPost Taper Lux  
PF-180 Intro Kit (722-3035).....Ea. \$235.40  
Refill Package 5/Pkg.....Ea. \$55.10

**BUY 3, GET 1 FREE**

BUY 3 ParaPost Fiber Lux or  
ParaPost Taper Lux Refills,  
GET 1 of the same FREE!

Taper Lux

Fiber Lux



ParaPost Fiber Lux  
PF-170 Intro Kit (722-3308).....Ea. \$235.40  
Refill Package 5/Pkg.....Ea. \$55.10

## JET BITE

- Smooth, mousse like consistency
- No resistance to closure
- High dimensional stability/no distortion
- Not brittle - easy to grind and trim

**C6400 (718-6299) \$45.85**  
**50ml Single Pack**  
Includes: 2 cartridges,  
12 mixing tips, 12 spreader tips



## BioSonic® Solutions

Convenient, Concentrated MeterDose® Dispensing

- Super concentrated
- Super effective
- Super convenient



UC30 General Purpose Concentrated Solution  
(718-5002)  
16 oz. \$37.65 6 x 16 oz. (Case) \$225.90

UC32 Enzymatic Ultrasonic Cleaning Concentrate  
(718-5028)  
8 oz. \$36.20 6 x 8 oz. (Case) \$217.20

**BUY 6, GET 2 OF THE SAME FREE**



- A ideal offering of Silicone or Diamond Impregnated Polishers for dental or laboratory practice
- Polisher designed to handle the unique characteristics of porcelain, veneers and inlays and onlays, chairside and in the lab
- Long lasting shape integrity

**NOW AVAILABLE!**

**BUY 3, GET 1 FREE**  
of equal or lesser value / 5 per pack



**COMPOSITE**

	MINIPOINT	POINT	CUP	SMALL LENS	MINIPOINT	POINT	CUP	SMALL LENS
<b>Cat No.</b>								
<b>1 @ \$11.65</b>	RP230055RA	RP230059RA	RP230063RA	RP230069RA	RP230072RA	RP230076RA	RP230080RA	RP230086RA
<b>3 @ \$11.20</b>	748-0015	748-0023	748-0031	748-0049	748-0056	748-0064	748-0072	748-0080
<b>Spk</b>	<b>FINE</b>				<b>EXTRA FINE</b>			
	STEP 1: SHAPING AND SMOOTHING				STEP 2: HIGH LUSTER POLISHING			

**COMPOSITE PLUS**

	MINIPOINT	POINT	CUP	SMALL LENS	MINIPOINT	POINT	CUP	SMALL LENS	MINIPOINT	POINT	CUP	SMALL LENS
<b>Cat No.</b>												
<b>1 @ \$55.90</b>	RP230264RA	RP230265RA	RP230267RA	RP230271RA	RP230291RA	RP230293RA	RP230295RA	RP230297RA	RP230299RA	RP230301RA	RP230303RA	RP230305RA
<b>3 @ \$53.75</b>	748-0247	748-0254	748-0262	748-0270	748-0163	748-0171	748-0189	748-0197	748-0205	748-0213	748-0221	748-0239
<b>Spk</b>	<b>MEDIUM</b>				<b>FINE</b>				<b>ULTRA FINE</b>			
	ONE STEP: MIRROR FINISHING				STEP 1: SHAPING AND SMOOTHING				STEP 2: HIGH LUSTER POLISHING			

**CERAMIC**

	POINT	CUP	LENS	WHEEL	POINT	CUP	LENS	WHEEL
<b>Cat No.</b>								
<b>1 @ \$11.65</b>	RP230114RA	RP230122RA	RP230134UM	RP230136UM	RP230139RA	RP230147RA	RP230159UM	RP230161UM
<b>3 @ \$11.20</b>	748-0361	748-0379	748-0387	748-0395	748-0403	748-0411	748-0429	748-0437
<b>Spk</b>	<b>MEDIUM</b>				<b>FINE</b>			
	STEP 1: REMOVING AND CONTOURING				STEP 2: HIGH LUSTER POLISHING			

**CERAMIC PLUS**

	MINIPOINT	FLAME	CUP	SMALL LENS	WHEEL	DISC	LARGE LENS	MINIPOINT	FLAME	CUP	SMALL LENS	
<b>Cat No.</b>												
<b>1 @ \$55.90</b>	RP230313RA	RP230317RA	RP230319RA	RP230321RA	RP230325RA	RP230327HP	RP230329HP	RP230331RA	RP230335RA	RP230337RA	RP230339RA	
<b>3 @ \$53.75</b>	748-0510	748-0536	748-0544	748-0551	748-0569	748-0916	748-0817	748-0585	748-0601	748-0619	748-0627	
<b>1 @ \$55.90</b>		RP230315HP			RP230323HP	<b>1 @ \$188.70</b>	<b>1 @ \$163.10</b>		RP230333HP			
<b>3 @ \$53.75</b>		748-0528			748-0577	<b>3 @ \$179.20</b>	<b>3 @ \$156.80</b>		748-0593			
<b>Spk</b>	<b>MEDIUM</b>							<b>FINE</b>				
	STEP 1: REMOVING AND SHAPING							STEP 2: SMOOTHING				

**AVAILABLE GRITS**

**Composite**

- F = Fine : 45 µm
- EF = Extra Fine : 25 µm

**Composite PLUS**

- F = Fine : 45 µm
- UF = Ultra Fine : 15 µm
- M = Medium : 110 µm

**Ceramic**

- M = Medium : 110 µm
- F = Fine : 45 µm

**Ceramic PLUS**

- M = Medium : 110 µm
- F = Fine : 45 µm
- UF = Ultra Fine : 15 µm

**CAD/CAM Restoration PLUS**

- SC = Super Coarse : 150 µm
- M = Medium : 110 µm
- F = Fine : 45 µm


**Brush PLUS**

- F = Fine : 45 µm

**Universal**

- SC = Super Coarse : 150 µm
- F = Fine : 45 µm

**For more information or to place your order, please contact your Patterson representative, local branch or call 1-800-873-7683**  
To qualify for FREE goods on individual promotion, full quantity must be purchased at one time. For redemption of free goods send Patterson invoice to: Attn: Redemption Center, Coltene/Whaledent Inc., 235 Ascot Parkway, Cuyahoga Falls, OH 44223 or fax to 330-916-7096 or email: freegoods@coltenewhaledent.com. All invoices must be received by July 31, 2011. For any other information, please call 1-800-221-3046 or visit us at www.coltenewhaledent.com. \*All items on one order; cannot be combined with any free goods offers. One iPad per order.



	WHEEL	DISC	LARGE LENS	MINIPOINT	FLAME	CUP	SMALL LENS	WHEEL	DISC	LARGE LENS
Cat No.										
<b>1 @ \$55.90</b> <b>3 @ \$53.75</b>	RP230343RA 748-0643	RP230345HP 748-0924	RP230347HP 748-0825	RP230349RA 748-0650	RP230353RA 748-0676	RP230355RA 748-0684	RP230357RA 748-0692	RP230361RA 748-0718	RP230363HP 748-0932	RP230365HP 748-0833
<b>1 @ \$55.90</b> <b>3 @ \$53.75</b>	RP230341HP 748-0635	<b>1 @ \$188.70</b> <b>3 @ \$179.20</b>	<b>1 @ \$163.10</b> <b>3 @ \$156.80</b>		RP230351HP 748-0668			RP230359HP 748-0700	<b>1 @ \$188.70</b> <b>3 @ \$179.20</b>	<b>1 @ \$163.10</b> <b>3 @ \$156.80</b>
<b>5pk</b>	FINE STEP 2: SMOOTHING			ULTRA FINE STEP 3: HIGH LUSTER POLISHING						

**CAD/CAM RESTORATION PLUS**

	LENS	SMALL WHEEL	WHEEL	TAPER	CYLINDER	FLAME	CUP	WHEEL	LARGE LENS	DISC
Cat No.										
<b>1 @ \$69.70</b> <b>3 @ \$67.20</b>	RP230367HP 748-1260	RP230369HP 748-1278	RP230371HP 748-1286	RP230373HP 748-1294	RP230375HP 748-1302	RP230379RA 748-1328	RP230381RA 7481-336	RP230383HP 748-1344	RP230385HP 748-1468	RP230387HP 748-1567
<b>1 @ \$69.70</b> <b>3 @ \$67.20</b>						RP230377HP 748-1310			<b>1 @ \$163.10</b> <b>3 @ \$156.80</b>	<b>1 @ \$188.70</b> <b>3 @ \$179.20</b>
<b>5pk</b>	SUPER COARSE STEP 1: REMOVING AND SHAPING					MEDIUM STEP 2: SMOOTHING				

**MANDREL Handpiece**

	FLAME	CUP	WHEEL	LARGE LENS	DISC	WIDE CUP	CUP	FLAME	MANDREL Handpiece		
Cat No.											
<b>1 @ \$69.70</b> <b>3 @ \$67.20</b>	RP230391RA 748-1369	RP230393RA 748-1377	RP230395HP 748-1385	RP230397HP 748-1476	RP230399HP 748-1575	RP232002RA 748-1666	RP232004RA 748-1674	RP232006RA 748-1682	RP255047HP 000-0000	RP255049RA 000-0000	
<b>1 @ \$69.70</b> <b>3 @ \$67.20</b>	RP230389HP 748-1351			<b>1 @ \$163.10</b> <b>3 @ \$156.80</b>	<b>1 @ \$188.70</b> <b>3 @ \$179.20</b>			<b>1 @ \$23.30</b> <b>3 @ \$22.40</b>		<b>1 @ \$13.95</b> <b>3 @ \$13.45</b>	
<b>5pk</b>	FINE STEP 3: HIGH LUSTER POLISHING						POLISHES CERAMIC, COMPOSITES AND METALS			STAINLESS STEEL MANDRELS FOR UNMOUNTED POLISHERS	

**UNIVERSAL (COMPOSITE AND CERAMIC)**

	WIDE CUP	SMALL LENS	LENS	SMALL WHEEL	WHEEL	MINIPOINT	FLAME	BULLET	CUP	WIDE CUP	
Cat No.											
<b>1 @ \$7.00</b> <b>3 @ \$6.75</b>	RP230180UM 748-1013	RP230186UM 748-1021	RP230188UM 748-1039	RP230194UM 748-1047	RP230196UM 748-1054	RP230201RA 748-1112	RP230203RA 748-1138	RP230207RA 748-1153	RP230211RA 748-1179	RP230215RA 748-1187	
						<b>1 @ \$8.15</b> <b>3 @ \$7.85</b>		RP230205HP 748-1120	RP230209HP 748-1146	RP230213HP 748-1161	
<b>5pk</b>	SUPER COARSE ONE STEP: REMOVING, SHAPING AND POLISHING					<b>5pk</b>	FINE ONE STEP: REMOVING, SHAPING AND POLISHING				

**For more information or to place your order, please contact your Patterson representative, local branch or call 1-800-873-7683**  
 To qualify for FREE goods on individual promotion, full quantity must be purchased at one time. For redemption of free goods send Patterson invoice to: Attn: Redemption Center, Coltene/Whaledent Inc., 235 Ascot Parkway, Cuyahoga Falls, OH 44223 or fax to 330-916-7096 or email: freegoods@coltenewhaledent.com. All invoices must be received by July 31, 2011. For any other information, please call 1-800-221-3046 or visit us at www.coltenewhaledent.com. \*All items on one order; cannot be combined with any free goods offers. One iPad per order.



**3+1**



- Eliminates risk of joint breakage
- Cuts up to 4x faster
- Lasts longer

Operatives

**BUY 3 PACKS, GET 1 FREE**  
(10 per pack)

Operatives Bulk Packs

**BUY 3 PACKS, GET 1 FREE**  
(100 per pack)

Trimming & Finishing

**BUY 3 PACKS, GET 1 FREE**  
(5 per pack)



**3+1**



- Designed to rapidly cut Amalgam, Porcelain, Metal, Enamel and Dentin
- Reduces prep and chair time

Speedster

**BUY 3 PACKS, GET 1 FREE**  
(10 per pack)

Speedster Bulk Packs

**BUY 3 PACKS, GET 1 FREE**  
(100 per pack)



**3+1**



**Diamond-Shaped One-Piece Carbides for Rapid Crown and Bridge Preparation**

- Cuts 45% faster than diamonds, without clogging

Speedster C&B

**BUY 3 PACKS, GET 1 FREE**  
(5 per pack)



**5+2**

Multi-layer diamond application assures maximum, uniform crystal coverage for rapid, smooth and comfortable cutting.

Multi Use Diamonds

**BUY 5 PACKS, GET 2 FREE**  
(5 per pack)



**3+1**



**A Smart Combination of... Convenience, Efficiency and Infection Control**

- A sharp, pre-sterilized bur for every procedure, right out of the box
- Promotes Advanced Infection Control and patient safety
- Alpen x1 lets you keep your edge without emptying your wallet

**BUY 3, GET 1 FREE**  
(25 per pack)



**3+1**



**Ideal for Zirconia & Alumina Frameworks PLUS...**

- Porcelain
- IPS e.max - Lithium Disilicate
- many other hard substructures
- Perfect for lab & chairside adjustment of CAD/CAM frameworks and implant abutments
- Unique, synthetic diamond coating allows for chatter-free endo access
- Fine grit enables exact interproximal contact and occlusal adjustment

**BUY 3, GET 1 FREE**  
(5 per pack)